
## DIRECTORY

|                 |    |                |    |
|-----------------|----|----------------|----|
| English.....    | 1  | Norsk.....     | 29 |
| Français.....   | 3  | Suomi.....     | 30 |
| Deutsch.....    | 5  | Български..... | 32 |
| Español.....    | 7  | Română.....    | 34 |
| Italiano.....   | 9  | Eesti.....     | 36 |
| Nederlands..... | 11 | Lietuvių.....  | 38 |
| Dansk.....      | 13 | Latviešu.....  | 40 |
| Svenska.....    | 15 | Türkçe.....    | 42 |
| Ελληνικά.....   | 17 | Русский.....   | 44 |
| Português.....  | 19 | Srpski.....    | 46 |
| Česky.....      | 21 | 简体中文.....      | 48 |
| Magyar.....     | 23 | 한국어.....       | 50 |
| Polski.....     | 25 | 日本語.....       | 52 |
| Slovensky.....  | 27 |                |    |

## English

# HemoSphere Pressure-Out Cable

HMS-1020

Edwards, Edwards Lifesciences, the stylized E logo, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

**Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.** Refer to the HemoSphere advanced monitor operator's manual, available at [eifu.edwards.com](http://eifu.edwards.com), for comprehensive monitoring procedures.

## 1.0 Intended Use

The HemoSphere pressure-out cable is a reusable accessory that enables the user to output pressure monitored by the HemoSphere advanced monitoring platform to a compatible patient monitor for standard pressure monitoring. Up to three waveforms can be outputted: arterial pressure (AP, red), central venous pressure (CVP, blue), and pulmonary artery pressure (PAP, yellow).

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

## 2.0 Warnings


This product contains metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

## 3.0 Precautions

Inspect the pressure-out cable before each use. Do not use if the cable is damaged. Damage may include broken insulation, frays, deteriorated strain relief, bent/recessed connector pins or any signs that the cable may be compromised.

## 4.0 HemoSphere Pressure-Out Cable Setup

Refer to the numbered figure and to the HemoSphere advanced monitor operator's manual, available at [eifu.edwards.com](http://eifu.edwards.com), for comprehensive monitoring procedures with that monitor.

| Step | Procedure                                                                                                                                                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Plug the 18-pin connector ① of the pressure-out cable into the rear panel of the HemoSphere advanced monitor. The female port is denoted by an analog output symbol.<br> |
| 2    | Use the two jack screws ② to properly secure the pressure-out cable plug in place.                                                                                                                                                                          |

| Step | Procedure                                                                                                                                                                                                                                                                                                              |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 3    | Connect the desired pressure signal plug(s) into the compatible patient monitor:<br><br>③ arterial pressure (AP, red)<br>④ pulmonary artery pressure (PAP, yellow)<br>⑤ central venous pressure (CVP, blue)<br><br>Ensure that the selected connector is fully engaged. Refer to patient monitor instructions for use. |
| 4    | Refer to the HemoSphere advanced monitor operator's manual for instructions on configuring the desired pressure signal output and for proper zeroing instructions.                                                                                                                                                     |

## 5.0 MRI Information

**MR Unsafe**

Do not use the HemoSphere pressure-out cable in an MR environment. The cable is MR unsafe since the accessory contains metallic components, which can experience RF-induced heating in the MRI environment.

## 6.0 Operating Conditions

Temperature: 10 to 32.5 °C

Humidity range: 20 to 90% non-condensing

Altitude (Atmospheric Pressure): 0 m/0 ft (1013 hPa) to 3048 m/10,000 ft (697 hPa)

## 7.0 Maintenance

It is recommended to inspect the HemoSphere pressure-out cable periodically for signs of damage. Check the condition of the cable and replace if any of the following are present: exposed wire, splits, cracks, or signs of stress.

## 8.0 Cleaning

The HemoSphere pressure-out cable is a reusable accessory and can be surface cleaned per hospital aseptic procedures. Accessory has been shown to handle cloth wipe down using the following chemical agents: 70% isopropyl alcohol solution, 2% glutaraldehyde, 10% bleach solution, or quaternary ammonium solution.

---

## CAUTIONS:

**Do not use any other cleaning agents, spray, or pour cleaning solution directly on the cable.**

**Do not steam, radiate, or EtO sterilize the cable.**

**Do not immerse.**

## 9.0 Storage

Store in a cool, dry place in original packaging to prevent damage.

**CAUTION: Do not coil the cable against its natural orientation or coil the cable tightly.**

Temperature: -18 to 45 °C

Relative humidity: 90% non-condensing at 45 °C

Altitude: 0 to 20,000 feet (6096 m)

## 10.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): .....800.822.9837

Outside the U.S. and Canada (24 hours): ...949.250.2222

In the UK: .....0870 606 2040 - Option 4

In Ireland: ..... 01 8211012 - Option 4

## 11.0 Notices

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Prices, specifications and model availability are subject to change without notice.

## 12.0 Disposal

To avoid contaminating or infecting personnel, the environment or other equipment, make sure the HemoSphere pressure-out cable is disinfected and decontaminated appropriately in accordance with your country's laws for equipment containing electrical and electronic parts prior to disposal.

Dispose of in accordance with national and local legislation and requirements for equipment containing electrical and electronic parts.

**Refer to the symbol legend at the end of this document.**

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|                                                                                                                           | Latviešu                                                                                                                                | Türkçe                                                                                         | Русский                                                                                                                                           | Srpski                                                                                      |
|---------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|                                          | Kataloga numurs                                                                                                                         | Katalog numarası                                                                               | Номер по каталогу                                                                                                                                 | Kataloški broj                                                                              |
| <br>eifu.edwards.com<br>+ 1 888 570 4016 | Levērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē                                                                          | Web sitesinde bulunan kullanım talimatlarını takip edin                                        | Следуйте инструкциям по применению на веб-сайте                                                                                                   | Pridržavajte se uputstva za upotrebu sa veb-sajta                                           |
|                                          | Sērijas numurs                                                                                                                          | Seri Numarası                                                                                  | Серийный номер                                                                                                                                    | Serijski broj                                                                               |
|                                          | Ražotājs                                                                                                                                | Üretici                                                                                        | Производитель                                                                                                                                     | Proizvođač                                                                                  |
|                                          | Izgatavošanas datums                                                                                                                    | Üretim tarihi                                                                                  | Дата производства                                                                                                                                 | Datum proizvodnje                                                                           |
|                                          | Pilnvarotais pārstāvis Eiropas Kopienā                                                                                                  | Avrupa Topluluğu'nda Yetkili Temsilci                                                          | Официальный представитель в Европейском сообществе                                                                                                | Ovlašćeni predstavnik u Evropskoj zajednici                                                 |
|                                          | Nedrīkst lietot MR vidē                                                                                                                 | MR için Güvenli Değildir                                                                       | Опасно при проведении МРТ                                                                                                                         | Nije bezbedno za MR                                                                         |
|                                          | Mitrumizturīgs aprikojums                                                                                                               | Su sıçramasına karşı korumalı ekipman                                                          | Водостойкое оборудование                                                                                                                          | Oprema sa zaštitom od kapanja                                                               |
|                                          | Aizsardzības līmenis pret priekšmetu iekļūšanu                                                                                          | Nesne girişine karşı koruma derecesi                                                           | Степень пылезащиты                                                                                                                                | Stepen zaštite od prodora objekata                                                          |
|                                          | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.                                             | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.    | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.                                                       | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|                                          | Pret defibrilāciju noturīga CF tipa daļa, kas saskaras ar pacientu                                                                      | Defibrilasyon korumalı CF tipi hastayla temas eden parça                                       | Контактирующий с пациентом элемент типа CF с защитой от разряда дефибриллятора                                                                    | Primenjeni deo tipa CF otporan na defibrilaciju                                             |
|                                        | Bīstamo vielu ierobežošanas elektriskās un elektroniskās iekārtās (RoHS — Restriction of Hazardous Substances) atbilstība — tikai Ķīnai | Tehnikeli Maddelerin Sınırlandırılması (RoHS) uygunluğu - Yalnızca Çin için                    | Соответствие ограничениям по опасным веществам (RoHS) — только для Китая                                                                          | Usaglašenost sa Ograničenjem upotrebe opasnih supstanci (RoHS) — samo Kina                  |
|                                        | Federālā komunikāciju komisija                                                                                                          | Federal Haberleşme Komisyonu                                                                   | Федеральная комиссия по связи                                                                                                                     | Savezna komisija za komunikacije                                                            |
|                                        | Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES direktīvu 2012/19/ES                                                  | 2012/19/EU sayılı EC direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır | Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU                                              | Odvojeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU   |
|                                        | Temperatūras ierobežojums                                                                                                               | Sıcaklık sınırı                                                                                | Ограничение по температуре                                                                                                                        | Ograničenje temperature                                                                     |
|                                        | Mitruma ierobežojums                                                                                                                    | Nem sınırlaması                                                                                | Ограничение по влажности                                                                                                                          | Ograničenje vlažnosti                                                                       |
|                                        | Sargāt no saules gaismas                                                                                                                | Güneş ışığından uzak tutun                                                                     | Беречь от солнечных лучей                                                                                                                         | Držati zaklonjeno od sunčeve svetlosti                                                      |
|                                        | Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju                                                                     | Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun                              | Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией                                                              | Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu                    |
|                                        | Conformité Européenne (CE markējums), ko piešķir TÜV SÜD Product Service GmbH (paziņotā struktūra)                                      | TÜV SÜD Product Service GmbH (onaylı kuruluş) Conformité Européenne (CE İşareti)               | Соответствие требованиям директив Европейского Союза (маркировка CE) подтверждено компанией TÜV SÜD Product Service GmbH (нотифицированный орган) | Conformité Européenne (CE znak) kompanije TÜV SÜD Product Service GmbH (ovlašćeno telo)     |
|                                        | Daudzums                                                                                                                                | Miktar                                                                                         | Количество                                                                                                                                        | Količina                                                                                    |

### Simbolu skaidrojums • Sembol Açıklamaları • Условные обозначения • Legenda sa simbolima

|                                                                                                                                                                                                                                                                                                                | Latviešu                                          | Türkçe                                                                               | Русский                                                | Srpski                                                                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------|--------------------------------------------------------------------------------------|--------------------------------------------------------|-------------------------------------------------------------------------|
|                                                                                                                                                                                                                               | Nejonizējošs elektromagnētiskais starojums        | İyonlaştırıcı olmayan elektromanyetik radyasyon                                      | Неионизирующее электромагнитное излучение              | Nejonizujuće elektromagnetno zračenje                                   |
|                                                                                                                                                                                                                               | Conformité Européenne (CE marķējums)              | Conformité Européenne (CE İşareti)                                                   | Conformité Européenne (маркировка CE)                  | Conformité Européenne (CE znak)                                         |
|                                                                                                                                                                                                                               | Nesatur svinu                                     | Kurşun içermez                                                                       | Не содержит свинец                                     | Bez olova                                                               |
|                                                                                                                                                                                                                               | UL apstiprinātas detaļas marķējums Kanādai un ASV | Kanada ve Amerika Birleşik Devletleri için UL Tarafından Kabul Gören Bileşen İşareti | Знак сертифицированного UL компонента для Канады и США | UL oznaka za odobrene komponente za Kanadu i Sjedinjene Američke Države |
|                                                                                                                                                                                                                               | Pārstrādājami litija joni                         | Geri Dönüştürülebilir Lityum İyon                                                    | Перерабатываемые литий-ионные компоненты               | Litijum-jonska baterija koju je moguće reciklirati                      |
|                                                                                                                                                                                                                               | Neizjaukt!                                        | Sökmeyin                                                                             | Разбирать запрещено                                    | Nemojte rastavljati                                                     |
|                                                                                                                                                                                                                               | Nedezzināt!                                       | Yakmayın                                                                             | При утилизации не сжигать                              | Nemojte paliti                                                          |
|                                                                                                                                                                                                                               | Medicīnas ierīce                                  | Tıbbi cihaz                                                                          | Медицинское устройство                                 | Medicinsko sredstvo                                                     |
|                                                                                                                                                                                                                              | Importētājs Eiropas Savienībā                     | Avrupa Birliği İthalatçısı                                                           | Импортер Европейского союза                            | Uvoznik za Evropsku uniju                                               |
| <b>Piezīme:</b> Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. • <b>Not:</b> Bu ürünün etiketinde tüm semboller yer almamış olabilir. • <b>Примечание:</b> На этикетках данного изделия могут указываться не все символы. • <b>Napomena:</b> Na ambalaži proizvoda se možda ne nalaze svi simboli. |                                                   |                                                                                      |                                                        |                                                                         |

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05/21  
10047195001 A / DOC-0174746 A  
© Copyright 2021, Edwards Lifesciences LLC  
All rights reserved.

**Edwards Lifesciences LLC**  
One Edwards Way  
Irvine, CA 92614 USA  
Made in China

**Edwards Lifesciences B.V.**  
Verlengde Poolseweg 16  
4818 CL Breda, Netherlands

|           |              |         |
|-----------|--------------|---------|
| Telephone | 949.250.2500 | Web IFU |
|           | 800.424.3278 |         |
| FAX       | 949.250.2525 |         |

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**Edwards Lifesciences Services GmbH**