DIRECTORY

English.....	1	Português.....	10	Eesti.....	19
Français.....	2	Česky.....	11	Lietuvių.....	20
Deutsch.....	3	Magyar.....	12	Latviešu.....	21
Español.....	4	Polski.....	13	Türkçe.....	22
Italiano.....	5	Slovensky.....	14	Русский.....	23
Nederlands.....	6	Norsk.....	15	Srpski.....	24
Dansk.....	7	Suomi.....	16	简体中文.....	25
Svenska.....	8	Български.....	17	한국어.....	26
Ελληνικά.....	9	Română.....	18	日本語.....	27

English

Roll Stand Assembly and Operation

REF: HEMRLSTD1000

For figures, refer to Figure 1 on page 29 through Figure 8 on page 30.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 INTENDED USE

This roll stand is intended to be used with the HemoSphere advanced monitor. See HemoSphere advanced monitor operator's manual for more details.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

Installation Note:

Separate installation kits and accessories may be included with the roll stand. Read all separate installation guides before assembling the roll stand or mounting instruments or accessories.

Tools Required:

- Phillips screwdriver (not provided)
- 5/32" [4 mm] hex wrench (provided)
- 2.5 mm hex wrench (provided)
- 1/2" [13 mm] wrench (not provided)

2.0 Attaching Roll Stand Post to the Base

Step	Procedure
1	Insert the post in the base and lay assembly on its side for access to bottom of the base (see Figure 1 on page 29).

Step	Procedure
2	Using a 1/2" [13 mm] wrench, attach the post to the base with one (1) 5/16-18 x 1" HHCS, 5/16 split lock washer, and 5/16" flat washer.

3.0 Attaching the Basket to the Post

Step	Procedure
1	Install both nylon clips at an angle, in the center of the basket. Rotate clips until both are horizontal (see Figure 2 on page 29).
2	Slide the clips to the top and bottom of basket. The top clip should be touching the top horizontal wire, and the bottom clip should be touching the lower reinforcing wire (see Figure 3 on page 29).
3	Slide the clips over the top of the post (see Figure 4 on page 29). Position the basket at the desired height. Tighten the clips using the supplied #8-32 machine screws and hex nuts (see Figure 5 on page 30).

4.0 Attaching Mounting Plate Assembly

Step	Procedure
1	Insert mounting plate assembly in top of post (see Figure 6 on page 30). Press downward while twisting until the assembly snaps into place.
2	Rotate the mounting plate until the spring-loaded plunger (indicated by ① in Figure 6 on page 30) is toward the front of the roll stand.

5.0 Attaching Roll Stand Handle

Fasten handle to roll stand post with two (2) #10-32 x 9/16" PHMS (see Figure 7 on page 30). Alternately tighten handle screws until handle does not rotate around post.

Note: Handle position shown is typical but handle may be attached at any desired location on the post.

6.0 Periodic Maintenance and Cleaning

All fasteners associated with the mounting system should be inspected periodically and tightened as necessary. Roll stand may be cleaned with any non-abrasive solution as approved by the facility in which the product will be used.

7.0 Warnings

- DO NOT use this product without first reading and understanding the instructions contained in this booklet. If you are unable to understand the Warnings, Caution or Instructions, contact a healthcare professional, dealer or technical personnel before use – otherwise serious bodily injury or product damage may occur.
- All wheels must be in contact with the floor at ALL TIMES.
- DO NOT use the roll stand while walking backwards down gradients or climbing stairs, curbs or to go over obstacles. Serious risk of fall or injury may occur.
- Engage all wheel locks when roll stand is not in motion.
- Do not place items in the utility basket exceeding 4.5 kg (10 lbs).
- Use care when moving roll stand with a patient. Avoid wet or slick surfaces, wires, cables, or tubes and any uneven surfaces that may cause roll stand to tip or fall.
- Only use accessories and spare parts authorized by GCX Corporation.
- DO NOT attach roll stand to other objects.
- DO NOT use in MRI areas.
- DO NOT use the roll stand for support or as a walking aid. Do not attempt to support your weight or use as a stand assist.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

8.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837
Outside the U.S. and Canada (24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications and model availability are subject to change without notice.

9.0 Disposal

Dispose of in accordance with hospital policy and local regulations.

Refer to the symbol legend at the end of this document.

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© Copyright 2021, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC
Edwards Lifesciences B.V.