
DIRECTORY

English.....	1	Česky	31	Latviešu	61
Français	4	Magyar	34	Türkçe	64
Deutsch	7	Polski	37	Русский	67
Español	10	Slovensky	40	Srpski	70
Italiano	13	Norsk	43	简体中文	73
Nederlands	16	Suomi	46	繁體中文 (台灣)	76
Dansk	19	Български	49	한국어	79
Svenska	22	Română	52	Hrvatski	82
Ελληνικά	25	Eesti	55		
Português	28	Lietuvių	58		

English

FloTrac Pressure Cable—BeneVision

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The FloTrac Pressure Cable—BeneVision is intended to be used by qualified personnel or trained clinicians in a critical care environment in a hospital setting. The FloTrac Pressure Cable—BeneVision is intended for use with compatible monitors, and Edwards pressure monitoring sensors and compatible catheters.

The FloTrac Pressure Cable—BeneVision is a reusable device that connects with a compatible monitor on one end ④ and an Edwards pressure monitoring sensor on the other end ① to continuously measure hemodynamic parameters. Refer to the monitor operator's manual for a list of supported hemodynamic parameters and applicable Edwards pressure monitoring sensor(s).

Note: The expected useful life of the FloTrac Pressure Cable—BeneVision is 5 years from date of purchase. If your equipment experiences a malfunction, please contact Technical Support or your local Edwards representative for further assistance.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established

Instructions for Use. Refer to the monitor operator's manual for comprehensive monitoring procedures and monitoring parameters.

2.0 Compatible Sensors

The FloTrac Pressure Cable—BeneVision is intended to be used in combination with a compatible monitor, an Edwards pressure monitoring sensor and compatible pressure monitoring catheter to measure intravascular pressure. Refer to the monitor operator's manual for a list of supported Edwards-compatible pressure monitoring sensors and catheters.

3.0 Intended Use/Purpose

The FloTrac Pressure Cable—BeneVision when used with an Edwards' pressure monitoring sensor and connected to a compatible monitor, offers continuous assessment of hemodynamic parameters. It may be used for monitoring hemodynamic parameters in conjunction with a perioperative goal directed therapy protocol in a hospital environment.

4.0 Indications for Use

The FloTrac Pressure Cable—BeneVision is indicated for use in adult critical care patients in which the balance between cardiac function, fluid status, vascular resistance and pressure needs continuous assessment.

5.0 Contraindications

There are no known contraindications.

6.0 Warnings

Compliance to IEC 60601-1 is only maintained when the FloTrac Pressure Cable—BeneVision (applied part accessory, defibrillation proof) is connected to a compatible monitoring platform. Connecting external equipment or configuring the system in a way not described in these instructions will not meet this standard. Failure to use the device as instructed may increase the risk of electrical shock to the patient/operator.

Do not modify, service or alter the product in any way. Servicing, alteration or modification may affect patient/operator safety and/or product performance.

7.0 Precautions

Do not use if the package or the FloTrac Pressure Cable—BeneVision is damaged. Damage may include cracks, scratches, dents or any signs that the housing may be compromised.

Risks of the FloTrac Pressure Cable—BeneVision include patient or clinician burns or electrical shock, inappropriate/unintended treatment, and/or delay in treatment.

8.0 FloTrac Pressure Cable—BeneVision Setup

Refer to the compatible monitor operator's manual for comprehensive monitoring procedures.

Step	Procedure
1	Connect the FloTrac Pressure Cable—BeneVision to a compatible monitor at the end indicated by ④ in Figure 1.
2	Setup and initiate the compatible monitor.
3	Insert the appropriate color insert at ② to indicate the monitored pressure type. The colors are as follows: • Red for arterial pressure (AP) • Blue for central venous pressure (CVP) • Yellow for pulmonary artery pressure (PAP) • Green for cardiac output (CO)
4	Connect the pressure monitoring sensor to the FloTrac Pressure Cable—BeneVision. ① The FloTrac Pressure Cable—BeneVision LED that surrounds the zero button at ③ will flash green indicating that the pressure monitoring sensor is detected. A yellow light indicates a sensor fault condition. If this occurs, refer to the monitor screen for specific fault condition details.
5	Follow all instructions found in pressure monitoring catheter IFU for catheter preparation and insertion procedures. The FloTrac Pressure Cable—BeneVision must be zeroed before monitoring.

9.0 Zero FloTrac Pressure Cable—BeneVision

Step	Procedure
1	To initiate the zeroing process, press and hold the zero button indicated by ③ in the figure or initiate zeroing process on the connected monitor screen. Refer to the compatible monitor's operator manual for on-screen zeroing instructions. CAUTION: To prevent cable damage, do not apply excessive force to the FloTrac Pressure Cable—BeneVision zero button.
2	Level the sensor stopcock valve (vent port) to the patient's phlebostatic axis position according to the instructions for use.
3	Open the stopcock valve to measure atmospheric conditions.

Step	Procedure
4	Press and hold the FloTrac Pressure Cable—BeneVision zero button or the on-screen button to complete zeroing process. The zero button LED will emit a solid green light when properly zeroed.
5	Confirm stable zero pressure and turn stopcock such that sensor is reading patient intravascular pressure.

10.0 MRI Information

MR Unsafe

Do not use the FloTrac Pressure Cable—BeneVision in an MR environment. The FloTrac Pressure Cable—BeneVision is MR unsafe since the device contains metallic components, which can experience RF-induced heating in the MRI environment.

11.0 Specifications

Physical Specifications

Cable Length: 3 m (10 ft) and 4.6 m (15 ft)

Operating Conditions

Temperature: 10 to 32.5 °C

Humidity range: 20 to 90% non-condensing

Altitude (Atmospheric Pressure): 0 m/0 ft (1013 hPa) to

3048 m/10,000 ft (697 hPa)

IP Rating: IPX4

12.0 Maintenance

When unplugging the FloTrac Pressure Cable—BeneVision from a compatible monitor or pressure monitoring sensor from the FloTrac Pressure Cable—BeneVision, always pull at the connection site. Do not pull from cables or use tools to disconnect.

CAUTION: Excessive dropping of the FloTrac Pressure Cable—BeneVision may result in cable damage and/or malfunction.

It is recommended to inspect the FloTrac Pressure Cable—BeneVision periodically for signs of damage. Check the condition of the cable and replace if any of the following are present: exposed wire, splits, cracks, or signs of stress.

The FloTrac Pressure Cable—BeneVision is a reusable device and can be surface cleaned per hospital aseptic procedures. Device has been shown to handle cloth wipe down using the following chemical agents: 70% isopropyl, 2% glutaraldehyde, 10% bleach solution, or quaternary ammonium solution. Disconnect the FloTrac Pressure Cable—BeneVision from the monitor to air dry the transducer connector. To blow dry the transducer connector, use clean, dry wall air, canned air, or CO₂ aerosol for at least two minutes. If left to dry under room conditions, allow the connector to dry for two days before using.

CAUTIONS: Do not use any other cleaning agents, spray, or pour cleaning solution directly on the FloTrac Pressure Cable—BeneVision.

Do not steam, radiate, or EtO sterilize the FloTrac Pressure Cable—BeneVision.

Device contains electronics. Handle with care.

Do not immerse.

Dispose of in accordance with hospital policy and local regulations.

13.0 Storage Conditions

Store in a cool, dry place.

Store in original packaging to prevent damage. Keep away from sunlight.

CAUTION: Do not coil the cable against its natural orientation or coil the cable tightly.

14.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837

Outside the U.S. and Canada (24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications, and model availability are subject to change without notice.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.


15.0

Table 1: FloTrac Pressure Cable--BeneVision Parameters:

Parameters ¹	Description	Sub-System Technology Used	Patient Population	Hospital Environment
CO	Continuous Cardiac Output	FloTrac Pressure Cable—BeneVision	Adult only	Operating room, intensive care unit, emergency room
CI	Continuous Cardiac Index			
CVP	Central Venous Pressure			
DIA_art	Diastolic Blood Pressure			
MAP	Mean Arterial Pressure			
MPAP	Mean Pulmonary Arterial Pressure			
PPV	Pulse Pressure Variation			
SV	Stroke Volume			
SVI	Stroke Volume Index			
SVR	Systemic Vascular Resistance			
SVRI	Systemic Vascular Resistance Index			
SVV	Stroke Volume Variation			
SYS_art	Systolic Blood Pressure			
PR	Pulse Rate			

¹Available parameters will vary depending on the pressure monitoring sensor that is connected, a list of supported parameters will be provided in the monitor's operations manual.

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© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences GmbH
Edwards Lifesciences LLC
Edwards Lifesciences B.V.