Edwards

Fogarty Venous Thrombectomy Catheter

Directory	
English (en).....	1
Français (fr).....	4
Deutsch (de).....	6
Español (es).....	8
Italiano (it).....	10
Nederlands (nl).....	12
Dansk (da).....	14
Svenska (sv).....	16
Ελληνικά (el).....	18
Português (pt).....	20
Česky (cs).....	22
Magyar (hu).....	24
Polski (pl).....	26
Slovensky (sk).....	28
Norsk (no).....	30
Suomi (fi).....	32
Български (bg).....	34
Română (ro).....	36
Eesti (et).....	38
Lietuvių (lt).....	40
Latviešu (lv).....	42
Türkçe (tr).....	44
Русский (ru).....	46
Srpski (sr).....	48
Hrvatski (hr).....	50
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	52

English

Fogarty Venous Thrombectomy Catheter

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

Fogarty venous thrombectomy catheters have a long, soft tip which permits advancement past the venous valves without undue trauma. The flexibility required for the venous procedure is supplied by the spring-wound body. The spring is covered with a braiding for additional strength. The gate valve is designed for single-handed operation and to minimize the possibility of leakage. The gate valve has an arrow to indicate the "open" and "closed" positions. There are three available catheter sizes (the 6 French (2.0 mm), 8 French (2.7 mm), and 8/10 French (2.7/3.3 mm)). Each may be quickly identified by the color-coded body. The size and filling capacity are printed on each catheter.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

This device aids in clot management and restoration of circulatory blood flow to veins of the non-central circulatory system and temporary occlusion of blood vessels.

The device is intended to be used by medical professionals who are trained in the safe use and clinical usage of endovascular catheter technologies as part of their respective institutional guidelines.

2.0 Intended Use/Purpose

Fogarty venous thrombectomy catheters are intended for use in the performance of venous thrombectomy and for the temporary occlusion of blood vessels.

3.0 Indications

Fogarty venous thrombectomy catheters are specifically designed for the performance of venous thrombectomy in adult patients. The Fogarty venous thrombectomy catheter takes the path of least resistance through the obstructing material. Once properly positioned, the balloon is inflated with a sterile medium. By the mechanism of fluid displacement, the balloon maintains uniform contact with the vessel wall as the catheter and occluding material are withdrawn. Additionally, it may be used for the temporary occlusion of blood vessels.

Fogarty venous thrombectomy catheters are designed for use in the veins of the non-central circulatory system only.

4.0 Contraindications

Venous thrombectomy catheters are not designed to withstand the additional strain encountered during endarterectomy procedures.

5.0 Warnings

- **Air should not be used for inflation in instances where balloon rupture could produce a dangerous air embolus.**
- **The maximum recommended volume should not be exceeded, as overinflation increases the possibility of balloon rupture. (See specification table).**
- **Undue manipulation or overdistension of the balloon can result in significant damage to the vessel lining, increasing the potential for recurrent thrombosis.**
- **Experience, along with the realization of the limitations of the instrument, have been the most significant factors in reducing the incidence of complications.**

6.0 Instructions

6.1 Balloon Inflation

The balloon should be inflated with a sterile, blood-compatible fluid which may include radiopaque solutions if they are highly dilute and nonparticulate.

Carbon dioxide should be substituted in cases where gas is used as the inflation fluid. However, it should be recognized that CO₂ will readily diffuse through the balloon wall, and corrections for gas loss must be made.

Use of a highly viscous or highly particulate contrast medium is not recommended for balloon inflation as the inflation passage may become occluded. Prior to use, the catheter, with the balloon inflated to its maximum recommended volume, should be carefully inspected by the surgeon.

The smallest syringe which will hold the stated maximum fluid capacity should be used, as the smaller diameter facilitates control of the withdrawal force because greater displacement of the plunger is required to make a given change in a balloon diameter.

The amount of fluid in the syringe should be checked before each inflation. If it is more than the maximum capacity, the syringe should be removed from the catheter and the barrel refilled to the proper volume, taking into account the amount of fluid remaining in the lumen of the catheter.

Inflation of the balloon is associated with a feeling of resistance. When no resistance is encountered, it should be assumed that the balloon has ruptured. Inflation should be discontinued at once, and the catheter withdrawn.

7.0 MRI Information

This product has not been tested for MRI compatibility.

8.0 Complications

As with all catheterization procedures, complications may occur. Incidents such as distal embolization, sepsis/infection, bleeding, blood loss, local hematoma, tissue/vascular trauma, soft tissue trauma/injury, air emboli, adverse reaction to device materials, ischemia and arteriovenous fistula have been reported with the use of thrombectomy catheters.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

9.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

STERILE	EO
---------	----

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

	Suomi	Български	Română	Eesti	Lietuvių
	Käyttöpituus	Използваема дължина	Lungimea utilă	Kasutatav pikkus	Naudingasis ilgis
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks	Vienkartinio naudojimo
QTY	Määrä	Количество	Cantitate	Kogus	Kiekis
LOT	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
EC REP	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană / Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Viimeinen käyttöpäivä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Naudoti iki nurodytos datos
STERILEO	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
BC	Pallon tilavuus	Вместимост на балона	Capacitatea balonului	Ballooni mahutavus	Balionėlio talpa
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
LATEX	Sisältää luonnokumilateksia tai sen jäämiä	Съдържа или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural	Sisaldab naturaalsed kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen
 ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli
 ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbole seletus
 ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба в уебсайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege kasutusjuhiseid veebisaidilt	Žr. naudojimo instrukcijas interneto svetainėje
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Koko	Размер	Mărime	Suurus	Dydis
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhiseiga	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinaï nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusi priemonės identifikatorius
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Mallinumero	Номер на модела	Număr model	Mudeli number	Modelio numeris

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebljiva dužina	Upotrebljiva dužina
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
	Daudzums	Miktar	Количество	Količina	Količina
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen
 ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli
 ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus
 ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sterilizēts ar etilēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidadom	Sterilizirano etilen-oksidadom
	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na internet stranici	Pogledajte upute za upotrebu na internetskoj stranici
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Izmērs	Boyut	Размер	Veličina	Veličina
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykiety niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei vältämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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05/23
10054348001 A / DOC-0208128 A
© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH