Edwards

DIRECTORY

English	1	Português	34	Eesti	66
Français	4	Česky	38	Lietuvių	70
Deutsch	8	Magyar	41	Latviešu	74
Español	12	Polski	45	Türkçe	77
Italiano	16	Slovensky	48	Русский	80
Nederlands	19	Norsk	52	Srpski	84
Dansk	23	Suomi	55	Hrvatski	88
Svenska	27	Български	59		
Ελληνικά	30	Română	63		

English

Swan-Ganz Pacing-TD Catheter D200F7, D205F7

D205F7 is not available in EU

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

For figures, see Figure 1 on page 92 through Figure 7 on page 97.

1.0 Description

The pacing-TD catheter is available in two models: the standard Model D200F7 and the Model D205F7. For additional clinical flexibility in the smaller anatomy, the electrodes have been moved distally on the Model D205F7. The pacing TD catheter has three atrial and two ventricular electrodes for atrial and ventricular pacing and atrioventricular (A-V) sequential pacing.

The Swan-Ganz pacing catheters are recommended for use *in situ* for up to 72 hours.

As part of the insertion procedure this product is used for ECG detection during placement, but is not intended for ECG monitoring.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The intravascular catheter is inserted through the central vein to connect to the right side of the heart and advance towards the pulmonary artery. Route of insertion can be internal jugular vein, antecubital or brachial vein. The cardiac structures in contact are the right atria, right ventricle, and pulmonary artery as well as the circulatory system.

The Swan-Ganz pacing-TD catheter provides temporary electrical cardiac stimulation to facilitate return to and/or continuance of normal cardiac rhythm.

Edwards, Edwards Lifesciences, the stylized E logo, CO-Set, CO-Set+, Swan, and Swan-Ganz are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

This device has not yet been tested in the pediatric population or in pregnant or lactating women.

2.0 Intended Use/Purpose

Swan-Ganz pacing-TD catheters are pulmonary artery catheters intended to serve as diagnostic and therapeutic tools for short-term use in the central circulatory system for patients who require intracardiac hemodynamic monitoring, temporary cardiac pacing, blood sampling and infusing solutions. Swan-Ganz pacing-TD catheters also provide electrocardiographic detection during catheter placement. When used with a compatible monitoring platform and accessories Swan-Ganz pacing-TD catheters offer a comprehensive hemodynamic profile to help clinicians assess cardiovascular function and guide treatment decisions.

3.0 Indications

The Swan-Ganz pacing-TD catheters are indicated for use for temporary, atrial, ventricular, or A-V sequential cardiac pacing of adult surgical and/or critically ill patients. This includes, but is not limited to, patients with arrhythmias, acute myocardial infarction, cardiac surgery, and coronary angiography.

In addition, these catheters are indicated for use in hemodynamic monitoring of patients experiencing medical conditions such as, but not limited to, post-major surgical recovery, trauma, sepsis, burns, pulmonary disease, pulmonary failure, stroke, and cardiac disease including heart failure.

4.0 Contraindications

Endocardial pacing catheters are contraindicated in patients with a mechanical valve. If the pacing lead is placed through a mechanical valve, the valve leaflet will be held permanently open and result in significant permanent valve regurgitation.

Relative contraindications may include patients with either recurrent sepsis, or with hypercoagulable state where the catheter could serve as a focus for septic or bland thrombus formation.

These products contain metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

5.0 Warnings

Air should never be used for balloon inflation in any situation where it may enter the arterial circulation, e.g. in all pediatric patients and adults with suspected right to left intracardiac intrapulmonary shunts.

Bacteria-filtered carbon dioxide is the recommended inflation medium because of its rapid absorption into the blood in the event of balloon rupture within the circulation. Carbon dioxide diffuses through the latex balloon, diminishing the balloon's flow-directed capability after 2 to 3 minutes of inflation.

Do not leave the catheter in a permanent wedge position. Furthermore, avoid lengthy balloon inflation while the catheter is

in a wedge position; this occlusive maneuver may result in pulmonary infarction.

This device is designed, intended and distributed for **SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE** this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Do not modify or alter the product in any way. Alteration or modification may affect patient/operator safety or product performance.

Cleaning and resterilization will damage the integrity of the latex balloon. Damage may not be obvious during routine inspection.

As part of the insertion procedure this product is used for ECG detection during placement, but is not intended for ECG monitoring.

Use of the catheter on patients with tricuspid tissue-based valve prosthesis can increase regurgitation because when the device lays across the tissue valve leaflets it may cause regurgitation. It should be noted that the lead should be placed into one of the commissures so not to touch the leaflet and cause regurgitation.

The pacing-TD catheter is not recommended for use in pacemaker-dependent patients. No absolute contraindications for the use of flow directed pulmonary artery catheters exist. However, a patient with a left bundle branch block may develop a right bundle branch block during catheter insertion, resulting in complete heart block. In such patients, temporary pacing modes should be immediately available.

Electrocardiographic monitoring during catheter passage is encouraged and is particularly important in the presence of either of the following conditions:

- Complete left bundle branch block, in which the risk of complete heart block is somewhat increased.
- Wolff-Parkinson-White syndrome and Ebstein's malformation, in which the risk of tachyarrhythmias is present.

Some models may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0. Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

6.0 Precautions

Clinicians using the device should be familiar with the device and understand its applications prior to use.

When handling indwelling leads, the terminal pins or exposed metal (on the product) are not to be touched nor be allowed to contact electrically conductive or wet surfaces.

Failure of a balloon flotation catheter to enter the right ventricle or pulmonary artery is rare, but may occur in patients with an enlarged right atrium or ventricle particularly if the cardiac output is low or in the presence of tricuspid or pulmonic incompetence or pulmonary hypertension. Deep inspiration by the patient during advancement may also facilitate passage.

7.0 Insertion

Swan-Ganz catheters can be inserted at the patient's bedside, without the aid of fluoroscopy, guided by continuous pressure monitoring.

Precaution: Pacing TD catheters are not suitable for placement through the inferior vena cava.

8.0 Equipment

WARNING: Compliance to IEC 60601-1 is only maintained when the catheter or probe (Type CF applied part, defibrillation proof) is connected to a patient monitor or equipment that has a Type CF defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the catheter or probe. Failure to ensure monitor or equipment compliance to IEC 60601-1 and catheter or probe compatibility may increase the risk of electrical shock to the patient/operator.

- Swan-Ganz pacing-TD catheter
- Percutaneous sheath introducer and contamination shield
- Compatible cardiac output computer, injectate probe, and appropriate connecting cable
- External pulse generator (ventricular demand or A-V sequential)
- External pulse generator cable adapters
- Sterile flush system and pressure transducers
- Bedside pressure monitor system

In addition, the following items should be immediately available if complications arise during catheter insertion: antiarrhythmic drugs, defibrillator, respiratory assist equipment, and temporary pacing equipment.

9.0 Preparation

Use aseptic technique

WARNING: This catheter requires special techniques for insertion and removal. Electrode dislodgment may result from pulling the catheter out through the percutaneous sheath.

Precaution: Avoid forceful wiping or stretching of the catheter during testing and cleaning as not to break the electrode or thermistor wire circuitry.

1. Flush catheter lumens with a sterile solution to ensure patency and to remove air.
2. Check balloon integrity. Inflate the balloon to the recommended volume and check for major asymmetry and for leaks by submerging in sterile saline or water.
3. Connect the catheter's injectate and pressure monitoring lumens to the flush system and pressure transducers. Ensure that the lines and transducers are free of air.
4. Test the thermistor's electrical continuity before insertion. Connect the catheter to the cardiac output computer and check for a "CATHETER FAULT".

10.0 Procedure

Visually inspect for breaches of packaging integrity prior to use.

The following two procedures are offered as an aid to the physician. The first procedure requires fluoroscopy and the second procedure uses pressure monitoring.

Precaution: In a minority of patients, the possibility exists that adequate thresholds will not be obtained or that capture will be lost. If either occurs, the use of a conventional pacing catheter should be considered.

Note: Due to the unique bend at the location of the ventricular electrodes, the catheter cannot be inserted from the inferior vena cava approach (see Specifications). Insertion must be accomplished through the superior vena cava.

10.1 Insertion Under Fluoroscopy

1. Introduce the catheter into the vein through a sheath introducer using percutaneous insertion using modified Seldinger technique.
2. Gently advance the catheter into the right atrium.

Note: When the catheter is near the junction of the right atrium and the superior vena cava of the typical adult patient, the tip has been advanced approximately 40 cm from the right or 50 cm from the left antecubital fossa, or 15 to 20 cm from the subclavian vein.

3. Should the catheter require stiffening during insertion, slowly perfuse the catheter with 5 to 10 ml cold sterile saline or 5% dextrose as the catheter is advanced through a peripheral vessel.
4. Using the syringe provided, inflate the balloon with CO₂ or air to the recommended volume (1.5 ml) printed on the catheter shaft (**Do not use liquid**). Note that an offset arrow on the gate valve indicates the "closed" position.

WARNING: Pulmonary complications may result from improper inflation technique. To avoid damage to the pulmonary artery and possible balloon rupture, do not inflate above the recommended volume. Use the volume-limited syringe provided in the catheter package.

Before reinflation with CO₂ or air, completely deflate the balloon by removing the syringe and opening the gate valve. Do not forcefully aspirate with the syringe, as this may damage the balloon. After deflation, reattach the syringe to the gate valve.

Precaution: It is recommended that the provided syringe be re-attached to the gate valve after balloon deflation to prevent inadvertent injection of liquids into the balloon lumen.

5. Advance the catheter in the usual manner until the inflated balloon is wedged in a central pulmonary artery (Figure 3 on page 95). Deflate the balloon and then pull the catheter back a few centimeters to remove any slack. Avoid lengthy balloon inflation while the catheter is in a wedge position; this is an occlusive maneuver and may result in pulmonary infarction.

Precaution: Catheter looping may occur when excessive length has been inserted, which could result in kinking or knotting (see Complications). If the right ventricle is not entered after advancing the catheter 15 cm beyond entry into the right atrium, the catheter may be looping, or the tip may be engaged in the neck of the vein with only the proximal shaft advancing into the heart. Deflate the balloon and withdraw the catheter until the 20 cm mark is visible. Re-inflate the balloon and advance the catheter.

6. Reinflate the balloon and advance the catheter until the ventricular electrodes contact the right ventricular wall, usually in the inflow tract. Deflate the balloon. Figure 4 on page 96 shows the final catheter position.
7. Determine pacing thresholds. If necessary, manipulate the catheter slightly to obtain good ventricular thresholds (1 or 2 mA). Check for stability of pacing with respiration and adjust catheter position if necessary.
8. Attempt atrial, ventricular, or A-V sequential pacing as required.

10.2 Insertion Using Pressure Monitoring

1. Insert the catheter and advance to the "wedge" position in the usual manner under continuous pressure monitoring (see steps 1 and 2 in **Insertion Under Fluoroscopy**). Figure 2 on page 94 illustrates the characteristic pressure waveforms.

Note: Use of a protective catheter sheath is recommended.

2. Deflate the balloon to verify the presence of a normal pulmonary artery pressure tracing.
3. Reinflate the balloon to determine the minimum inflation volume necessary to obtain a wedge tracing. If a wedge is obtained with less than the maximum recommended volume (see specifications table for balloon inflation capacity), the catheter must be withdrawn to a position where full inflation volume produces a wedge tracing.

CAUTION: Do not pull the catheter across the pulmonic valve while the balloon is inflated to avoid damage to the valve.

Note: When the catheter is near the junction of the right atrium and the superior or inferior vena cava of the typical adult patient, the tip has been advanced approximately 40 cm from the right or 50 cm from the left antecubital fossa, 15 to 20 cm

from the jugular vein, 10 to 15 cm from the subclavian vein, or about 30 cm from the femoral vein.

Note: After deflation, the catheter tip may tend to recoil towards the pulmonic valve and slip back into the right ventricle, requiring that the catheter be repositioned.

10.2.1 Atrial Pacing

1. Connect the "distal atrial" electrode (#3) to the negative pulse generator terminal. Connect the "central atrial" electrode (#4) to the positive terminal.

Note: To facilitate a connection between the catheter and pulse generator, a cable adapter may be required.

2. Adjust the pulse generator output to 0.1 milliamps and the rate to 15% above the patient's heart rate, or to a physiological rate. With the pulse generator turned on, slowly increase the output until atrial pacing occurs (see Figure 6 on page 97). A typical atrial pacing current threshold is 5 milliamps. If atrial pacing does not occur, slowly advance or withdraw the catheter 0.5 cm at a time with the pulse generator operating at 5 milliamps and the balloon deflated.

Note: Diaphragmatic pacing may occasionally occur; it can usually be alleviated by advancing the catheter 0.5 to 1 cm.

3. Check the catheter position once again to make sure it has not been advanced to a permanent wedge position.

Precaution: If atrial and ventricular pacing occur only with the catheter in a permanent wedge position, withdraw the catheter to a pulmonary artery position. In this situation, pacing should not be performed because of the likelihood of pulmonary infarction resulting from the catheter being permanently wedged (see Complications).

If using a contamination shield, extend the distal end towards the introducer valve. Extend the proximal end of the catheter contamination shield to desired length, and secure.

10.2.2 Ventricular Pacing

1. Connect the "distal ventricular" electrode (#1) to the negative pulse generator terminal. Connect the "proximal ventricular" electrode (#2) to the positive terminal.

Note: To facilitate a connection between the catheter and pulse generator, a cable adapter may be required.

2. Adjust the external generator R-wave sensitivity to approximately 3 millivolts to avoid competition between the pulse generator's rate and the patient's rate.
3. Repeat steps 2 and 3 in **Atrial Pacing** to achieve ventricular pacing. Ventricular stimulation at 3 milliamps or less is generally obtainable. Check for adequate sensitivity once pacing is established. Successful ventricular capture is shown in Figure 5 on page 96.
4. Unused electrode connectors must be capped to prevent their contact with a faulty ground.

10.2.3 A-V Sequential Pacing

1. After achieving atrial and ventricular pacing, connect the two ventricular electrodes to the A-V sequential pulse generator; attempt pacing. Successful A-V sequential pacing is shown in Figure 7 on page 97.

Note: To facilitate a connection between the catheter and pulse generator, a cable adapter may be required.

11.0 Maintenance and Use *in situ*

1. Keep the catheter tip centrally located in a main branch of the pulmonary artery. Ideally, the catheter tip should be located near the hilum of the lungs. The tip migrates towards the periphery of the lungs during balloon inflation. Therefore, a central location before inflation is important. Keep the tip in a position where a full or near full (1.0 to 1.5 ml) inflation volume is necessary to produce a "wedge" tracing.
2. Anticipate spontaneous catheter tip migration towards the periphery of the pulmonary bed. To avoid possible damage to the pulmonary artery, continuously monitor the catheter tip pressure while the catheter is in place. If a wedge pressure tracing is observed when the balloon is deflated, the catheter should be pulled back to a central pulmonary artery position.

Spontaneous catheter tip migration towards the periphery of the lung occurs during cardiopulmonary bypass. Partial catheter withdrawal (3 to 5 cm) just before bypass should be considered, as it may help reduce distal migration and prevent permanent catheter wedging postbypass.

After termination of bypass, the catheter may require repositioning. Check the distal pulmonary artery tracing before inflating the balloon.

Precaution: Over a period of time, the catheter tip may migrate towards the periphery of the pulmonary bed and lodge in a small vessel. Damage may occur either by prolonged occlusion or by over-distension of the vessel upon reinflation of the balloon (see Complications).

3. **Reinflation of the balloon while in the pulmonary artery should be performed gradually and with caution while observing the pressure tracing.**

Note: Inflation is usually associated with a feeling of resistance. On release, the syringe plunger should usually spring back. If no resistance to inflation is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may continue to be used for hemodynamic monitoring. However, be sure to take precautions to prevent infusion of air or liquid into the balloon lumen.

4. Inflate the balloon slowly to the **minimum** volume needed to obtain PAOP (never exceeding the recommended volume). If PAOP is obtained at volumes less than 1.0 ml, pull the catheter back to a position at which the full or near full inflation volume (1.0 to 1.5 ml) produces a wedge pressure tracing.
5. Measure pulmonary artery occlusion pressure only when necessary and keep the number of wedge pressure measurements and wedge time to a minimum (two respiratory cycles or 10 to 15 seconds), especially in patients with pulmonary hypertension.

In some patients, pulmonary arterial end-diastolic pressure can often be substituted for pulmonary artery occlusion pressure if the pressures are nearly identical, obviating the need for repeated balloon inflation.

Note: Avoid prolonged maneuvers to obtain wedge pressure. If difficulties are encountered, give up the "wedge".

6. Keep the inflation syringe attached to the gate valve to prevent inadvertent injection of liquid into the balloon inflation lumen.
7. Keep pressure monitoring lumens patent by intermittent flush or continuous slow infusion with heparinized saline solution.

WARNING: To avoid pulmonary artery rupture, never flush the catheter when the balloon is wedged in the pulmonary artery.

8. Periodically check IV lines, pressure lines, and transducer domes to keep them free of air. Also ensure that connecting lines and stopcocks remain tightly fitted.
9. Infusion of viscous solutions (e.g., whole blood or albumin) is not recommended, as they flow too slowly and may occlude the catheter lumen.
10. The catheter should remain indwelling only as long as is required by the patient's condition.

Precaution: The incidence of complications increases significantly with indwelling periods longer than 72 hours.

12.0 Cardiac Output Determination

To determine cardiac output by thermodilution, a known amount of sterile solution of known temperature is injected into the right atrium or vena cava, and the resultant change in blood temperature is measured in the pulmonary artery by the catheter thermistor. Cardiac output is inversely proportional to the integrated area under the resulting curve. This method has been shown to provide good correlation with the direct Fick method and dye dilution technique for cardiac output determination.

Consult the references on the use of iced versus room temperature injectate or open versus closed injectate delivery systems.

Refer to the appropriate cardiac output computer manual for specific instructions in the use of thermodilution catheters for cardiac output determination. Correction factors or computation constants needed to correct for indicator heat transfer are given in the specifications.

Edwards cardiac output computers require that a computation constant be used to correct for injectate temperature rise as it passes through the catheter. The computation constant is a function of injectate volume, temperature, and catheter dimensions. The computation constants listed in the specifications have been determined *in vitro*.

13.0 MRI Information

MR Unsafe

The Swan-Ganz device is MR unsafe as the result of the device containing metallic components, which experience RF-induced heating in the MRI environment; therefore the device poses hazards in all MRI environments.

14.0 Complications

All invasive procedures inherently involve some patient risks. Although serious complications associated with pulmonary artery catheters are relatively uncommon, the physician is advised before deciding to use the catheter to consider and weigh the potential benefits and risks associated with the use of the catheter against alternative procedures. Strict adherence to the foregoing instructions and the awareness of possible complications have been the most significant factors in reducing the incidence of complications.

14.1 Perforation of the Pulmonary Artery

Factors associated with the development of fatal pulmonary artery rupture during the use of flow-directed balloon-tipped catheters are pulmonary hypertension, advanced age, cardiac surgery with hypothermia and anti-coagulation, and distal catheter tip migration.

Extreme care should therefore be exercised during the measurement of pulmonary artery occlusion pressure in patients with pulmonary hypertension. The period of time during which the balloon remains inflated and wedged in these patients should be minimal and limited to two respiratory cycles, or 10 to 15 seconds.

Spontaneous catheter tip migration towards the periphery of the lung occurs during cardiopulmonary bypass. Partial catheter withdrawal just before bypass should be considered, as it may help reduce distal migration and prevent permanent catheter wedging postbypass. After termination of bypass, the catheter may require repositioning. Check the distal pulmonary artery tracing before inflating the balloon. **A central location of the catheter tip near the hilum of the lung may prevent pulmonary artery perforation.**

14.2 Pulmonary Infarction

Tip migration with spontaneous wedging, air embolism, and thromboembolism can lead to this complication.

14.3 Cardiac Arrhythmias

Although usually transient and self-limiting, arrhythmias may occur during insertion, removal, or following displacement of the tip from the pulmonary artery into the right ventricle. Whereas premature ventricular contractions are the most commonly encountered arrhythmias, ventricular tachycardia and atrial and ventricular fibrillation have also been reported. ECG monitoring and the immediate availability of antiarrhythmic drugs and defibrillating equipment are recommended. Prophylactic lidocaine may be helpful in decreasing the incidence of ventricular arrhythmias during catheterization.

14.4 Knotting

Flexible catheters have been reported to form knots, most often as a result of looping in the right ventricle. Sometimes the knot can be resolved by insertion of a suitable guidewire and manipulation of the catheter under fluoroscopy. If the knot does not include any intracardiac structures, the knot may be gently tightened and the catheter withdrawn through the site of entry.

14.5 Sepsis/Infection

Positive catheter-tip cultures resulting from contamination and colonization have been reported, as well as incidences of septic and aseptic vegetation in the right heart. Increased risks of septicemia and bacteremia have been associated with blood sampling, the infusing of fluids, and catheter-related thrombosis. Preventive measures should be taken to guard against infection.

14.6 Other Complications

Pulmonary artery catheters have also been associated with right bundle branch block and complete heart block, tricuspid and pulmonic valve damage, pneumothorax, thrombosis, blood loss, cardiac structure/wall injury or damage, hematoma, embolism, anaphylaxis, cardiac tissue/artery burn, thrombophlebitis, perforation of the right ventricle, and inadvertent stimulation/pacing.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance for this medical device. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

15.0 How Supplied

Swan-Ganz catheters are supplied sterile, unless otherwise stated. Do not use if package has been previously opened or damaged.

Catheters are for single use only. Do not clean or resterilize a used catheter.

The packaging is designed to avoid crushing of the catheter and to protect the balloon from exposure to the atmosphere. It is therefore recommended that the catheter remain inside the package until use.

16.0 Storage

Store in a cool, dry place.

17.0 Operating Conditions

Intended to operate under physiological conditions of the human body.

18.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the recommended time may result in balloon deterioration, since the natural latex rubber in the balloon is acted upon and deteriorated by the atmosphere and could lead to illness or an adverse event as the device may not function as originally intended.

19.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): . . . 800.822.9837

Outside the U.S. and Canada (24 hours): . . . 949.250.2222

In the UK: . . . 0870 606 2040 - Option 4

Ireland: . . . 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

20.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the latest version of the monitoring system operator's manual for more information.

Refer to the symbol legend at the end of this document.

Computation Constants			
Model Number	Temp (°C)	Volume (ml)	D200F7, D205F7 Computation Constants (cc)*
	0 - 5	10	0.542
		5	0.247
	19 - 22	3	0.132
		10	0.578
		5	0.274
		3	0.154
Computation Constants for CO-Set+			
Cold Injectate		10	0.561
10 ml: 6-12 °C		5	0.259
5 ml: 8-16 °C		3	---
Room Temperature		10	0.608
Injectate		5	0.301
5 or 10 ml: 18-25 °C		3	---
*CC = (1.08)C _T (60)V _I			
Specifications			
Function	Atrio-Ventricular Pacing and Thermodilution		
Model Number COC Compatibility	D200F7 Edwards	D205F7 Edwards	
Usable Length (cm)	110	110	
Catheter Body French Size	7F (2.3 mm)	7F (2.3 mm)	
Required Introducer Size	8F (2.7 mm)	8F (2.7 mm)	
Body Color	Yellow	Yellow	
Diameter of Inflated Balloon (mm)	13	13	
Balloon Inflation Capacity (ml)	1.5	1.5	
Injectate Port (cm from tip)	30	29.5	
Injectate Lumen Volume (ml)	1.20	0.93	
Electrode Connector	0.40" x 0.080" (10.16 mm x 2.03 mm) dia pin plugs		0.40" x 0.080" (10.16 mm x 2.03 mm) dia pin plugs
Atrial Electrodes (cm from tip)	28.5, 31.0, 33.5		25.5, 28.0, 30.5
Ventricular Electrodes (cm from tip)	18.5, 19.5		16.5, 17.5
Compatible Guidewire Diameter (in)	0.020 (0.51 mm)		0.020 (0.51 mm)
Frequency Response Distortion at 10 Hz	< 3 dB		< 3 dB
Thermistor Location (cm from tip)	4		4
All specifications given are nominal values.			

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
Figure 5: Surface ECG showing ventricular pacing

- Figure 5 : ECG de surface indiquant la stimulation ventriculaire
- Abbildung 5: Oberflächen-EKG bei Kammerstimulation
- Figura 5: ECG de superficie que muestra la estimulación ventricular
- Figura 5: ECG superficiale che mostra la stimolazione ventricolare
- Afbeelding 5: Oppervlakte-ECG met het verloop van ventriculaire stimulatie
- Figur 5: Overflate-EKG, der viser ventrikulær pacing
- Figur 5: Yt-EKG som visar kammarstimulering
- Εικόνα 5: ΗΚΓ επιφανείας που δείχνει κοιλιακή βηματοδότηση
- Figura 5: ECG de superficie a apresentar estimulação ventricular
- Obrázek 5: Povrchové EKG ukazující komorovou stimulaci
- 5. ábra: Kamrai ingerlést mutatós felszíni EKG
- Rysunek 5: Powierzchniowe EKG przedstawiające stymulację komorową
- Obrázok 5: Povrchové EKG zobrazujúce komorovú stimuláciu
- Figur 5: Overflate-EKG viser ventrikulær pacing
- Kuva 5: Pinta-EKG, jossa näkyy kammion tahdistus
- Фигура 5: Поверхностна ЕКГ, показваща камерно пейсирание
- Figura 5: ECG de suprafață care prezintă stimularea ventriculară
- Joonis 5. Pinna-EKG, mis näitab vatsakese stimulaatsiooni
- 5 pav. Paviršiaus EKG, rodanti skilvelių stimuliavimą
- 5. attēls. Virsmas EKG uzņēmumi, kuros attēlota ventrikulārā kardiostimulācija
- Şekil 5: Ventriküler pacing gösteren yüzey EKG görüntüsü
- Рис. 5. Кривая ЭКГ во время желудочковой стимуляции
- Slika 5: Površinski EKG koji prikazuje ventrikularni pejsing
- Slika 5: Površinski EKG koji prikazuje elektrostimulaciju ventrikula

DC-2026-7


Figure 6: Surface ECG showing atrial pacing

- Figure 6 : ECG de surface indiquant la stimulation auriculaire
- Abbildung 6: Oberflächen-EKG bei Vorhofstimulation
- Figura 6: ECG de superficie que muestra la estimulación auricular
- Figura 6: ECG superficiale che mostra la stimolazione atriale
- Afbeelding 6: Oppervlakte-EKG met het verloop van atriale stimulatie
- Figur 6: Overflade-EKG, der viser atriel pacing
- Figur 6: Yt-EKG som visar förmaksstimulering
- Εικόνα 6: ΗΚΓ επιφανείας που δείχνει κοιλική βηματοδότηση
- Figura 6: ECG de superficie a apresentar estimulação auricular
- Obrázek 6: Povrchové EKG ukazující síňovou stimulaci
- 6. ábra: Pitvari ingerlést mutató felszíni EKG
- Rysunek 6: Powierzchniowe EKG przedstawiające stymulację przedsionkową
- Obrázok 6: Povrchové EKG zobrazujúce predsieňovú stimuláciu
- Figur 6: Overflate-EKG viser atrial pacing
- Kuva 6: Pinta-EKG, jossa näkyy eteisen tahdistus
- Фигура 6: Поверхностна ЕКГ, показваща предсърдно пейсиране
- Figura 6: ECG de suprafață care prezintă stimularea atrială
- Joonis 6. Pinna-EKG, mis näitab koja stimulatsiooni
- 6 pav. Paviršiaus EKG, rodanti prieširdžių stimuliavimą
- 6. attēls. Virsmas EKG uzņēmumi, kuros attēlota atriālā kardiostimulācija
- Şekil 6: Atrial pacing gösteren yüzey EKG görüntüsü
- Рис. 6. Кривая ЭКГ во время предсердной стимуляции
- Slika 6: Površinski EKG koji prikazuje atrijalni pejsing
- Slika 6: Površinski EKG koji prikazuje elektrostimulaciju atrijske

DC-2026-8


Figure 7: Surface ECG showing A-V sequential pacing

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	English	Français	Deutsch	Español	Italiano	Nederlands
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello	Modelnummer
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad
	Minimum introducer size	Taille minimale de l'introducteur	Mindestgröße der Einführhilfe	Tamaño de introductor mínimo	Misura minima dell'introduttore	Minimum introducergrootte
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum
	Balloon capacity	Capacité du ballonnet	Ballonkapazität	Capacidad del balón	Capacità del palloncino	Balloncapaciteit
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo	Unieke instrumentidentificatie
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum
	Size	Taille	Größe	Tamaño	Misura	Afmetingen
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig
	Contains or presence of Phthalates: bis (2-ethylhexyl) phthalate (DEHP)	Traces ou présence de phthalates : phthalate de bis (2-éthylhexyle) (DEHP)	Enthält Phthalate oder kann Spuren von Phthalaten enthalten: Bis(2-ethylhexyl) phthalat (DEHP)	Contiene o hay presencia de ftalatos: Di[2-etilhexil] ftalato (DEHP)	Contiene o sono presenti ftalati: di (2-etilesil) ftalato (DEHP)	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl) ftalaat (DEHP)

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	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
#	Mudeli number	Modelio numeris	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
⌀	Vālislābimōot	Išorinis skersmuo	Ārējais diametrs	Dış çap	Внешний диаметр	Spoljašnji prečnik	Vanjski promjer
— cm —	Kasutatav pikkus	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebljiva dužina	Upotrebljiva dužina
GW	Juhtetraadi soovituslik suurus	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs	Kılavuz Tel İçin Önerilen Boyut	Рекомендованный размер проводника	Preporučena veličina žice vodice	Preporučena veličina žice vodilice
I	Sisesti väikseim suurus	Minimalus įvediklio dydis	Minimālais ievadītāja izmērs	Minimum introdüser boyutu	Минимальный размер интродьюсера	Minimalna veličina uvodnika	Minimalna veličina uvodnog instrumenta
⚠	Ettevaatust	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
⊗	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не использовать повторно	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
QTY	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina	Količina
LOT	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
⌚	Kōlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
BC	Balloon mahutavus	Balionių talpa	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
☀☔	Sāilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
STERILE EO	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidadom	Sterilizirano etilen-oksidadom
UDI	Seadme kordumatu identiftseerimistunnus	Unikalusis priemonės identifikatorius	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator proizvoda	Jedinstveni identifikator proizvoda
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
🏭	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvodač	Proizvodač
🏭	Tootmiskuupäev	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
SZ	Suurus	Dydis	Izmērs	Boyut	Размер	Veličina	Veličina
LATEX	Sisaldab naturaalsset kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
PHT DEHP	Sisaldab ftalaate: bis-(2-etüülheksüül)ftalaat (DEHP)	Sudėtyje yra ftalatų: bis-(2-etilheksil) ftalatas (DEHP)	Vērojama ftalātu klātbūtne vai satur ftalātus: bis-(2-etilheksil)ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭГФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)

	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Palun lugege kasutusjuhiseid veebisaidilt	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na internetskoj stranici
	Järgige kasutusjuhendit, mille leiata veebisaidilt	Vadovautis naudojimo instrukcijomis, pateiktomis svetainėje	Levėrojiēt lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Ohtlik magnetresonants-tomograafias	MR nesaugus	Nedrikst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR	Nije sigurno kod pregleda MR-om
	Vaadake kasutusjuhendit	Žr. naudojimo instrukcijas	Skatīt lietošanas instrukciju	Kullanım talimatlarına başvurun	См. инструкции по применению	Pogledajte uputstvo za upotrebu	Pogledajte upute za upotrebu
	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis	Levėrojiēt lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu
	Ärge resteriliseerige	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Mittepürogeenne	Nepirogeniškas	Nepirogēna	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno
	B-tüüpi kohaldatav osa	B tipo liečiamoji dalis	B tipa lietotā daļa	B Tipi uygulanan parça	Рабочая часть типа В	Primenjeni deo tipa B	Primjenjen dio tipa B
	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Defibrillatsioonikindel CF-tüüpi rakendusosa	Defibriliavimo poveikiui atspari CF tipo darbinė dalis	Pret defibrilāciju noturīga CF tipa daļa, kas saskaras ar pacientu	Defibrilasyon Korumalı CF Tipi hastayla temas eden parça	Контактирующий с пациентом элемент типа CF с защитой от разряда дефибриллятора	Primenjeni deo tipa CF otporan na defibrilaciju	Primjenjen dio tipa CF otporan na defibrilaciju
	Sisaldab ohtlikke aineid	Sudėtyje yra pavojingų medžiagų	Satur bīstamas vielas	Tehlikeli maddeler içerir	Содержит вредные вещества	Sadrži opasne supstance	Sadrži opasne tvari
	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Importija	Importuotojas	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Meditsiiniseade	Medicinos priemonė	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Kahekordne steriilne kattesüsteem	Dvigubo sterilaus barjero sistema	Divkārsa sterilitātes aizsargslāņa sistēma	Çift steril bariyer sistemi	Двойная барьерная система для стерилизации	Sistem sa dve sterilne barijere	Sustav dvostruke sterilne barijere
Märkus. Kõik sümbolid ei pruugi antud toote etiketil esineda. • Pastaba. Šio gaminio etiketėje gali būti pateikti ne visi simboliai. • Piezīme. Ši izstrādājuma marķējumā var nebūt iekļauti visi simboli. • Not: Bu ürünün etiketinde tüm semboller yer almayış olabılır. • Примечание. На этикетках данного изделия могут указываться не все символы. • Napomena: Na ambalaži proizvoda se možda ne nalaze svi simboli. • Napomena: Na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.							

This page intentionally left blank.
Cette page est intentionnellement blanche.
Diese Seite wurde absichtlich frei gelassen.
Pàgina en blanco.
Questa pagina è stata lasciata intenzionalmente in bianco.
Deze pagina is opzettelijk leeg gelaten.
Denne side skal være tom.
Denna sida har avsiktligt lämnats tom.
Η σελίδα αυτή αφέθηκε σκόπιμα κενή.
Esta página foi intencionalmente deixada em branco.
Tato strana je záměrně ponechána prázdná.
Ez az oldal szándékosan maradt üresen.
Niniejsza strona pozostaje pusta.
Táto stránka je zámerne ponechaná prázdna.
Denne siden er tom med hensikt.
Tämä sivu on jätetty tarkoituksella tyhjäksi.
Тази страница е умишлено оставена празна.
Această pagină a fost lăsată albă în mod intenționat.
See leht on tahtlikult tühjaks jäetud.
Šis puslapis tikslingai paliktas tuščias.
Ši lappuse atstāta tukša ar nolūku.
Bu sayfa özellikle boş bırakılmıştır.
Эта страница оставлена пустой намеренно.
Ova stranica je namerno ostavljena prazna.
Ova stranica namjerno je ostavljena praznom

This page intentionally left blank.
Cette page est intentionnellement blanche.
Diese Seite wurde absichtlich frei gelassen.
Pàgina en blanco.
Questa pagina è stata lasciata intenzionalmente in bianco.
Deze pagina is opzettelijk leeg gelaten.
Denne side skal være tom.
Denna sida har avsiktligt lämnats tom.
Η σελίδα αυτή αφέθηκε σκόπιμα κενή.
Esta página foi intencionalmente deixada em branco.
Tato strana je záměrně ponechána prázdná.
Ez az oldal szándékosan maradt üresen.
Niniejsza strona pozostaje pusta.
Táto stránka je zámerne ponechaná prázdna.
Denne siden er tom med hensikt.
Tämä sivu on jätetty tarkoituksella tyhjäksi.
Тази страница е умишлено оставена празна.
Această pagină a fost lăsată albă în mod intenționat.
See leht on tahtlikult tühjaks jäetud.
Šis puslapis tikslingai paliktas tuščias.
Ši lappuse atstāta tukša ar nolūku.
Bu sayfa özellikle boş bırakılmıştır.
Эта страница оставлена пустой намеренно.
Ova stranica je namerno ostavljena prazna.
Ova stranica namjerno je ostavljena praznom

EC REP

Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

01/25
10007454004 B / DOC-0210800 B
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU