Edwards

DIRECTORY

English (en)	1	Português (pt)	46	Română (ro)	86
Français (fr)	6	Česky (cs)	51	Eesti (et)	91
Deutsch (de)	11	Magyar (hu)	56	Lietuvių (lt)	96
Español (es)	16	Polski (pl)	61	Latviešu (lv)	101
Italiano (it)	21	Slovensky (sk)	66	Türkçe (tr)	106
Nederlands (nl)	26	Norsk (no)	71	Русский (ru)	111
Dansk (da)	31	Suomi (fi)	76	Srpski (sr)	116
Svenska (sv)	36	Български (bg)	81	Hrvatski (hr)	121
Ελληνικά (el)	41				

Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai
■ Attēli ■ Şekiller ■ Рисунок ■ Slike ■ Slike

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos
■ Legenda se symboly ■ Jelmagyarozat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sūbolite seletus ■ Simbolių paaiškinimas
■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola


English

Swan-Ganz

Hi-Shore True Size Thermodilution Catheters: 131F7, 131F7P, 141F7
Thermodilution Catheter (with contamination shield): 131VF7P
True Size Thermodilution Catheter (CardioCath): 143TF7
Thermodilution Catheter (S-Tip): 151F7
VIP Thermodilution Catheters: 831F75, 831F75P
Venous Infusion Port (VIP) Thermodilution Catheter (with contamination shield): 831VF75P
VIP+ Tri-Lumen Infusion Thermodilution Catheters: 834F75, 834F75P
True Size Thermodilution Catheters: 096F6, 096F6P, TS105F5, 132F5
Polymer Blend True Size C-Tip ControlCath Thermodilution Catheters: C144F7, C145F6N
Polymer Blend True Size S-Tip ControlCath Thermodilution Catheter: S144F7
Synthetic ControlCath Thermodilution Catheter: C146F7
Catheter/Introducer System Catheter: S9FC146F7
Polymer Blend Torque Support True Size T-Tip TD: T173F6

TS105F5, 096F6, C145F6N, and T173F6 are not available in EU.

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific country.

Swan-Ganz VIP+ Tri-Lumen Infusion Thermodilution Catheter (Model 834F75)

Model 834F75 is illustrated. The models listed above contain some but not all of the features shown. Model C146F7 and S9FC146F7 are illustrated – No components of these models or the products they contain are made from natural rubber latex or dry natural rubber.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

CAUTION: This product contains natural rubber latex which may cause allergic reactions (with the exception of models C146F7 and S9C146F7).

For Single Use Only

Please refer to Figure 1 on page 126 thru Figure 3 on page 128.

1.0 Description

The device is used by medical professionals who have been trained in safe use of invasive hemodynamic technologies and clinical usage of pulmonary artery catheters as part of their respective institutional guidelines.

The Swan-Ganz catheters are flow-directed pulmonary artery catheters used to monitor hemodynamic pressures. The distal (pulmonary artery) port on Swan-Ganz catheters allow sampling of mixed venous blood for the assessment of oxygen transport balance and the calculation of derived parameters such as oxygen consumption, oxygen utilization coefficient, and intrapulmonary shunt fraction when used with compatible hemodynamic monitor and accessories.

The Swan-Ganz Hi-Shore, CardioCath, and "S-Tip" thermodilution catheters (Models 141F7, 143TF7 and 151F7 respectively) have the same functions as the standard Swan-Ganz thermodilution catheter (Models 131F7 and 131VF7). The "S-Tip" catheter (Model 151F7) is designed for femoral vein insertion. The Swan-Ganz Hi-Shore thermodilution catheter (Model 141F7), is slightly stiffer than the standard Swan-Ganz catheter, and may be used when more torque control and maneuverability is needed (i.e., from the femoral approach). The CardioCath thermodilution catheter (Model 143TF7) is made from the same material as the Swan-Ganz Hi-Shore TD catheter but is even stiffer due to the stiffening rod in one lumen. CardioCath catheter has a "C" tip bend. In addition to the standard distal (pulmonary artery) and injectate lumens, the Swan-Ganz VIP thermodilution catheter (Models 831F75, 831F75P and 831VF75P) has an additional lumen that provides direct access to the right atrium. The Swan-Ganz VIP+ tri-lumen infusion thermodilution catheter (Models 834F75 and 834F75P) is equipped with a right atrium lumen and an additional lumen. The catheter's right ventricular (RV) lumen terminates 19 cm from the catheter tip and the right atrial (RA) lumen terminates at 31 cm. The VIP lumen allows for continuous infusion, even during cardiac output determinations.

ControlCath thermodilution catheters are made from a proprietary polymer blend and offer improved torque control and maneuverability over the standard Swan-Ganz thermodilution catheter. ControlCath thermodilution catheters feature two tip configurations designed for flexibility in insertion. The "S-Tip" configuration is designed specifically for femoral vein insertion,

Edwards, Edwards Lifesciences, the stylized E logo, CardioCath, ControlCath, CO-Set, CO-Set+, Hi-Shore, Swan, Swan-Ganz, True Size, VIP, and VIP+ are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

while the “C-Tip” configuration allows for both femoral or SVC approaches. ControlCath catheters also offer improved column strength for easy passage through the introducer sheath, and will not lose their torque properties at body temperature as much as other thermodilution catheters. This results in excellent handling characteristics and consistent waveform fidelity throughout the procedure.

Torque support thermodilution catheters are made from a proprietary polymer blend and offer improved torque control and maneuverability over standard Swan-Ganz thermodilution catheters. Torque support thermodilution catheters feature three tip configurations designed for flexibility in insertion. The “S-Tip” configuration is designed specifically for femoral vein insertion, while the “C-Tip” and “T-Tip” configurations allow for both femoral or SVC approaches. Torque support catheters also offer improved column strength for easy passage through the introducer sheath, and will not lose their torque properties at body temperature as much as other thermodilution catheters. This results in excellent handling characteristics and consistent waveform fidelity throughout the procedure.

The intravascular catheter is inserted through the central vein into the right side of the heart and is advanced towards the pulmonary artery. Route of insertion can be internal jugular, femoral, antecubital and brachial veins. The body parts in contact are the atrium, ventricles, pulmonary artery and circulatory system.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended for use in adult critically ill or surgical patient populations. The device has not yet been tested in the pediatric population or in pregnant or lactating women.

2.0 Intended Use / Purpose

Swan-Ganz catheters are pulmonary artery catheters intended for short-term use in the central circulatory system for patients who require intracardiac hemodynamic monitoring, blood sampling and infusing solutions. When used with a compatible monitoring platform and accessories, Swan-Ganz catheters offer a comprehensive hemodynamic profile to help clinicians assess cardiovascular function and guide treatment decisions.

3.0 Indications for Use

Swan-Ganz catheters are diagnostic and monitoring tools used for hemodynamic monitoring of adult critically ill patients including but not limited to post major surgical recovery, trauma, sepsis, burns, pulmonary disease, pulmonary failure, cardiac disease including heart failure.

4.0 Contraindications

These products contain metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

Patients with either recurrent sepsis, or with hypercoagulopathy, in which the catheter could serve as a focal point for septic or bland thrombus formation, should not be considered candidates for a balloon flotation catheter.

5.0 Warnings

No absolute contraindications to the use of flow-directed pulmonary artery catheters exist. However, a patient with a left bundle branch block may develop a right bundle branch block during catheter insertion, resulting in complete heart block. In such patients, temporary pacing modes should be immediately available. Electrocardiographic monitoring during catheter passage is encouraged and is particularly important in the presence of either of the following conditions:

- Complete left bundle branch block, in which the risk of complete heart block is somewhat increased.
- Wolff-Parkinson-White syndrome and Ebstein’s malformation, in which the risk of tachyarrhythmias is present.

Do not modify or alter the product in any way. Alteration or modification may affect product performance.

Air should never be used for balloon inflation in any situation where it may enter the arterial circulation, e.g., in all pediatric patients and in adults with suspected right to left intracardiac or intrapulmonary shunts. Bacteria-filtered carbon dioxide is the recommended inflation medium because of its rapid absorption into the blood in the event of balloon rupture within the circulation. Carbon dioxide diffuses through the latex balloon, diminishing the balloon’s flow-directed capability after 2 to 3 minutes of inflation.

Do not leave the catheter in a permanent wedge position. Furthermore, avoid lengthy balloon inflation while the catheter is in a wedge position; this occlusive maneuver may result in pulmonary infarction.

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data

to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Cleaning and resterilization will damage the integrity of the latex balloon. Damage may not be obvious during routine inspection.

6.0 Precautions

Failure of a balloon flotation catheter to enter the right ventricle or pulmonary artery is rare, but may occur in patients with an enlarged right atrium or ventricle particularly if the cardiac output is low or in the presence of tricuspid or pulmonic incompetence or pulmonary hypertension. Use of a Swan-Ganz Hi-Shore thermodilution catheter (Model 141F7) may be helpful in these patients. Deep inspiration by the patient during advancement may also facilitate passage.

Clinicians using the device should be familiar with the device and understand its applications prior to use.

7.0 Recommended Equipment

WARNING: Compliance to IEC 60601-1 is only maintained when the catheter or probe (CF applied part, defibrillation proof) is connected to a patient monitor or equipment that has a Type CF defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the catheter or probe. Failure to ensure monitor or equipment compliance to IEC 60601-1 and catheter or probe compatibility may increase the risk of electrical shock to the patient/operator.

- Swan-Ganz catheter; ControlCath thermodilution catheter or torque support thermodilution catheter.
- Percutaneous sheath introducer and contamination shield.
- Compatible cardiac output computer, compatible injectate probe, and connecting cable or compatible computer.
- Sterile flush system and pressure transducers.
- Bedside ECG and pressure monitor system.

In addition, the following items should be immediately available if complications arise during catheter insertion: antiarrhythmic drugs, defibrillator, respiratory assist equipment, and a means for temporary pacing.

8.0 Catheter Preparation

Use aseptic technique.

Visually inspect for breaches of packaging integrity prior to use.

Note: Use of a protective catheter sheath is recommended.

Precaution: Avoid forceful wiping or stretching of the catheter during testing and cleaning as not to break the thermistor wire circuitry if present.

Step	Procedure
1	Flush catheter lumens with a sterile solution to ensure patency and to remove air.
2	Check balloon integrity by inflating it to the recommended volume. Check for major asymmetry and for leaks by submerging in sterile saline or water. Deflate balloon before insertion.
3	Connect the catheter’s injectate and pressure monitoring lumens to the flush system and pressure transducers. Ensure that the lines and transducers are free of air.
4	Test the thermistor’s electrical continuity before insertion (refer to computer operations manual for detailed information).

9.0 Insertion Procedure

Swan-Ganz catheters can be inserted at the patient’s bedside without the aid of fluoroscopy, guided by continuous pressure monitoring. Simultaneous pressure monitoring from the distal lumen is recommended. Fluoroscopy is recommended for femoral vein insertion.

Note: The catheter should pass easily through the right ventricle and pulmonary artery and into a wedge position in less than a minute.

Note: Should the catheter require stiffening during insertion, slowly perfuse the catheter with 5 ml to 10 ml of cold sterile saline or 5% dextrose as the catheter is advanced through a peripheral vessel.

Although a variety of techniques can be used for insertion, the following guidelines are provided as an aid to the physician:

Step	Procedure
1	Introduce the catheter into the vein through a sheath introducer using percutaneous insertion using modified Seldinger technique.

Step	Procedure
2	Under continuous pressure monitoring, with or without the aid of fluoroscopy, gently advance the catheter into the right atrium. Entry of the catheter tip into the thorax is signaled by an increased respiratory fluctuation in pressure. Figure 3 on page 128 shows the characteristic intracardiac and pulmonary pressure waveforms. Note: When the catheter is near the junction of the right atrium and the superior or inferior vena cava of the typical adult patient, the tip has been advanced approximately 40 cm from the right or 50 cm from the left antecubital fossa, 15 to 20 cm from the jugular vein, 10 to 15 cm from the subclavian vein, or about 30 cm from the femoral vein.
3	Using the syringe provided, inflate the balloon with CO₂ or air to the maximum recommended volume. Do not use liquid. Note that an offset arrow on the gate valve indicates the “closed” position. Note: Inflation is usually associated with a feeling of resistance. On release, the syringe plunger should usually spring back. If no resistance to inflation is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may continue to be used for hemodynamic monitoring. However, be sure to take precautions to prevent infusion of air or liquid into the balloon lumen. WARNING: Pulmonary complications may result from improper inflation technique. To avoid damage to the pulmonary artery and possible balloon rupture, do not inflate above the recommended volume.
4	Advance the catheter until pulmonary artery occlusion pressure (PAOP) is obtained, then passively deflate the balloon by removing the syringe from the gate valve. Do not forcefully aspirate as this may damage the balloon. After deflation, re-attach the syringe. Note: Avoid prolonged maneuvers to obtain wedge pressure. If difficulties are encountered, give up the “wedge”. Note: Before re-inflation with CO₂ or air, completely deflate the balloon by removing the syringe and opening the gate valve. Precaution: It is recommended that the provided syringe be re-attached to the gate valve after balloon deflation to prevent inadvertent injection of liquids into the balloon lumen. Precaution: If a right ventricular pressure tracing is still observed after advancing the catheter several centimeters beyond the point where the initial right ventricular pressure tracing was observed, the catheter may be looping in the right ventricle which can result in kinking or knotting of the catheter (see Complications). Deflate the balloon and withdraw the catheter into the right atrium. Reinflate the balloon and readvance the catheter to a pulmonary artery wedge position, then deflate the balloon. Precaution: Catheter looping may occur when excessive length has been inserted, which could result in kinking or knotting (see Complications). If the right ventricle is not entered after advancing the catheter 15 cm beyond entry into the right atrium, the catheter may be looping, or the tip may be engaged in a neck vein with only the proximal shaft advancing into the heart. Deflate the balloon and withdraw the catheter until the 20 cm mark is visible. Reinflate the balloon and advance the catheter.
5	Reduce or remove any excessive length or loop in the right atrium or ventricle by slowly pulling the catheter back approximately 2 to 3 cm. Precaution: Do not pull the catheter across the pulmonic valve while the balloon is inflated to avoid damage to the valve.

Step	Procedure
6	Reinflate the balloon to determine the minimum inflation volume necessary to obtain a wedge tracing. If a wedge is obtained with less than the maximum recommended volume (see specifications table for balloon inflation capacity), the catheter must be withdrawn to a position where full inflation volume produces a wedge tracing. Precaution: Overtightening the proximal Tuohy-Borst adaptor of the contamination shield may impair catheter function by potentially compressing and occluding the lumens.
7	Confirm final catheter tip position with chest X-ray. Note: If using a contamination shield, extend the distal end towards the introducer valve. Extend the proximal end of the catheter contamination shield to desired length, and secure.

Note: After deflation, the catheter tip may tend to recoil towards the pulmonic valve and slip back into the right ventricle, requiring that the catheter be repositioned.

10.0 Guidelines for Femoral Insertion

Fluoroscopy is recommended for femoral vein insertion.

Note: The “S-Tip” catheter is designed for femoral vein insertion only.

Precaution: Femoral insertion may lead to a redundancy of the catheter length in the right atrium and difficulties in obtaining a pulmonary artery wedge (occlusion) position.

Precaution: With femoral insertion, it is possible to transfix the femoral artery in some situations during percutaneous entry into the vein. Proper femoral vein puncture technique should be followed, including removal of the innermost occluding stylet when the insertion set needle is advanced toward the vein.

- When advancing the catheter into the inferior vena cava, the catheter may slip into the opposite iliac vein. Pull the catheter back into the ipsilateral iliac vein, inflate the balloon, and let the bloodstream carry the balloon into the inferior vena cava.
- If the catheter does not pass from the right atrium into the right ventricle, it may be necessary to change the orientation of the tip. Gently rotate the catheter and simultaneously withdraw it several centimeters. Care must be exercised so that the catheter is not kinked as it is rotated.
- If difficulty is encountered in positioning the catheter, a suitable sized guidewire may be inserted to stiffen the catheter.

Precaution: To avoid damage to intracardiac structures, do not advance the guidewire beyond the catheter tip. The tendency for thrombus formation will increase with the duration of guidewire use. Keep the period of time that the guidewire is used to a minimum; aspirate 2 to 3 ml from the catheter lumen and flush twice after guidewire removal.

11.0 Maintenance and Use *in situ*

The catheter should remain indwelling only as long as is required by the patient's condition.

Precaution: The incidence of complications increases significantly with indwelling periods longer than 72 hours.

11.1 Catheter Tip Position

Keep catheter tip centrally located in a main branch of the pulmonary artery near the hilum of the lungs. Do not advance tip too far peripherally. Tip should be kept where full or near full inflation volume is required to produce a wedge tracing. The tip migrates toward periphery during balloon inflation.

After deflation, the catheter tip may tend to recoil towards the pulmonic valve and may slip back into the right ventricle, requiring that the catheter be repositioned.

11.2 Catheter Tip Migration

Anticipate spontaneous catheter tip migration towards periphery of pulmonary bed. Continuously monitor distal lumen pressure to verify tip position. If wedge tracing is observed when balloon is deflated, pull catheter back. Damage may be caused by prolonged occlusion or over-distension of vessel upon re-inflation of the balloon.

Spontaneous catheter tip migration towards the periphery of the lung occurs during cardiopulmonary bypass. Partial catheter withdrawal (3 to 5 cm) just before bypass should be considered, as it may help reduce distal migration and prevent permanent catheter wedging postbypass. After termination of bypass, the catheter may require repositioning. Check the distal pulmonary artery tracing before inflating the balloon.

Precaution: Over a period of time, the catheter tip may migrate towards the periphery of the pulmonary bed and lodge in a small vessel. Damage may occur either by prolonged occlusion or by

over-distension of the vessel upon reinflation of the balloon (see Complications).

PA pressures should be continuously monitored with the alarm parameter set to detect physiologic changes as well as spontaneous wedge.

11.3 Balloon Inflation and Wedge Pressure Measurement

Re-inflation of the balloon should be performed gradually while monitoring pressures. Inflation is usually associated with a feeling of resistance. If no resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may still be used for hemodynamic monitoring, however, take precautions against infusion of air or liquids into the balloon lumen. During normal catheter use, keep inflation syringe attached to gate valve to prevent inadvertent injection of liquid into the balloon inflation lumen.

Measure wedge pressure only when necessary and only when tip is properly positioned (see above). Avoid prolonged maneuvers to obtain wedge pressure and keep wedge time to a minimum (two respiratory cycles or 10 - 15 seconds), especially in patients with pulmonary hypertension. If difficulties are encountered, discontinue wedge measurements. In some patients, pulmonary arterial end-diastolic pressure can often be substituted for pulmonary artery wedge pressure if the pressures are nearly identical, obviating the need for repeated balloon inflation.

11.4 Spontaneous Tip Wedging

The catheter may migrate into the distal pulmonary artery and spontaneous tip wedging may occur. To avoid this complication, pulmonary artery pressure should be monitored continuously with a pressure transducer and display monitor.

Forward advancement should never be forced if resistance is encountered.

11.5 Catheter Contamination Shield

Secure distal Tuohy-Borst adapter to the catheter.

Extend the proximal end of the catheter contamination shield to desired length, and secure proximal Tuohy-Borst adapter to the catheter.

11.6 Patency

All pressure monitoring lumens should be filled with a sterile, heparinized saline solution (e.g., 500 I.U. heparin in 500 ml saline) and flushed at least once each half hour or by continuous slow infusion. If loss of patency occurs and cannot be corrected by flushing, the catheter should be removed.

11.7 General

Keep pressure monitoring lumens patent by intermittent flush or continuous slow infusion with heparinized saline solution or use of a heparin lock using the provided injection caps with heparinized saline solution. Infusion of viscous solutions (e.g., whole blood or albumin) is not recommended, as they flow too slowly and may occlude the catheter lumen.

To use injection caps:

Step	Procedure
1	Disinfect injection caps before entry with syringe needle (see Complications).
2	Use a small bore needle (22 gauge (0.7 mm) or smaller) to puncture and inject through the injection caps.

WARNING: To avoid pulmonary artery rupture never flush the catheter when the balloon is wedged in the pulmonary artery.

Periodically check IV lines, pressure lines, and transducers to keep them free of air. Also ensure that connecting lines and stopcocks remain tightly fitted.

12.0 Cardiac Output Determination

To determine cardiac output by thermodilution, a known amount of sterile solution of known temperature is injected into the right atrium or vena cava, and the resultant change in blood temperature is measured in the pulmonary artery by the catheter thermistor. Cardiac output is inversely proportional to the integrated area under the resulting curve. This method has been shown to provide good correlation with the direct Fick method and dye dilution technique for cardiac output determination.

Consult the references on the use of iced versus room temperature injectate or open versus closed injectate delivery systems.

Refer to the appropriate cardiac output computer manual for specific instructions in the use of thermodilution catheters for cardiac output determination. Correction factors or computation constants needed to correct for indicator heat transfer are given in the specifications.

Edwards cardiac output computers require that a computation constant be used to correct for injectate temperature rise as it passes through the catheter. The computation constant is a function of injectate volume, temperature, and catheter dimensions. The computation constants listed in the specifications have been determined *in vitro*.

13.0 MRI Information

The Swan-Ganz device is MR unsafe as the result of the device containing metallic components, which experience RF-induced heating in the MRI environment; therefore the device poses hazards in all MRI environments.

14.0 Complications

Invasive procedures involve some patient risks. Although serious complications are relatively uncommon, the physician is advised, before deciding to insert or use the catheter, to consider the potential benefits in relation to the possible complications. The techniques for insertion, methods of using the catheter to obtain patient data information, and the occurrence of complications is well described in the literature.

The general risks and complications associated with indwelling catheters are described in the literature.

Strict adherence to these instructions and awareness of risks reduces the incidence of complications. Several known complications include:

14.1 Perforation of the Pulmonary Artery

Factors which are associated with the development of fatal pulmonary artery rupture include pulmonary hypertension, advanced age, cardiac surgery with hypothermia and anticoagulation, distal catheter tip migration, and arteriovenous fistula formation and other vascular traumas.

Extreme care should therefore be exercised during the measurement of pulmonary artery wedge pressure in patients with pulmonary hypertension.

In all patients, balloon inflation should be limited to two respiratory cycles, or 10 to 15 seconds.

A central location of the catheter tip near the hilum of the lung may prevent pulmonary artery perforation.

14.2 Pulmonary Infarction

Tip migration with spontaneous wedging, air embolism, and thromboembolism can lead to pulmonary artery infarction.

14.3 Cardiac Arrhythmias

Although usually transient and self-limited, arrhythmias may occur during insertion, withdrawal, and repositioning of the tip from the pulmonary artery into the right ventricle. Whereas premature ventricular contractions are the most commonly encountered arrhythmias, ventricular tachycardia and atrial and ventricular fibrillation have also been reported. ECG monitoring and the immediate availability of antiarrhythmic drugs and defibrillating equipment are recommended. Use of prophylactic lidocaine should be considered to decrease the incidence of ventricular arrhythmias during catheterization.

14.4 Knotting

Flexible catheters have been reported to form knots, most often as a result of looping in the right ventricle. Sometimes the knot can be resolved by insertion of a suitable guidewire and manipulation of the catheter under fluoroscopy. If the knot does not include any intracardiac structures, the knot may be gently tightened and the catheter withdrawn through the site of entry.

14.5 Sepsis/Infection

Positive catheter-tip cultures resulting from contamination and colonization have been reported, as well as incidences of septic and aseptic vegetation in the right heart. Increased risks of septicemia and bacteremia have been associated with blood sampling, the infusing of fluids, and catheter-related thrombosis. Preventive measures should be taken to guard against infection.

14.6 Other Complications

Other complications include right bundle branch block and complete heart block, tricuspid and pulmonic valve damage, pneumothorax, thrombosis, blood loss, cardiac structure/wall injury or damage, hematoma, embolism, anaphylaxis, and cardiac tissue/artery burn.

In addition, allergic reactions to latex have been reported. Physicians should identify latex sensitive patients and be prepared to treat allergic reactions promptly.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance for this medical device. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

14.7 Long Term Monitoring

The duration of catheterization should be the minimum required by the patient's clinical state since the risk of thromboembolic and infection complications increases with time. The incidence of complications increases

significantly with indwelling periods longer than 72 hours. Prophylactic systemic anticoagulation and antibiotic protection should be considered when long-term catheterization (i.e. over 48 hours) is required, as well as in cases involving increased risk of clotting or infection.

15.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

The packaging is designed to avoid crushing of the catheter and to protect the balloon from exposure to the atmosphere. It is therefore recommended that the catheter remain inside the package until use.

16.0 Storage

Store in a cool, dry place.

17.0 Operating Conditions/Use Environment

Intended to operate under physiological conditions of the human body in a controlled clinical environment.

18.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

19.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

20.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Computation Constants									
Model		131F7, 131F7P, 131VF7P, 141F7, 151F7	143TF7	831F75, 831F75P, and 831VF75P	834F75 and 834F75P	096F6 and 096F6P	132F5	TS105F5	C144F7, S144F7, C145F6N, C146F7, S9FC146F7, T173F6
Injectate Temp (°C)	Injectate Volume (ml)	Computation Constants (CC)***							
0-5	10	0.542	0.554	0.564	0.564	0.547	--	0.542	0.547
	5	0.247	0.259	0.262	0.257	0.259	0.270	0.255	0.252
	3	0.132	--	0.139	0.143	0.144	0.154	0.143	0.134
	1	--	--	--	--	--	0.037	--	--
19-22	10	0.578	0.587	0.588	0.582	0.582	--	--	0.579
	5	0.274	0.286	0.283	0.277	0.280	0.292	--	0.277
	3	0.154	--	0.158	0.156	0.161	0.170	--	0.156
	1	--	--	--	--	--	0.048	--	--
23-25	10	0.595	0.599	0.612	0.607	0.608	--	0.605	0.597
	5	0.287	0.291	0.301	0.294	0.305	0.307	0.297	0.295
	3	0.165	--	0.177	0.170	0.180	0.181	0.175	0.169
	1	--	--	--	--	--	0.055	--	--
Computation Constants for C0-Set+ closed injectate delivery system									
Cold Injectate									
6 °C - 12 °C	10	0.561	0.569	0.578	0.574	0.558	--	0.552	0.570
8 °C - 12 °C	5	0.259	0.266	0.272	0.287	0.277	0.285	0.265	0.271
Room Temperature									
18 °C - 25 °C	10	0.608	0.589	0.592	0.595	0.607	--	0.589	0.585
18 °C - 25 °C	5	0.301	0.287	0.290	0.298	0.301	0.307	0.289	0.287
***CC = (1.08)C _T (60)(V _I)									

Specifications

Function	Thermodilution catheter				Thermodilution and Infusion catheter		True Size catheter			Polymer Blend True Size catheter			
Model Number	131F7, 131F7P, 131VF7P (Standard)	141F7 (Hi-Shore catheter)	143TF7 (CardioCath catheter)	151F7 ("S" Tip)	831F75, 831F75P, 831VF75P (VIP catheter)	834F75, 834F75P (VIP+ catheter)	096F6, 096F6P	132F5**	TS105F5	C144F7, S144F7	C145F6N	C146F7, S9FC146F7	T173F6
Usable Length (cm)	110	110	110	110	110	110	110	75	100	110	110	110	100
Catheter Body French Size	7F (2.3 mm)	7F (2.3 mm)	7F (2.3 mm)	7F (2.3 mm)	7.5F (2.5 mm)	7.5F (2.5 mm)	6F (2.0 mm)	5F (1.7 mm)	5F (1.7 mm)	7F (2.3 mm)	6F (2.0 mm)	7F (2.3 mm)	6F (2.0 mm)
Body Color	-	-	-	-	-	-	-	-	-	White	White	White	White
Minimum Recommended Introducer Size	8F (2.7 mm)	8F (2.7 mm)	8F (2.7 mm)	8F (2.7 mm)	8.5F (2.8 mm)	8.5F (2.8 mm)	7F (2.3 mm)	6F (2.0 mm)	6F (2.0 mm)	8F (2.7 mm)	7F (2.3 mm)	8.5F (2.8 mm)	7F (2.3 mm)
Diameter of Inflated Balloon (mm)	13	13	13	13	13	13	11	8	8	13	11	13	11
Balloon Inflation Capacity (ml)	1.5	1.5	1.5	1.5	1.5	1.5	1.0	0.7	0.8	1.5	1.0	1.5	1.0
Distance from Tip (cm)													
Injectate Port	30	30	30	30	30	30	30	15	24	30	30	30	23
Infusion VIP Port	-	-	-	-	31	31	-	-	-	-	-	-	-
RV Infusion VIP+ Port	-	-	-	-	-	19	-	-	-	-	-	-	-
Thermistor	4	4	4	4	4	4	4	1.5	-	4	4	4	4
Lumen Volume (ml)													
Distal Lumen	1.02	1.01	0.89	1.02	0.86	0.89	0.92	0.64	0.67	1.01	0.94	1.01	0.90
Injectate Lumen	0.81	0.81	0.73	0.81	0.86	0.75	0.74	0.57	0.58	0.76	0.74	0.76	0.74
VIP/RA Lumen	-	-	-	-	0.87	0.97	-	-	-	-	-	-	-
RV Lumen	-	-	-	-	-	0.93	-	-	-	-	-	-	-
Infusion Rate* (ml/min)													
Infusion VIP Lumen	-	-	-	-	13	14	5 (distal)	-	-	-	-	-	-
RV Infusion VIP+ Lumen	-	-	-	-	-	7	7 (injectate)	-	-	-	-	-	-
Distal Lumen	-	-	-	-	-	5	-	-	1.7	7	6	7	6
Proximal Injectate Lumen	-	-	-	-	-	9	-	-	2.6	10	7	10	8
Maximum Guidewire Diameter													
Distal Lumen	0.025 in (0.64 mm)	0.025 in (0.64 mm)	0.025 in (0.64 mm)	0.025 in (0.64 mm)	0.025 in (0.64 mm)	0.025 in (0.64 mm)	0.025 in (0.64 mm)	0.020 in (0.51 mm)	0.018 in (0.46 mm)	0.028 in (0.71 mm)	0.025 in (0.64 mm)	0.028 in (0.71 mm)	0.025 in (0.64 mm)
Frequency Response Distortion at 10 Hz													
Distal Lumen	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB	< 3 dB

All specifications given are nominal values. A syringe is provided with each catheter.

* Using normal saline at room temperature, 1 m above insertion site, gravity drip; rates represent average values.

**In situations where the insertion site or patient physiology requires a larger insertion distance, a longer catheter model or larger introducer size should be selected.

†A "C" in the model number denotes the "C-Tip" configuration; an "S" in the model number denotes the "S-Tip" configuration; a "T" in the model number denotes the "T-Tip" configuration.

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

Figure 1: Swan-Ganz ControlCath Thermodilution Catheter or Swan-Ganz Torque Support Thermodilution Catheter

- Figure 1 : Cathéter de thermodilution ControlCath Swan-Ganz ou cathéter de thermodilution à contrôle de torque amélioré Swan-Ganz
- Abbildung 1: Swan-Ganz ControlCath Thermodilutionskatheter oder Swan-Ganz Thermodilutionskatheter mit Drehmomentkontrolle
- Figura 1: Cateter de termodilución Swan-Ganz ControlCath o Cateter de termodilución con control del par de torsión Swan-Ganz
- Figura 1: Catetere per termodiluzione ControlCath Swan-Ganz o catetere per termodiluzione Swan-Ganz con supporto di torsione
- Afbeelding 1: Swan-Ganz ControlCath -thermodilutiekatheter of Swan-Ganz -thermodilutiekatheter met draaicontrol
- Figur 1: Swan-Ganz ControlCath termodilutionskateter eller Swan-Ganz termodilutionskateter med torsjonsstøtte
- Figur 1: Swan-Ganz ControlCath termodilutionskateter eller Swan-Ganz termodilutionskateter med vridstöd
- Εικόνα 1: Καθετήρας θερμοαραίωσης Swan-Ganz ControlCath ή καθετήρας θερμοαραίωσης Swan-Ganz με υποστήριξη ροπής στρέψης
- Figura 1: Cateter de termodiluição Swan-Ganz ControlCath ou cateter de termodiluição Swan-Ganz com suporte de torque
- Obrázek 1: Termodiluční katétr Swan-Ganz ControlCath nebo termodiluční katétr Swan-Ganz s podporou točivého momentu
- 1. ábra: Swan-Ganz ControlCath termodilúciós katéter vagy Swan-Ganz elforgatáskönnyített termodilúciós katéter
- Rysunek 1: Cewnik do termodylucji Swan-Ganz ControlCath lub cewnik do termodylucji Swan-Ganz ze wspornikiem momentu obrotowego
- Obrázok 1: Termodilučný katéter Swan-Ganz ControlCath alebo termodilučný katéter Swan-Ganz s podporou krútenia
- Figur 1: Swan-Ganz ControlCath -termodilusionskateter eller Swan-Ganz -termodilusionskateter med torsjonsstøtte
- Kuva 1: Swan-Ganz ControlCath -termodilutiokatetri tai taivutettava Swan-Ganz -termodilutiokatetri
- Фигура 1: Катетър Swan-Ganz ControlCath за термодилуция или катетър Swan-Ganz за термодилуция с контрол на въртящия момент
- Figura 1: Cateter de termodilufiie Swan-Ganz ControlCath sau cateter de termodilufiie cu suport pentru torsiune Swan-Ganz
- Joonis 1. Termidilutsiooni kateeter Swan-Ganz ControlCath või pööramistoega termidilutsiooni kateeter Swan-Ganz
- 1 pav. „Swan-Ganz ControlCath“ termidilucijos kateteriai arba „Swan-Ganz“ sukima palaikantys termidilucijos kateteriai
- 1. attēls. Swan-Ganz ControlCath termidilucijas katetrs vai Swan-Ganz termidilucijas katetrs ar vērpes atbalstu
- Şekil 1: Swan-Ganz ControlCath Termidilüsyon Kateteri veya Swan-Ganz Tork Destekli Termidilüsyon Kateteri
- Рис. 1. Термодилуционный катетер Swan-Ganz ControlCath или термодилуционный катетер Swan-Ganz с передачей вращения
- Slika 1: Swan-Ganz ControlCath termidilucioni kateter ili Swan-Ganz termidilucioni kateter sa torzionom podrškom
- Slika 1: Kateter za termodiluciju Swan-Ganz ControlCath ili kateter za termodiluciju s olakšanim zakretanjem Swan-Ganz

1. Proximal Injectate Port @ 23 cm ■ Orifice proximal d'injection à 23 cm ■ Proximale Injektatöffnung bei 23 cm ■ Puerto del injectable proximal a 23 cm ■ Porta di iniezione prossimale a 23 cm ■ Proximale injectaatpoort @ 23 cm ■ Proksimal injektatport @ 23 cm ■ Proximal injektatport vid 23 cm ■ Εγγύς θύρα εγχέομενου διαλύματος στα 23 cm ■ Porta proximal do injetado @ 23 cm ■ Proximální port pro vstřikovaný roztok ve 23 cm ■ Proximális injektátumnyílás 23 cm -nél ■ Proksymalny port do iniekcji w odległości 23 cm ■ Proximálny port na injektát pri 23 cm ■ Proksimal injektatport ved 23 cm ■ Proksimaalinen injektioportti 23 cm:n kohdalla ■ Проксимален порт за инжектирано вещество при 23 cm ■ Port proximal de injectare la 23 cm ■ Proksimaalne süstelahuse ava 23 cm juures ■ Proksimalinė leidžiama tirpalo anga ties 23 cm ■ Proksimālā injektāta pieslēgvieta pie 23 cm atzīmes ■ Proksimal Enjektat Portu (23 cm uzaklıkta) ■ Проксимальный порт для введения раствора на расстоянии 23 cm ■ Proksimalni otvor za injektat na 23 cm ■ Proksimalni priključak za injektat @ 23 cm
2. Thermistor ■ Thermistance ■ Thermistor ■ Termistor ■ Termistore ■ Thermistor ■ Termistor ■ Termistor ■ Θερμική αντίσταση ■ Termistor ■ Termistor ■ Termisztor ■ Termistor ■ Termistor ■ Termistor ■ Termistori ■ Термистор ■ Termistor ■ Termistor ■ Termistorius ■ Termistors ■ Termistör ■ Термистор ■ Termistor ■ Termistor
3. Balloon ■ Ballonnet ■ Ballon ■ Balón ■ Palloncino ■ Ballon ■ Ballon ■ Ballong ■ Μπαλόνι ■ Balão ■ Balónek ■ Ballon ■ Balon ■ Balónik ■ Ballong ■ Pallo ■ Балон ■ Balon ■ Balloon ■ Balionelis ■ Balons ■ Balon ■ Баллон ■ Balon ■ Balon
4. Distal Lumen Port ■ Orifice de la lumière distale ■ Distale Lumenöffnung ■ Puerto de la luz distal ■ Porta del lume distale ■ Distale lumenpoort ■ Distal lumenport ■ Distal lumenport ■ Θύρα περιφερικού αυλού ■ Porta do lúmen distal ■ Port distálného lumina ■ Disztális lumen nyílása ■ Port kanálu dystalnego ■ Port distálneho lúmenu ■ Distal lumenport ■ Distaalisen luumenin portti ■ Порт на дистален лумен ■ Port lumen distal ■ Distaalse valendiku port ■ Distalinio spindžio anga ■ Distālā lūmena pieslēgvieta ■ Distal Lūmen Portu ■ Порт дистального просвета ■ Otvor distalnog lumena ■ Priključak distalnog lumena
5. Proximal Injectate Port @ 23 cm ■ Orifice proximal d'injection à 23 cm ■ Proximale Injektatöffnung bei 23 cm ■ Puerto del injectable proximal a 23 cm ■ Porta di iniezione prossimale a 23 cm ■ Proximale injectaatpoort @ 23 cm ■ Proksimal injektatport @ 23 cm ■ Proximal injektatport vid 23 cm ■ Εγγύς θύρα εγχέομενου διαλύματος στα 23 cm ■ Porta proximal do injetado @ 23 cm ■ Proximální port pro vstřikovaný roztok ve 23 cm ■ Proximális injektátumnyílás 23 cm -nél ■ Proksymalny port do iniekcji w odległości 23 cm ■ Proximálny port na injektát pri 23 cm ■ Proksimal injektatport ved 23 cm ■ Proksimaalinen injektioportti 23 cm:n kohdalla ■ Проксимален порт за инжектирано вещество при 23 cm ■ Port proximal de injectare la 23 cm ■ Proksimaalne süstelahuse ava 23 cm juures ■ Proksimalinė leidžiama tirpalo anga ties 23 cm ■ Proksimālā injektāta pieslēgvieta pie 23 cm atzīmes ■ Proksimal Enjektat Portu (23 cm uzaklıkta) ■ Проксимальный порт для введения раствора на расстоянии 23 cm ■ Proksimalni otvor za injektat na 23 cm ■ Proksimalni priključak za injektat @ 23 cm
6. Thermistor ■ Thermistance ■ Thermistor ■ Termistor ■ Termistore ■ Thermistor ■ Termistor ■ Termistor ■ Θερμική αντίσταση ■ Termistor ■ Termistor ■ Termisztor ■ Termistor ■ Termistor ■ Termistor ■ Termistori ■ Термистор ■ Termistor ■ Termistor ■ Termistorius ■ Termistors ■ Termistör ■ Термистор ■ Termistor ■ Termistor
7. Balloon ■ Ballonnet ■ Ballon ■ Balón ■ Palloncino ■ Ballon ■ Ballon ■ Ballong ■ Μπαλόνι ■ Balão ■ Balónek ■ Ballon ■ Balon ■ Balónik ■ Ballong ■ Pallo ■ Балон ■ Balon ■ Balloon ■ Balionelis ■ Balons ■ Balon ■ Баллон ■ Balon ■ Balon
8. Distal Lumen Port ■ Orifice de la lumière distale ■ Distale Lumenöffnung ■ Puerto de la luz distal ■ Porta del lume distale ■ Distale lumenpoort ■ Distal lumenport ■ Distal lumenport ■ Θύρα περιφερικού αυλού ■ Porta do lúmen distal ■ Port distálného lumina ■ Disztális lumen nyílása ■ Port kanálu dystalnego ■ Port distálneho lúmenu ■ Distal lumenport ■ Distaalisen luumenin portti ■ Порт на дистален лумен ■ Port lumen distal ■ Distaalse valendiku port ■ Distalinio spindžio anga ■ Distālā lūmena pieslēgvieta ■ Distal Lūmen Portu ■ Порт дистального просвета ■ Otvor distalnog lumena ■ Priključak distalnog lumena
9. Proximal Injectate Port @ 23 cm ■ Orifice proximal d'injection à 23 cm ■ Proximale Injektatöffnung bei 23 cm ■ Puerto del injectable proximal a 23 cm ■ Porta di iniezione prossimale a 23 cm ■ Proximale injectaatpoort @ 23 cm ■ Proksimal injektatport @ 23 cm ■ Proximal injektatport vid 23 cm ■ Εγγύς θύρα εγχέομενου διαλύματος στα 23 cm ■ Porta proximal do injetado @ 23 cm ■ Proximální port pro vstřikovaný roztok ve 23 cm ■ Proximális injektátumnyílás 23 cm -nél ■ Proksymalny port do iniekcji w odległości 23 cm ■ Proximálny port na injektát pri 23 cm ■ Proksimal injektatport ved 23 cm ■ Proksimaalinen injektioportti 23 cm:n kohdalla ■ Проксимален порт за инжектирано вещество при 23 cm ■ Port proximal de injectare la 23 cm ■ Proksimaalne süstelahuse ava 23 cm juures ■ Proksimalinė leidžiama tirpalo anga ties 23 cm ■ Proksimālā injektāta pieslēgvieta pie 23 cm atzīmes ■ Proksimal Enjektat Portu (23 cm uzaklıkta) ■ Проксимальный порт для введения раствора на расстоянии 23 cm ■ Proksimalni otvor za injektat na 23 cm ■ Proksimalni priključak za injektat @ 23 cm
10. Thermistor ■ Thermistance ■ Thermistor ■ Termistor ■ Termistore ■ Thermistor ■ Termistor ■ Termistor ■ Θερμική αντίσταση ■ Termistor ■ Termistor ■ Termisztor ■ Termistor ■ Termistor ■ Termistor ■ Termistori ■ Термистор ■ Termistor ■ Termistor ■ Termistorius ■ Termistors ■ Termistör ■ Термистор ■ Termistor ■ Termistor
11. Distal Lumen Port ■ Orifice de la lumière distale ■ Distale Lumenöffnung ■ Puerto de la luz distal ■ Porta del lume distale ■ Distale lumenpoort ■ Distal lumenport ■ Distal lumenport ■ Θύρα περιφερικού αυλού ■ Porta do lúmen distal ■ Port distálného lumina ■ Disztális lumen nyílása ■ Port kanálu dystalnego ■ Port distálneho lúmenu ■ Distal lumenport ■ Distaalisen luumenin portti ■ Порт на дистален лумен ■ Port lumen distal ■ Distaalse valendiku port ■ Distalinio spindžio anga ■ Distālā lūmena pieslēgvieta ■ Distal Lūmen Portu ■ Порт дистального просвета ■ Otvor distalnog lumena ■ Priključak distalnog lumena
12. Balloon ■ Ballonnet ■ Ballon ■ Balón ■ Palloncino ■ Ballon ■ Ballon ■ Ballong ■ Μπαλόνι ■ Balão ■ Balónek ■ Ballon ■ Balon ■ Balónik ■ Ballong ■ Pallo ■ Балон ■ Balon ■ Balloon ■ Balionelis ■ Balons ■ Balon ■ Баллон ■ Balon ■ Balon

Figure 2: Swan-Ganz ControlCath Thermodilution Catheter tip configurations ■ Figure 2 : Configurations des extrémités de cathéter de thermodilution ControlCath Swan-Ganz ■ Abbildung 2: Spitzenvarianten beim Swan-Ganz ControlCath Thermodilutionskatheter ■ Figura 2: Configuración de la punta del catéter de termodilución ControlCath Swan-Ganz ■ Figura 2: Configurazioni della punta del catetere per termodiluzione ControlCath Swan-Ganz ■ Afbeelding 2: Tipconfiguraties Swan-Ganz ControlCath -thermodilutiekatheter ■ Figur 2: Konfiguration af Swan-Ganz ControlCath thermodilutionskaterspids ■ Figur 2: Konfigurationer av Swan-Ganz ControlCath thermodilutionskaterspets ■ Εικόνα 2: Διαμορφώσεις άκρου καθετήρα θερμοδiluσιώνς Swan-Ganz ControlCath ■ Figura 2: Configurações das pontas do cateter de termodiluição Swan-Ganz ControlCath ■ Obrázek 2: Termodiluční katétr Swan-Ganz ControlCath – konfigurace hrotu ■ 2. ábra: A Swan-Ganz ControlCath termodilúciós katétervég konfigurációi ■ Rysunek 2: Konfiguracje końcówki cewnika do termodylucji Swan-Ganz ControlCath ■ Obrázok 2: Konfigurácie hrotov termodilučných katétrov Swan-Ganz ControlCath ■ Figur 2: Swan-Ganz ControlCath termodilusionskaterspiss-konfigurasjoner ■ Kuva 2: Swan-Ganz ControlCath -termodilutiokatetrin T-kärvivaihtoehdot ■ Фигура 2: Конфигурации на върх на катетър Swan-Ganz ControlCath за термодилуция ■ Figura 2: Configurații pentru vârful cateterului de termodilufie Swan-Ganz ControlCath ■ Joonis 2. Swan-Ganz ControlCath termodilutsiooni kateetri otsa konfiguratsioonid ■ 2 pav. „Swan-Ganz ControlCath“ termodilucijos kateterių galiukų konfigūracijos ■ 2. attēls. Swan-Ganz ControlCath termodilucijas katetra gala konfigurācijas ■ Šekil 2: Swan-Ganz ControlCath Termodilüsyon Kateteri uç yapılandırılmaları ■ Рис. 2. Конфигурации кончиков термодилуционного катетера Swan-Ganz ControlCath ■ Slika 2: Konfiguracije vrha Swan-Ganz ControlCath termodilucionog katetera ■ Slika 2: Konfiguracije vrha katetera za termodiluciju Swan-Ganz ControlCath

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	English	Français	Deutsch	Español	Italiano	Nederlands
	Number of lumens	Nombre de lumières	Anzahl der Lumen	Número de luces	Numero di lumi	Aantal lumina
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad
	Lumen size	Taille de la lumière	Lumendurchmesser	Tamaño de la luz	Dimensioni del lume	Lumengrootte
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello	Modelnummer
	Minimum introducer size	Taille minimale de l'introducteur	Mindestgröße der Einführhilfe	Tamaño de introductor mínimo	Misura minima dell'introdotore	Minimum introduceergrootte
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum
	Inner diameter	Diamètre interne	Innendurchmesser	Diámetro interior	Diametro interno	Binnendiameter
	Balloon capacity	Capacité du ballonnet	Ballonkapazität	Capacidad del balón	Capacità del palloncino	Balloncapaciteit
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto	Droog houden
	Sterile	Stérile	Steril	Estéril	Sterile	Steriel
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni	Gesteriliseerd met behulp van straling
	Sterilized using steam or dry heat	Stérilisé à la vapeur d'eau ou à la chaleur sèche	Mit Wasserdampf oder trockener Hitze sterilisiert	Esterilizado con vapor o calor seco	Sterilizzato con vapore o calore secco	Gesteriliseerd met behulp van stoom of droge warmte
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig
	Not made with natural rubber latex	Non fabriqué à partir de latex de caoutchouc naturel	Nicht aus Latex hergestellt	No contiene látex de caucho natural.	Non realizzato in lattice di gomma naturale	Bevat geen natuurlijk rubberlatex
	Size	Taille	Größe	Tamaño	Misura	Afmetingen
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Consultare le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website

	English	Français	Deutsch	Español	Italiano	Nederlands
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM	Non compatibile con RM	MRI-onveilig
	MR Safe	Aucun risque en milieu RM	MR-sicher	Seguro para RM	Compatibile con RM	MRI-veilig
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata	MRI-veilig onder bepaalde voorwaarden
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Consultare le istruzioni per l'uso	Volg de gebruiksaanwijzing
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare	Niet opnieuw steriliseren
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen
	Type B applied part	Pièce appliquée de type B	Anwendungsteil vom Typ B	Pieza aplicada tipo B	Parte applicata di tipo B	Toegepast onderdeel van type B
	Type CF applied part	Pièce appliquée de type CF	Anwendungsteil vom Typ CF	Pieza aplicada tipo CF	Parte applicata di tipo CF	Toegepast onderdeel van type CF
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing
	Open	Ouvrir	Offen	Abrir	Aperto	Open
	Aspirate Balloon -0.5 cc Before Introduction or Withdrawal.	Aspirer -0,5 cm ³ dans le ballonnet avant introduction ou retrait.	Vor der Einführung oder dem Entfernen mit der Spritze ein Vakuum (-0,5 ml) im Ballon erzeugen.	Aspire el balón -0,5 cm ³ antes de la introducción o la extracción.	Aspirare nel palloncino -0,5 cc prima dell'introduzione o del ritiro.	Zuig 0,5 cc uit de ballon alvorens hem in te brengen of te verwijderen.
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/ l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)	Conformité Européenne (CE-markering)
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel
	Importer	Importateur	Importeur	Importador	Importatore	Importeur
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren
	Double sterile barrier system	Système de barrière stérile double	Doppeltes Sterilbarriersystem	Un sistema doble de barrera estéril	Sistema a doppia barriera sterile	Dubbel steriel barrièresysteem
	Fragile, handle with care	Fragile, manipuler avec précaution	Vorsicht! Zerbrechlich!	Frágil; manipular con cuidado	Fragile, maneggiare con cura	Breekbaar, voorzichtigheid geboden
	This way up	Flèches indiquant le haut du carton	Oben	Este lado arriba	Alto	Deze kant boven
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo	Unieke instrumentidentificatie

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen.

**Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων
■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat**

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	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Valendike arv	Spindžių skaičius	Lūmenu skaits	Lūmen sayısı	Количество просветов	Broj lumena	Broj lumena
	Vālislābimōot	Išorinis skersmuo	Ārējais diametrs	Dış çap	Внешний диаметр	Spoljašnji prečnik	Vanjski promjer
	Kasutata pikkus	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebljiva dužina	Upotrebljiva dužina
	Juhtetraadi soovitusslik suurus	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs	Kılavuz Tel İçin Önerilen Boyut	Рекомендованный размер проводника	Preporučena veličina žice vodice	Preporučena veličina žice vodilice
	Valendiku suurus	Spindžio dydis	Lūmena izmērs	Lūmen boyutu	Размер просвета	Veličina lumena	Veličina lumena
	Mudeli number	Modelio numeris	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Sisesti väikeim suurus	Minimalus įvediklio dydis	Minimālais ievadītāja izmērs	Minimum introdüser boyutu	Минимальный размер интродюсера	Minimalna veličina uvodnika	Minimalna veličina uvodnog instrumenta
	Ettevaatust	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не использовать повторно	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina	Količina
	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Kõlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sisemine läbimōot	Vidinis skersmuo	Iekšējais diametrs	İç çap	Внутренний диаметр	Unutrašnji prečnik	Unutrašnji promjer
	Balloon mahutavus	Balionių talpa	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Hoida kuivalt	Laikykite sausoje vietoje	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
	Steriilne	Sterilus	Sterils	Steril	Стерильно	Sterilno	Sterilno
	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Steriliseeritud kiirgust kasutades	Sterilizuota švitinant	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Steriliseeritud aur või kuiva kuumusega	Sterilizuota naudojant garus arba sausą karštį	Sterilizēts, izmantojot tvaiku vai sausu karstumu	Buhar veya kuru ısı kullanılarak sterilize edilmiştir	Стерилизовано паром или сухим жаром	Sterilizovano korišćenjem pare ili suvom toplotom	Sterilizirano parom ili vrućim, suhim zrakom
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Tootmiskupäev	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Sisaldab naturaalseid kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Pole valmistatud looduslikust kummilateksist.	Pagaminta be natūralios gumos latekso	Izgatavošanā nav izmantots dabiskā kaučuka latekss	Doğal kauçuk lateksten yapılmamıştır	Не содержит натурального каучукового латекса	Ne sadrži lateks od prirodne gume	Nije izrađeno od prirodnog gumenog lateksa
	Suurus	Dydis	Izmērs	Boyut	Размер	Veličina	Veličina
	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas interneto svetainėje	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na internet stranici	Pogledajte upute za upotrebu na internetskoj stranici
	Järgige kasutusjuhendit, mille leiata veebisaidilt	Vadovaukites naudojimo instrukcijomis, pateiktomis svetainėje	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici

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Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

01/24
10007450008 A / DOC-0210899 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone	949.250.2500
	800.424.3278
FAX	949.250.2525

Web IFU