Edwards Lifesciences

Directory

English (en)	1	Magyar (hu)	45	Română (ro)	69
Français (fr)	5	Polski (pl)	49	Eesti (et)	73
Deutsch (de)	9	Slovensky (sk)	53	Lietuvių (lt)	77
Español (es)	13	Norsk (no)	57	Latviešu (lv)	81
Italiano (it)	17	Suomi (fi)	61	Türkçe (tr)	85
Nederlands (nl)	21	Български (bg)	65	Русский (ru)	89
Dansk (da)	25			Srpski (sr)	93
Svenska (sv)	29			Hrvatski (hr)	97
Ελληνικά (el)	33				
Português (pt)	37				
Česky (cs)	41				
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Φιγυρι ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike					101
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola					115

English

Model 93600 CO-Set+ Closed Injectate Delivery System for Cold Injectate

Model 93610 CO-Set+ Closed Injectate Delivery System for Room Temperature Injectate

Model 93520 CO-Set+ Cooling Container


Model 93505 Flow-Through Housing

Model 93650 CO-Set+ Syringe

Model 93522 CO-Set+ In-Line Injectate Temperature Probe

Model 9850A Injectate Temperature Probe

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Model 93600 CO-Set+ Closed Injectate Delivery System, Cold Injectate

- | | |
|---|--|
| 1. Swan-Ganz Thermodilution Catheter | 10. Snap Clamp |
| 2. Catheter Connecting Cable | 11. Cooling Container |
| 3. Patient Monitor | 12. Tubing Connector |
| 4. Sterile Injectate Solution (User Supplied) | 13. Three Way Stopcock |
| 5. Non-Vented IV Spike | 14. Injectate Temperature Probe |
| 6. Tubing Coil | 15. Flow-Through Check Valve |
| 7. Lid | 16. 10 ml syringe with jacket and contamination shield |
| 8. Thermistor Connector | 17. Cooling Container Bracket |
| 9. Injectate Probe Connector | |

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For single use only

1.0 Description

CO-Set+ is a closed injectate delivery system. It is an accessory to be used with Swan-Ganz catheters that intermittently measure cardiac output (CO). It allows for delivery of injectate for the measurement of Cardiac Output by thermodilution (TD). The unit is designed for use with a patient monitor and a Swan-Ganz thermodilution catheter. These devices are intended to be used by licensed physicians, surgeons, including cardiac, and vascular surgeons.

CO-Set+ closed injectate delivery system for cold injection (Figure 1 on page 101) and for room temperature injectate (Figure 2 on page 103) is connected between a suitable supply of sterile injectate (D5W) and a Swan-Ganz thermodilution catheter.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The measurement of cardiac outputs may improve data driven clinical decision making for medically necessary intervention and/or clinical re-evaluation.

2.0 Intended Use/Purpose

The CO-Set+ closed injectate delivery system is intended to be a pathway to deliver thermodilution injectate (used to make cardiac output determinations) into a thermodilution catheter.

3.0 Indications for Use

The CO-Set+ closed injectate delivery systems used for cold injectate (Model 93600) or room temperature injectate (Model 93610), provide simple and convenient delivery of cold or room temperature injectate for the measurement of cardiac output by thermodilution. The devices are intended for use in adult critically ill patients.

4.0 Contraindications

There are no known contraindications.

5.0 CO-Set+ closed injectate delivery system for cold injectate (Model 93600)

After the system is primed and the cooling coil is immersed in an ice bath, the injectate will reach operating temperature within 5 minutes.

The covered ice bath will remain usable for up to 6 hours; for best results, the ice bath should be checked periodically to ensure that ice is still present.

6.0 Flow-Through Housing (Model 93505)

When used in conjunction with the Injectate Temperature Probe, the Flow-Through Housing (Model 93505) allows direct temperature measurement of the thermal indicator at the proximal lumen infusion site. Measurement of the injectate temperature at the infusion site helps reduce errors frequently encountered in standard injection techniques (for example, injectate warming caused by handling the syringe). The temperature of the injectate is fed directly by the injectate temperature probe into a patient monitor for curve analysis and determination of cardiac output. An appropriate computation constant must be entered into the computer to determine the correct cardiac output. A 3-way stopcock should be placed proximal to the flow-through housing (Figure 11 on page 112) for standard catheter maintenance and fluid withdrawal.

7.0 CO-Set+ Syringe (Model 93650)

The CO-Set+ syringe is a 10 ml syringe for use with the Model 93600 (Cold) or Model 93610 (Room Temperature) closed injectate delivery systems. This individually packaged and sterilized syringe is intended to replace the originally supplied syringe if the need arises.

Additionally, the CO-Set+ syringe may be used with the Model 93505 Flow-Through Housing in place of the user supplied 10 ml syringe described in its Instruction for Use.

8.0 CO-Set+ In-Line Injectate Temperature Probe (Model 93522)

The model 93522 CO-Set+ in-line injectate temperature probe is used in conjunction with injectate delivery systems: model 93600 (Cold) CO-Set+ or model 93610 (Room Temperature) CO-Set+ or model 93505 Flow-Through Housing. The CO-Set+ in-line injectate temperature probe is a three element network consisting of two resistors and one thermistor configured to produce a linear resistance versus temperature curve. The nominal measurement accuracy over the range of 0 ° to 25 °C is ± 0.3 °C. The CO-Set+ in-line injectate temperature probe is stationed at the point where the syringe injects the thermal indicator into the catheter's proximal lumen hub. The probe is able to precisely measure the injectate temperature as it is being introduced into the catheter.

9.0 Injectate Temperature Probe (Model 9850A)

The Model 9850A Injectate Temperature Probe is a dual thermistor composite that produces a linearized response to temperature change. The nominal measurement accuracy over the range of 0 ° to 25 °C is ± 0.3 °C. The probe's thermistors are electrically insulated in a tubular stainless steel probe specifically designed for liquid immersion. The probe should not be immersed in liquid past the point where the stainless steel tube attaches to the spring-cord assembly.

The injectate temperature probe attaches to the properly marked branch of the Edwards connecting cable for the

Edwards Vigilance monitors. The model 9850A Injectate Temperature Probe is an in-line bath probe used for bolus cardiac output. The probe is 3 1/2 in (8.9 cm) in length and has a 10 ft (3 m) long cable.

10.0 Instructions for Use

Equipment Needed for Thermodilution

- Swan-Ganz thermodilution catheter
- Flow-Through Injectate Temperature Probe
- Sterile D5W in IV container (bag recommended) for use as injectate (500 ml minimum)
- Patient Monitor
- Connecting cable for patient monitor
- Empty flush container for collecting discarded injectate
- Continuous flush device

Equipment Needed for cold injectate:

- CO-Set+ closed injectate delivery system for cold injectate (Model 93600)
- Cooling Container (Model 93520 or equivalent)
- Crushed ice
- Iced water
- Cooling Container IV Pole Mounting Bracket (Model 93521) (optional)

Equipment Needed for room temperature injectate:

- CO-Set+ closed injectate delivery system for room injectate (Model 93610)

10.1 Connecting Injectate Container to Closed System (Figure 4)

Use aseptic technique.

Step	Procedure
1	Suspend injectate container (IV bag or bottle) at the desired height; do not place it higher than three feet (90 cm) above the catheter insertion site.
2	Remove CO-Set+ closed injectate delivery system for cold injection from its package.
3	Close the snap clamp.
4	Remove the protective cap from the IV spike and insert the nonvented spike into the injectate container. Note: If a non-vented glass bottle is used, an airway needle or vented extension set is required to vent the bottle.
5	For cold injectate (Model 93600) only: Pull apart the TUBING COIL (not cooling coil) to obtain the desired tubing length between the injectate container and cooling container. WARNING: Do not pull apart the Cooling Coil. Injectate may not cool to appropriate temperature if Cooling Coil is uncoiled. If additional tubing length is needed, pull apart the tubing coil.

10.2 Connecting Injectate Container to Closed System (Figure 5)

Use Aseptic Technique

Step	Procedure
1	Suspend injectate container (IV bag or bottle) at the desired height; do not place it higher than 3 feet (90 cm) above the catheter insertion site.

Step	Procedure
2	Remove a CO-Set+ system from its package.
3	Close the snap clamp.
4	Remove the protective cap from the IV spike and insert the non-vented spike into the injectate container.

Note: If a non-vented glass bottle is used, an airway needle or vented extension set is required to vent the bottle.

10.3 Flushing Air from Fluid Path

Use aseptic technique.

Step	Procedure
1	Connect the 10 ml syringe to the flow-through/check valve (Figure 6 on page 108). CAUTION: Do not connect the syringe to the catheter side of the flow-through/check valve. Damage to the flow-through/check valve can occur if the flow-through/check valve is flushed or pressurized from the wrong side. Disabling back-pressurization of the flow-through/check valve can also occur if the 3-way stopcock between the flow-through/check valve and the catheter is not opened to the catheter before injecting. If the flow-through/check valve becomes damaged, then it is possible that blood could be aspirated back into the CO-Set+ system syringe. Should this occur, use of the CO-Set+ system (syringe and flow-through/check valve) should be discontinued and the system replaced.
2	Ensure that all the components are securely connected.
3	Open the snap clamp to allow solution to flow from the injectate container. To collect the discarded solution, hold the flow-through/check valve over an empty container.
4	Slowly pull the syringe plunger to prime the system; then push the plunger in. Repeat 5 or 6 times, or until the system is completely free of air. The system may also be filled by squeezing the IV bag or plastic bottle. Remove the syringe from the flow-through/check valve when filling the system by this method.
5	Return the syringe plunger to the fully depressed position.
6	Close the snap clamp.

10.4 Preparing Cooling Container (Model 93520) (See Fig. 7)

Step	Procedure
1	Pack the container with crushed ice up to the inside ridge in the internal ribs.
2	Seat the cooling coil on the first ridge. For maximum cooling, position the coil with the outlet at the bottom as shown.

Step	Procedure
3	Place the inlet and outlet tubing in the slots provided at the top of the container. Be sure that the white plastic tubing connector is inside the container as shown.
4	Completely cover the cooling coil with crushed ice.
5	Add iced water until it is visible above the ice level.
6	Place any unused length of tubing coil into the ice container to provide additional cold solution.
7	Lock the lid on the cooling container.
8	Let sit for 5 minutes.
9	Flush system once to ensure cold fluid is used for injection (Figure 8 on page 110).

Note: Use the spigot to drain excess water before refilling with more ice. There is no need to disassemble the system.

Note: An optional bracket (Model 93521) is available to mount the container on an IV pole with a diameter between 1/2 and 1 1/2 inches (1.27 and 3.81 cm).

10.5 For Models 93600 and 93610: Connecting System to Catheter and Patient Monitor

Use aseptic technique.

Step	Procedure
1	Verify the absence of air in the fluid path. Attach the flow-through/check valve to the catheter injectate lumen stopcock (Figure 9 on page 110).
2	Insert the injectate temperature probe into the flow-through/check valve until it snaps into place (See Figure 10 on page 111). Secure the system to minimize movement with respect to the patient.
3	Connect the injectate temperature probe cable to the "Injectate Probe" lead of the monitor cable.

10.6 For Models 93505: Connecting System Housing to Catheter and Patient Monitor

Use aseptic technique.

Step	Procedure
1	Verify that the fluid path is completely free of air. Attach the flow-through housing to the catheter injectate lumen stopcock. (Figure 11 on page 112)
2	Insert the reusable injectate temperature probe into the flow-through housing (Figure 12 on page 113). Secure the system to minimize movement with respect to the patient.
3	Connect the injectate temperature probe cable to the "Injectate Probe" lead of the computer cable.

10.7 Prepare Patient Monitor for Operation

Step	Procedure
1	Select injectate volume and catheter size on monitor to determine computation constant or enter the computation constant corresponding to the injectate volume and catheter model being used (see package insert accompanying the catheter for computation constant).
2	Prepare the monitor for operation according to the directions in the appropriate monitor operations manual. CO-Set+ closed injectate delivery system is now ready for use.

10.8 Cardiac Output Measurement

Step	Procedure
1	Open the snap clamp.
2	Turn the stopcock at the catheter injectate hub to close the IV flush and to open the fluid path between the syringe and catheter.
3	Withdraw desired amount (5 to 10 ml) of injectate into syringe.
4	Press start on patient monitor.
5	Immediately perform injections at a rate of at least 2.5 ml/sec. Note: The syringe plunger and the finger grips on the syringe jacket have been designed to help accommodate a wide range of user hand sizes for efficient and comfortable injection of fluid. To ensure uniform fluid injections, care should be taken to avoid pushing in the syringe plunger at an angle to the syringe barrel during injections.
6	Perform cardiac output determination according to instructions in the appropriate patient monitor operations manual. Note: To improve cardiac output reproducibility, space injections approximately one minute apart and increase injectate volume to 10 ml if not already doing so.
7	At the end of the measurements, return the stopcock to its original (flush) position (assuring that the syringe is empty), and close the snap clamp.

10.9 Injectate Temperature

Lower or higher temperatures will affect cardiac output accuracy when using the computation constants provided.

Ensure that the injectate temperature falls within the following ranges:

Model 93600:

Injectate Volume	Temperature Range
10 ml	6 °C - 12 °C
5 ml	8 °C - 16 °C

Model 93610:

Ensure that the injectate temperature falls within the range of 18 °C to 25 °C for 5 or 10 ml injectate.

Model 93505:

Ensure that the injectate temperature falls within 6 °C to 16 °C for cold injectate and 18 °C and 25 °C for room temperature injectate.

Models 93522 and 93650:

Please refer to package insert provided with Models 93600, 93610, or 93505 for detailed instructions for use.

10.10 Removing the Injectate Temperature Probe (See Fig. 13)

Grasp the flow-through/check valve, and carefully withdraw (unsnap) the Injectate Temperature Probe. Do not squeeze or pinch the area around the probe connector as this will make probe removal from the flow-through/check valve difficult.

10.11 Model 93505 only: To Remove Flow-Through Temperature Probe (Figure 12)

Grasp the flow-through housing, gently press both tabs, and withdraw the temperature probe.

11.0 Warnings

With the exception of the temperature probes, these devices are designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

12.0 Precautions

Device risks include infection, peripheral air embolism, adverse reaction to device materials, burns or shock, arrhythmia, inappropriate/unintended treatment, and/or delay in monitoring and/or treatment.

Care should be taken to avoid air entrapment during setup and use. The flow-through/check valve and syringe should be secured to minimize movement with respect to the patient.

Models 93522 and 9850A:

Do not immerse any electrical connectors of the cable or accessory equipment in liquids.

The recommended Use By Date imprinted on the flag label on the cable represents the recommended date the cable should be replaced or returned to Edwards Lifesciences for evaluation.

Please contact Edwards Technical Support or your local Edwards representative for further assistance.

13.0 How Supplied

Models 93600, 93610, 93650, and 93505:

Contents sterile and nonpyrogenic if package is undamaged or unopened. Do not use if package is opened or damaged. Do not resterilize. Visually inspect for breaches of packaging integrity prior to use.

The CO-Set+ closed injectate delivery system is recommended for single patient use only. It is recommended that CO-Set+ closed injectate delivery system be changed (replaced) for each patient according to hospital procedures for changing tubing sets containing dextrose solutions.

Models 93522 and 9850A are reusable and shipped non-sterile. Clean the probe by wiping it with a cleaning solution such as 70% isopropyl alcohol or 10% Bleach Solution.

14.0 Storage

Store in a cool, dry place.

15.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

16.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: This device is designed, intended, and distributed for single use only.

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

17.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Computation Constants for CO-Set+ devices - see package insert accompanying the catheter for computation constant.

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	Suomi	Български	Română	Eesti	Lietuvių
#	Mallinnumero	Номер на модела	Număr model	Mudeli number	Modelio numeris
QTY	Määrä	Количество	Cantitate	Kogus	Kiekis
!	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
2	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks.	Vienkartinio naudojimo
!	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
STERILIZE	Älä steriloι uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
☀ ☔	Säilytettävä kuivassa ja viileässä.	Да се съхранява на хладно, сухо място.	Stocați într-un loc rece și uscat.	Säilitada jahedas ja kuivas kohas.	Laiykite vėsioje, sausoje vietoje.
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
CE 0123	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
MD	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
!	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
! eifu.edwards.com +1 888 570 4016	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
!	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Основные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
 eifu.edwards.com + 1 888 570 4016	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiata veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/ Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/ Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijose / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Viimeinen käyttöpäivä	Срок на годност	Data expirării	Kõlblik kuni	Naudoti iki nurodytos datos
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Epästeriili	Нестерилно	Nesteril	Mittesteriilne	Nesterilus
	Sisältää luonnokumilateksia tai sen jäämiä	Съдържание или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural	Sisaldab naturaalseid kummlatekseid	Turi savo sudėtyje ar yra natūralios gumos latekso
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна барьерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Основные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
#	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
QTY	Daudzums	Miktar	Количество	Količina	Količina
!	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
2	Vienreizējai lietošanai	Tek Kullanımlık	Только для одноразового использования	Za jednokratnu upotrebu	Za jednokratnu upotrebu
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
2 STERILIZE	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Uzglabāt vēsā, sausā vietā.	Serin, kuru yerde saklayın.	Хранить в прохладном, сухом месте.	Čuvati na hladnom i suvom mestu.	Pohranite na hladnom, suhom mjestu.
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
CE 0123	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
MD	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
eifu.edwards.com + 1 888 570 4016	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na web-stranici
	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символах ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Основные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
 eifu.edwards.com + 1 888 570 4016	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/ Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Nesterils	Steril Değildir	Нестерильно	Nesterilno	Nesterilno
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilenoksidom	Sterilizirano etilenoksidom
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykietcie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei vältämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminių etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almayış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

EC REP

Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

Made in Dominican Republic
CCT Critical Care Technologies, S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

12/2024
10056131001 A / DOC-0216441 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU