Edwards

FloTrac Pressure Cable—Philips

Directory	
English (en).....	1
Français (fr).....	4
Deutsch (de).....	7
Español (es).....	10
Italiano (it).....	13
Nederlands (nl).....	16
Dansk (da).....	19
Svenska (sv).....	22
Ελληνικά (el).....	25
Português (pt).....	29
Česky (cs).....	32
Magyar (hu).....	35
Polski (pl).....	38
Slovensky (sk).....	41
Norsk (no).....	44
Suomi (fi).....	47
Български (bg).....	50
Română (ro).....	54
Eesti (et).....	57
Lietuvių (lt).....	60
Latviešu (lv).....	63
Türkçe (tr).....	66
Русский (ru).....	69
Srpski (sr).....	73

Hrvatski (hr).....	76
简体中文 (zh-cn).....	79
繁體中文 (zh-tw).....	82
한국어 (ko).....	85
Slovenščina (sl).....	88
Symbol Legend ■ Légende des symboles	
■ Zeichenerklärung ■ Significado de los símbolos	
■ Legenda dei simboli ■ Lijst met symbolen	
■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα	
συμβόλων ■ Legenda de símbolos ■ Legenda se symboly	
■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky	
k symbolom ■ Symbolforklaring ■ Merkkien selitykset	
■ Легенда на символите ■ Legendă de simboluri	
■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu	
skaidrojums ■ Sembol Açıklamaları ■ Условные	
обозначения ■ Legenda sa simbolima ■ Legenda simbola	
■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov	
..... 91	

English

FloTrac Pressure Cable—Philips


Figure 1: FloTrac Pressure Cable—Philips

④ Connector may vary depending on the applicable compatible Philips monitor.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The FloTrac Pressure Cable—Philips is intended to be used by qualified personnel or trained clinicians in a critical care environment in a hospital setting. The FloTrac Pressure

Edwards, Edwards Lifesciences, the stylized E logo, and FloTrac are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Cable—Philips is intended for use with compatible Philips monitors, and Edwards pressure monitoring sensors and compatible catheters.

The FloTrac Pressure Cable—Philips is a reusable device that connects with a compatible Philips monitor on one end ④ and an Edwards pressure monitoring sensor on the other end ① to continuously measure hemodynamic parameters. Refer to the Philips monitor's Instructions for Use for a list of supported hemodynamic parameters and applicable Edwards pressure monitoring sensor(s).

Note: The expected useful life of the FloTrac Pressure Cable—Philips is 5 years from date of purchase. If your equipment experiences a malfunction, please contact Technical Support or your local Edwards representative for further assistance.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions for Use. Refer to the Philips monitor's Instructions for Use for comprehensive monitoring procedures and monitoring parameters.

2.0 Compatible Sensors

The FloTrac Pressure Cable—Philips is intended to be used in combination with a compatible Philips monitor, an Edwards pressure monitoring sensor and compatible pressure monitoring catheter to measure intravascular pressure. Refer to the Philips monitor's Instructions for Use for a list of supported Edwards-compatible pressure monitoring sensors and catheters.

3.0 Intended Use/Purpose

The FloTrac Pressure Cable—Philips when used with an Edwards' pressure monitoring sensor and connected to a compatible monitor, offers continuous assessment of hemodynamic parameters. It may be used for monitoring hemodynamic parameters in conjunction with a perioperative goal directed therapy protocol in a hospital environment.

4.0 Indications for Use

The FloTrac Pressure Cable—Philips is indicated for use in adult critical care patients in which the balance between cardiac function, fluid status, vascular resistance and pressure needs continuous assessment.

5.0 Contraindications

There are no known contraindications.

6.0 Warnings

Compliance to IEC 60601-1 is only maintained when the FloTrac Pressure Cable—Philips (applied part accessory, defibrillation proof) is connected to a compatible Philips monitoring platform. Connecting external equipment or configuring the system in a way not described in these instructions will not meet this standard. Failure to use the device as instructed may increase the risk of electrical shock to the patient/operator.

Do not modify, service or alter the product in any way. Servicing, alteration or modification may affect patient/operator safety and/or product performance.

7.0 Precautions

Do not use if the package or the FloTrac Pressure Cable—Philips is damaged. Damage may include cracks, scratches, dents or any signs that the housing may be compromised.

Risks of the FloTrac Pressure Cable—Philips include patient or clinician burns or electrical shock, inappropriate/unintended treatment, and/or delay in treatment.

8.0 FloTrac Pressure Cable—Philips Setup

Refer to the compatible Philips monitor's instruction for use comprehensive monitoring procedures.

Step	Procedure
1	Connect the FloTrac Pressure Cable—Philips to a compatible Philips monitor at the end indicated by ④ in Figure 1.
2	Setup and initiate the compatible Philips monitor.
3	Insert the appropriate color insert at ② to indicate the monitored pressure type. The colors are as follows: - Red for arterial pressure (AP) - Blue for central venous pressure (CVP) - Yellow for pulmonary artery pressure (PAP) - Green for cardiac output (CO)
4	Connect the pressure monitoring sensor to the FloTrac Pressure Cable—Philips. ① The FloTrac Pressure Cable—Philips LED that surrounds the zero button at ③ will flash green indicating that the pressure monitoring sensor is detected. A yellow light indicates a sensor fault condition. If this occurs, refer to the monitor screen for specific fault condition details.
5	Follow all instructions found in pressure monitoring catheter IFU for catheter preparation and insertion procedures. The FloTrac Pressure Cable—Philips must be zeroed before monitoring.

9.0 Zero FloTrac Pressure Cable—Philips

Step	Procedure
1	To initiate the zeroing process, press and hold the zero button indicated by ③ in the figure or initiate zeroing process on the connected monitor screen. Refer to the compatible Philips monitor's Instruction for Use for on-screen zeroing instructions. CAUTION: To prevent cable damage, do not apply excessive force to the FloTrac Pressure Cable—Philips zero button.
2	Level the sensor stopcock valve (vent port) to the patient's phlebostatic axis position according to the instructions for use.
3	Open the stopcock valve to measure atmospheric conditions.
4	Press and hold the FloTrac Pressure Cable—Philips zero button or the on-screen button to complete zeroing process. The zero button LED will turn off when properly zeroed.
5	Confirm stable zero pressure and turn stopcock such that sensor is reading patient intravascular pressure.

10.0 MRI Information

MR Unsafe

Do not use the FloTrac Pressure Cable—Philips in an MR environment. The FloTrac Pressure Cable—Philips is MR unsafe since the device contains metallic components, which can experience RF-induced heating in the MRI environment.

11.0 Specifications

Physical Specifications

Cable Length: 3 m (10 ft)

Operating Conditions

Temperature: 10 to 32.5 °C

Humidity range: 20 to 90% non-condensing

Altitude (Atmospheric Pressure): 0 m/0 ft (1013 hPa) to

3048 m/10,000 ft (697 hPa)

IP Rating: IPX4

12.0 Maintenance

When unplugging the FloTrac Pressure Cable—Philips from a compatible Philips monitor or pressure monitoring sensor from the FloTrac Pressure Cable—Philips, always pull at the connection site. Do not pull from cables or use tools to disconnect.

CAUTION: Excessive dropping of the FloTrac Pressure Cable—Philips may result in cable damage and/or malfunction.

It is recommended to inspect the FloTrac Pressure Cable—Philips periodically for signs of damage. Check the condition of the cable and replace if any of the following are present: exposed wire, splits, cracks, or signs of stress.

The FloTrac Pressure Cable—Philips is a reusable device and can be surface cleaned per hospital aseptic procedures. Device has been shown to handle cloth wipe down using the following chemical agents: 70% isopropyl, 2% glutaraldehyde, 10% bleach solution, or quaternary ammonium solution. Disconnect the FloTrac Pressure Cable—Philips from the monitor to air dry the transducer connector. To blow dry the transducer connector, use clean, dry wall air, canned air, or CO₂ aerosol for at least two minutes. If left to dry under room conditions, allow the connector to dry for two days before using.

CAUTION: Do not use any other cleaning agents, spray, or pour cleaning solution directly on the FloTrac Pressure Cable—Philips.

Do not steam, radiate, or EtO sterilize the FloTrac Pressure Cable—Philips.

Device contains electronics. Handle with care.

Do not immerse.

Dispose of in accordance with hospital policy and local regulations.

13.0 Storage Conditions

Store in a cool, dry place.

Store in original packaging to prevent damage. Keep away from sunlight.

CAUTION: Do not coil the cable against its natural orientation or coil the cable tightly.

14.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications, and model availability are subject to change without notice.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

15.0

Table 1: FloTrac Pressure Cable—Philips Parameters:

Parameters ¹	Description	Sub-System Technology Used	Patient Population	Hospital Environment
CO	Continuous Cardiac Output	FloTrac Pressure Cable—Philips	Adult only	Operating room, intensive care unit, emergency room
CI	Continuous Cardiac Index			
CVP	Central Venous Pressure			
DIA_art	Diastolic Blood Pressure			
MAP	Mean Arterial Pressure			
MPAP	Mean Pulmonary Arterial Pressure			
PPV	Pulse Pressure Variation			
SV	Stroke Volume			
SVI	Stroke Volume Index			
SVR	Systemic Vascular Resistance			
SVRI	Systemic Vascular Resistance Index			
SVV	Stroke Volume Variation			
SYS_art	Systolic Blood Pressure			
PR	Pulse Rate			

¹Available parameters will vary depending on the pressure monitoring sensor that is connected, a list of supported parameters will be provided in the Philips monitor's Instructions for Use.


Figure 1 : Câble de pression FloTrac — Philips

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


	English	Français	Deutsch	Español	Italiano
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Consultare le istruzioni per l'uso sul sito Web
	Serial Number	Numéro de série	Seriennummer	Número de serie	Numero di serie
	Reference Number	Numéro de référence	Referenznummer	Número de referencia	Numero di riferimento
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea
	Service number for Philips	Numéro d'assistance Philips	Servicenummer für Philips	Número de servicio de Philips	Numero dell'assistenza per Philips
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM	Non compatibile con RM
	Provides protection against water splashing in any direction to IPX4 standard	Assure une protection contre les projections d'eau multidirectionnelles conformément à la norme IPX4	Schutz gegen Spritzwasser aus allen Richtungen gemäß IPX4	Proporciona protección contra las salpicaduras de agua en cualquier dirección hasta el estándar IPX4	Fornisce protezione contro gli spruzzi d'acqua in qualsiasi direzione secondo lo standard IPX4
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Defibrillator Proof Type CF applied part	Partie appliquée de type CF résistante à la défibrillation	Defibrillationsgeschütztes Anwendungsteil des Typs CF	Pieza-aplicada de tipo CF a prueba de desfibrilación	Parte applicata di tipo CF a prova di defibrillazione
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Consultare le istruzioni per l'uso

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov

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


	Nederlands	Dansk	Svenska	Ελληνικά	Português
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Defibrillatiebestendig toegepast onderdeel van type CF	Defibrillatorsikker type CF-anvendt del	Defibrilleringssäker patientansluten del typ CF	Εφαρμοζόμενο εξάρτημα τύπου CF ανθεκτικό σε απινίδωση	Peça aplicada de tipo CF à prova de desfibrilhação
	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen	Ακολουθείτε τις οδηγίες χρήσης	Seguir as instruções de utilização
UDI	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Milieuvriendelijke gebruikperiode (EFUP) - Alleen China	Periode for miljøvenlig brug (EFUP, Environmentally friendly use period) – kun i Kina	Miljövänlig användningsperiod (EFUP, Environmentally Friendly Use Period) – endast Kina	Περίοδος χρήσης φιλική προς το περιβάλλον (EFUP) - Μόνο για την Κίνα	Período de utilização seguro para o ambiente (EFUP, Environmentally Friendly Use Period) – Apenas China
FC	Federal Communications Commission	Den Føderale Kommunikationskommission	Federal Communications Commission	Ομοσπονδιακή Επιτροπή Επικοινωνιών	Comissão Federal de Comunicações
	Gescheiden inzameling voor elektrische en elektronische apparatuur in overeenstemming met EG-richtlijn 2012/19/EU.	Separat indsamling af elektrisk og elektronisk udstyr iht. EF-direktivet 2012/19/EU.	Separat insamling för elektrisk och elektronisk utrustning i enlighet med EU-direktiv 2012/19/EU.	Ξεχωριστή συλλογή για ηλεκτρικό και ηλεκτρονικό εξοπλισμό σύμφωνα με την Οδηγία του Ευρωπαϊκού Κοινοβουλίου 2012/19/ΕΕ.	Recolha em separado para o equipamento elétrico e eletrônico, de acordo com a Diretiva CE 2012/19/UE.
	Uit het zonlicht houden	Beskyttes mod sollys	Skyddas mot solljus	Φυλάσσετε μακριά από την άμεση ηλιακή ακτινοβολία	Manter afastado da luz solar
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
CE 0123	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marcação CE)
QTY	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Markering voor UL-erkende component voor Canada en de Verenigde Staten	UL-godkendelsesmærke for Canada og USA	UL-godkännandemärke för komponentmärke för Kanada och USA	Σήμα εξαρτήματος με πιστοποίηση UL για τον Καναδά και τις Ηνωμένες Πολιτείες	Marca de componente reconhecido UL para o Canadá e os Estados Unidos
	Medisch hulpmiddel	Medicinsk udstyr	Medicinsk utrustning	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Breekbaar, voorzichtigheid geboden	Skørbelig, skal håndteres forsigtigt	Ömtåligt, hanteras försiktigt	Εύθραστο: χειριστείτε με προσοχή	Frágil, manusear com cuidado
	Deze kant boven	Denne side opad	Denna sida upp	Αυτή η πλευρά προς τα πάνω	Este lado para cima

	Česky	Magyar	Polski	Slovensky	Norsk
	Číslo modelu	Típuszám	Numer modelu	Číslo modelu	Modellnummer
	Dodržujte návod k použití uvedený na webových stránkách.	Kövesse a következő honlapon található használati utasítást	Postępować zgodnie z instrukcją stosowania na stronie internetowej	Postupujte podľa návodu na použitie na webovej stránke	Følg bruksanvisningen på nettsiden
	Sériové číslo	Sorozatszám	Numer seryjny	Sériové číslo	Serienummer
	Referenční číslo	Referenciaszám	Numer referencyjny	Referenčné číslo	Referansenummer
	Výrobce	Gyártó	Producent	Výrobca	Produsent
	Datum výroby	Gyártás ideje	Data produkcji	Dátum výroby	Produksjonsdato
	Autorizovaný zástupce v Evropském společenství / Evropské unii	Hivatalos képviselő az Európai Közösségekben/ Európai Unióban	Autoryzowany przedstawiciel we Wspólnocie Europejskiej/Unii Europejskiej	Autorizovaný zástupca v krajinách Európskeho spoločenstva/ Európskej únie	Autorisert representant i Det europeiske fellesskap / EU

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov

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	Suomi	Български	Română	Eesti	Lietuvių
	Viitenumero	Референтен номер	Număr de referință	Viitenumber	Nuorodinis numeris
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/ Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Philips-huoltonumero	Сервизен номер за Philips	Număr operațiuni de service pentru Philips	Philips hooldusnumber	„Philips“ skirtas techninės priežiūros numeris
	Ei sovi magneettikuvaukseen	Небезопасно при MR	Incompatibil RM	Ohtlik magnetresonantsto mograafias	MR nesaugus
	Kotelo-intiloukka IPX4 (suojaus roiskuvasta vedeltä)	Предоставя защита срещу пръски вода във всички посоки според стандарт IPX4	Asigură protecție la standardul IPX4 împotriva apei împrorscate în orice direcție	Kaitseb mis tahes suunas langevate veepritsmete eest standardi IPX4 järgi.	Apsauga nuo aptaškymo vandeniu bet kuria kryptimi pagal IPX4 standartą
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Defibrillointia kestävä CF-tyypin liityntäosa	Устойчива на дефибрилации приложна част тип CF	Componentă aplicată, tip CF, rezistentă la defibrilare	Defibrillaatorikindel CF-tüüpi kohaldatav osa	Defibriliavimo poveikiui atspari CF tipo darbinė dalis
	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Ympäristöystävällinen käyttöaika (EFUP) – vain Kiina	Екологичен период на използване (EFUP) – само за Китай	Perioadă de utilizare ecologică (EFUP) – Numai pentru China	Keskkonnahoidlik kasutusaeg (EFUP) – ainult Hiina	Aplinkai nekenksmingo naudojimo laikotarpis (EFUP) – tik Kinijoje




Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov

	Suomi	Български	Română	Eesti	Lietuvių
	Yhdysvaltojen telehallintovirasto	Федерална комисия за комуникации	Comisia Federală pentru Comunicații	Föderaalne Sidekomisjon	Federalinė ryšių komisija
	Eröttele sähköiset ja elektroniset laitteet EY-direktiivin 2012/19/EU mukaisesti.	Разделно събиране на електрическо и електронно оборудване според Директивата на ЕС 2012/19/ЕС.	Colectarea separată a echipamentelor electrice și electronice, în conformitate cu directiva CE 2012/19/UE.	Elektri- ja elektroonikaseadmet e eraldi kogumine kooskõlas EÜ direktiiviga 2012/19/EL.	Atskiras elektros ir elektroninės įrangos atliekų surinkimas pagal EB direktyvą 2012/19/ES.
	Säilytä auringonvalolta suojattuna	Да се пази от слънчева светлина	A se feri de lumina solară	Hoida eemal päikesevalgusest	Saugoti nuo saulės šviesos
	Älä käyttää, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Määrä	Количество	Cantitate	Kogus	Kiekis
	UL Recognized Component Mark -laatuimerkki Kanadassa ja Yhdysvalloissa	Маркировка за разпознат UL компонент за Канада и САЩ	Marcă de componentă recunoscută de UL pentru Canada și Statele Unite	UL-i tunnustatud komponendi märgis Kanadas ja Ameerika Ühendriikides	UL pripažįstamo komponento ženklas Kanadoje ir Jungtinėse Amerikos Valstijose
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Särkyvä, käsiteltävä varoen	Чупливо, да се борава внимателно	Fragil, a se manipula cu atenție	Õrn, käsitseda ettevaatlikult	Trapus, elgtis atsargiai
	Tämä puoli ylöspäin	С тази страна нагоре	Cu partea aceasta în sus	See pool üles	Šia puse į viršų

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	levērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Sērijas numurs	Seri Numarası	Серийный номер	Serijski broj	Serijski broj
	Atsauces numurs	Referans Numarası	Ссылочный номер	Broj reference	Referentni broj
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Pilnvarotais pārstāvis Eiropas Kopienā/ Eiropas Savienībā	Avrupa Birliği'nde/ Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Philips apkopes numurs	Philips için servis numarası	Сервисный номер Philips	Broj servisa kompanije Philips	Broj servisa za Philips
	Nedrīkst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR	Nije sigurno kod pregleda MR-om
	Nodrošina aizsardzību pret ūdens šļakstīšanos jebkurā virzienā atbilstoši IPX4 standartam.	IPX4 standardına uygun olarak herhangi bir yönden sıçrayan suya karşı koruma sağlar	Обеспечивает защиту от брызг воды (в любом направлении) по стандарту IPX4	Obezbeđuje zaštitu od prskanja vode u bilo kom pravcu prema IPX4 standardu	Osigurava zaštitu od prskanja vode u bilo kojem smjeru u skladu sa standardom IPX4
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Pret defibrilāciju noturīga CF tipa daļa, kas saskaras ar pacientu	Defibrilasyona Dayanıklı CF Tipi hastaya temas eden parça	Контактирующий с пациентом элемент типа CF с защитой от разрядов дефибриллятора	Primenjeni deo tipa CF otporan na defibrilaciju	Primjenjen dio tipa CF otporan na defibrilaciju
	Levērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Vidi saudzējošas lietošanas periods (EFUP) — tikai Ķīna	Çevre dostu kullanım süresi (EFUP) — Yalnızca Çin	Срок экологически благоприятного использования (EFUP) — только для Китая	Period ekološki prihvatljive upotrebe (EFUP) – samo za Kinu	Ekološki prihvatljivo razdoblje upotrebe (Environmentally friendly use period, EFUP) – samo za Kinu
	Federālā komunikāciju komisija	Federal Haberleşme Komisyonu	Федеральная комиссия по связи	Savezna komisija za komunikacije	Savezna komisija za komunikacije
	Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar EK Direktīvu 2012/19/ES.	2012/19/EU sayılı AB direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır.	Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU.	Odvojeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU.	Odvojeno prikupljanje električne i elektroničke opreme u skladu s direktivom EZ-a 2012/19/EU.
	Sargāt no saules gaismas	Güneş ışığından uzak tutun	Беречь от солнечных лучей	Držati zaklonjeno od sunčeve svjetlosti	Držite dalje od sunčeve svjetlosti
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Daudzums	Miktar	Количество	Količina	Količina
	UL apstiprinātas detaļas marķējums Kanādai un Amerikas Savienotajām Valstīm	Kanada ve Amerika Birleşik Devletleri için UL Tarafından Kabul Gören Bileşen İşareti	Знак сертифицированного о UL компонента для Канады и США	UL oznaka za odobrene komponente za Kanadu i Sjedinjene Američke Države	Oznaka komponente koju je priznala organizacija UL za Kanadu i Sjedinjene Američke Države
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 符號圖例 ■ 기호 범례 ■ Legenda simbolov

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
	Ar šo pusi uz augšu	Bu taraf yukarı	Верхняя сторона	Ova strana mora da bude prema gore	Ova strana prema gore

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10/2023
10057683001 A / DOC-0224646 A
© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH