Directory

English (en).....	1
Français (fr).....	3
Deutsch (de).....	4
Español (es).....	5
Italiano (it).....	6
Nederlands (nl).....	7
Dansk (da).....	8
Svenska (sv).....	9
Ελληνικά (el).....	10
Português (pt).....	11
Česky (cs).....	12
Magyar (hu).....	13
Polski (pl).....	14
Slovensky (sk).....	15
Norsk (no).....	16
Suomi (fi).....	17
Български (bg).....	18
Română (ro).....	19
Eesti (et).....	20
Lietuvių (lt).....	21
Latviešu (lv).....	22
Türkçe (tr).....	23
Русский (ru).....	24
Srpski (sr).....	25
Hrvatski (hr).....	26
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	27

English

23001A and 23001B CO-Set+ injectate temperature probes

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Indications

This CO-Set+ injectate temperature probe is designed for Philips Medical Systems Cardiac Output Modules, for measurement of iced or room temperature injectate with CO-Set+ system during cardiac output procedure.

Accuracy

Accuracy is +/- 0.2 °C from 0 to 26 °C.

Warnings

- All wire-lead patient-connected sensor assemblies (including temperature probes) are subject to reading error, local heating and possible damage from high intensity sources of RF energy. Electrosurgical (ESU) equipment represents one such source. Capacitively coupled currents may seek alternate paths to ground

through probe cables and associated instruments. Patient burns may result.

- If at all possible, remove the probe from patient contact before activating the surgical unit or other RF source. If probes must be used simultaneously with electrosurgical apparatus, the instruments attached to the probes should be checked for adequate isolation from RF grounds. Hazards can be reduced by selecting a temperature monitoring point located away from the expected RF current path between the active probe and the return path.
- Do not use these temperature probes on patients who are undergoing Magnetic Resonance Imaging (MRI) procedures.

2.0 Handling Precautions

Be sure to observe the following:

- Never boil or autoclave.
- Handle gently – do not pull on the junction between the wire and the plug, or the junction between the wire and the sensing element (white tip, flat disk, or metal probe).
- Hold the lead wire (not the sensing element) while wiping the lead wire.
- Do not tie the probe cables in knots or bunch tightly.
- Do not wrap the probe cables around the instrument.
- Avoid soaking in alcohols, strong detergents, or highly alkaline solutions.
- Avoid contact with aromatic or chlorinated ketones, phenols, ether or ester solvents.
- Do not use gas disinfection.

Prolonged exposure to the above substances will cause the vinyl to lose flexibility and eventually crack. The user must determine that a probe is suitable and sufficiently flexible for use.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

Note: The expected useful life of the CO-Set+ injectate temperature probe is 3.5 years from date of purchase, at which point the probe should be replaced and returned to Edwards Lifesciences. Please contact Technical Support or your local Edwards representative for further assistance.

Refer to Instructions For Use provided with the corresponding injectate delivery system (model 93600 or model 93610 or model 93505 Flow-Through Housing).

3.0 Cleaning

Clean the probe by wiping it with a cleaning solution such as 70% isopropyl alcohol or 10% Bleach Solution.

4.0 Precaution

Do not immerse any electrical connectors of the probe or accessory equipment in liquids.

Please contact Edwards Technical Support or your local Edwards representative for further assistance.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

5.0 Sterilization

These devices are reusable and shipped non-sterile.

Edwards does not recommend ethylene oxide sterilization or autoclaving of these devices.

6.0 Warranty

Edwards Lifesciences warrants that the equipment is fit for the purposes and indications described in the labeling for a period of one (1) year after the date of purchase when used in accordance with the directions for use. Unless the equipment is used in accordance with such instructions, this warranty is void and of no effect. No other express or implied warranty exists, including any warranty of merchantability or fitness for a particular purpose. Edwards’ sole obligation and purchaser’s exclusive remedy for breach of warranty shall be limited to repair or replacement at Edwards’ option. Edwards shall not be liable for proximate, incidental, or consequential damages.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications, and model availability are subject to change without notice.

7.0 Storage

Probes and leads should be formed in loose loops to prevent damage during storage.

8.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

Refer to the symbol legend at the end of this document.

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

	English	Français	Deutsch	Español	Italiano
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea
	Separate collection for electrical and electronic equipment in accordance with EC Directive 2012/19/EU	Tri sélectif des équipements électriques et électroniques conformément à la Directive CE 2012/19/UE	Getrennte Sammlung von Elektro- und Elektronikgeräten gemäß EU-Richtlinie 2012/19/EU	Recogida por separado de equipos eléctricos y electrónicos según la directiva de la CE 2012/19/UE	Raccolta separata per attrezzature elettriche ed elettroniche in conformità alla Direttiva CE 2012/19/UE
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Consultare le istruzioni per l'uso
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Importer	Importateur	Importeur	Importador	Importatore
	Conformité Européenne (CE Mark) of TÜV SÜD Product Service GmbH (notified body)	Conformité européenne (marquage CE) de TÜV SÜD Product Service GmbH (organisme notifié)	Conformité Européenne (CE-Kennzeichnung) der TÜV SÜD Product Service GmbH (Benannte Stelle)	Conformité Européenne (Marca CE) de TÜV SÜD Product Service GmbH (organismo notificado)	Conformité Européenne (Marchio CE) di TÜV SÜD Product Service GmbH (ente notificato)
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto
	Quantity	Quantité	Menge	Cantidad	Quantità
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen
 ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli
 ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių
 paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español	Italiano
	Not made with natural rubber latex	Non fabriqué à partir de latex de caoutchouc naturel	Nicht aus Latex hergestellt	No contiene látex de caucho natural.	Non realizzato in lattice di gomma naturale
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM	Non compatibile con RM

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Lotnummer	Partinummer	Lotnummer	Αριθμός παρτίδας	Número de lote
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Gescheiden inzameling voor elektrische en elektronische apparatuur in overeenstemming met EG-richtlijn 2012/19/EU	Separat indsamling af elektrisk og elektronisk udstyr iht. EF-direktivet 2012/19/EU	Separat insamling för elektrisk och elektronisk utrustning i enlighet med EU-direktiv 2012/19/EU	Ξεχωριστή συλλογή για ηλεκτρικό και ηλεκτρονικό εξοπλισμό σύμφωνα με την Οδηγία του Ευρωπαϊκού Κοινοβουλίου 2012/19/ΕΕ	Recolha em separado para o equipamento elétrico e eletrônico, de acordo com a Diretiva CE 2012/19/UE.
	Niet-steriel	Ikke-steril	Icke-steril	Μη αποστειρωμένο	Não esterilizado
	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen	Ακολουθείτε τις οδηγίες χρήσης	Seguir as instruções de utilização
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Conformité Européenne (CE-markeering) van TÜV SÜD Product Service GmbH (aangemelde instantie)	Conformité Européenne (CE-mærke) fra TÜV SÜD Product Service GmbH (bemyndiget organ)	Conformité Européenne (CE-märke) för TÜV SÜD Product Service GmbH (anmält organ)	Conformité Européenne (σήμανση CE) της TÜV SÜD Product Service GmbH (κοινοποιημένος οργανισμός)	Conformité Européenne (Marcação CE) da TÜV SÜD Product Service GmbH (organismo notificado)
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco



Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Μοναδικό αναγνωριστικό συσκευής	Identificador único de dispositivo
	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen	Ακολουθείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Seguir as instruções de utilização no site
	Bevat geen natuurlijk rubberlatex	Ikke fremstillet med naturligt gummilatex	Innehåller inte naturgummilatex	Δεν κατασκευάζεται από φυσικό λάτεξ	Não fabricado com látex de borracha natural
	MRI-onveilig	MR-usikker	MR-farlig	Μη ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)	Utilização não segura em ambiente de RM

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
12/2023
10057449001 A / DOC-0224871 A
© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in Malaysia

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH