
## Directory

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |    |                      |    |                     |     |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----|----------------------|----|---------------------|-----|
| English (en) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 01 | Português (pt) ..... | 44 | Română (ro) .....   | 81  |
| Français (fr) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 06 | Česky (cs) .....     | 49 | Eesti (et) .....    | 86  |
| Deutsch (de) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 11 | Magyar (hu) .....    | 53 | Lietuvių (lt) ..... | 90  |
| Español (es) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 16 | Polski (pl) .....    | 58 | Latviešu (lv) ..... | 94  |
| Italiano (it) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 21 | Slovensky (sk) ..... | 63 | Türkçe (tr) .....   | 99  |
| Nederlands (nl) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 26 | Norsk (no) .....     | 68 | Русский (ru) .....  | 103 |
| Dansk (da) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 31 | Suomi (fi) .....     | 72 | Srpski (sr) .....   | 108 |
| Svenska (sv) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 35 | Български (bg) ..... | 76 | Hrvatski (hr) ..... | 113 |
| Ελληνικά (el) .....                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 39 |                      |    |                     |     |
| Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Φιγυρι ■ Figuri ■ Joonised ■ Paveikslėliai<br>■ Attēli ■ Šekiller ■ Рису́нки ■ Slike ■ Slike .....                                                                                                                                                                                                                                                                                     |    |                      |    |                     |     |
| Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos<br>■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Šūbolite seletus ■ Simbolių paaiškinimas<br>■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ..... |    |                      |    |                     |     |

## English

### Swan-Ganz IQ Pulmonary Arterial Catheter SvO<sub>2</sub>/VIP Catheter: AIQSGF8

Swan-Ganz IQ Pulmonary Arterial Catheter (Model AIQSGF8)

1. EEPROM Connector
2. Thermistor Connector
3. Balloon Inflation Valve
4. PA Distal Lumen
5. Balloon

6. Thermistor @ 4 cm
7. RV Port @ 12.7 cm
8. Optical Module Connector
9. RV Port Hub
10. PA Distal Lumen Hub

11. Infusion (VIP) Lumen Hub
12. Proximal Injectate Lumen Hub
13. Infusion (VIP) Port @ 30 cm
14. Proximal Injectate Port @ 26 cm

### For use with a compatible cardiac output computer†

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

**CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.**

Edwards, Edwards Lifesciences, the stylized E logo, CO-Set, CO-Set+, Swan, Swan-Ganz, Swan-Ganz IQ, and VIP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

† Measurement capabilities (i.e. CCO, CCO/SvO<sub>2</sub>) of the compatible cardiac output computer vary by model number. Ensure that the monitor being used is able to measure the desired parameters.

### For Single Use Only

For figures, see Figure 1 on page 118 and Figure 2 on page 119.

#### 1.0 Description

The device is used by medical professionals who have been trained in safe use of invasive hemodynamic technologies and clinical usage of pulmonary artery catheters as part of their respective institutional guidelines.

The Swan-Ganz IQ pulmonary arterial catheter is a flow-directed pulmonary artery catheter designed to enable the monitoring of hemodynamic pressures. The pulmonary artery (PA) distal lumen terminates at the distal tip. The proximal injectate lumen terminates at a port located 26 cm from the distal tip. When the distal tip is located in the pulmonary artery, the proximal injectate port will reside in the right atrium or vena cava, allowing for bolus cardiac output injections, right atrial pressure and right ventricular pressure monitoring, blood sampling, or infusion of solutions.

When used with a compatible cardiac output computer†, the Swan-Ganz IQ SvO<sub>2</sub>/VIP catheter allows for intermittent calculation and display of cardiac output and mixed venous oxygen saturation. The oximetry lumen (optical module connector) terminates at the distal tip. This lumen contains the fibers that transmit the light to the pulmonary artery for measurement of mixed venous oxygen saturation. Mixed venous oxygen saturation is monitored by fiberoptic reflectance spectrophotometry. The amount of light absorbed, refracted, and reflected depends on the relative amounts of oxygenated and deoxygenated hemoglobin in the blood.

The SvO<sub>2</sub>/VIP catheter provides an additional VIP lumen that allows for continuous infusion. The VIP catheter lumen (proximal infusion lumen) terminates at a port located 30 cm from the distal tip. This port allows for infusion of solutions, pressure monitoring or blood sampling.

The intravascular catheter is inserted through the central vein into the right side of the heart and is advanced towards the pulmonary

artery. Route of insertion can be internal jugular, femoral, antecubital and brachial veins. The body parts in contact are the atrium, ventricles, pulmonary artery and circulatory system.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended for use in adult critically ill or surgical patient populations. The device has not yet been tested in pediatric population or in pregnant or lactating women.

2.0 Intended Use / Purpose

The Swan-Ganz IQ catheter (model AIQSGF8) is a pulmonary artery catheter intended for short-term use in the central circulatory system for patients who require intracardiac hemodynamic monitoring, blood sampling and infusing solutions. When used with a compatible monitoring platform and accessories, the Swan-Ganz IQ catheter offers a comprehensive hemodynamic profile to help clinicians assess cardiovascular function and guide treatment decisions.

3.0 Indications for Use

Swan-Ganz catheters are diagnostic and monitoring tools used for hemodynamic monitoring of adult critically ill patients including but not limited to post major surgical recovery, trauma, sepsis, burns, pulmonary disease, pulmonary failure, and cardiac disease including heart failure.

4.0 Contraindications

Patients with either recurrent sepsis, or with hypercoagulopathy, in which the catheter could serve as a focal point for septic or bland thrombus formation, should not be considered candidates for a balloon flotation catheter.

These products contain metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

5.0 Warnings

No absolute contraindications to the use of flow-directed pulmonary artery catheters exist. However, a patient with a left bundle branch block may develop a right bundle branch block during catheter insertion, resulting in complete heart block. In such patients, temporary pacing modes should be immediately available.

Electrocardiographic monitoring during catheter passage is encouraged and is particularly important in the presence of either of the following conditions:

- Complete left bundle branch block, in which the risk of complete heart block is somewhat increased.
- Wolff-Parkinson-White syndrome and Ebstein's malformation, in which the risk of tachyarrhythmias is present.

Do not modify or alter the product in any way. Alteration or modification may affect product performance.

Air should never be used for balloon inflation in any situation where it may enter the arterial circulation, e.g., in all pediatric patients and in adults with suspected right to left intracardiac or intra-pulmonary shunts. Bacteria-filtered carbon dioxide is the recommended inflation medium because of its rapid absorption into the blood in the event of balloon rupture within the circulation. Carbon dioxide diffuses through the latex balloon, diminishing the balloon's flow-directed capability after 2 to 3 minutes of inflation.

Do not leave the catheter in a permanent wedge position. Furthermore, avoid lengthy balloon inflation while the catheter is in a wedge position; this occlusive maneuver may result in pulmonary infarction.

While the RV Port can be used for infusion of solutions such as normal saline or Ringer's lactate, infusion of

medication may lead to complications including Cardiac Arrhythmias.

This device is designed, intended, and distributed for SINGLE USE ONLY. DO NOT RESTERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity, and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Cleaning and resterilization will damage the integrity of the latex balloon. Damage may not be obvious during routine inspection.

6.0 Precautions

Failure of a balloon flotation catheter to enter the right ventricle or pulmonary artery is rare, but may occur in patients with an enlarged right atrium or ventricle particularly if the cardiac output is low or in the presence of tricuspid or pulmonic incompetence or pulmonary hypertension. Placement of the catheter may be difficult in patients with a history of annuloplasty rings and the use of fluoroscopy during insertion is recommended. Deep inspiration by the patient during advancement may also facilitate passage.

Clinicians using the device should be familiar with the device and understand its applications prior to use.

7.0 Recommended Equipment

WARNING: Compliance to IEC 60601-1 is only maintained when the catheter or probe (Type CF applied part, defibrillation proof) is connected to a patient monitor or equipment that has a Type CF defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the catheter or probe. Failure to ensure monitor or equipment compliance to IEC 60601-1 and catheter or probe compatibility may increase the risk of electrical shock to the patient/operator.

- Swan-Ganz catheter
- Percutaneous sheath introducer and contamination shield
- Compatible cardiac output platform for intermittent cardiac output and mixed venous oxygen saturation
- Connecting cables
- Model OM2 or OM2E Optical Module (Model AIQSGF8)
- Sterile flush system and pressure transducers
- Bedside ECG and pressure monitor system

In addition, the following items should be immediately available if complications arise during catheter insertion: antiarrhythmic drugs, defibrillator, respiratory assist equipment, and a means for temporary pacing.

8.0 Monitor Set-Up and Calibration for Mixed Venous Oxygen Saturation Monitoring

This section only applies to the AIQSGF8 model with SvO<sub>2</sub> capability.

The compatible cardiac output computer can be calibrated prior to catheter insertion by performing an *in vitro* calibration. When performing an *in vitro* calibration, do so before preparing the catheter (i.e. flushing the lumens). The catheter tip must not get wet before an *in vitro* calibration is performed. An *in vivo* calibration is required if an *in vitro* calibration is not done. *In vivo* calibration may be used to periodically recalibrate the monitor. Refer to the monitor operator's manual for detailed calibration instructions.

9.0 Catheter Preparation

Visually inspect for breaches of packaging integrity prior to use.

Use aseptic technique.

Note: Use of a protective catheter sheath is recommended.

Precaution: Avoid forceful wiping or stretching of the catheter during testing and cleaning as not to break the thermistor wire circuitry.

Before the catheter is inserted, the following preparation procedure should be followed:

| Step | Procedure                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Perform <i>in vitro</i> calibration (when measuring mixed venous oxygen saturation).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 2    | To open the catheter for calibration, peel back the top left perforated portion of the lid and fold diagonally according to the dotted fold line. In order to gain access to the catheter lumens for flushing, peel lid from the top right corner and fold according to the dotted fold line; otherwise peel off the entire tray lid by peeling back the lid from the bottom right corner. Flush lumens with sterile saline or dextrose solution to ensure patency and to remove air.                                                                                                                                                                             |
| 3    | Gently lift the catheter up and remove it from the silicone gripper (see Figure 1 on page 118, Step 1).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 4    | Once the catheter has cleared the silicone gripper, pull the balloon out of the calibrator cup and remove the catheter from the tray (see Figure 1 on page 118, Step 2).<br><b>Note: To avoid damaging the balloon, do not pull the balloon through the silicone gripper.</b>                                                                                                                                                                                                                                                                                                                                                                                     |
| 5    | Check balloon integrity by inflating it to the recommended volume. Check for major asymmetry and for leaks by submerging in sterile saline or water. Deflate balloon before insertion.                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 6    | Connect the catheter's injectate and pressure monitoring lumens to the flush system and pressure transducers. Ensure that the lines and transducers are free of air.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 7    | Test the thermistor's electrical continuity before insertion. Connect the thermistor to the monitor and confirm that no fault messages appear.<br><b>Note: The AIQSGF8 catheter is not capable of CCO monitoring. The EEPROM connector is not required for intermittent cardiac monitoring. The AIQSGF8 catheter contains no thermal filament. On some monitors when the thermistor and EEPROM connector are connected together, the message "Connect thermal filament for CO monitoring" or "Fault: CO - Check Thermal Filament Connection" will display and an alarm may sound. As such, leave unconnected and intermittent cardiac monitoring may proceed.</b> |
| 8    | If using a compatible cardiac output computer for intermittent cardiac output measurement, connect the thermistor to the monitor.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

10.0 Insertion Procedure

Swan-Ganz catheters can be inserted at the patient's bedside without the aid of fluoroscopy, guided by continuous pressure monitoring.

Simultaneous pressure monitoring from the distal lumen is recommended. Fluoroscopy is recommended for femoral vein insertion.

**Note:** Should the catheter require stiffening during insertion, slowly perfuse the catheter with 5 ml to 10 ml of cold sterile saline or 5% dextrose as the catheter is advanced through a peripheral vessel.

**Note:** The catheter should pass easily through the right ventricle and pulmonary artery and into a wedge position in less than a minute.

Although a variety of techniques can be used for insertion, the following guidelines are provided as an aid to the physician:

| Step | Procedure                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Introduce the catheter into the vein through a sheath introducer using percutaneous insertion using modified Seldinger technique.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| 2    | <p>Under continuous pressure monitoring, with or without the aid of fluoroscopy, gently advance the catheter into the right atrium. Entry of the catheter tip into the thorax is signaled by an increased respiratory fluctuation in pressure. Figure 2 on page 119 shows the characteristic intracardiac and pulmonary pressure waveforms.</p> <p><b>Note:</b> When the catheter is near the junction of the right atrium and the superior or inferior vena cava of the typical adult patient, the tip has been advanced approximately 40 cm from the right or 50 cm from the left antecubital fossa, 15 to 20 cm from the jugular vein, 10 to 15 cm from the subclavian vein, or about 30 cm from the femoral vein.</p>                                                                                                                                                                     |
| 3    | <p>Using the syringe provided, inflate the balloon with CO<sub>2</sub> or air to the maximum recommended volume. <b>Do not use liquid.</b> Note that an offset arrow on the gate valve indicates the "closed" position.</p> <p><b>Note:</b> Inflation is usually associated with a feeling of resistance. On release, the syringe plunger should usually spring back. If no resistance to inflation is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may continue to be used for hemodynamic monitoring. However, be sure to take precautions to prevent infusion of air or liquid into the balloon lumen.</p> <p><b>WARNING:</b> Pulmonary complications may result from improper inflation technique. To avoid damage to the pulmonary artery and possible balloon rupture, do not inflate above the recommended volume.</p> |

| Step | Procedure                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 4    | <p>Advance the catheter until pulmonary artery occlusion pressure (PAOP) is obtained, then passively deflate the balloon by removing the syringe from the gate valve. Do not forcefully aspirate as this may damage the balloon. After deflation, re-attach the syringe.</p> <p><b>Note:</b> Avoid prolonged maneuvers to obtain wedge pressure. If difficulties are encountered, give up the "wedge".</p> <p><b>Note:</b> Before reinflation with CO<sub>2</sub> or air, completely deflate the balloon by removing the syringe and opening the gate valve.</p> <p><b>Precaution:</b> It is recommended that the provided syringe be re-attached to the gate valve after balloon deflation to prevent inadvertent injection of liquids into the balloon lumen.</p> <p><b>Precaution:</b> If a right ventricular pressure tracing is still observed after advancing the catheter several centimeters beyond the point where the initial right ventricular pressure tracing was observed, the catheter may be looping in the right ventricle which can result in kinking or knotting of the catheter (see Complications). Deflate the balloon and withdraw the catheter into the right atrium. Reinflate the balloon and readvance the catheter to a pulmonary artery wedge position, then deflate the balloon.</p> <p><b>Precaution:</b> Catheter looping may occur when excessive length has been inserted, which could result in kinking or knotting (see Complications). If the right ventricle is not entered after advancing the catheter 15 cm beyond entry into the right atrium, the catheter may be looping, or the tip may be engaged in a neck vein with only the proximal shaft advancing into the heart. Deflate the balloon and withdraw the catheter until the 20 cm mark is visible. Reinflate the balloon and advance the catheter.</p> |
| 5    | <p>Reduce or remove any excessive length or loop in the right atrium or ventricle by slowly pulling the catheter back approximately 2 to 3 cm.</p> <p><b>Precaution:</b> Do not pull the catheter across the pulmonic valve while the balloon is inflated to avoid damage to the valve.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 6    | <p>Reinflate the balloon to determine the minimum inflation volume necessary to obtain a wedge tracing. If a wedge is obtained with less than the maximum recommended volume (see specifications table for balloon inflation capacity), the catheter must be withdrawn to a position where full inflation volume produces a wedge tracing.</p> <p><b>Precaution:</b> Overtightening the proximal Tuohy-Borst adapter of the contamination shield may impair catheter function.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| 7    | Confirm final catheter tip position with chest X-ray.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |

**Note:** After deflation, the catheter tip may tend to recoil towards the pulmonic valve and slip back into the right ventricle, requiring that the catheter be repositioned.

**Note:** If using a contamination shield, extend the distal end towards the introducer valve. Extend the proximal end of the catheter contamination shield to desired length, and secure.

11.0 Bolus Thermodilution Method

Bolus TD cardiac output measurement is made by injecting an exact amount of physiological solution (saline or dextrose) of known temperature into the right atrium or superior vena cava and by using the thermistor in the pulmonary artery to detect the resultant change in blood temperature. Cardiac output is inversely proportional to the area under the temperature-versus-time curve. The accuracy of this method depends on the accuracy with which the quantity and temperature of the injectate are known. The accuracy of the thermodilution method correlates well with the dye dilution technique and with the direct Fick method.

Edwards Lifesciences cardiac output computers require that a computation constant be used to correct for injectate temperature rise as it passes through the catheter. The computation constant is a function of injectate volume, temperature, and catheter dimensions. The computation constants listed in the specifications have been determined *in vitro*.

12.0 Maintenance and Use *in situ*

The catheter should remain indwelling only as long as is required by the patient's condition.

**Precaution:** The incidence of complications increases significantly with indwelling periods longer than 72 hours.

12.1 Catheter Tip Position

Keep catheter tip centrally located in a main branch of the pulmonary artery near the hilum of the lungs. Do not advance tip too far peripherally. Tip should be kept where full or near full inflation volume is required to produce a wedge tracing. The tip migrates toward periphery during balloon inflation.

12.2 Catheter Tip Migration

Anticipate spontaneous catheter tip migration towards periphery of pulmonary bed. Continuously monitor distal lumen pressure to verify tip position. If wedge tracing is observed when balloon is deflated, pull catheter back. Damage may be caused by prolonged occlusion or over-distension of vessel upon re-inflation of the balloon.

Spontaneous catheter tip migration towards the periphery of the lung occurs during cardiopulmonary bypass. Partial catheter withdrawal (3 to 5 cm) just before bypass should be considered, as it may help reduce distal migration and prevent permanent catheter wedging postbypass. After termination of bypass, the catheter may require repositioning. Check the distal pulmonary artery tracing before inflating the balloon.

**Precaution:** Over a period of time, the catheter tip may migrate towards the periphery of the pulmonary bed and lodge in a small vessel. Damage may occur either by prolonged occlusion or by over-distension of the vessel upon reinflation of the balloon (see Complications).

PA pressures should be continuously monitored with the alarm parameter set to detect physiologic changes as well as spontaneous wedge.

12.3 Balloon Inflation and Wedge Pressure Measurement

Re-inflation of the balloon should be performed gradually while monitoring pressures. Inflation is usually associated with a feeling of resistance. If no resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may still be used for hemodynamic monitoring, however, take precautions against infusion of air or liquids into the balloon lumen. During normal catheter use, keep inflation syringe attached to gate valve to prevent inadvertent injection of liquid into the balloon inflation lumen.

Measure wedge pressure only when necessary and only when tip is properly positioned (see above). Avoid prolonged maneuvers to

obtain wedge pressure and keep wedge time to a minimum (two respiratory cycles or 10 - 15 seconds), especially in patients with pulmonary hypertension. If difficulties are encountered, discontinue wedge measurements. In some patients, pulmonary arterial end-diastolic pressure can often be substituted for pulmonary artery wedge pressure if the pressures are nearly identical, obviating the need for repeated balloon inflation.

12.4 Spontaneous Tip Wedging

The catheter may migrate into the distal pulmonary artery and spontaneous tip wedging may occur. To avoid this complication, pulmonary artery pressure should be monitored continuously with a pressure transducer and display monitor.

Forward advancement should never be forced if resistance is encountered.

12.5 Patency

All pressure monitoring lumens should be filled with a sterile, heparinized saline solution (e.g., 500 I.U. heparin in 500 ml saline) and flushed at least once each half hour or by continuous slow infusion. If loss of patency occurs and cannot be corrected by flushing, the catheter should be removed.

12.6 General

Keep pressure monitoring lumens patent by intermittent flush or continuous slow infusion with heparinized saline solution. Infusion of viscous solutions (e.g., whole blood or albumin) is not recommended, as they flow too slowly and may occlude the catheter lumen.

**WARNING: To avoid pulmonary artery rupture, never flush the catheter when the balloon is wedged in the pulmonary artery.**

Periodically check IV lines, pressure lines, and transducers to keep them free of air. Also ensure that connecting lines and stopcocks remain tightly fitted.

13.0 MRI Information

The Swan-Ganz device is MR unsafe as the result of the device containing metallic components, which experience RF-induced heating in the MRI environment; therefore the device poses hazards in all MRI environments.

14.0 Complications

Invasive procedures involve some patient risks. Although serious complications are relatively uncommon, the physician is advised, before deciding to insert or use the catheter, to consider the potential benefits in relation to the possible complications. The techniques for insertion, methods of using the catheter to obtain patient data information, and the occurrence of complications is well described in the literature.

Strict adherence to these instructions and awareness of risks reduce the incidence of complications. Several known complications include:

14.1 Perforation of the Pulmonary Artery

Factors associated with fatal pulmonary artery rupture include pulmonary hypertension, advanced age, cardiac surgery with hypothermia and anticoagulation, distal catheter tip migration, arteriovenous fistula formation and other vascular traumas.

Extreme care should be used during the measurement of pulmonary artery wedge pressure in patients with pulmonary artery hypertension. In all patients, balloon inflation should be limited to two respiratory cycles, or 10 to 15 seconds.

A central location of the catheter tip near the hilum of the lung may prevent pulmonary artery perforation.

14.2 Pulmonary Infarction

Tip migration with spontaneous wedging, air embolism, and thromboembolism can lead to pulmonary artery infarction.

14.3 Cardiac Arrhythmias

Cardiac arrhythmias may occur during insertion, withdrawal, and repositioning, but are usually transient and self-limited. Premature ventricular contractions are the most commonly observed arrhythmia. Ventricular tachycardia and atrial tachycardia have been reported. Use of prophylactic lidocaine should be considered to decrease the incidence of ventricular arrhythmias during catheterization. ECG monitoring and immediate availability of antiarrhythmic drugs and defibrillator equipment is recommended.

14.4 Knotting

Flexible catheters have been reported to form knots, most often as a result of looping in the right ventricle. Sometimes the knot can be resolved by insertion of a suitable guidewire and manipulation of the catheter under fluoroscopy. If the knot does not include any intracardiac structures, the knot may be gently tightened and the catheter withdrawn through the site of entry.

14.5 Sepsis/Infection

Positive catheter-tip cultures resulting from contamination and colonization have been reported, as well as incidences of septic and aseptic vegetation in the right heart. Increased risks of septicemia and bacteremia have been associated with blood sampling, the infusing of fluids, and catheter-related thrombosis. Preventive measures should be taken to guard against infection.

14.6 Other Complications

Other complications include right bundle branch block, complete heart block, tricuspid and pulmonic prosthetic mechanical valve damage, blood loss, cardiac structure/wall/vessel injury or damage, hematoma, embolism, anaphylaxis, pneumothorax, additional surgical procedure due to catheter entanglement with other intra-vascular devices, necrosis and thrombosis. In addition, allergic reactions to latex have been reported. Physicians should identify latex sensitive patients and be prepared to treat allergic reactions promptly.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance for this medical device. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

15.0 Long-Term Monitoring

The duration of catheterization should be the minimum required by the patient's clinical state since the risk of thromboembolic and infectious complications increases with time. The incidence of complications increases significantly with indwelling periods longer than 72 hours. Prophylactic systemic anticoagulation and antibiotic protection should be considered when long-term catheterization (i.e., over 48 hours) is required, as well as in cases involving increased risk of clotting or infection.

16.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

The packaging within is designed to avoid crushing of the catheter and to protect the balloon from exposure to the atmosphere. It is therefore recommended that the catheter remain inside the package until use.

17.0 Storage

Store in a cool, dry place.

18.0 Operating Conditions/Use Environment

Intended to operate under physiological conditions of the human body in a controlled clinical environment.

19.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the recommended expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

**Note: Resterilization will not extend the shelf life.**

20.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers: Inside the U.S. and Canada (24 hours): ..... 800.822.9837  
Outside the U.S. and Canada (24 hours): ..... 949.250.2222  
In the UK: ..... 0870 606 2040 - Option 4  
In Ireland: ..... 01 8211012 - Option 4

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**

21.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

**Refer to the symbol legend at the end of this document.**

Specifications:

|                                   | Model<br>AIQSGF8      |
|-----------------------------------|-----------------------|
| Body Color                        | yellow                |
| Usable Length (cm)                | 110                   |
| Catheter Body                     | 8F (2.7 mm)           |
| Diameter of Inflated Balloon (mm) | 13                    |
| Required Introducer Size          | 9F (3.0 mm)           |
| Balloon Inflation Capacity (ml)   | 1.5                   |
| Distance from Distal Tip (cm)     |                       |
| Thermistor                        | 4                     |
| RV Port                           | 12.7                  |
| Injectate Port                    | 26                    |
| Infusion (VIP) Port               | 30                    |
| Distance Between Markings (cm)    | 10                    |
| Lumen Volumes (ml)                |                       |
| Distal Lumen                      | 0.90                  |
| RV Lumen                          | 0.78                  |
| Injectate Lumen                   | 0.85                  |
| Infusion (VIP) Lumen              | 1.10                  |
| Infusion Rate* (ml/min)           |                       |
| Distal Lumen                      | 4                     |
| RV Lumen                          | 5                     |
| Injectate Lumen                   | 9                     |
| Infusion (VIP) Lumen              | 16                    |
| Compatible Guidewire Diameter     | 0.018 in<br>(0.46 mm) |
| Frequency Response                |                       |
| Distortion at 10 Hz               |                       |
| Distal Lumen                      | < 3 dB                |

|                           | <b>Model</b>               |
|---------------------------|----------------------------|
|                           | <b>AIQSGF8</b>             |
| RV Lumen                  | < 3 dB                     |
| Balloon Inflation Syringe | 3 ml,<br>limited to 1.5 ml |

All specifications given are nominal values.

\* Using normal saline at room temperature, 1 m above insertion site, gravity drip.

#### Catheter Model and Functions:

|                  | <b>AIQSGF8</b> |
|------------------|----------------|
| IC0              | X              |
| SvO <sub>2</sub> | X              |
| Infusion (VIP)   | X              |

#### Computation Constants

##### For use with bath temperature probes


| <b>Model</b>               | <b>AIQSGF8</b>               |                                     |
|----------------------------|------------------------------|-------------------------------------|
| <b>Injectate Temp (°C)</b> | <b>Injectate Volume (ml)</b> | <b>Computation Constants (CC)**</b> |
| 0 - 5                      | 10                           | 0.550                               |
|                            | 5                            | 0.256                               |
|                            | 3                            | —                                   |
| 19 - 22                    | 10                           | 0.585                               |
|                            | 5                            | 0.282                               |
|                            | 3                            | —                                   |
| 23 - 25                    | 10                           | 0.600                               |
|                            | 5                            | 0.292                               |

##### Computation Constants for CO-Set+ Delivery System

|               |    |       |
|---------------|----|-------|
| 6 °C - 12 °C  | 10 | 0.559 |
| 8 °C - 16 °C  | 5  | 0.263 |
| 18 °C - 25 °C | 10 | 0.602 |
|               | 5  | 0.295 |

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
|                | English                                                                                     | Français                                                                                    | Deutsch                                                                                     | Español                                                                                     | Italiano                                                                                    |
|----------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|                | Number of lumens                                                                            | Nombre de lumières                                                                          | Anzahl der Lumen                                                                            | Número de luces                                                                             | Numero di lumi                                                                              |
|                | Exterior diameter                                                                           | Diamètre externe                                                                            | Außendurchmesser                                                                            | Diámetro exterior                                                                           | Diametro esterno                                                                            |
|                | Usable length                                                                               | Longueur utile                                                                              | Nutzlänge                                                                                   | Longitud útil                                                                               | Lunghezza utile                                                                             |
|                | Recommended guidewire size                                                                  | Taille de fil-guide recommandée                                                             | Empfohlene Führungsdrahtgröße                                                               | Tamaño de guía recomendado                                                                  | Misura filo guida consigliata                                                               |
|                | Lumen size                                                                                  | Taille de la lumière                                                                        | Lumendurchmesser                                                                            | Tamaño de la luz                                                                            | Dimensioni del lume                                                                         |
|                | Model Number                                                                                | Référence                                                                                   | Modellnummer                                                                                | Número de modelo                                                                            | Numero di modello                                                                           |
|                | Minimum introducer size                                                                     | Taille minimale de l'introducteur                                                           | Mindestgröße der Einführhilfe                                                               | Tamaño de introductor mínimo                                                                | Misura minima dell'introdotto                                                               |
|                | Caution                                                                                     | Avertissement                                                                               | Vorsicht                                                                                    | Aviso                                                                                       | Attenzione                                                                                  |
|                | Do not re-use                                                                               | Ne pas réutiliser                                                                           | Nicht wiederverwenden                                                                       | No reutilizar                                                                               | Non riutilizzare                                                                            |
|                | Quantity                                                                                    | Quantité                                                                                    | Menge                                                                                       | Cantidad                                                                                    | Quantità                                                                                    |
|                | Lot Number                                                                                  | N° du lot                                                                                   | Chargenbezeichnung                                                                          | Número de lote                                                                              | Numero di lotto                                                                             |
|                | Use-by date                                                                                 | Date d'expiration                                                                           | Verwendbar bis                                                                              | Fecha de caducidad                                                                          | Utilizzare entro                                                                            |
|                | Inner diameter                                                                              | Diamètre interne                                                                            | Innendurchmesser                                                                            | Diámetro interior                                                                           | Diametro interno                                                                            |
|                | Balloon capacity                                                                            | Capacité du ballonnet                                                                       | Ballonkapazität                                                                             | Capacidad del balón                                                                         | Capacità del palloncino                                                                     |
|                | Sterile                                                                                     | Stérile                                                                                     | Steril                                                                                      | Estéril                                                                                     | Sterile                                                                                     |
|                | Sterilized using ethylene oxide                                                             | Stérilisé à l'oxyde d'éthylène                                                              | Mit Ethylenoxid sterilisiert                                                                | Esterilizado con óxido de etileno                                                           | Sterilizzato con ossido di etilene                                                          |
|                | Sterilized using irradiation                                                                | Stérilisé par irradiation                                                                   | Durch Bestrahlung sterilisiert                                                              | Esterilizado con radiación                                                                  | Sterilizzato mediante radiazioni                                                            |
|                | Sterilized using steam or dry heat                                                          | Stérilisé à la vapeur d'eau ou à la chaleur sèche                                           | Mit Wasserdampf oder trockener Hitze sterilisiert                                           | Esterilizado con vapor o calor seco                                                         | Sterilizzato con vapore o calore secco                                                      |
| <b>Rx only</b> | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|                | Manufacturer                                                                                | Fabricant                                                                                   | Hersteller                                                                                  | Fabricante                                                                                  | Produttore                                                                                  |
|                | Date of manufacture                                                                         | Date de fabrication                                                                         | Herstellungsdatum                                                                           | Fecha de fabricación                                                                        | Data di produzione                                                                          |
|                | Contains or presence of natural rubber latex                                                | Contient ou présence de latex de caoutchouc naturel                                         | Enthält Naturlatex                                                                          | Contiene o hay presencia de látex de caucho natural                                         | Contiene o è presente lattice di gomma naturale                                             |
|                | Size                                                                                        | Taille                                                                                      | Größe                                                                                       | Tamaño                                                                                      | Misura                                                                                      |
|                | Consult instructions for use on the website                                                 | Consulter le mode d'emploi sur notre site Web                                               | Gebrauchsanweisung auf der Website beachten                                                 | Consulte las instrucciones de uso en el sitio web                                           | Consultare le istruzioni per l'uso sul sito Web                                             |

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|                | Nederlands                                                                                  | Dansk                                                                                       | Svenska                                                                                     | Ελληνικά                                                                                    | Português                                                                                   |
|----------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|                | Aantal lumina                                                                               | Antal lumens                                                                                | Antal lumen                                                                                 | Αριθμός αυλών                                                                               | Números de lúmenes                                                                          |
|                | Buitendiameter                                                                              | Udvendig diameter                                                                           | Ytterdiameter                                                                               | Εξωτερική διάμετρος                                                                         | Diâmetro exterior                                                                           |
|                | Bruikbare lengte                                                                            | Anvendelig længde                                                                           | Brukbar längd                                                                               | Ωφέλιμο μήκος                                                                               | Comprimento útil                                                                            |
|                | Aanbevolen maat voerdraad                                                                   | Anbefalet guidewirestørrelse                                                                | Rekommenderad ledarstorlek                                                                  | Συνιστώμενο μέγεθος οδηγού σύρματος                                                         | Tamanho recomendado do fio-guia                                                             |
|                | Lumengrootte                                                                                | Lumenstørrelse                                                                              | Lumenstorlek                                                                                | Μέγεθος αυλού                                                                               | Tamanho do lúmen                                                                            |
|                | Modelnummer                                                                                 | Modelnummer                                                                                 | Modellnummer                                                                                | Αριθμός μοντέλου                                                                            | Número do modelo                                                                            |
|                | Minimum introducergrootte                                                                   | Minimum introducerstørrelse                                                                 | Minsta storlek på införsare                                                                 | Ελάχιστο μέγεθος εισαγωγέα                                                                  | Tamanho mínimo do introdutor                                                                |
|                | Let op                                                                                      | Forsigtig                                                                                   | Var försiktig                                                                               | Προσοχή                                                                                     | Aviso                                                                                       |
|                | Niet hergebruiken                                                                           | Må ikke genanvendes                                                                         | Får inte återanvändas                                                                       | Μην επαναχρησιμοποιείτε                                                                     | Não reutilizar                                                                              |
|                | Hoeveelheid                                                                                 | Mængde                                                                                      | Antal                                                                                       | Ποσότητα                                                                                    | Quantidade                                                                                  |
|                | Lotnummer                                                                                   | Partinummer                                                                                 | Lotnummer                                                                                   | Αριθμός παρτίδας                                                                            | Número de lote                                                                              |
|                | Vervaldatum                                                                                 | Anvendes inden                                                                              | Sista förbrukningsdag                                                                       | Ημερομηνία λήξης                                                                            | Data de vencimento                                                                          |
|                | Binnendiameter                                                                              | Indre diameter                                                                              | Invändig diameter                                                                           | Εσωτερική διάμετρος                                                                         | Diâmetro interno                                                                            |
|                | Balloncapaciteit                                                                            | Ballonkapacitet                                                                             | Ballongkapacitet                                                                            | Χωρητικότητα μπαλονιού                                                                      | Capacidade do balão                                                                         |
|                | Steriel                                                                                     | Steril                                                                                      | Steril                                                                                      | Αποστειρωμένο                                                                               | Esterilizado                                                                                |
|                | Gesteriliseerd met behulp van ethyleenoxide                                                 | Steriliseret ved brug af ethylenoxid                                                        | Steriliserad med etylenoxid                                                                 | Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου                                              | Esterilizado com óxido de etileno                                                           |
|                | Gesteriliseerd met behulp van straling                                                      | Steriliseret ved brug af bestråling                                                         | Steriliserad med strålning                                                                  | Αποστειρωμένο με χρήση ακτινοβολίας                                                         | Esterilizado por irradiação                                                                 |
|                | Gesteriliseerd met behulp van stoom of droge warmte                                         | Dampsteriliseret eller tørsteriliseret                                                      | Steriliserad med ånga eller torrvarme                                                       | Αποστειρωμένο με χρήση ατμού ή ξηρής θερμότητας                                             | Esterilizado por vapor ou calor seco                                                        |
| <b>Rx only</b> | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|                | Fabrikant                                                                                   | Producent                                                                                   | Tillverkare                                                                                 | Κατασκευαστής                                                                               | Fabricante                                                                                  |
|                | Fabricagedatum                                                                              | Fremstillingsdato                                                                           | Tillverkningsdatum                                                                          | Ημερομηνία κατασκευής                                                                       | Data de fabrico                                                                             |
|                | Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig                          | Indeholder eller har spor af naturlig gummilatex                                            | Innehåller eller har spår av naturgummilatex                                                | Περιέχει λάτεξ φυσικού καουτσούκ                                                            | Contém ou está presente látex de borracha natural                                           |
|                | Afmetingen                                                                                  | Størrelse                                                                                   | Storlek                                                                                     | Μέγεθος                                                                                     | Tamanho                                                                                     |
|                | Raadpleeg de gebruiksaanwijzing op de website                                               | Se brugsanvisningen på webstedet                                                            | Se bruksanvisningen på webbplatsen                                                          | Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web                        | Consultar as instruções de utilização no site                                               |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                   | Nederlands                                                                        | Dansk                                                                     | Svenska                                                                  | Ελληνικά                                                                                    | Português                                                                                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------|--------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
|  eifu.edwards.com<br>+ 1 888 570 4016                                                                                                                                                                                                                                                                                                                             | Volg de gebruiksaanwijzing op de website                                          | Følg brugsanvisningen på hjemmesiden                                      | Följ bruksanvisningen på webbplatsen                                     | Ακολουθείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web                           | Seguir as instruções de utilização no site                                                 |
|                                                                                                                                                                                                                                                                                                                                                                  | MRI-onveilig                                                                      | MR-usikker                                                                | MR-farlig                                                                | Μη ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)                                        | Utilização não segura em ambiente de RM                                                    |
|                                                                                                                                                                                                                                                                                                                                                                  | MRI-veilig                                                                        | MR-sikker                                                                 | MR-säker                                                                 | Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)                                           | Utilização segura em ambiente de RM                                                        |
|                                                                                                                                                                                                                                                                                                                                                                  | MRI-veilig onder bepaalde voorwaarden                                             | MR-betinget                                                               | MR-villkorlig                                                            | Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις                          | Utilização condicionada em ambiente de RM                                                  |
|                                                                                                                                                                                                                                                                                                                                                                  | Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing            | Se brugsanvisningen, eller se elektronisk brugsanvisning                  | Se bruksanvisningen eller den elektroniska bruksanvisningen              | Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης          | Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas |
|                                                                                                                                                                                                                                                                                                                                                                  | Volg de gebruiksaanwijzing                                                        | Følg brugsanvisningen                                                     | Följ bruksanvisningen                                                    | Ακολουθείτε τις Οδηγίες χρήσης                                                              | Seguir as instruções de utilização                                                         |
|                                                                                                                                                                                                                                                                                                                                                                  | Niet opnieuw steriliseren                                                         | Må ikke resteriliseres                                                    | Får inte omsteriliseras                                                  | Μην επαναποστειρώνετε                                                                       | Não voltar a esterilizar                                                                   |
|                                                                                                                                                                                                                                                                                                                                                                  | Niet-pyrogeen                                                                     | Ikke-pyrogen                                                              | Icke-pyrogen                                                             | Μη πυρετογόνο                                                                               | Não pirogénico                                                                             |
|                                                                                                                                                                                                                                                                                                                                                                  | Toegepast onderdeel van type B                                                    | Type B anvendt del                                                        | Patientansluten del av B-typ                                             | Εφαρμοζόμενο εξάρτημα τύπου B                                                               | Peça aplicada de tipo B                                                                    |
|                                                                                                                                                                                                                                                                                                                                                                 | Toegepast onderdeel van type CF                                                   | Anvendt del af typen CF                                                   | Patientansluten del av CF-typ                                            | Εφαρμοζόμενο εξάρτημα τύπου CF                                                              | Peça aplicada de tipo CF                                                                   |
|                                                                                                                                                                                                                                                                                                                                                                | Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing | Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen   | Använd inte om förpackningen är skadad och se bruksanvisningen           | Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης | Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização     |
|                                                                                                                                                                                                                                                                                                                                                                | Open                                                                              | Åben                                                                      | Öppen                                                                    | Ανοικτό                                                                                     | Abrir                                                                                      |
|                                                                                                                                                                                                                                                                                                                                                                 | Zuig 0,5 cc uit de ballon alvorens hem in te brengen of te verwijderen.           | Aspirer ballonen -0,5 cm <sup>3</sup> før indføring eller tilbagetrækning | Aspirera ballongen -0,5 cm <sup>3</sup> före införsel eller utdragning.  | Αναρροφήστε από το μπαλόνι -0,5 cm <sup>3</sup> πριν την εισαγωγή ή την απόσυρση            | Aspirar o balão -0,5 cm <sup>3</sup> antes da introdução ou da retirada.                   |
|                                                                                                                                                                                                                                                                                                                                                                 | Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/Europese Unie         | Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union | Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen | Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/Ευρωπαϊκή Ένωση                      | Representante autorizado na Comunidade Europeia/União Europeia                             |
|                                                                                                                                                                                                                                                                                                                                                                | Conformité Européenne (CE-markering)                                              | Conformité Européenne (CE-mærke)                                          | Conformité Européenne (CE-märke)                                         | Conformité Européenne (Σήμανση CE)                                                          | Conformité Européenne (Marcação CE)                                                        |
|                                                                                                                                                                                                                                                                                                                                                                | Medisch hulpmiddel                                                                | Medicinsk udstyr                                                          | Medicinteknisk produkt                                                   | Ιατροτεχνολογικό προϊόν                                                                     | Dispositivo médico                                                                         |
|                                                                                                                                                                                                                                                                                                                                                                | Importeur                                                                         | Importør                                                                  | Importör                                                                 | Εισαγωγέας                                                                                  | Importador                                                                                 |
|                                                                                                                                                                                                                                                                                                                                                                | Op een koele en droge plaats bewaren                                              | Skal opbevares køligt og tørt                                             | Förvara svalt och torrt.                                                 | Φυλάσσετε σε δροσερό και ξηρό χώρο                                                          | Guardar num local fresco e seco                                                            |
|                                                                                                                                                                                                                                                                                                                                                                | Dubbel steriel barrièresysteem                                                    | Dobbelt sterilt barriersystem                                             | Dubbelt sterilt barriärsystem                                            | Σύστημα διπλού αποστειρωμένου φραγμού                                                       | Sistema de barreira dupla esterilizada                                                     |
|                                                                                                                                                                                                                                                                                                                                                                | Unieke instrumentidentificatie                                                    | Unik udstyrsidentifikation                                                | Unik produktidentifiering                                                | Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος                                      | Identificador único de dispositivo                                                         |
|                                                                                                                                                                                                                                                                                                                                                                | Deze kant boven                                                                   | Denne side opad                                                           | Denna sida upp                                                           | Αυτή η πλευρά προς τα πάνω                                                                  | Este lado para cima                                                                        |
|                                                                                                                                                                                                                                                                                                                                                                | Breekbaar, voorzichtigheid geboden                                                | Skrøbelig, skal håndteres forsigtigt                                      | Ömtåligt, hanteras försiktigt                                            | Εύθραστο: χειριστείτε με προσοχή                                                            | Frágil, manusear com cuidado                                                               |
| <b>Opmerking:</b> Het label van dit product bevat misschien niet alle symbolen. ■ <b>Bemærk:</b> Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ <b>Obs!</b> Alla symboler kanske inte används för märkning av denna produkt. ■ <b>Σημείωση:</b> Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ <b>Nota:</b> poderão não estar incluídos todos os símbolos na rotulagem deste produto. |                                                                                   |                                                                           |                                                                          |                                                                                             |                                                                                            |

Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring

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|                | Suomi                                                                                       | Български                                                                                   | Română                                                                                      | Eesti                                                                                       | Lietuvių                                                                                    |
|----------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|                | Luumenien määrä                                                                             | Брой лумени                                                                                 | Număr de Lumenuri                                                                           | Valendike arv                                                                               | Spindžių skaičius                                                                           |
|                | Ulkohalkaisija                                                                              | Външен диаметър                                                                             | Diametru exterior                                                                           | Välisläbimõõt                                                                               | Išorinis skersmuo                                                                           |
|                | Käyttöpituus                                                                                | Използваема дължина                                                                         | Lungimea utilă                                                                              | Kasutatav pikkus                                                                            | Naudingasis ilgis                                                                           |
|                | Ohjainlangan suositeltu koko                                                                | Препоръчителен размер телен водач                                                           | Mărimea recomandată a firului de ghidaj                                                     | Juhtetraadi soovituslik suurus                                                              | Rekomenduojamas kreipiamosios vielos dydis                                                  |
|                | Luumenin koko                                                                               | Размер на лумена                                                                            | Mărimea lumenului                                                                           | Valendiku suurus                                                                            | Spindžio dydis                                                                              |
|                | Mallinnumero                                                                                | Номер на модела                                                                             | Număr model                                                                                 | Mudeli number                                                                               | Modelio numeris                                                                             |
|                | Sisäänviejän vähimmäiskoko                                                                  | Минимален размер на интродюсера                                                             | Mărimea minimă a dispozitivului de introduce                                                | Sisesti väikseim suurus                                                                     | Minimalus įvediklio dydis                                                                   |
|                | Tärkeä huomautus                                                                            | Внимание                                                                                    | Atenție                                                                                     | Ettevaatust                                                                                 | Perspėjimas                                                                                 |
|                | Ei uudelleenkäytettävä                                                                      | Да не се използва повторно                                                                  | A nu se reutiliza                                                                           | Mitte korduskasutada                                                                        | Nenaudoti pakartotinai                                                                      |
|                | Määrä                                                                                       | Количество                                                                                  | Cantitate                                                                                   | Kogus                                                                                       | Kiekis                                                                                      |
|                | Eränumero                                                                                   | Партиден номер                                                                              | Număr de lot                                                                                | Partii number                                                                               | Partijos numeris                                                                            |
|                | Viimeinen käyttöpäivä                                                                       | Срок на годност                                                                             | A se utiliza până la data de                                                                | Kõlblik kuni                                                                                | Naudoti iki nurodytos datos                                                                 |
|                | Sisäläpimita                                                                                | Вътрешен диаметър                                                                           | Diametru interior                                                                           | Sisemine läbimõõt                                                                           | Vidinis skersmuo                                                                            |
|                | Pallon tilavuus                                                                             | Вместимост на балона                                                                        | Capacitatea balonului                                                                       | Ballooni mahutavus                                                                          | Balionių talpa                                                                              |
|                | Steriili                                                                                    | Стерилно                                                                                    | Steril                                                                                      | Steriilne                                                                                   | Sterilus                                                                                    |
|                | Steriloitu etyleenioksidilla                                                                | Стерилизирано с етиленов оксид                                                              | Sterilizat cu oxid de etilenă                                                               | Steriliseeritud etüleenoksiidi kasutades                                                    | Sterilizuota etileno oksidu                                                                 |
|                | Steriloitu säteilyttämällä                                                                  | Стерилизирано с облъчване                                                                   | Sterilizat prin iradiere                                                                    | Steriliseeritud kiiritamist kasutades                                                       | Sterilizuota švitinant                                                                      |
|                | Lämpö- tai höyrysteriloitu                                                                  | Стерилизирано чрез автоклавиране или сух стерилизатор                                       | Sterilizat cu abur sau prin căldură uscată                                                  | Steriliseeritud auru või kuiva kuumusega                                                    | Sterilizuota naudojant garus arba sausą karštį                                              |
| <b>Rx only</b> | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|                | Valmistaja                                                                                  | Производител                                                                                | Producător                                                                                  | Tootja                                                                                      | Gamintojas                                                                                  |
|                | Valmistuspäivämäärä                                                                         | Дата на производство                                                                        | Data fabricației                                                                            | Tootmiskuupäev                                                                              | Pagaminimo data                                                                             |
|                | Sisältää luonnokumilateksia tai sen jäämiä                                                  | Съдържание или наличие на естествен каучуков латекс                                         | Conține sau prezintă latex din cauciuc natural                                              | Sisaldab naturaalsed kummilateksit                                                          | Turi savo sudėtyje ar yra natūralios gumos latekso                                          |
|                | Koko                                                                                        | Размер                                                                                      | Mărime                                                                                      | Suurus                                                                                      | Dydis                                                                                       |
|                | Katso käyttöohjeet verkkosivustolta                                                         | Направете справка с инструкциите за употреба, намиращи се на уеб сайта                      | Consultați instrucțiunile de utilizare de pe site-ul web                                    | Palun lugege veebisaidil olevaid kasutusjuhiseid                                            | Žr. naudojimo instrukcijas, pateiktas interneto svetainėje                                  |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Suomi                                                           | Български                                                                                    | Română                                                                                                 | Eesti                                                                        | Lietuvių                                                                |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------|----------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------|-------------------------------------------------------------------------|
|  eifu.edwards.com<br>+ 1 888 570 4016                                                                                                                                                                                                                                                                                                                                | Noudata käyttöohjeita, jotka ovat verkkosivustolla              | Спазвайте инструкциите за употреба, дадени на уебсайта                                       | Respectați instrucțiunile de utilizare disponibile pe site-ul web                                      | Järgige kasutusjuhendit, mille leiате veebisaidilt                           | Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje           |
|                                                                                                                                                                                                                                                                                                                                                                     | Ei sovi magneettikuvaukseen                                     | Небезопасно при MR                                                                           | Incompatibil RM                                                                                        | Ohtlik magnetresonantstomograafias                                           | MR nesaugus                                                             |
|                                                                                                                                                                                                                                                                                                                                                                     | Sopii magneettikuvaukseen                                       | Безопасно при MR                                                                             | Sigur în utilizarea la RM                                                                              | Ohutu magnetresonantstomograafias                                            | MR saugus                                                               |
|                                                                                                                                                                                                                                                                                                                                                                     | Ehdollisesti turvallinen magneettikuvauksessa                   | Безопасно при MR при определени условия                                                      | Condiționat RM                                                                                         | Ohutu magnetresonantstomograafias teatud tingimustel                         | Sąlyginis MR                                                            |
|                                                                                                                                                                                                                                                                                                                                                                     | Katso käyttöohjeet tai sähköiset käyttöohjeet                   | Направте справка с инструкциите за употреба или електронните инструкции за употреба          | Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic | Tutvuge paberandjal või elektroonilise kasutusjuhisega                       | Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas     |
|                                                                                                                                                                                                                                                                                                                                                                     | Noudata käyttöohjeita                                           | Спазвайте инструкциите за употреба                                                           | Respectați instrucțiunile de utilizare                                                                 | Järgige kasutusjuhendit                                                      | Vadovautis naudojimo instrukcijomis                                     |
|                                                                                                                                                                                                                                                                                                                                                                     | Älä steriloi uudelleen                                          | Не стерилизирайте повторно                                                                   | A nu se resteriliza                                                                                    | Ärge resteriliseerige                                                        | Kartotinai nesterilizuokite                                             |
|                                                                                                                                                                                                                                                                                                                                                                     | Pyrogeeniton                                                    | Непирогенно                                                                                  | Non pirogen                                                                                            | Mittepürogeenne                                                              | Nepirogeniškas                                                          |
|                                                                                                                                                                                                                                                                                                                                                                     | B-tyypin liityntäosa                                            | Приложна част тип B                                                                          | Componentă aplicată de tip B                                                                           | B-tüüpi kohaldatav osa                                                       | B tipo liečiamoji dalis                                                 |
|                                                                                                                                                                                                                                                                                                                                                                   | CF-tyypin liityntäosa                                           | Приложна част тип CF                                                                         | Componentă aplicată de tip CF                                                                          | CF-tüüpi kohaldatav osa                                                      | CF tipo liečiamoji dalis                                                |
|                                                                                                                                                                                                                                                                                                                                                                   | Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet | Да не се използва, ако опаковката е повредена, и направте справка с инструкциите за употреба | A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare                 | Ärge kasutage, kui pakend on kahjustatud, ja lugege kasutusjuhiseid          | Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas          |
|                                                                                                                                                                                                                                                                                                                                                                   | Auki                                                            | Отворено                                                                                     | Deschidere                                                                                             | Ava                                                                          | Atidarymas                                                              |
|                                                                                                                                                                                                                                                                                                                                                                    | Aspiroi pallosta -0,5 ml ennen sisäänvientiä tai ulosvetämistä. | Аспирирайте балона с -0,5 cc преди въвеждане или изтегляне                                   | Aspirați din balon 0,5 cc înainte de inserție sau de retragere                                         | Enne sisestamist või eemaldamist aspiireerige balloon -0,5 cm <sup>3</sup> . | Prieš įvesdami arba ištraukdami, iš balionėlio išsiurbkite -0,5 ml oro. |
|                                                                                                                                                                                                                                                                                                                                                                    | Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa   | Упълномощен представител в Европейската общност / Европейския съюз                           | Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană                                      | Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus                        | Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje             |
|                                                                                                                                                                                                                                                                                                                                                                   | Conformité Européenne (CE-merkintä)                             | Conformité Européenne (CE маркировка)                                                        | Conformité Européenne (marcay CE)                                                                      | Conformité Européenne (CE-märk)                                              | Conformité Européenne (CE žymė)                                         |
|                                                                                                                                                                                                                                                                                                                                                                   | Lääkinnällinen laite                                            | Медицинско изделие                                                                           | Dispozitiv medical                                                                                     | Meditsiiniseade                                                              | Medicinos priemonė                                                      |
|                                                                                                                                                                                                                                                                                                                                                                   | Maahantuoja                                                     | Вносител                                                                                     | Importator                                                                                             | Importija                                                                    | Importuotojas                                                           |
|                                                                                                                                                                                                                                                                                                                                                                   | Säilytettävä kuivassa ja viileässä                              | Да се съхранява на хладно, сухо място                                                        | Stocați într-un loc rece și uscat                                                                      | Säilitage jahedas ja kuivas kohas.                                           | Laikyti vėsioje, sausoje vietoje                                        |
|                                                                                                                                                                                                                                                                                                                                                                   | Kaksinkertainen steriili sulkujärjestelmä                       | Двойна стерилна бариерна система                                                             | Sistem cu barieră sterilă dublă                                                                        | Kahekordne steriilne kattesüsteem                                            | Dvigubo sterilaus barjero sistema                                       |
|                                                                                                                                                                                                                                                                                                                                                                   | Laitteen yksilöivä tunnus                                       | Уникален идентификатор на изделието                                                          | Identificator unic al dispozitivului                                                                   | Seadme kordumatu identifitseerimistunnus                                     | Unikalasis priemonės identifikatorius                                   |
|                                                                                                                                                                                                                                                                                                                                                                   | Tämä puoli ylöspäin                                             | С тази страна нагоре                                                                         | Cu partea aceasta în sus                                                                               | See pool üles                                                                | Šia puse į viršų                                                        |
|                                                                                                                                                                                                                                                                                                                                                                   | Särkyvä, käsiteltävä varoen                                     | Чупливо, да се борава внимателно                                                             | Fragil, a se manipula cu atenție                                                                       | Ūrn, käsitseda ettevaatlikult                                                | Trapus, elgtis atsargiai                                                |
| <b>Huomautus:</b> kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ <b>Забележка:</b> Възможно е не всички символи да са включени в етикета на този продукт.<br><b>Notă:</b> este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ <b>Märkus:</b> Koik simboliid ei pruugi antud toote etiketil esineda. ■ <b>Pastaba:</b> Šio gaminio etiketėje gali būti pateikti ne visi simboliai. |                                                                 |                                                                                              |                                                                                                        |                                                                              |                                                                         |

**Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola**

|                | Latviešu                                                                                    | Türkçe                                                                                      | Русский                                                                                     | Srpski                                                                                      | Hrvatski                                                                                    |
|----------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|                | Lūmenu skaits                                                                               | Lümen sayısı                                                                                | Количество просветов                                                                        | Broj lumena                                                                                 | Broj lumena                                                                                 |
|                | Ārējais diametrs                                                                            | Dış çap                                                                                     | Внешний диаметр                                                                             | Spoljašnji prečnik                                                                          | Vanjski promjer                                                                             |
|                | Izmantojamais garums                                                                        | Kullanılabilir uzunluk                                                                      | Рабочая длина                                                                               | Upotrebljiva dužina                                                                         | Upotrebljiva dužina                                                                         |
|                | Ieteicamais vadītājstīgas izmērs                                                            | Kılavuz Tel İçin Önerilen Boyut                                                             | Рекомендованный размер проводника                                                           | Preporučena veličina žice vodice                                                            | Preporučena veličina žice vodilice                                                          |
|                | Lūmena izmērs                                                                               | Lümen boyutu                                                                                | Размер просвета                                                                             | Veličina lumena                                                                             | Veličina lumena                                                                             |
|                | Modeļa numurs                                                                               | Model Numarası                                                                              | Номер модели                                                                                | Broj modela                                                                                 | Broj modela                                                                                 |
|                | Minimālais ievadītāja izmērs                                                                | Minimum introdüser boyutu                                                                   | Минимальный размер интродьюсера                                                             | Minimalna veličina uvodnika                                                                 | Minimalna veličina uvodnog instrumenta                                                      |
|                | Uzmanību!                                                                                   | Dikkat                                                                                      | Предостережение                                                                             | Oprez                                                                                       | Oprez                                                                                       |
|                | Nelietot atkārtoti                                                                          | Yeniden kullanmayın                                                                         | Не подлежит повторному использованию                                                        | Ne koristiti ponovo                                                                         | Nemojte ponovno upotrebljavati                                                              |
|                | Daudzums                                                                                    | Miktar                                                                                      | Количество                                                                                  | Količina                                                                                    | Količina                                                                                    |
|                | Partijas numurs                                                                             | Lot Numarası                                                                                | Номер партии                                                                                | Broj partije                                                                                | Broj serije                                                                                 |
|                | Derīguma termiņš                                                                            | Son Kullanma Tarihi                                                                         | Использовать до                                                                             | Rok korišćenja                                                                              | Rok upotrebe                                                                                |
|                | Iekšējais diametrs                                                                          | İç çap                                                                                      | Внутренний диаметр                                                                          | Unutrašnji prečnik                                                                          | Unutrašnji promjer                                                                          |
|                | Balona ietilpība                                                                            | Balon Kapasitesi                                                                            | Рабочий объем баллона                                                                       | Kapacitet balona                                                                            | Kapacitet balona                                                                            |
|                | Sterils                                                                                     | Steril                                                                                      | Стерильно                                                                                   | Sterilno                                                                                    | Sterilno                                                                                    |
|                | Sterilizēts ar etilēna oksīdu                                                               | Etilen oksit kullanılarak sterilize edilmiştir                                              | Стерилизовано этиленоксидом                                                                 | Sterilizovano etilen-oksidom                                                                | Sterilizirano etilen-oksidom                                                                |
|                | Sterilizēts apstarojot                                                                      | Işınlama yoluyla sterilize edilmiştir                                                       | Стерилизовано излучением                                                                    | Sterilizovano zračenjem                                                                     | Sterilizirano zračenjem                                                                     |
|                | Sterilizēts, izmantojot tvaiku vai sausu karstumu                                           | Buhar veya kuru ısı kullanılarak sterilize edilmiştir                                       | Стерилизовано паром или сухим жаром                                                         | Sterilizovano korišćenjem pare ili suvom toplotom                                           | Sterilizirano parom ili vrućim, suhim zrakom                                                |
| <b>Rx only</b> | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|                | Ražotājs                                                                                    | Üretici                                                                                     | Производитель                                                                               | Proizvođač                                                                                  | Proizvođač                                                                                  |
|                | Izgatavošanas datums                                                                        | Üretim tarihi                                                                               | Дата производства                                                                           | Datum proizvodnje                                                                           | Datum proizvodnje                                                                           |
|                | Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā                                       | Doğal kauçuk lateks içerir veya mevcuttur                                                   | Содержит натуральный латекс или изготовлено из натурального латекса                         | Sadrži ili ima prisustva prirodne lateks gume                                               | Sadrži prirodni gumeni lateks ili je prisutan                                               |
|                | Izmērs                                                                                      | Boyut                                                                                       | Размер                                                                                      | Veličina                                                                                    | Veličina                                                                                    |
|                | Skatīt lietošanas instrukciju timeklā vietnē                                                | Web sitesindeki kullanım talimatlarına başvurun                                             | См. инструкции по применению на веб-сайте                                                   | Pogledajte uputstva za upotrebu na veb-sajtu                                                | Pogledajte upute za upotrebu na web-stranici                                                |

|                                                                                                                                                                                                                                                                                                                                                                                                 | Latviešu                                                                                | Türkçe                                                                            | Русский                                                                                    | Srpski                                                                         | Hrvatski                                                                              |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  eifu.edwards.com<br>+ 1 888 570 4016                                                                                                                                                                                                                                                                           | Ievērojiet lietošanas instrukciju,<br>kas pieejama tīmekļa vietnē                       | Web sitesinde bulunan kullanım<br>talimatlarını takip edin                        | Следуйте инструкциям по<br>применению на веб-сайте                                         | Pridržavajte se uputstva za<br>upotrebu sa veb-sajta                           | Pridržavajte se uputa za upotrebu<br>na internetskoj stranici                         |
|                                                                                                                                                                                                                                                                                                                | Nedrīkst lietot MR vidē                                                                 | MR için Güvenli Değildir                                                          | Опасно при проведении МРТ                                                                  | Nije bezbedno za MR                                                            | Nije sigurno kod pregleda MR-om                                                       |
|                                                                                                                                                                                                                                                                                                                | Drošs lietošanai MR vidē                                                                | MR için Güvenli                                                                   | Безопасно при проведении МРТ                                                               | Bezbedno za MR                                                                 | Sigurno kod pregleda MR-om                                                            |
|                                                                                                                                                                                                                                                                                                                | Izmantojams MR vidē, ievērojot<br>noteiktus nosacījumus                                 | MR Koşullu                                                                        | Условно безопасно при<br>проведении МРТ                                                    | Uslovno bezbedno za MR                                                         | Uvjetno siguran kod pregleda<br>MR-om                                                 |
|                                                                                                                                                                                                                                                                                                                | Skatiet lietošanas instrukcijas vai<br>elektroniskās lietošanas<br>instrukcijas         | Kullanım talimatlarına veya<br>elektronik kullanım talimatlarına<br>başvurun      | См. инструкции по применению<br>либо электронные инструкции<br>по применению               | Pogledajte uputstvo za upotrebu ili<br>elektronsko uputstvo za upotrebu        | Pogledajte upute za upotrebu ili<br>elektroničke upute za upotrebu                    |
|                                                                                                                                                                                                                                                                                                                | Ievērojiet lietošanas instrukciju                                                       | Kullanım talimatlarını takip edin                                                 | Следуйте инструкциям по<br>применению                                                      | Pridržavajte se uputstva za<br>upotrebu                                        | Pridržavajte se uputa za upotrebu                                                     |
|                                                                                                                                                                                                                                                                                                                | Nesterilizēt atkārtoti                                                                  | Yeniden sterilize etmeyin                                                         | Не стерилизовать повторно                                                                  | Nemojte ponovo sterilisati                                                     | Nemojte ponovno sterilizirati                                                         |
|                                                                                                                                                                                                                                                                                                                | Nepirogēns                                                                              | Nonpirogeniktir                                                                   | Апирогенно                                                                                 | Apirogeno                                                                      | Nepirogeno                                                                            |
|                                                                                                                                                                                                                                                                                                                | B tipa lietotā daļa                                                                     | B Tipi uygulanan parça                                                            | Рабочая часть типа B                                                                       | Primenjeni deo tipa B                                                          | Primjenjen dio tipa B                                                                 |
|                                                                                                                                                                                                                                                                                                               | CF tipa kontaktējošā daļa                                                               | CF Tipi uygulanan parça                                                           | Рабочая часть типа CF                                                                      | Primenjeni deo tipa CF                                                         | Primijenjen dio vrste CF                                                              |
|                                                                                                                                                                                                                                                                                                              | Nelietot, ja iepakojums ir bojāts,<br>un skatīt lietošanas instrukciju                  | Ambalaj hasarlıysa kullanmayın ve<br>kullanım talimatlarına başvurun              | Не использовать, если упаковка<br>повреждена; использовать в<br>соответствии с инструкцией | Ne koristite ako je pakovanje<br>oštećeno i pogledajte uputstvo za<br>upotrebu | Nemojte upotrebljavati ako je<br>pakiranje oštećeno i pogledajte<br>upute za upotrebu |
|                                                                                                                                                                                                                                                                                                              | Atvērt                                                                                  | Açık                                                                              | Открыть                                                                                    | Otvoreno                                                                       | Otvoreno                                                                              |
|                                                                                                                                                                                                                                                                                                              | Pirms balona ievadišanas vai<br>izņemšanas atsūknējiet balonu<br>-0,5 cm <sup>3</sup> . | Uygulama veya Geri Çekme<br>Öncesinde Balonu -0,5 cm <sup>3</sup> Aspire<br>Edin. | Перед введением или<br>извлечением удалите из<br>баллона 0,5 куб. см<br>содержимого.       | Aspirirajte balon -0,5 ml pre<br>uvođenja ili izvlačenja                       | Aspirirajte balon -0,5 cm <sup>3</sup> prije<br>uvođenja ili izvlačenja               |
|                                                                                                                                                                                                                                                                                                              | Pilnvarotais pārstāvis Eiropas<br>Kopienā/Eiropas Savienībā                             | Avrupa Birliği'nde/Avrupa<br>Topluluğu'nda Yetkili Temsilci                       | Официальный представитель в<br>Европейском сообществе/<br>Европейском союзе                | Ovlašćeni predstavnik u Evropskoj<br>zajednici / Evropskoj Uniji               | Ovlašteni predstavnik u Evropskoj<br>zajednici / Europskoj uniji                      |
|                                                                                                                                                                                                                                                                                                              | Conformité Européenne (CE<br>marķējums)                                                 | Conformité Européenne (CE İşareti)                                                | Conformité Européenne<br>(маркировка CE)                                                   | Conformité Européenne (CE znak)                                                | Conformité Européenne<br>(oznaka CE)                                                  |
|                                                                                                                                                                                                                                                                                                              | Medicīniska ierīce                                                                      | Tıbbi cihaz                                                                       | Медицинское устройство                                                                     | Medicinsko sredstvo                                                            | Medicinski uređaj                                                                     |
|                                                                                                                                                                                                                                                                                                              | Importētājs                                                                             | İthalatçı                                                                         | Импортер                                                                                   | Uvoznik                                                                        | Uvoznik                                                                               |
|                                                                                                                                                                                                                                                                                                              | Uzglabāt vēsā, sausā vietā                                                              | Serin, kuru yerde saklayın                                                        | Хранить в прохладном, сухом<br>месте                                                       | Čuvati na hladnom i suvom mestu                                                | Pohranite na hladnom, suhom<br>mjestu                                                 |
|                                                                                                                                                                                                                                                                                                              | Divkārsa sterilitātes aizsargslāņa<br>sistēma                                           | Çift steril bariyer sistemi                                                       | Двойная барьерная система для<br>стерилизации                                              | Sistem sa dve sterilne barijere                                                | Sustav dvostruke sterilne barijere                                                    |
|                                                                                                                                                                                                                                                                                                              | Unikāls ierīces identifikators                                                          | Benzersiz cihaz tanımlayıcısı                                                     | Уникальный идентификатор<br>устройства                                                     | Jedinstveni identifikator sredstva                                             | Jedinstveni identifikator proizvoda                                                   |
|                                                                                                                                                                                                                                                                                                              | Ar šo pusi uz augšu                                                                     | Bu taraf yukarı                                                                   | Верхняя сторона                                                                            | Ova strana mora da bude<br>prema gore                                          | Ova strana prema gore                                                                 |
|                                                                                                                                                                                                                                                                                                              | Trausls, rīkoties piesardzīgi                                                           | Kırılabilir, dikkatli taşıyın                                                     | Хрупкое изделие, обращаться с<br>осторожностью                                             | Lomljivo, rukovati sa oprezom                                                  | Lomljivo, pažljivo rukovati                                                           |
| <b>Piezīme.</b> Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ <b>Not:</b> Bu ürünün etiketinde tüm semboller yer almayış olabilir. ■ <b>Примечание.</b> На этикетках данного изделия могут указываться не все символы. ■ <b>Napomena:</b> Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ <b>Napomena:</b> na oznakama ovog proizvoda ne moraju se nalaziti svi simboli. |                                                                                         |                                                                                   |                                                                                            |                                                                                |                                                                                       |

**Edwards Lifesciences GmbH**  
Parking 30  
85748 Garching bei München  
Germany

10/24  
10048653001 A / DOC-0230850 A  
© Copyright 2024, Edwards Lifesciences LLC  
All rights reserved.

  
**Edwards Lifesciences LLC**  
One Edwards Way  
Irvine, CA 92614 USA  
Made in USA

 **Edwards Lifesciences B.V.**  
Verlengde Poolseweg 16  
4818 CL Breda, Netherlands

Telephone 949.250.2500  
800.424.3278  
FAX 949.250.2525

Web IFU