Edwards


Directory

English (en).....	1	Slovensky (sk).....	28
Français (fr).....	4	Norsk (no).....	30
Deutsch (de).....	6	Suomi (fi).....	32
Español (es).....	8	Български (bg).....	34
Italiano (it).....	10	Română (ro).....	36
Nederlands (nl).....	12	Eesti (et).....	38
Dansk (da).....	14	Lietuvių (lt).....	40
Svenska (sv).....	16	Latviešu (lv).....	42
Ελληνικά (el).....	18	Türkçe (tr).....	44
Português (pt).....	20	Русский (ru).....	46
Česky (cs).....	22	Srpski (sr).....	48
Magyar (hu).....	24	Hrvatski (hr).....	50
Polski (pl).....	26		

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola52

English

HemoSphere Pressure-Out Cable

Edwards, Edwards Lifesciences, the stylized E logo, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device. Refer to the HemoSphere advanced monitor operator's manual, available at eifu.edwards.com, for comprehensive monitoring procedures.

1.0 Intended Use

The HemoSphere pressure-out cable is a reusable accessory that enables the user to output pressure monitored by the HemoSphere advanced monitoring platform to a compatible patient monitor for standard pressure monitoring. Up to three waveforms can be outputted: arterial pressure (AP, red), central venous pressure (CVP, blue), and pulmonary artery pressure (PAP, yellow).

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

2.0 Warnings


This product contains metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

3.0 Precautions

Inspect the pressure-out cable before each use. Do not use if the cable is damaged. Damage may include broken insulation, frays, deteriorated strain relief, bent/recessed connector pins or any signs that the cable may be compromised.

4.0 HemoSphere Pressure-Out Cable Setup

Refer to the numbered figure and to the HemoSphere advanced monitor operator's manual, available at eifu.edwards.com, for comprehensive monitoring procedures with that monitor.

Step	Procedure
1	Plug the 18-pin connector ① of the pressure-out cable into the rear panel of the HemoSphere advanced monitor. The female port is denoted by an analog output symbol. 

Step	Procedure
2	Use the two jack screws ② to properly secure the pressure-out cable plug in place.
3	Connect the desired pressure signal plug(s) into the compatible patient monitor: ③ arterial pressure (AP, red) ④ pulmonary artery pressure (PAP, yellow) ⑤ central venous pressure (CVP, blue) Ensure that the selected connector is fully engaged. Refer to patient monitor instructions for use.
4	Refer to the HemoSphere advanced monitor operator's manual for instructions on configuring the desired pressure signal output and for proper zeroing instructions.

5.0 MRI Information

MR Unsafe

Do not use the HemoSphere pressure-out cable in an MR environment. The cable is MR unsafe since the accessory contains metallic components, which can experience RF-induced heating in the MRI environment.

6.0 Operating Conditions

Temperature: 10 to 32.5 °C

Humidity range: 20 to 90% non-condensing

Altitude (Atmospheric Pressure): 0 m/0 ft (1013 hPa) to 3048 m/10,000 ft (697 hPa)

7.0 Maintenance

It is recommended to inspect the HemoSphere pressure-out cable periodically for signs of damage. Check the condition of the cable and replace if any of the following are present: exposed wire, splits, cracks, or signs of stress.

8.0 Cleaning

The HemoSphere pressure-out cable is a reusable accessory and can be surface cleaned per hospital aseptic procedures. Accessory has been shown to handle cloth wipe down using the following chemical agents: 70% isopropyl alcohol solution, 2% glutaraldehyde, 10% bleach solution, or quaternary ammonium solution.

CAUTIONS:

Do not use any other cleaning agents, spray, or pour cleaning solution directly on the cable.

Do not steam, radiate, or EtO sterilize the cable.

Do not immerse.

9.0 Storage

Store in a cool, dry place in original packaging to prevent damage.

CAUTION: Do not coil the cable against its natural orientation or coil the cable tightly.

10.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours):800.822.9837

Outside the U.S. and Canada (24 hours): ...949.250.2222

In the UK:0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

11.0 Notices

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Prices, specifications and model availability are subject to change without notice.

12.0 Disposal

To avoid contaminating or infecting personnel, the environment or other equipment, make sure the HemoSphere pressure-out cable is disinfected and decontaminated appropriately in accordance with your country's laws for equipment containing electrical and electronic parts prior to disposal.

Dispose of in accordance with national and local legislation and requirements for equipment containing electrical and electronic parts.

Refer to the symbol legend at the end of this document.

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	English	Français	Deutsch	Español
 eifu.edwards.com + 1 888 570 4016	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web
	Serial Number	Numéro de série	Seriennummer	Número de serie
	Manufacturer	Fabricant	Hersteller	Fabricante
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación
	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM
IPX1	Drip proof equipment	Équipement anti-gouttes	Tropfwassergeschütztes Gerät	Equipo a prueba de goteo
IPX4	Extent of protection against ingress of objects	Extension de la protection contre l'entrée d'objets	Erweiterter Schutz gegen feste Fremdkörper	Grado de protección contra la entrada de objetos
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Defibrillation Proof Type CF applied part	Partie appliquée de type CF résistante à la défibrillation	Defibrillationsgeschütztes Anwendungsteil vom Typ CF	Pieza aplicada de tipo CF a prueba de desfibrilación
	Do not disassemble	Ne pas démonter	Nicht zerlegen	No desmontar
	Lead-free	Sans plomb	Bleifrei	Sin plomo
	Importer	Importateur	Importeur	Importador
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario
	Restriction of hazardous substances (RoHS) – China only	RoHS (Restriction sur l'usage de certaines substances dangereuses) – Chine uniquement	Vermeidung gefährlicher Substanzen (RoHS) – nur China	Restricción de Sustancias Peligrosas (RoHS); solo en China

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos

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	Italiano	Nederlands	Dansk	Svenska
 eifu.edwards.com + 1 888 570 4016	Seguire le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen
	Numero di serie	Serienummer	Serienummer	Serienummer
	Produttore	Fabrikant	Producent	Tillverkare
	Data di produzione	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum
	Rappresentante autorizzato nella Comunità Europea/ Unione Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/ Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen
	Non compatibile con RM	MRI-onveilig	MR-usikker	MR-farlig
IPX1	Attrezzatura a prova di gocciolamento	Druipwaterdichte apparatuur	Drypsikkert udstyr	Droppskyddad utrustning
IPX4	Protezione dall'ingresso di oggetti	Mate van bescherming tegen het binnendringen van deeltjes	Skal beskyttes mod indtrængen af genstande	Graden av skydd mot intrång från objekt
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Parte applicata di tipo CF a prova di defibrillazione	Defibrillatiebestendig toegepast onderdeel van type CF	Defibrilleringssikker type CF-anvendt del	Defibrilleringssäker typ CF tillämpad del
	Non smontare	Niet uit elkaar nemen	Må ikke adskilles	Får inte monteras isär
	Senza piombo	Loodvrij	Blyfri	Blyfri
	Importatore	Importeur	Importør	Importör
	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt
	Limitazioni sulle sostanze pericolose (RoHS) - solo Cina	Beperking op het gebruik van gevaarlijke stoffen (RoHS) - Alleen China	RoHS (Restriction of Hazardous Substances) – kun Kina	Restriction of Hazardous Substances (RoHS) – endast Kina

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	Български	Română	Eesti	Lietuvių
 eifu.edwards.com + 1 888 570 4016	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje
	Сериен номер	Număr de serie	Seerianumber	Serijos numeris
	Производител	Producător	Tootja	Gamintojas
	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Упълномощен представител в Европейската общност/ Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/ Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Небезопасно при MR	Incompatibil RM	Ohtlik magnetresonantstomograafias	MR nesaugus
IPX1	Защитено от капки оборудване	Echipament protejat împotriva picăturilor de apă	Tilkumiskindel seade	Įranga su apsauga nuo pūslių
IPX4	Степен на защита срещу проникване на предмети	Grad de protecție împotriva pătrunderii obiectelor	Kaitsmise määr objektide sissetungi eest	Objektų patekimo į prietaisą apsaugos lygis
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Приложна част от тип CF със защита от дефибрилация	Piesă aplicată de tip CF, rezistentă la defibrilare	Defibrillatsioonikindel CF-tüüpi rakendusosa	Defibriliavimo poveikiui atspari CF tipo darbinė dalis
	Не разглобявайте	A nu se demonta	Ärge demonteerige	Neardyti
	Не съдържа олово	Nu conține plumb	Pliivaba	Be švino
	Вносител	Importator	Importija	Importuotojas
MD	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Ограничение за употребата на определени опасни вещества (RoHS) – само за Китай	Restricții de utilizare a substanțelor periculoase (RoHS) – numai China	Vastab ohtlike ainete kasutamist piiravale direktiivile Restriction of Hazardous Substance (RoHS) – ainult Hiina	Atitiktis pavojingų medžiagų naudojimo apribojimams (RoHS) – tik Kinijoje

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03/25
10020061004 A / DOC-0232370 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in China

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Conformité Européenne (oznaka CE) tijela TÜV SÜD Product Service GmbH (prijavljeno tijelo)
Edwards Lifesciences GmbH