Edwards

Directory

English (en).....	1
Français (fr).....	3
Deutsch (de).....	5
Español (es).....	7
Italiano (it).....	9
Nederlands (nl).....	11
Dansk (da).....	13
Svenska (sv).....	15
Ελληνικά (el).....	17
Português (pt).....	19
Česky (cs).....	21
Magyar (hu).....	23
Polski (pl).....	25
Slovensky (sk).....	27
Norsk (no).....	29
Suomi (fi).....	31
Български (bg).....	33
Română (ro).....	35
Eesti (et).....	37
Lietuvių (lt).....	38
Latviešu (lv).....	40
Türkçe (tr).....	42
Русский (ru).....	44
Srpski (sr).....	46
Hrvatski (hr).....	48
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike.....	50
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola	54

English

Roll Stand Assembly and Operation

REF: HEMRLSTD1000

For figures, refer to Figure 1 on page 50 through Figure 8 on page 53.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 INTENDED USE

This roll stand is intended to be used with the HemoSphere advanced monitor. See HemoSphere advanced monitor operator's manual for more details.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

Installation Note:

Separate installation kits and accessories may be included with the roll stand. Read all separate installation guides before assembling the roll stand or mounting instruments or accessories.

Tools Required:

- Phillips screwdriver (not provided)
- 5/32" [4 mm] hex wrench (provided)
- 2.5 mm hex wrench (provided)
- 1/2" [13 mm] wrench (not provided)

2.0 Attaching Roll Stand Post to the Base

Step	Procedure
1	Insert the post in the base and lay assembly on its side for access to bottom of the base (see Figure 1 on page 50).
2	Using a 1/2" [13 mm] wrench, attach the post to the base with one (1) 5/16-18 x 1" HHCS, 5/16 split lock washer, and 5/16" flat washer.

3.0 Attaching the Basket to the Post

Step	Procedure
1	Install both nylon clips at an angle, in the center of the basket. Rotate clips until both are horizontal (see Figure 2 on page 50).
2	Slide the clips to the top and bottom of basket. The top clip should be touching the top horizontal wire, and the bottom clip should be touching the lower reinforcing wire (see Figure 3 on page 51).
3	Slide the clips over the top of the post (see Figure 4 on page 51). Position the basket at the desired height. Tighten the clips using the supplied #8-32 machine screws and hex nuts (see Figure 5 on page 52).

4.0 Attaching Mounting Plate Assembly

Step	Procedure
1	Insert mounting plate assembly in top of post (see Figure 6 on page 52). Press downward while twisting until the assembly snaps into place.
2	Rotate the mounting plate until the spring-loaded plunger (indicated by ① in Figure 6 on page 52) is toward the front of the roll stand.

5.0 Attaching Roll Stand Handle

Fasten handle to roll stand post with two (2) #10-32 x 9/16" PHMS (see Figure 7 on page 53). Alternately tighten handle screws until handle does not rotate around post.

Note: Handle position shown is typical but handle may be attached at any desired location on the post.

6.0 Periodic Maintenance and Cleaning

All fasteners associated with the mounting system should be inspected periodically and tightened as necessary. Roll stand may be cleaned with any non-abrasive solution as approved by the facility in which the product will be used.

7.0 Warnings

- DO NOT use this product without first reading and understanding the instructions contained in this booklet. If you are unable to understand the Warnings, Caution or Instructions, contact a healthcare professional, dealer or technical personnel before use – otherwise serious bodily injury or product damage may occur.
- All wheels must be in contact with the floor at ALL TIMES.
- DO NOT use the roll stand while walking backwards down gradients or climbing stairs, curbs or to go over obstacles. Serious risk of fall or injury may occur.
- Engage all wheel locks when roll stand is not in motion.
- Do not place items in the utility basket exceeding 4.5 kg (10 lbs).
- Use care when moving roll stand with a patient. Avoid wet or slick surfaces, wires, cables, or tubes and any uneven surfaces that may cause roll stand to tip or fall.
- Only use accessories and spare parts authorized by GCX Corporation.
- DO NOT attach roll stand to other objects.
- DO NOT use in MRI areas.
- DO NOT use the roll stand for support or as a walking aid. Do not attempt to support your weight or use as a stand assist.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

8.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 821 1012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications and model availability are subject to change without notice.

9.0 Disposal

Dispose of in accordance with hospital policy and local regulations.

Refer to the symbol legend at the end of this document.

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
	Nederlands	Dansk	Svenska	Ελληνικά	Português
QTY	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
SN	Serienummer	Serienummer	Serienummer	Αριθμός σειράς	Número de série
	Massa van rolstandaard	Rullestanders vægt	Rullstativets massa	Βάρος τροχήλατης βάσης	Massa do suporte de rolo
	Gescheiden inzameling voor elektrische en elektronische apparatuur in overeenstemming met EG-richtlijn 2012/19/EU	Separat indsamling af elektrisk og elektronisk udstyr iht. EF-direktivet 2012/19/EU	Separat insamling för elektrisk och elektronisk utrustning i enlighet med EU-direktiv 2012/19/EU	Ξεχωριστή συλλογή για ηλεκτρικό και ηλεκτρονικό εξοπλισμό σύμφωνα με την Οδηγία του Ευρωπαϊκού Κοινοβουλίου 2012/19/ΕΕ	Recolha em separado para o equipamento elétrico e eletrônico, de acordo com a Diretiva CE 2012/19/UE.
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
 eifu.edwards.com +1 888 570 4016	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen	Ακολουθείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Seguir as instruções de utilização no site
EC REP	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/ Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/União Europeia
	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen	Ακολουθείτε τις Οδηγίες χρήσης	Seguir as instruções de utilização
 + 	Massa van rolstandaard + veilige werklust	Rullestanders vægt + sikker arbejdsbelastning	Rullstativets massa + säker belastning	Βάρος τροχήλατης βάσης + ασφαλές φορτίο εργασίας	Massa do suporte de rolo + carga de trabalho segura
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Importeur	Importør	Importör	Εισαγωγέας	Importador

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola



	Nederlands	Dansk	Svenska	Ελληνικά	Português
MD	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
#	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
UDI	Unieke instrument-identificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Μοναδικό αναγνωριστικό συσκευής	Identificador único de dispositivo

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
	Latviešu	Türkçe	Русский	Srpski	Hrvatski
 eifu.edwards.com + 1 888 570 4016	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na web-stranici
QTY	Daudzums	Miktar	Количество	Količina	Količina
SN	Sērijas numurs	Seri Numarası	Серийный номер	Serijski broj	Serijski broj
	Statīva uz ritenīšiem svārs	Yuvarlak sehpanın kütlesi	Масса передвижной стойки	Masa stalka sa točkičima	Masa držača role
	Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES Direktīvu 2012/19/EU	2012/19/EU sayılı EC direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır	Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU	Odvojeno sakupljanje električne i elektronske opreme u skladu s direktivom EZ 2012/19/EU	Odvojeno prikupljanje električne i elektroničke opreme u skladu s direktivom EZ-a 2012/19/EU
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
 eifu.edwards.com + 1 888 570 4016	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
EC REP	Pilnvarotais pārstāvis Eiropas Kopienā/ Eiropas Savienībā	Avrupa Birliği'nde/ Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu
	Statīva uz ritenīšiem svārs + droša darba slodze	Yuvarlak sehpanın kütlesi + güvenli iş yükü	Масса передвижной стойки + безопасная рабочая нагрузка	Masa stalka sa točkičima + bezbedno opterećenje pri radu	Masa držača role + sigurno radno opterećenje
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykiety niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei vältämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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02/25
10047205002 A / DOC-0233092 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH