Edwards

DIRECTORY

English	1
Français	4
Deutsch	7
Español	10
Italiano	13
Nederlands	16
Dansk	19
Svenska	22
Ελληνικά	25
Português	28
Česky	31

Magyar	34
Polski	37
Slovensky	40
Norsk	43
Suomi	46
Български	49
Română	52
Eesti	55
Lietuvių	58
Latviešu	61
Türkçe	64

Русский	67
Srpski	71
Hrvatski	74


Symbol legend ■ Légende des symboles ■ Zeichenerklärung
■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst
met symbolen ■ Symbolforklaring ■ Symbolförklaring ■
Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se
symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky
k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■
Легенда на символите ■ Legendă de simboluri ■ Sümbolite
seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■
Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa
simbolima ■ Legenda simbola

English

FloTrac Jr sensor

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

T4076

FloTrac Jr sensor

1. To Patient
2. Sensor Vent Port
3. Snap-Tab Device
4. To IV Set
5. To Pressure monitor (red)
6. To FloTrac Jr sensor compatible monitor (green)

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

These are general instructions for setting up a pressure monitoring system and/or Edwards compatible hardware with the FloTrac Jr sensor. Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

CAUTION: The use of lipids with the FloTrac Jr sensor may compromise product integrity.

1.0 Concept/Description

The FloTrac Jr sensor is a sterile, single use kit that monitors pressures when attached to pressure monitoring catheters. When connected to a compatible monitor, the FloTrac Jr sensor minimally-invasively measures cardiac output and key hemodynamic parameters, which assist the clinician in assessing the patient's physiologic status and support clinical decisions related to hemodynamic optimization. They are intended to be used in the surgical and critical care environments by medical professionals who have been trained in the safe use of hemodynamic technologies and clinical usage of arterial line pressure monitoring. The disposable sterile cable with a red connector interfaces exclusively with an Edwards cable that is specifically wired for the pressure monitor being used. The disposable sterile cable with a green connector interfaces exclusively with the Edwards cables for use with the Edwards arterial pressure based cardiac output monitoring devices or hardware.

The FloTrac Jr sensor has a straight, flow-through design across the pressure sensors with an integral flush device.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

When used with a compatible monitoring platform, the FloTrac Jr sensor provides information regarding the hemodynamic status of the patient, which may lead to improved data driven clinical decision making for medically necessary intervention and/or clinical re-evaluation.

2.0 Intended Use/Purpose

The sensors measure an electrical resistance to provide blood pressure and arterial pressure based cardiac output.

3.0 Indications

The FloTrac Jr sensor is indicated for use in intravascular pressure monitoring. It is also indicated for use with the Edwards arterial pressure based cardiac output monitoring devices or hardware to measure cardiac output. The FloTrac Jr sensor is indicated for use in pediatric patients ≥ 12 years of age.

4.0 Contraindications

There are no absolute contraindications for using the FloTrac Jr sensor in patients requiring invasive pressure monitoring.

5.0 Warnings

- **Do not allow air bubbles to enter the setup. See complications section of this IFU regarding air emboli and abnormal pressure readings.**
- **Do not use the flush device during intracranial pressure monitoring.**

6.0 Specifications*

FloTrac Jr sensor	
Operating Pressure Range	-50 to +300 mmHg
Operating Temperature Range	15 to 40 °C
Sensitivity	5.0 μ V/V/mmHg $\pm$ 1%
Nonlinearity and Hysteresis	$\pm$ 1.5% of reading or $\pm$ 1 mmHg, whichever is greater
Excitation Frequency	dc to 5,000 Hz
Excitation Impedance	350 ohms $\pm$ 10% with typical Edwards Monitor Cable attached
Phase Shift	<5°
Signal Impedance	300 ohms $\pm$ 5%
Zero Offset	$\leq \pm$ 25 mmHg
Zero Thermal Drift	$\leq \pm$ 0.3 mmHg/°C
Output Drift	$\pm$ 1 mmHg per 8 hours after 20 second warm-up
Sensitivity Thermal Drift	$\leq \pm$ 0.1%/°C
Natural Frequency	40 Hz nominal for a standard kit (48"/12") (122 cm/30.5 cm); >200 Hz for sensor alone
Defibrillator Challenge	withstands 5 repeated discharges of 360 Joules within 5 minutes delivered into a 50 ohm load
Leakage Current	<2 μ amps at 120V RMS 60 Hz
Overpressure Tolerance	-500 to +5000 mmHg
Shock Resistance	withstands 3 drops from 1 meter
Light Sensitivity	< 1 mmHg at 6 volts excitation when exposed to a 3400 °K tungsten light source at 3000 foot candles (32,293 lumen/m ²)
Volumetric Displacement	$\leq$ 0.03 mm ³ /100 mmHg for transducer without flush device
Flow rate across flush device with IV bag pressurized to 300 mmHg Blue Snap-Tab device	3 $\pm$ 1 ml/hr

*at 6.00 VDC and 25 °C unless otherwise stated.

All specifications meet or exceed the AAMI BP22 Standard for performance interchangeability of resistance bridge type blood pressure transducers.

7.0 Procedure

Step	Procedure
1	Ensure that the cables are compatible with the monitors being used. Connect the reusable cables to the monitors, and turn the monitors on to allow the electronics to warm up.
2	Using aseptic technique, remove the sensor and kit from the sterile packaging. CAUTION: Modifying any FloTrac Jr kit may reduce dynamic response resulting in compromised hemodynamic monitoring performance.
3	Ensure that all connections are secure but are not overtightened.
4	Connect the sensor cable with the red shield to the reusable cable appropriate for the pressure monitor.
5	Connect the sensor cable with the green shield to the Edwards' FloTrac Jr connecting cable.

- **High pressures, which may be generated by an infusion pump at certain flow rates, may override the flush device restriction, resulting in fast flushing at the rate set by the pump.**
- **Avoid contact with any topical cream or ointment that attacks polymeric materials. May affect product integrity.**
- **Do not expose electrical connections to fluid contact. May result in electrical shock to user or patient, and/or arrhythmia.**
- **Do not autoclave the reusable cable as it may affect product integrity.**
- **This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.**
- **Compliance to IEC 60601-1 is only maintained when the FloTrac Jr sensor is connected to a patient monitor or equipment that has a Type CF defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the FloTrac Jr sensor. Failure to ensure monitor or equipment compliance to IEC 60601-1 and FloTrac Jr sensor compatibility may increase the risk of electrical shock to the patient/operator.**
- **Do not modify or alter the product in any way. Alteration or modification may affect patient/operator safety or product performance.**
- **For a patient undergoing an MRI examination, please refer to the MRI Safety Information section for specific conditions to ensure patient safety.**

Step	Procedure
6	Remove all air from the IV flush solution bag. CAUTION: If all air is not removed from the bag, air may be forced into the patient's vascular system when the solution is exhausted. See complications section of this IFU regarding air emboli.
7	Heparin anticoagulation therapy to be administered per hospital protocol.
8	Close the roller clamp on the IV set and connect the IV set to the IV flush bag. Hang the bag approximately 2 feet (60 cm) above the patient. This height will provide approximately 45 mmHg of pressure to prime the setup.
9	Fill the drip chamber halfway with flush solution by squeezing the drip chamber. Open the roller clamp.
10	Flow is provided by pulling on the Snap-Tab device and discontinued by releasing the Snap-Tab device.
11	Prime system using gravity only (do not pressurize bag) to decrease fluid turbulence and mitigation of bubbles.

Step	Procedure
12	For kits with IV sets attached, open the sensor vent port by turning the stopcock handle. Deliver flush solution first through the sensor and out through the vent port, then through the remaining kit by turning the appropriate stopcocks. Remove all air bubbles. CAUTION: Significant distortion of the pressure waveform or air emboli can result from air bubbles in the setup. See complications section of this IFU regarding air emboli.
13	Replace all vented caps on sideports of the stopcocks with nonvented caps.
14	Mount the sensor on an IV pole using the appropriate clamp and holder. Snap sensor into place in holder.
15	Pressurize IV flush solution bag after initial gravity prime. Flow rate will vary with pressure across the flush device. The flow rate with the IV bag pressurized to 300 mmHg: 3 ± 1 ml/hr.
16	Connect pressure tubing to the catheter per manufacturer's instructions.
17	Flush system per hospital policy. CAUTION: After each fast-flush operation, observe the drip chamber to verify that the continuous flush rate is as desired. See complications section of this IFU regarding clotted catheter and bleed-back, and overinfusion.

8.0 Zeroing and Calibration

Step	Procedure
1	Adjust the level of the sensor vent port (the fluid-air interface) to correspond to the chamber where pressure is being measured. For example, in cardiac monitoring, zero at level of the right atrium. This is at the phlebostatic axis, determined by the intersection of the midaxillary line and the fourth intercostal space.
2	Remove the non-vented cap from the stopcock above the transducer and open the port to the atmosphere while maintaining sterility of both the inside of the cap and port.
3	Zero both the bedside monitor and the Edwards monitoring system to the atmospheric air per monitor's instructions.
4	Close the vent port to the atmosphere and then replace the non-vented cap.
5	The system is ready to begin monitoring.

9.0 Testing Dynamic Response

The assembly may be tested for dynamic response by observing the pressure waveform on an oscilloscope or monitor. Bedside determination of the dynamic response of the catheter, monitor, kit, and sensor system is done after the system is flushed, attached to the patient, zeroed, and calibrated. A square-wave test may be performed by pulling the Snap-Tab device and releasing quickly.

Note: Poor dynamic response can be caused by air bubbles, clotting, excessive tubing length, excessively compliant pressure tubing, small bore tubing, loose connections, or leaks.

10.0 Routine Maintenance

Follow hospital policies and procedures for frequency of zeroing the sensor and pressure monitor and for replacing and maintaining pressure monitoring lines. The FloTrac Jr sensor is precalibrated and has a negligible drift with time (see **Specifications**).

Step	Procedure
1	Adjust zero pressure reference each time level of the patient is changed. CAUTION: When rechecking zero or verifying accuracy, ensure that the non-vented cap is removed before opening the sensor vent port to the atmosphere. See complications section of this IFU regarding abnormal pressure readings.
2	Periodically check fluid path for air bubbles. Ensure that connecting lines and stopcocks remain tightly fitted.
3	Periodically observe the drip chamber to verify that the continuous flush rate is as desired.
4	The Centers for Disease Control recommends replacing disposable or reusable transducers at 96-hour intervals. Replace other components of the system, including the tubing, continuous-flush device, and flush solution, at the time the sensor is replaced.

11.0 MRI Safety Information

MR Conditional

The following device was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International Designation: F2503-13,

Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

11.1 FloTrac Jr Sensors

- Non-clinical testing demonstrated that the FloTrac Jr sensor is MR Conditional according to the following conditions:
- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m) or less
- This device and the associated cable are not intended for use inside of the bore of the MR system and should not be in direct contact with the patient; failure to follow this guideline may result in a serious patient injury
- This device and the associated cable may be in the MR system room but not in operation or connected to a pressure monitoring system during an MRI examination
- The pressure tubing connected to the FloTrac Jr sensor is MR Safe and can be placed inside of the bore of the MR system during an MRI examination

12.0 Complications

Device risks may include hyper/hypovolemia, blood loss, and/or patient burn.

12.1 Sepsis/Infection

Positive cultures can result from contamination of the pressure setup. Increased risks of septicemia and bacteremia have been associated with blood sampling, infusing fluids, and catheter related thrombosis.

12.2 Air Emboli

Air can enter the patient through stopcocks that are inadvertently left open, from accidental disconnection of the pressure setup, or from flushing residual air bubbles into the patient.

12.3 Clotted Catheter and Bleed-Back

If the flush system is not adequately pressurized relative to the patient's blood pressure, blood bleed-back and catheter clotting may occur.

12.4 Overinfusion

Excessive flow rates may result from pressures greater than 300 mmHg. This may result in a potentially harmful increase in blood pressure and fluid overdose.

12.5 Abnormal Pressure Readings

Pressure readings can change quickly and dramatically because of loss of proper calibration, loose connection, or air in the system.

CAUTION: Abnormal pressure readings should correlate with the patient's clinical manifestations. Verify system integrity before instituting therapy.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

13.0 How Supplied

The FloTrac Jr sensor is supplied sterile in preconnected monitoring kits. Contents are sterile and the fluid path is nonpyrogenic if package is undamaged or unopened. Do not use if package is opened or damaged. Do not resterilize. Visually inspect for breaches of packaging integrity prior to use.

14.0 Storage

Store in a cool, dry place.

15.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

16.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

17.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers: Inside the U.S. and Canada (24 hours): 800.822.9837
Outside the U.S. and Canada (24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications, and model availability are subject to change without notice.

Refer to the latest version of the monitoring system operator's manual for more information.

Refer to the symbol legend at the end of this document.

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	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello	Modelnummer	Modelnummer
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op	Forsigtig
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken	Må ikke genanvendes
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid	Mængde
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer	Partinummer
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum	Anvendes inden
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant	Producent
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum	Fremstillingsdato
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare	Niet opnieuw steriliseren	Må ikke resteriliseres
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen	Ikke-pyrogen
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata	MRI-veilig onder bepaalde voorwaarden	MR-betinget
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Importer	Importateur	Importeur	Importador	Importatore	Importeur	Importør
	Contents	Contenu	Inhalt	Contenido	Contenuti	Inhoud	Indhold
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarriersystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barriersystem
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo	Unieke instrumentidentificatie	Unik udstyrsidentifikation
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Consultare le istruzioni per l'uso	Volg de gebruiksaanwijzing	Følg brugsanvisningen
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/ Den Europæiske Union
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr

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	Slovensky	Norsk	Suomi	Български	Română	Eesti
#	Číslo modelu	Modellnummer	Mallinúmero	Номер на модела	Număr model	Mudeli number
!	Upozornenie	Forsiktig	Tärkeä huomautus	Внимание	Atenție	Ettevaatust
ⓧ	Nepoužívajte opakovane	Må ikke gjenbrukes	Ei uudelleenkäytettävä	Да не се използва повторно	A nu se reutiliza	Mitte korduskasutada
QTY	Množstvo	Antall	Määrä	Количество	Cantitate	Kogus
LOT	Číslo šarže	Lotnummer	Eränumero	Партиден номер	Număr de lot	Partii number
⌚	Dátum spotreby	Utløpsdato	Viimeinen käyttöpäivä	Срок на годност	A se utiliza până la data de	Kõlblik kuni
STERILE EO	Sterilizované pomocou etylénoksidu	Sterilisert med etylenoksid	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades
ⓧ	Nepoužívajte, ak je obal poškodený, a prečítajte si návod na použitie	Skal ikke brukes hvis pakningen er skadet. Se bruksanvisningen	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid
🏭	Výrobca	Produsent	Valmistaja	Производител	Producător	Tootja
🏭	Dátum výroby	Produksjonsdato	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev
ⓧ	Opakovane nesterilizujte	Må ikke resteriliseres	Älä steriloι uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige
eifu.edwards.com + 1 888 570 4016	Postupujte podľa návodu na použitie na webovej stránke	Følg bruksanvisningen på nettsiden	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiata veebisaidilt
ⓧ	Nepyrogeenne	Ikke-pyrogen	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne
MR	Podmienene bezpečné v prostredí MR	MR-sikker under spesifiserte forhold	Ehdollisesti turvallinen magneettikuvauksessa	Безопасно при MR при определени условия	Condiționat RM	Ohutu magnetresonantsomograafias teatud tingimustel
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
🌐	Dovozca	Importør	Maahantuoja	Вносител	Importator	Importija
📖	Obsah	Innhold	Sisältö	Съдържание	Cuprins	Sisu
🛡️	Systém jednej sterilnej bariéry	Enkelt, sterilt barriersystem	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem
UDI	Jedinečný identifikátor pomôcky	Unik enhetsidentifikator	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus
👤	Postupujte podľa návodu na použitie	Følg bruksanvisningen	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit
☀️☔	Skladujte na chladnom a suchom mieste	Oppbevares tørt og kjølig	Säilytettyä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.
EC REP	Autorizovaný zástupca v krajinách Európskeho spoločenstva/Európskej únie	Autorisert representant i Det europeiske fellesskap / EU	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/ Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus
CE 0123	Conformité Européenne (značka CE)	Conformité Européenne (CE-merke)	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)
MD	Zdravotnícka pomôcka	Medisinsk utstyr	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символах ■ Legenda de simboluri ■ Sümboleite seletus ■ Simbolijy paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Основные обозначения ■ Legenda sa simbolima ■ Legenda simbola

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Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

Made in Dominican Republic
CCT Critical Care Technologies
S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

06/24
10059781001 A / DOC-0236906 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda,
Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU