DIRECTORY


English	1	Português	37	Eesti	74
Français	5	Česky	41	Lietuvių	78
Deutsch	9	Magyar	45	Latviešu	82
Español	13	Polski	49	Türkçe	86
Italiano	17	Slovensky	53	Русский	90
Nederlands	21	Norsk	57	Srpski	95
Dansk	25	Suomi	61	Hrvatski	99
Svenska	29	Български	65		
Ελληνικά	33	Română	70		

English

Swan-Ganz

CCO Thermodilution Catheter: 139F75
 CCOmbo CCO/SvO₂ Thermodilution Catheter: 744F75
 CCOmbo CCO/SvO₂/VIP Thermodilution Catheter: 746F8
 CCOmbo V CCO/CEDV Thermodilution Catheter: 177F75N
 CCOmbo V (CCO/SvO₂/CEDV) Thermodilution Catheter: 774F75
 CCOmbo V (CCO/SvO₂/CEDV/VIP) Thermodilution Catheter: 777F8

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

Please refer to Figure 1 and Figure 2 on page 103.

Edwards, Edwards Lifesciences, the stylized E logo, CCOmbo, CCOmbo V, CO-Set, CO-Set+, Swan, Swan-Ganz, Vigilance, and VIP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

1.0 Description

The device is used by medical professionals who have been trained in safe use of invasive hemodynamic technologies and clinical usage of pulmonary artery catheters as part of their respective institutional guidelines.

The Swan-Ganz continuous cardiac output thermodilution catheters are flow-directed pulmonary artery catheters designed to enable the monitoring of hemodynamic pressures and to provide continuous cardiac output. The pulmonary artery (PA) distal lumen terminates at the distal tip. The proximal injectate lumen terminates at a port located 26 cm from the distal tip. When the distal tip is located in the pulmonary artery, the proximal injectate port will reside in the right atrium or vena cava, allowing for bolus cardiac output injections, right atrial pressure monitoring, blood sampling, or infusion of solutions.

When used with a compatible cardiac output computer¹, the Swan-Ganz continuous cardiac output (CCO) thermodilution catheter (Models 139F75, and 177F75N) allows for continuous calculation and display of cardiac output. To measure cardiac output continuously, a compatible cardiac output computer uses thermal energy emitted by the thermal filament located on the catheter to calculate cardiac output using thermodilution principles. Alternatively,

cardiac output can be measured by the traditional bolus thermodilution method.

When used with a compatible cardiac output computer¹, the Swan-Ganz CCO/SvO₂ catheter (Models 744F75, 774F75) allows for continuous calculation and display of cardiac output and mixed venous oxygen saturation. The oximetry lumen (optical module connector) terminates at the distal tip. This lumen contains the fibers that transmit the light to the pulmonary artery for measurement of mixed venous oxygen saturation. Mixed venous oxygen saturation is monitored by fiberoptic reflectance spectrophotometry. The amount of light absorbed, refracted, and reflected depends on the relative amounts of oxygenated and deoxygenated hemoglobin in the blood.

When used with a compatible cardiac output computer¹ the CCO/SvO₂/VIP catheter (Models 746F8, 777F8) provides an additional VIP lumen that allows for continuous infusion. The VIP lumen (proximal infusion lumen) terminates at a port located 30 cm from the distal tip. This port allows for infusion of solutions, pressure monitoring or blood sampling.

In addition to the capabilities listed above, models 177F75N, 774F75, 777F8 are designed to enable continuous end diastolic volume (CEDV) when used with the compatible cardiac output computer. To measure end diastolic volume continuously, the compatible cardiac output computer uses thermal energy emitted by the thermal filament located on the catheter, and an ECG monitor signal (preferably in lead II configuration) "slaved" into the compatible cardiac output computer (refer to the appropriate operator's manual for information regarding "slaving" techniques) to calculate ejection fraction using thermodilution principles. CEDV is then derived from the ejection fraction and cardiac output measurements.

The intravascular catheter is inserted through the central vein into the right side of the heart and is advanced towards the pulmonary artery. Route of insertion can be internal jugular, femoral, antecubital and brachial veins. The body parts in contact are the atrium, ventricles, pulmonary artery and circulatory system.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended for use in adult critically ill or surgical patient populations. The device has not yet been tested in pediatric population or in pregnant or lactating women.

2.0 Intended Use / Purpose

Swan-Ganz catheters are pulmonary artery catheters intended for short-term use in the central circulatory system for patients who require intracardiac hemodynamic monitoring, blood sampling and infusing solutions. When used with a compatible monitoring platform and accessories, Swan-Ganz catheters offer a comprehensive hemodynamic profile to help clinicians assess cardiovascular function and guide treatment decisions.

¹ Measurement capabilities (i.e. CCO, CCO/SvO₂ or CCO/SvO₂/CEDV) of the compatible cardiac output computer vary by model number. Ensure that the monitor being used is able to measure the desired parameters.

3.0 Indications

Swan-Ganz catheters are diagnostic and monitoring tools used for hemodynamic monitoring of adult critically ill patients including but not limited to post major surgical recovery, trauma, sepsis, burns, pulmonary disease, pulmonary failure, cardiac disease including heart failure.

4.0 Contraindications

Patients with either recurrent sepsis, or with hypercoagulopathy, in which the catheter could serve as a focal point for septic or bland thrombus formation, should not be considered candidates for a balloon flotation catheter.

These products contain metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

5.0 Warnings

No absolute contraindications to the use of flow-directed pulmonary artery catheters exist. However, a patient with a left bundle branch block may develop a right bundle branch block during catheter insertion, resulting in complete heart block. In such patients, temporary pacing modes should be immediately available.

Electrocardiographic monitoring during catheter passage is encouraged and is particularly important in the presence of either of the following conditions:

- Complete left bundle branch block, in which the risk of complete heart block is somewhat increased.
- Wolff-Parkinson-White syndrome and Ebstein's malformation, in which the risk of tachyarrhythmias is present.

Do not modify or alter the product in any way. Alteration or modification may affect product performance.

Air should never be used for balloon inflation in any situation where it may enter the arterial circulation, e.g., in all pediatric patients and in adults with suspected right to left intracardiac or intra-pulmonary shunts. Bacteria-filtered carbon dioxide is the recommended inflation medium because of its rapid absorption into the blood in the event of balloon rupture within the circulation. Carbon dioxide diffuses through the latex balloon, diminishing the balloon's flow-directed capability after 2 to 3 minutes of inflation.

Do not leave the catheter in a permanent wedge position. Furthermore, avoid lengthy balloon inflation while the catheter is in a wedge position; this occlusive maneuver may result in pulmonary infarction.

CCO monitoring should always be discontinued when blood flow around the thermal filament is stopped to avoid thermal tissue injury. Clinical situations where CCO monitoring should be discontinued include, but are not limited to:

- Time periods when a patient is on cardiopulmonary bypass,
- Partial withdrawal of the catheter so that the thermistor is not in the pulmonary artery, or
- Removal of the catheter from the patient.

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Cleaning and resterilization will damage the integrity of the latex balloon. Damage may not be obvious during routine inspection.

6.0 Precautions

Failure of a balloon flotation catheter to enter the right ventricle or pulmonary artery is rare, but may occur in patients with an enlarged right atrium or ventricle particularly if the cardiac output is low or in the presence of tricuspid or pulmonic incompetence or pulmonary hypertension. Deep inspiration by the patient during advancement may also facilitate passage.

Clinicians using the device should be familiar with the device and understand its applications prior to use.

7.0 Recommended Equipment

WARNING: Compliance to IEC 60601-1 is only maintained when the catheter or probe (Type CF applied part, defibrillation proof) is connected to a patient monitor or equipment that has a Type CF

defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the catheter or probe. Failure to ensure monitor or equipment compliance to IEC 60601-1 and catheter or probe compatibility may increase the risk of electrical shock to the patient/operator.

- Swan-Ganz catheter
- Percutaneous sheath introducer and contamination shield
- Vigilance monitor for continuous cardiac output, mixed venous oxygen saturation and continuous end diastolic volume measurement (or other compatible cardiac output computer for measuring cardiac output by the bolus thermodilution method)
- Injectate temperature sensing probe (if performing bolus thermodilution measurements)
- Connecting cables
- Model OM2 or OMZE Optical Module (Models 744F75, 746F8, 774F75, and 777F8)
- Sterile flush system and pressure transducers
- Bedside ECG and pressure monitor system
- Appropriate ECG "slave" cables for CEDV models (177F75N, 774F75, and 777F8)

In addition, the following items should be immediately available if complications arise during catheter insertion: antiarrhythmic drugs, defibrillator, respiratory assist equipment, and a means for temporary pacing.

Monitor Set-Up and Calibration for Mixed Venous Oxygen Saturation Monitoring

This section only applies to the following models which have SvO₂ capability: 744F75, 746F8, 774F75, and 777F8.

The compatible cardiac output computer can be calibrated prior to catheter insertion by performing an *in vitro* calibration. When performing an *in vitro* calibration, do so before preparing the catheter (i.e. flushing the lumens). **The catheter tip must not get wet before an *in vitro* calibration is performed.** An *in vivo* calibration is required if an *in vitro* calibration is not done. *In vivo* calibration may be used to periodically recalibrate the monitor. Refer to the monitor operator's manual for detailed calibration instructions.

8.0 Catheter Preparation

Visually inspect for breaches of packaging integrity prior to use.

Use aseptic technique.

Note: Use of a protective catheter sheath is recommended.

Precaution: Avoid forceful wiping or stretching of the catheter during testing and cleaning as not to break the thermistor wire circuitry or detach thermal filament leads from other circuit components if present.

Before the catheter is inserted, the following preparation procedure should be followed:

1. Perform *in vitro* calibration (when measuring mixed venous oxygen saturation).
2. To open the catheter for calibration, peel back the top left perforated portion of the lid and fold diagonally according to the dotted fold line. In order to gain access to the catheter lumens for flushing, peel lid from the top right corner and fold according to the dotted fold line; otherwise peel off the entire tray lid by peeling back the lid from the bottom right corner. Flush lumens with sterile saline or dextrose solution to ensure patency and to remove air.
3. Gently lift the catheter up and remove it from the silicone gripper (see Figure 1 on page 103, Step 1).
4. Once the catheter has cleared the silicone gripper, pull the balloon out of the calibrator cup and remove the catheter from the tray (see Figure 1 on page 103, Step 2).

Note: To avoid damaging the balloon, do not pull the balloon through the silicone gripper.

5. Check balloon integrity by inflating it to the recommended volume. Check for major asymmetry and for leaks by submerging in sterile saline or water. Deflate balloon before insertion.

6. Connect the catheter's injectate and pressure monitoring lumens to the flush system and pressure transducers. Ensure that the lines and transducers are free of air.
7. Test the thermistor's electrical continuity before insertion. Connect the thermistor to the monitor and confirm that no fault messages appear.
8. If using a compatible cardiac output computer for continuous cardiac output measurement, connect the thermistor and thermal filament to the monitor and observe the following message: "Press START to begin CCO monitoring."
9. If using a compatible cardiac output computer for CEDV measurement, connect the thermistor and thermal filament to the monitor and observe the following message: "Press START to begin CCO monitoring."

Note: If proper ECG "slave" is implemented, CEDV measurement will begin.

9.0 Insertion Procedure

Swan-Ganz catheters can be inserted at the patient's bedside without the aid of fluoroscopy, guided by continuous pressure monitoring.

Simultaneous pressure monitoring from the distal lumen is recommended. Fluoroscopy is recommended for femoral vein insertion.

Note: Should the catheter require stiffening during insertion, slowly perfuse the catheter with 5 ml to 10 ml of cold sterile saline or 5% dextrose as the catheter is advanced through a peripheral vessel.

Note: The catheter should pass easily through the right ventricle and pulmonary artery and into a wedge position in less than a minute.

Although a variety of techniques can be used for insertion, the following guidelines are provided as an aid to the physician:

1. Introduce the catheter into the vein through a sheath introducer using percutaneous insertion using modified Seldinger technique.
2. Under continuous pressure monitoring, with or without the aid of fluoroscopy, gently advance the catheter into the right atrium. Entry of the catheter tip into the thorax is signaled by an increased respiratory fluctuation in pressure. Figure 2 on page 103 shows the characteristic intracardiac and pulmonary pressure waveforms.
3. Using the syringe provided, inflate the balloon with CO₂ or air to the maximum recommended volume. **Do not use liquid.** Note that an offset arrow on the gate valve indicates the "closed" position.

Note: Inflation is usually associated with a feeling of resistance. On release, the syringe plunger should usually spring back. If no resistance to inflation is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may continue to be used for hemodynamic monitoring. However, be sure to take precautions to prevent infusion of air or liquid into the balloon lumen.

WARNING: Pulmonary complications may result from improper inflation technique. To avoid damage to the pulmonary artery and possible balloon rupture, do not inflate above the recommended volume.

4. Advance the catheter until pulmonary artery occlusion pressure (PAOP) is obtained, then passively deflate the balloon by removing the syringe from the gate valve. Do not forcefully aspirate as this may damage the balloon. After deflation, re-attach the syringe.

Note: Avoid prolonged maneuvers to obtain wedge pressure. If difficulties are encountered, give up the "wedge".

Note: Before reinflation with CO₂ or air, completely deflate the balloon by removing the syringe and opening the gate valve.

Precaution: It is recommended that the provided syringe be re-attached to the gate valve after balloon deflation to prevent inadvertent injection of liquids into the balloon lumen.

Precaution: If a right ventricular pressure tracing is still observed after advancing the catheter several centimeters

beyond the point where the initial right ventricular pressure tracing was observed, the catheter may be looping in the right ventricle which can result in kinking or knotting of the catheter (see Complications). Deflate the balloon and withdraw the catheter into the right atrium. Reinflate the balloon and readvance the catheter to a pulmonary artery wedge position, then deflate the balloon.

Precaution: Catheter looping may occur when excessive length has been inserted, which could result in kinking or knotting (see Complications). If the right ventricle is not entered after advancing the catheter 15 cm beyond entry into the right atrium, the catheter may be looping, or the tip may be engaged in a neck vein with only the proximal shaft advancing into the heart. Deflate the balloon and withdraw the catheter until the 20 cm mark is visible. Reinflate the balloon and advance the catheter.

5. Reduce or remove any excessive length or loop in the right atrium or ventricle by slowly pulling the catheter back approximately 2 to 3 cm.

Precaution: Do not pull the catheter across the pulmonic valve while the balloon is inflated to avoid damage to the valve.

6. Reinflate the balloon to determine the minimum inflation volume necessary to obtain a wedge tracing. If a wedge is obtained with less than the maximum recommended volume (see specifications table for balloon inflation capacity), the catheter must be withdrawn to a position where full inflation volume produces a wedge tracing.

Precaution: Overtightening the proximal Tuohy-Borst adapter of the contamination shield may impair catheter function.

7. Confirm final catheter tip position with chest X-ray.

Note: After deflation, the catheter tip may tend to recoil towards the pulmonic valve and slip back into the right ventricle, requiring that the catheter be repositioned.

Note: If using a contamination shield, extend the distal end towards the introducer valve. Extend the proximal end of the catheter contamination shield to desired length, and secure.

10.0 Maintenance and Use *in situ*

The catheter should remain indwelling only as long as is required by the patient's condition.

Precaution: The incidence of complications increases significantly with indwelling periods longer than 72 hours.

10.1 Catheter Tip Position

Keep catheter tip centrally located in a main branch of the pulmonary artery near the hilum of the lungs. Do not advance tip too far peripherally. Tip should be kept where full or near full inflation volume is required to produce a wedge tracing. The tip migrates toward periphery during balloon inflation.

Insertion of the thermal filament beyond the pulmonic valve may result in erroneous continuous cardiac output measurements.

10.2 Catheter Tip Migration

Anticipate spontaneous catheter tip migration towards periphery of pulmonary bed. Continuously monitor distal lumen pressure to verify tip position. If wedge tracing is observed when balloon is deflated, pull catheter back. Damage may be caused by prolonged occlusion or over-distension of vessel upon re-inflation of the balloon.

Spontaneous catheter tip migration towards the periphery of the lung occurs during cardiopulmonary bypass. Partial catheter withdrawal (3 to 5 cm) just before bypass should be considered, as it may help reduce distal migration and prevent permanent catheter wedging postbypass. After termination of bypass, the catheter may require repositioning. Check the distal pulmonary artery tracing before inflating the balloon.

Precaution: Over a period of time, the catheter tip may migrate towards the periphery of the pulmonary bed and lodge in a small vessel. Damage may occur either by prolonged occlusion or by over-distension of the vessel upon reinflation of the balloon (see Complications).

PA pressures should be continuously monitored with the alarm parameter set to detect physiologic changes as well as spontaneous wedge.

10.3 Balloon Inflation and Wedge Pressure Measurement

Re-inflation of the balloon should be performed gradually while monitoring pressures. Inflation is usually associated with a feeling of resistance. If no

resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may still be used for hemodynamic monitoring, however, take precautions against infusion of air or liquids into the balloon lumen. During normal catheter use, keep inflation syringe attached to gate valve to prevent inadvertent injection of liquid into the balloon inflation lumen.

Measure wedge pressure only when necessary and only when tip is properly positioned (see above). Avoid prolonged maneuvers to obtain wedge pressure and keep wedge time to a minimum (two respiratory cycles or 10 - 15 seconds), especially in patients with pulmonary hypertension. If difficulties are encountered, discontinue wedge measurements. In some patients, pulmonary arterial end-diastolic pressure can often be substituted for pulmonary artery wedge pressure if the pressures are nearly identical, obviating the need for repeated balloon inflation.

10.4 Spontaneous Tip Wedging

The catheter may migrate into the distal pulmonary artery and spontaneous tip wedging may occur. To avoid this complication, pulmonary artery pressure should be monitored continuously with a pressure transducer and display monitor.

Forward advancement should never be forced if resistance is encountered.

10.5 Patency

All pressure monitoring lumens should be filled with a sterile, heparinized saline solution (e.g., 500 I.U. heparin in 500 ml saline) and flushed at least once each half hour or by continuous slow infusion. If loss of patency occurs and cannot be corrected by flushing, the catheter should be removed.

10.6 General

Keep pressure monitoring lumens patent by intermittent flush or continuous slow infusion with heparinized saline solution. Infusion of viscous solutions (e.g., whole blood or albumin) is not recommended, as they flow too slowly and may occlude the catheter lumen.

WARNING: To avoid pulmonary artery rupture, never flush the catheter when the balloon is wedged in the pulmonary artery.

Periodically check IV lines, pressure lines, and transducers to keep them free of air. Also ensure that connecting lines and stopcocks remain tightly fitted.

11.0 Cardiac Output Measurement

11.1 Continuous

Continuous cardiac output measurement is made by periodically warming the blood in the right atrium or ventricle with a known quantity of heat. The catheter thermistor detects the small change in blood temperature downstream, and the compatible cardiac output computer computes a dilution curve via a modified Stewart-Hamilton indicator dilution equation. This measurement technique is conducted without additional instrument calibration, material preparation, or operator intervention. If confirmation of the displayed continuous cardiac output value is deemed necessary, performing a bolus TD cardiac output measurement is recommended. Refer to the compatible cardiac output computer operator's manual for more information.

11.2 Bolus Thermodilution Method

Bolus TD cardiac output measurement is made by injecting an exact amount of physiological solution (saline or dextrose) of known temperature into the right atrium or superior vena cava and by using the thermistor in the pulmonary artery to detect the resultant change in blood temperature. Cardiac output is inversely proportional to the area under the temperature-versus-time curve. The accuracy of this method depends on the accuracy with which the quantity and temperature of the injectate are known. The accuracy of the thermodilution method correlates well with the dye dilution technique and with the direct Fick method.

Edwards Lifesciences cardiac output computers require that a computation constant be used to correct for injectate temperature rise as it passes through the catheter. The computation constant is a function of injectate volume, temperature, and catheter dimensions. The computation constants listed in the specifications have been determined *in vitro*.

12.0 Continuous Volumetric Measurement

Continuous volumetric measurements are made by periodically warming the blood in the right atrium or ventricle with a known quantity of heat, and by sensing heart rate from a "slaved" ECG signal (refer to the appropriate operator's manual for information regarding "slaving" techniques). The catheter thermistor detects the small change in blood temperature downstream, and the compatible cardiac output computer computes an ejection fraction based on thermodilution principles. Subsequently, continuous stroke, end systolic and end diastolic volume measurements are

derived from the ejection fraction and cardiac output measurements as follows:

CSV = CCO/HR

CEdV = CSV/CEF

CESV = CEdV - CSV

Where:

CSV = Continuous Stroke Volume

CCO = Continuous Cardiac Output

HR = Heart Rate

CEF = Continuous Ejection Fraction

CEdV = Continuous End Diastolic Volume

CESV = Continuous End Systolic Volume

This measurement technique is conducted without additional instrument calibration, material preparation or operator intervention. Refer to the compatible cardiac output computer operator's manual for more information.

13.0 MRI Information

MR Unsafe

The Swan-Ganz device is MR unsafe as the result of the device containing metallic components, which experience RF-induced heating in the MRI environment; therefore the device poses hazards in all MRI environments.

14.0 Complications

Invasive procedures involve some patient risks. Although serious complications are relatively uncommon, the physician is advised, before deciding to insert or use the catheter, to consider the potential benefits in relation to the possible complications. The techniques for insertion, methods of using the catheter to obtain patient data information, and the occurrence of complications is well described in the literature.

Strict adherence to these instructions and awareness of risks reduces the incidence of complications. Several known complications include:

14.1 Perforation of the Pulmonary Artery

Factors associated with fatal pulmonary artery rupture include pulmonary hypertension, advanced age, cardiac surgery with hypothermia and anticoagulation, distal catheter tip migration, arteriovenous fistula formation and other vascular traumas.

Extreme care should be used during the measurement of pulmonary artery wedge pressure in patients with pulmonary artery hypertension. In all patients, balloon inflation should be limited to two respiratory cycles, or 10 to 15 seconds.

A central location of the catheter tip near the hilum of the lung may prevent pulmonary artery perforation.

14.2 Pulmonary Infarction

Tip migration with spontaneous wedging, air embolism, and thromboembolism can lead to pulmonary artery infarction.

14.3 Cardiac Arrhythmias

Cardiac arrhythmias may occur during insertion, withdrawal, and repositioning, but are usually transient and self-limited. Premature ventricular contractions are the most commonly observed arrhythmia. Ventricular tachycardia and atrial tachycardia have been reported. Use of prophylactic lidocaine should be considered to decrease the incidence of ventricular arrhythmias during catheterization. ECG monitoring and immediate availability of antiarrhythmic drugs and defibrillator equipment is recommended.

14.4 Knotting

Flexible catheters have been reported to form knots, most often as a result of looping in the right ventricle. Sometimes the knot can be resolved by insertion of a suitable guidewire and manipulation of the catheter under fluoroscopy. If the knot does not include any intracardiac structures, the knot may be gently tightened and the catheter withdrawn through the site of entry.

14.5 Sepsis/Infection

Positive catheter-tip cultures resulting from contamination and colonization have been reported, as well as incidences of septic and aseptic vegetation in the right heart. Increased risks of septicemia and bacteremia have been associated with blood sampling, the infusing of fluids, and catheter-related thrombosis. Preventive measures should be taken to guard against infection.

14.6 Other Complications

Other complications include right bundle branch block, complete heart block, tricuspid and pulmonic valve damage, pneumothorax, thrombosis, blood loss, cardiac structure/wall injury or damage, hematoma, embolism, anaphylaxis, cardiac tissue/artery burn. In addition, allergic reactions to latex have been reported. Physicians should identify latex sensitive patients and be prepared to treat allergic reactions promptly.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance for this medical device. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

15.0 Long-Term Monitoring

The duration of catheterization should be the minimum required by the patient's clinical state since the risk of thromboembolic and infectious complications increases with time. The incidence of complications increases significantly with indwelling periods longer than 72 hours. Prophylactic systemic anticoagulation and antibiotic protection should be considered when long-term catheterization (i.e., over 48 hours) is required, as well as in cases involving increased risk of clotting or infection.

16.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

The packaging within is designed to avoid crushing of the catheter and to protect the balloon from exposure to the atmosphere. It is therefore recommended that the catheter remain inside the package until use.

17.0 Storage

Store in a cool, dry place.

18.0 Operating Conditions / Use Environment

Intended to operate under physiological conditions of the human body in a controlled clinical environment.

19.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the recommended expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

20.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837
Outside the U.S. and Canada (24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

21.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Catheter Models and Functions:

	139F75	177F75N	744F75	774F75	746F8	777F8
ICO	X	X	X	X	X	X
CCO	X	X	X	X	X	X
SvO ₂			X	X	X	X
VIP	X	X			X	X
CEDV		X		X		X

Specifications:

	Models 139F75, 177F75N	Models 744F75, 774F75	Models 746F8, 777F8
Body Color	yellow	yellow	yellow
Usable Length (cm)	110	110	110
Catheter Body	7.5F (2.5 mm)	7.5F (2.5 mm)	8F (2.7 mm)
Diameter of Inflated Balloon (mm)	13	13	13
Required Introducer Size	8.5F (2.8 mm)	8.5F (2.8 mm) or 9F (3.0 mm)	9F (3.0 mm)
Balloon Inflation Capacity (ml)	1.5	1.5	1.5
Distance from Distal Tip (cm)			
Thermistor	4	4	4
Thermal Filament	14-25	14-25	14-25
Injectate Port	26	26	26
VIP Port	30	-	30
Distance Between Markings (cm)	10	10	10
Lumen Volumes (ml)			
Distal Lumen	0.96	0.96	0.90
Injectate Lumen	0.8	0.95	0.85
Infusion Lumen	0.95	-	1.10
Infusion Rate* (ml/min)			
Distal Lumen	6	6	4
Injectate Lumen	9	14	9
VIP Lumen	16	-	16
Compatible Guidewire Diameter	0.025 in (0.64 mm)	0.018 in (0.46 mm)	0.018 in (0.46 mm)

	Models 139F75, 177F75N	Models 744F75, 774F75	Models 746F8, 777F8
Frequency Response			
Distortion at 10 Hz			
Distal Lumen	< 3 dB	< 3 dB	< 3 dB
Balloon Inflation Syringe	3 ml, limited to 1.5 ml	3 ml, limited to 1.5 ml	3 ml, limited to 1.5 ml

All specifications given are nominal values.

* Using normal saline at room temperature, 1 m above insertion site, gravity drip.

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
	English	Français	Deutsch	Español	Italiano	Nederlands
	Number of lumens	Nombre de lumières	Anzahl der Lumen	Número de luces	Numero di lumi	Aantal lumina
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad
	Lumen size	Taille de la lumière	Lumendurchmesser	Tamaño de la luz	Dimensioni del lume	Lumengrootte
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello	Modelnummer
	Minimum introducer size	Taille minimale de l'introducteur	Mindestgröße der Einführhilfe	Tamaño de introductor mínimo	Misura minima dell'introduttore	Minimum introducergrootte
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum
	Inner diameter	Diamètre interne	Innendurchmesser	Diámetro interior	Diametro interno	Binnendiameter
	Balloon capacity	Capacité du ballonnet	Ballonkapazität	Capacidad del balón	Capacità del palloncino	Balloncapaciteit
	Sterile	Stérile	Steril	Estéril	Sterile	Steriel
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni	Gesteriliseerd met behulp van straling
	Sterilized using steam or dry heat	Stérilisé à la vapeur d'eau ou à la chaleur sèche	Mit Wasserdampf oder trockener Hitze sterilisiert	Esterilizado con vapor o calor seco	Sterilizzato con vapore o calore secco	Gesteriliseerd met behulp van stoom of droge warmte
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig
	Size	Taille	Größe	Tamaño	Misura	Afmetingen
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel
	Consult instructions for use on the website	Consulter le mode d'emploi sur le site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Consultare le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM	Non compatibile con RM	MRI-onveilig

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	Polski	Slovensky	Norsk	Suomi	Български	Română
	Liczba kanałów	Počet lúmenov	Antall lumen	Luumenien määrä	Брой лумени	Număr de Lumenuri
	Średnica zewnętrzna	Vonkajší priemer	Ytre diameter	Ulkohalkaisija	Външен диаметър	Diametru exterior
	Długość użytkowa	Použitelná dĺžka	Anvendelig lengde	Käyttöpituus	Използваема дължина	Lungimea utilă
	Zalecany rozmiar przewoźnika	Odporúčaná veľkosť vodiaceho drôtu	Anbefalt ledediameterstørrelse	Ohjainlangan suositeltu koko	Препоръчителен размер телен водач	Mărimea recomandată a firului de ghidaj
	Rozmiar kanału	Veľkosť lúmenu	Lumenstørrelse	Luumenin koko	Размер на лумена	Mărimea lumenului
	Numer modelu	Číslo modelu	Modellnummer	Mallinnumero	Номер на модела	Număr model
	Minimalny rozmiar introduktora	Minimálna veľkosť zavadzača	Minste tillatte innførerstørrelse	Sisäänviejän vähimmäiskoko	Минимален размер на интродюсера	Mărimea minimă a dispozitivului de introducere
	Przeostroga	Upozornenie	Forsiktig	Tärkeä huomautus	Внимание	Atenție
	Nie używać ponownie	Nepoužívajte opakovane	Må ikke gjenbrukes	Ei uudelleenkäytettävä	Да не се използва повторно	A nu se reutiliza
	Ilość	Množstvo	Antall	Määrä	Количество	Cantitate
	Numer serii	Číslo šarže	Lotnummer	Eränumero	Партиден номер	Număr de lot
	Data przydatności do użycia („Zużyć do”)	Dátum spotreby	Utløpsdato	Viimeinen käyttöpäivä	Срок на годност	Data expirării
	Średnica wewnętrzna	Vnútný priemer	Innvendig diameter	Sisäläpimita	Вътрешен диаметър	Diametru interior
	Pojemność balonu	Objem balónika	Ballongkapasitet	Pallon tilavuus	Вместимост на балона	Capacitatea balonului
	Jałowy	Sterilný	Steril	Steriili	Стерилно	Steril
	Wysterylizowano przy użyciu tlenku etylenu	Sterilizované pomocou etylénoxidu	Sterilisert med etylenoksid	Steriloitu etyleenioksidilla	Стерилизирано с этиленов оксид	Sterilizat cu oxid de etilenă
	Wysterylizowano przez napromienienie	Sterilizované pomocou ožarovania	Sterilisert med stråling	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere
	Wysterylizowano parą lub suchym, gorącym powietrzem	Sterilizované parou alebo suchým teplom	Sterilisert med stråling eller tørr varme	Lämpö- tai höyrysteriloitu	Стерилизирано чрез автоклавиране или сух стерилизатор	Sterilizat cu abur sau prin căldură uscată
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Producent	Výrobca	Produsent	Valmistaja	Производител	Producător
	Data produkcji	Dátum výroby	Produksjonsdato	Valmistuspäivämäärä	Дата на производство	Data fabricației
	Zawiera lateks pochodzący z naturalnego kauczuku lub wykazuje jego obecność	Obsahuje alebo je prítomný prírodný kaučukový latex	Inneholder eller tilstedeværelse av naturgummilateks	Sisältää luonnokumilateksia tai sen jäämiä	Съдържа или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural
	Rozmiar	Veľkosť	Størrelse	Koko	Размер	Mărime
	Wyrób medyczny	Zdravotnícka pomôcka	Medisinsk utstyr	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical
	Należy zapoznać się z instrukcją użycia na stronie internetowej	Prečítajte si návod na použitie na webovej stránke	Se bruksanvisningen på nettstedet	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба в уебсайта	Consultați instrucțiunile de utilizare de pe site-ul web
	Postępować zgodnie z instrukcją stosowania na stronie internetowej	Postupujte podľa návodu na použitie na webovej stránke	Følg bruksanvisningen på nettsiden	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web
	Produktu nie można bezpiecznie używać w środowisku badań metodą rezonansu magnetycznego	Nie je bezpečné v prostredí MR	MR-usikker	Ei sovi magneettikuvaukseen	Небезопасно при MR	Incompatibil RM

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	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Valendike arv	Spindžių skaičius	Lūmenu skaits	Lümen sayısı	Количество просветов	Broj lumena	Broj lumena
	Vālislābimōot	Išorinis skersmuo	Ārējais diametrs	Dış çap	Внешний диаметр	Spoljašnji prečnik	Vanjski promjer
	Kasutatav pikkus	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjljiva dužina	Upotrebjljiva dužina
	Juhtetraadi soovituslik suurus	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs	Kılavuz Tel İçin Önerilen Boyut	Рекомендованный размер проводника	Preporučena veličina žice vodice	Preporučena veličina žice vodilice
	Valendiku suurus	Spindžio dydis	Lūmena izmērs	Lümen boyutu	Размер просвета	Velicina lumena	Velicina lumena
	Mudeli number	Modelio numeris	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Sisesti väikseim suurus	Minimalus įvediklio dydis	Minimālais ievadītāja izmērs	Minimum introdüser boyutu	Минимальный размер интродьюсера	Minimalna veličina uvodnika	Minimalna veličina uvodnog instrumenta
	Ettevaatust	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не использовать повторно	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina	Količina
	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Kõlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sisemine läbimōot	Vidinis skersmuo	Iekšējais diametrs	İç çap	Внутренний диаметр	Unutrašnji prečnik	Unutrašnji promjer
	Ballooni mahutavus	Balionių talpa	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Sterilne	Sterilus	Sterils	Steril	Стерильно	Sterilno	Sterilno
	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Steriliseeritud auru või kuiva kuumusega	Sterilizuota naudojant garus arba sausą karštį	Sterilizēts, izmantojot tvaiku vai sausu karstumu	Buhar veya kuru ısı kullanılarak sterilize edilmiştir	Стерилизовано паром или сухим жаром	Sterilizovano korišćenjem pare ili suvom toplotom	Sterilizirano parom ili vrućim, suhim zrakom
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Tootmiskuupäev	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Sisaldab naturaalselt kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Suurus	Dydis	Izmērs	Boyut	Размер	Velicina	Velicina
	Meditsiiniseade	Medicinos priemonė	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Palun lugege kasutusjuhiseid veebisaidilt	Žr. naudojimo instrukcijas interneto svetainėje	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na internet stranici	Pogledajte upute za upotrebu na internetskoj stranici
	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovaukitės naudojimo instrukcijomis, pateiktomis svetainėje	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa web-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Ohtlik magnetresonantstomograafias	MR nesaugus	Nedrīkst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR	Nije sigurno kod pregleda MR-om

	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Ohutu magnetresonantstomograafias	MR saugus	Drošs lietošanai MR vidē	MR için Güvenli	Безопасно при проведении МРТ	Bezbedno za MR	Sigurno kod pregleda MR-om
	Ohutu magnetresonantstomograafias teatud tingimustel	Sąlyginis MR	Izmantojams MR vidē, ievērojot noteiktus nosacījumus	MR Koşullu	Условно безопасно при проведении МРТ	Uslovno bezbedno za MR	Uvjetno siguran kod pregleda MR-om
	Tutvuge paberkandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu
	Ärge resteriliseerige	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Mittepüroogenne	Nepirogeniškas	Nepirogēna	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno
	B-tüüpi kohaldatav osa	B tipo liečiamoji dalis	B tipa lietotā daļa	B Tipi uygulanan parça	Рабочая часть типа B	Primenjeni deo tipa B	Primjenjen dio tipa B
	CF-tüüpi kohaldatav osa	CF tipo liečiamoji dalis	CF tipa kontaktējošā daļa	CF Tipi uygulanan parça	Рабочая часть типа CF	Primenjeni deo tipa CF	Primjenjen dio vrste CF
	Ärge kasutage, kui pakend on kahjustatud, ja lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Ava	Atidarymas	Atvērt	Açık	Открыть	Otvoreno	Otvoreno
	Enne sisestamist või väljavõtmist aspireerige balloon –0,5 cm ³ ulatuses.	Prieš įvesdami arba ištraukdami, iš balionėlio išsiurbkite –0,5 ml oro.	Pirms balona ievadīšanas vai izņemšanas atsūknējiet balonu –0,5 cm ³ .	Uygulama veya Geri Çekme Öncesinde Balonu –0,5 cm ³ Aspire Edin.	Перед введением или извлечением удалите из баллона 0,5 куб. см содержимого.	Aspirirajte balon –0,5 ml pre uvođenja ili izvlačenja.	Aspirirajte balon –0,5 cm ³ prije uvođenja ili izvlačenja.
	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Įgaliojasis atstovas Europos Bendrijoje / Europos Sąjungoje	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašteni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Importija	Importuotojas	Importētājs	İthalatçı	Импортёр	Uvoznik	Uvoznik
	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Kahekordne steriilne kattesüsteem	Dvigubo sterilaus barjero sistema	Divkārsa sterilitātes aizsargslāņa sistēma	Çift steril bariyer sistemi	Двойная барьерная система для стерилизации	Sistem sa dve sterilne barijere	Sustav dvostruke sterilne barijere
	Seadme kordumatu identifitseerimistunnus	Unikalusis priemonės identifikatorius	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	See pool üles	Šia puse į viršų	Ar šo pusi uz augšu	Bu taraf yukarı	Верхняя сторона	Ova strana mora da bude prema gore	Ova strana prema gore
	Ūrņ, kārtsēda ettevaatlīkult	Trapus, elgtis atsargiai	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
Märkus. Kõik sümbolid ei pruugi antud toote etiketil esineda. • Pastaba. Šio gaminio etiketėje gali būti pateikti ne visi simboliai. • Piezime. Ši išstrādājuma marķējuma var nebūt iekļauti visi simboli. • Not: Bu üürünün etiketinde tüm semboller yer almamış olabilir. • Примечание. На этикетках данного изделия могут указываться не все символы. • Napomena: Na ambalaži proizvoda se možda ne nalaze svi simboli. • Napomena: na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.							

Edwards

Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

Made in Dominican Republic
CCT Critical Care Technologies, S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

06/24
10060790001 A / DOC-0243683 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU