Edwards

Directory

English (en)	1	Česky (cs)	18	Latviešu (lv)	35
Français (fr)	3	Magyar (hu)	20	Türkçe (tr)	36
Deutsch (de)	4	Polski (pl)	21	Русский (ru)	38
Español (es)	6	Slovensky (sk)	23	Srpski (sr)	40
Italiano (it)	8	Norsk (no)	25	Hrvatski (hr)	41
Nederlands (nl)	10	Suomi (fi)	26	简体中文 (zh-cn)	43
Dansk (da)	11	Български (bg)	28	한국어 (ko)	45
Svenska (sv)	13	Română (ro)	30	日本語 (ja)	47
Ελληνικά (el)	15	Eesti (et)	31		
Português (pt)	16	Lietuvių (lt)	33		

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring
 ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom
 ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas
 ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ■ 符号图例 ■ 기호 범례
 ■ 記号の凡例

49


English

ForeSight IQ

Large Sensor

≥ 40 kg - Adult

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The sensor, when used in combination with the ForeSight oximeter cable is a single use applied part that measures hemoglobin allowing the clinician to continuously and accurately determine absolute levels of blood oxygenation saturation in the tissue (StO₂). This device is supplied non-sterile.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions for Use.

The device is used by clinicians who have been trained in the use of tissue oximetry devices in accordance with their institutional policies.

The benefit of using ForeSight sensor is to non-invasively provide tissue oxygen saturation and hemoglobin values to allow clinicians to manage their patients. Potential risks include burns, electrical shock, tissue damage, transient hypoxia, adverse reaction to device materials, and/or inappropriate/unintended treatment.

2.0 Intended Use/Purpose

When used with the HemoSphere Alta advanced monitoring platform in combination with the ForeSight oximeter cable, the intended purpose is to monitor absolute regional hemoglobin oxygen saturation of blood under the sensor.

3.0 Indications for Use

When used in conjunction with the ForeSight oximeter cable:

The ForeSight IQ large sensor is indicated for monitoring of absolute regional hemoglobin oxygen saturation of blood under the sensor in individuals at risk for reduced flow or no-flow ischemic states. It is intended for use on adults and transitional adolescents ≥ 40 kg.

4.0 Contraindications

The sensor is contraindicated for use on patients:

- With a physical site area too limited for proper sensor placement
- With allergic reactions to sensor adhesive
- Undergoing an MRI scan because of associate risk of injury

5.0 Warnings

- Assess the sensor site at least every 12 hours, or more often as required by the institution's protocol.
- Remove the sensor if the circulatory condition or skin integrity has deteriorated.
- Do not lay patient on the sensor or cable.
- Do not use in an MR environment due to risk of burn as result of sensor heating.
- Do not attach the sensor to skin with unapproved devices, such as headbands, hats, etc.
- Do not place the sensor or accessories over eyes, nose, or mouth.
- Do not cut the sensor. Cutting the sensor can result in injury to the patient.
- Do not use in an MRI environment.
- When used in settings with LED lighting, sensors should be covered with a light blocker prior to connection to the preamp cable, as some high intensity systems can interfere with the sensor's near infrared light detection.
- Failure to remove the protective liner may cause erroneous StO₂ readings.

6.0 Cautions

- Avoid positioning the sensor over hair, air sinus, hematomas, or broken skin.
- Avoid attaching the sensor to sites with excess adipose, ascites, or edema.
- The materials used in the manufacture of the sensor were NOT designed for reuse. Reuse can cause the sensor not to perform as intended.

7.0 Instructions

Step	Procedure
1	Remove the sensor from the package. Carefully inspect the sensor for damage. Discard and replace if damage is found.
2	Select sensor location on the monitor.
3	Clean and dry the sensor site (1).
4	Remove protective liner from the sensor (2).
5	Apply the sensor to the patient: a) Cerebral Use (Fig. 3a): Select the site on the forehead well above the eyebrow and just below the hairline. Do not attach sensor over hair. b) Somatic Use (Fig. 3b, 3c, 3d): Select the site that provides ideal access to the desired skeletal muscle tissue, for example: • Biceps (upper arm), brachioradialis (forearm), or deltoid (shoulder). • Quadriceps (upper leg), or gastrocnemius and tibialis (calf). Apply the sensor with the connector towards the feet. • Latissimus dorsi (flank) or external oblique (abdomen). Note: Somatic monitoring only available with ForeSight oximeter cable software version 4.0.1 or above. Note: Do not attach sensor over hair. Note: You may use Tegaderm under the sensor in patients with delicate skin or edema.
6	Insert the sensor connector straight into the sensor cable connector until it snaps into place (Figure 4). Use the bedsheet clip to secure the cable and prevent pulling on the sensor.
7	If needed, fold the sensor flat cable to route it in the desired direction.

For use only with ForeSight oximeter cable software version 2.5.7 or above.

8.0 Disposal

Sensors are designed for single-patient use, and are not to be reprocessed. Re-used sensors present a risk of cross-contamination or infection. Use a new sensor for each patient and discard it after use. Disposal should follow in accordance with local hospital and institution policies.

9.0 Storage

Store in a cool, dry place.

10.0 Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration and malfunction.

11.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Refer to the latest version of the monitoring system operator's manual for more information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.


Refer to the symbol legend at the end of this document.

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
	English	Français	Deutsch	Español	Italiano	Nederlands
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile	Niet-steriel
	Not made with natural rubber latex	Non fabriqué à partir de latex de caoutchouc naturel	Nicht aus Latex hergestellt	No contiene látex de caucho natural.	Non realizzato in lattice di gomma naturale	Bevat geen natuurlijk rubberlatex
 eifu.edwards.com + 1 888 570 4016	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Importer	Importateur	Importeur	Importador	Importatore	Importeur
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo	Unieke instrumentidentificatie
Note: Not all symbols may be included in the labeling of this product. ■ Remarque : il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ Hinweis: Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ Nota: Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ Nota: alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ Opmerking: Het label van dit product bevat misschien niet alle symbolen.						

Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat

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	Polski	Slovensky	Norsk	Suomi	Български	Română
	Niejałowy	Nesterilné	Usteril	Epästeriili	Нестерилно	Nesteril
	Nie zawiera lateksu naturalnego	Neobsahuje prírodný kaučukový latex	Fremstilt uten naturgummilateks	Ei sisällä luonnonkumilateksia	Не е произведено с естествен каучуков латекс.	Nu este fabricat cu latex din cauciuc natural
	Postępować zgodnie z instrukcją stosowania na stronie internetowej	Postupujte podľa návodu na použitie na webovej stránke	Følg bruksanvisningen på nettsiden	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Importer	Dovozca	Importør	Maahantuoja	Вносител	Importator
	Należy przechowywać w chłodnym i suchym miejscu	Skladujte na chladnom a suchom mieste	Oppbevares tørt og kjølig	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat
UDI	Unikalny identyfikator wyrobu medycznego	Jedinečný identifikátor pomôcky	Unik enhetsidentifikator	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului
Uwaga: Na etykiecie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ Poznámka: Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ Merk: Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ Huomautus: kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ Забележка: Възможно е не всички символи да са включени в етикета на този продукт. ■ Notă: este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs.						

**Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları
■ Условные обозначения ■ Legenda sa simbolima**

	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski
#	Mudeli number	Modelio numeris	Modeļa numurs	Model Numarası	Номер модели	Broj modela
MD	Meditsiiniseade	Medicinos priemonė	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo
SN	Seerianumber	Serijos numeris	Sērijas numurs	Seri Numarası	Серийный номер	Serijski broj
LOT	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije
QTY	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina
CE 0123	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)
EC REP	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Įgaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Tootmiskoopäev	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje
	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač
	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu
	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo
	Kõlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok upotrebe
	Elektri- ja elektroonikaseadmete eraldi kogumise koostöös EÜ direktiiviga 2012/19/EL	Atskiras elektros ir elektroninės įrangos atliekų surinkimas pagal EB direktyvą 2012/19/ES	Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES direktīvu 2012/19/ES	2012/19/EU sayılı EC direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır	Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU	Odvojeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU
	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridrżavajte se uputstva za upotrebu
	Ohtlik magnetresonantstomograafias	MR nesaugus	Nedrīkst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR
	Õrn, käsitseada ettevaatlikult	Trapus, elgtis atsargiai	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom

Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları
■ Условные обозначения ■ Legenda sa simbolima

	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski
	Mittesteriilne	Nesterilus	Nesterils	Steril değildir	Нестерильно	Nesterilno
	Pole valmīstatud looduslikust kummilateksist.	Pagaminta be natūralios gumos latekso	Izgatavošanā nav izmantots dabiskā kaučuka latekss	Doğal kauçuk lateksten yapılmamıştır	Не содержит натурального каучукового латекса	Ne sadrži lateks od prirodne gume
 eifu.edwards.com +1 888 570 4016	Jārgīge kasutusjuhendit, mille leiate veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Importija	Importuotojas	Importētājs	İthalatçı	Импортер	Uvoznik
	Sāļītada jāhēdas ja kuīvas kohas.	Laikyti vėsioje, sausoje vietoje	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu
	Seadme kordumatu identititseerimistunnus	Unikalusi priemonės identifikatorius	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva
Mārkus. Kōīk sūmbolid ei pruugi antud toote etiketil esineda. ■ Pastaba. Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ Piezīme. Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ Not: Bu ūrūnūn etiketinē tūm semboller yer almamış olabilir. ■ Примечание. На этикетках данного изделия могут указываться не все символы. ■ Napomena: Na ambalaži proizvoda se možda ne nalaze svi simboli.						

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10/2024
10061272001 A / DOC-0247126 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in China

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH