
# Edwards

## English

### ForeSight Jr Small Sensor

< 8 kg (cerebral), < 5 kg (non-cerebral) – Pediatric

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

#### 1.0 Description

The sensor, when used in combination with the ForeSight oximeter cable or in combination with the FORE-SIGHT ELITE absolute tissue oximeter, is a single use applied part that measures hemoglobin allowing the clinician to continuously and accurately determine absolute levels of blood oxygenation saturation in the tissue (StO<sub>2</sub>). This device is supplied non-sterile.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions for Use.

The device is used by clinicians who have been trained in the use of tissue oximetry devices in accordance with their institutional policies.

The benefit of using ForeSight sensor is to non-invasively provide tissue oxygen saturation and hemoglobin values to allow clinicians to manage their patients. Potential risks include burns, electrical shock, tissue damage, transient hypoxia, adverse reaction to device materials, and/or inappropriate/unintended treatment.

#### 2.0 Intended Use/Purpose

When used with the HemoSphere advanced monitoring platform in combination with the HemoSphere technology module, the intended purpose is to monitor absolute regional hemoglobin oxygen saturation of blood under the sensor.

#### 3.0 Indications for Use

When used in conjunction with the FORE-SIGHT ELITE absolute tissue oximeter or ForeSight oximeter cable:

The small sensor is indicated for monitoring of absolute regional hemoglobin oxygen saturation of blood under the sensor in individuals at risk for reduced flow or no-flow ischemic states. It is intended for cerebral use on pediatric subjects < 8 kg and non-cerebral use on pediatric subjects < 5 kg.

#### 4.0 Contraindications

The sensor is contraindicated for use on patients:

- With a physical site area too limited for proper sensor placement
- With allergic reactions to sensor adhesive
- Undergoing an MRI scan because of associate risk of injury

#### 5.0 Warnings

- **Assess the sensor site at least every 12 hours, or more often as required by the institution's protocol.**
- **Remove the sensor if the circulatory condition or skin integrity has deteriorated.**
- **Do not attach the sensor to damaged or irritated skin.**

- Do not lay patient on the sensor or cable.
- Do not use in an MR environment due to risk of burn as result of sensor heating.
- Do not attach the sensor to skin with unapproved devices, such as headbands, hats, wraps, etc.
- Do not place the sensor or accessories over eyes, nose, or mouth.
- Do not cut the sensor. Cutting the sensor can result in injury to the patient.
- Do not use in an MRI environment.
- When used in settings with LED lighting, sensors should be covered with a light blocker prior to connection to the preamp cable, as some high intensity systems can interfere with the sensor's near infrared light detection.
- The use of barrier film dressings other than Tegaderm could affect the accuracy of the StO<sub>2</sub> readings.
- Failure to remove the sensor protective liner prior to monitoring may cause erroneous StO<sub>2</sub> readings.

## 6.0 Cautions

- Avoid positioning the sensor over hair, air sinus, hematoma, birthmark, or broken skin.
- Avoid attaching the sensor to sites with excess adipose, ascites, or edema.
- The materials used in the manufacture of the sensor were NOT designed for reuse. Reuse can cause the sensor not to perform as intended.

## 7.0 Instructions

| Step | Procedure                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Remove the sensor from the package. Carefully inspect the sensor for damage. Discard and replace if damage is found.                                                                                                                                                                                                                                                                                                                             |
| 2    | Select sensor location on the monitor.                                                                                                                                                                                                                                                                                                                                                                                                           |
| 3    | Clean and dry the sensor site (Figure 1).                                                                                                                                                                                                                                                                                                                                                                                                        |
| 4    | Remove protective liner from the sensor (Figure 2).                                                                                                                                                                                                                                                                                                                                                                                              |
| 5    | Apply the sensor to the patient:<br>a) <b>Cerebral Use</b> (Figure 3a): Select the site on the forehead well above the eyebrow and just below the hairline.<br>b) <b>Non-cerebral Use</b> (Figures 3b, 3c): Select the site that provides ideal access to the desired tissue, for example, latissimus dorsi (flank), external oblique or abdomen.<br><b>Note: You may use Tegaderm under the sensor in patients with delicate skin or edema.</b> |
| 6    | Insert the sensor connector straight into the sensor cable connector until it snaps into place (Figure 4). Use the bedsheet clip to secure the cable and prevent pulling on the sensor.                                                                                                                                                                                                                                                          |
| 7    | If needed, fold the sensor flat cable to route it in the desired direction (Figure 4).                                                                                                                                                                                                                                                                                                                                                           |

For use only with ForeSight oximeter cable software version 2.5.7 or above and with FORE-SIGHT ELITE absolute tissue oximeter software version 4.5.6 or above.

## 8.0 Disposal

Sensors are designed for single-patient use, and are not to be reprocessed. Re-used sensors present a risk of cross-contamination or infection. Use a new sensor for each patient and discard it after use. Disposal should follow in accordance with local hospital and institution policies.

## 9.0 Storage

Store in a cool, dry place.

## 10.0 Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration and malfunction.

## 11.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada  
(24 hours): ..... 800.822.9837  
Outside the U.S. and Canada  
(24 hours): ..... 949.250.2222  
In the UK: ..... 0870 606 2040 - Option 4  
In Ireland: ..... 01 8211012 - Option 4

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**

Refer to the latest version of the monitoring system operator's manual for more information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

**Refer to the symbol legend at the end of this document.**

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
|                                                                                     | Latviešu                                                                                   | Türkçe                                                                                         | Русский                                                                                              | Srpski                                                                                     | Hrvatski                                                                                     |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
|    | Ražotājs                                                                                   | Üretici                                                                                        | Производитель                                                                                        | Proizvođač                                                                                 | Proizvođač                                                                                   |
|    | Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju                        | Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun                              | Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией                 | Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu                   | Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu              |
|    | Nelietot atkārtoti                                                                         | Yeniden kullanmayın                                                                            | Не подлежит повторному использованию                                                                 | Ne koristiti ponovo                                                                        | Nemojte ponovno upotrebljavati                                                               |
|    | Derīguma termiņš                                                                           | Son kullanma tarihi                                                                            | Использовать до                                                                                      | Rok korišćenja                                                                             | Rok upotrebe                                                                                 |
|   | Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES direktīvu 2012/19/ES     | 2012/19/EU sayılı EC direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır | Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU | Odvojeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU  | Odvojeno prikupljanje električne i elektroničke opreme u skladu s direktivom EZ-a 2012/19/EU |
|  | Levērojiēt lietošanas instrukciju                                                          | Kullanım talimatlarını takip edin                                                              | Следуйте инструкциям по применению                                                                   | Pridržavajte se uputstva za upotrebu                                                       | Pridržavajte se uputa za upotrebu                                                            |
|  | Nedrīkst lietot MR vidē                                                                    | MR için Güvenli Değildir                                                                       | Опасно при проведении МРТ                                                                            | Nije bezbedno za MR                                                                        | Nije sigurno kod pregleda MR-om                                                              |
|  | Uzglabāt sausu                                                                             | Kuru halde tutun                                                                               | Беречь от влаги                                                                                      | Držite suvo                                                                                | Držati suhim                                                                                 |
|  | Trausls, rīkoties piesardzīgi                                                              | Kırılabilir, dikkatli taşıyın                                                                  | Хрупкое изделие, обращаться с осторожностью                                                          | Lomljivo, rukovati sa oprezom                                                              | Lomljivo, pažljivo rukovati                                                                  |
|  | Nesterils                                                                                  | Steril değildir                                                                                | Нестерильно                                                                                          | Nesterilno                                                                                 | Nesterilno                                                                                   |
|  | Izgatavošanā nav izmantots dabiskā kaučuka latekss                                         | Doğal kauçuk lateksten yapılmamıştır                                                           | Не содержит натурального каучукового латекса                                                         | Ne sadrži lateks od prirodne gume                                                          | Nije izrađeno od prirodnog gumenog lateksa                                                   |
|  | Levērojiēt lietošanas instrukciju, kas pieejama tīmekļa vietnē                             | Web sitesinde bulunan kullanım talimatlarını takip edin                                        | Следуйте инструкциям по применению на веб-сайте                                                      | Pridržavajte se uputstva za upotrebu sa veb-sajta                                          | Pridržavajte se uputa za upotrebu na internetskoj stranici                                   |
| <b>Rx only</b>                                                                      | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician     | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician           | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician   |

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

|                                                                                   | Latviešu                       | Türkçe                        | Русский                             | Srpski                             | Hrvatski                            |
|-----------------------------------------------------------------------------------|--------------------------------|-------------------------------|-------------------------------------|------------------------------------|-------------------------------------|
|  | Importētājs                    | İthalatçı                     | Импортер                            | Uvoznik                            | Uvoznik                             |
|  | Uzglabāt vēsā, sausā vietā     | Serin, kuru yerde saklayın    | Хранить в прохладном, сухом месте   | Čuvati na hladnom i suvom mestu    | Pohranite na hladnom, suhom mjestu  |
|  | Unikāls ierīces identifikators | Benzersiz cihaz tanımlayıcısı | Уникальный идентификатор устройства | Jedinstveni identifikator sredstva | Jedinstveni identifikator proizvoda |

**Note:** Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámká:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykietach niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámká:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забележка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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06/25  
10042624003 A / DOC-0249412 A  
© Copyright 2025, Edwards Lifesciences LLC  
All rights reserved

**Edwards Lifesciences LLC**  
One Edwards Way  
Irvine, CA 92614 USA  
Made in Malaysia

**Edwards Lifesciences B.V.**  
Verlengde Poolseweg 16  
4818 CL Breda, Netherlands

Telephone 949.250.2500  
800.424.3278  
FAX 949.250.2525

Web IFU

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**Edwards Lifesciences GmbH**