Edwards

Directory

English (en).....	1
Français (fr).....	4
Deutsch (de).....	7
Español (es).....	10
Italiano (it).....	13
Nederlands (nl).....	16
Dansk (da).....	19
Svenska (sv).....	21
Ελληνικά (el).....	23
Português (pt).....	26
Česky (cs).....	29
Magyar (hu).....	31
Polski (pl).....	34
Slovensky (sk).....	37
Norsk (no).....	40
Suomi (fi).....	42
Български (bg).....	44
Română (ro).....	47

Eesti (et).....	50
Lietuvių (lt).....	52
Latviešu (lv).....	54
Türkçe (tr).....	56
Русский (ru).....	58
Srpski (sr).....	61
Hrvatski (hr).....	63
Figures ■ Figures ■ Abbildungen	
■ Figuras ■ Figure ■ Afbeeldingen	
■ Figurer ■ Figurer ■ Εικόνες ■ Figuras	
■ Obrázky ■ Ábrák ■ Rysunki	
■ Obrázky ■ Figurer ■ Kuvat ■ Фигури	
■ Figuri ■ Joonised ■ Paveikslėliai	
■ Attēli ■ Şekiller ■ Рисунки ■ Slike	
■ Slike	65
Symbol Legend ■ Légende des	
symboles ■ Zeichenerklärung	
■ Significado de los símbolos	
■ Legenda dei simboli ■ Lijst met	


symbolen ■ Symbolforklaring	
■ Symbolforklaring ■ Υπόμνημα	
συμβόλων ■ Symbol Legend	
■ Legenda se symboly	
■ Jelmagyarázat ■ Legenda simboli	
■ Vysvětlivky k symbolom	
■ Symbolforklaring ■ Merkkien	
selitykset ■ Легенда на символите	
■ Legenda de simboluri ■ Sūmbolite	
seletus ■ Simbolių paaiškinimas	
■ Simbolu skaidrojums ■ Sembol	
Açıklamaları ■ Условные обозначения	
■ Legenda sa simbolima ■ Legenda	
simbola.....	67

English

Acumen IQ Finger Cuff Adult

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

AIQCA

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Acumen IQ finger cuff adult, when used with an appropriate Edwards monitoring system, provides continuous, noninvasive hemodynamic monitoring and the applicable derived parameters.

Edwards, Edwards Lifesciences, the stylized E logo, Acumen, Acumen IQ, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

The Acumen IQ finger cuff adult utilizes the volume-clamp method to measure blood pressure with an inflatable bladder wrapped around the middle phalanx of the index, middle or fourth/ring finger. This device is supplied non-sterile. Potential risks include inappropriate/unintended treatment, adverse reaction to device materials, patient or clinician burns or electrical shock, peripheral ischemia, and/or minor tissue damage.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

Device is used by clinicians who have been trained in safe use of noninvasive hemodynamic technologies in accordance with their institutional guidelines.

2.0 Intended Use/Purpose

The intended purpose is to noninvasively measure blood pressure and use the information to derive hemodynamic parameters when connected to a HemoSphere advanced monitoring platform.

3.0 Indications

The Acumen IQ finger cuff adult is indicated for patients over 18 years of age to noninvasively measure blood pressure and associated hemodynamic parameters when used with a HemoSphere advanced monitoring platform.

4.0 Contraindications

In some patients with extreme contraction of the smooth muscle in the arteries and arterioles in the lower arm and hand, such as may be present in patients with Raynaud's disease, blood pressure measurement can become impossible.

5.0 Single Patient Use

The Acumen IQ finger cuff adult is designed for single patient use. Upon starting measurement, the finger cuff can be used and re-applied for up to 72 hours on one patient. After 8 hours of continuous monitoring on a single finger, the finger cuff should be re-applied to another finger.

Do not attempt to clean and reuse the finger cuff on more than a single patient, such action could lead to illness or an adverse event as the device may not function as originally intended.

6.0 Warnings

- **Improper placement or alignment of the finger cuff can lead to inaccurate monitoring.**
- **Do not apply the finger cuff on injured skin as this may cause further injury.**
- **To reduce the risk of skin irritation and tissue damage, do not monitor longer than 8 hours continuously on a single finger. To continue to monitor beyond 8 hours, use an additional finger cuff on another finger or move the cuff in use to another finger.**
- **Do not use two finger cuffs simultaneously on the same finger.**
- **Measurement on one finger in contradiction with the instructions for use may affect patient comfort and/or lead to minor injuries.**
- **Do not apply finger cuff(s) on a hand or finger when external constriction (that may prevent circulation to the hand or finger) is present.**

- **Do not apply finger cuff(s) on a hand or finger when a second blood pressure measurement device is actively monitoring on the same arm (or hand or finger).**
- **Do not use the finger cuff with magnetic resonance imaging.**

7.0 Cautions

- **Do not use a damaged finger cuff. This may result in inaccurate measurements or may damage the Edwards monitoring system.**
- **Never bend a finger cuff to a flat shape as it will damage the finger cuff and affect measurement accuracy.**
- **Excessive ambient lighting may interfere with the finger cuff measurement. Avoid using the finger cuff under close, direct lighting.**
- **Always disconnect the finger cuff when it is not applied to a finger, to prevent damage by accidental overinflation.**
- **The effectiveness of the finger cuffs have not been established in pre-eclamptic patients.**

8.0 Instructions for Use

Refer to Figure 1 on page 65 through Figure 3 on page 66 for figures corresponding to the steps below.

8.1 Apply the Finger Cuff

Step	Procedure
1	Gently open the finger cuff and place the finger cuff on the middle phalanx of the index, middle or fourth/ring finger. Ensure the finger cuff is centered between the second and third knuckles, and the two green lines on the inside of the finger cuff (Figure 1). Note: Do not apply the finger cuff on the thumb, small finger, or previously fractured fingers.
2	Align the distal end of the finger with the center line of the alignment tab on the finger cuff (Figure 1).
3	Allow the finger cuff to close around the finger (Figure 2). Ensure that the finger cuff remains aligned and does not rotate.
4	Remove the adhesive backing from the finger cuff and secure around the finger (Figure 2-A). Remove the adhesive backing from the distal end of the alignment tab and secure to the top of the finger cuff (Figure 2-B).
5	Lead the finger cuff cable between two fingers to the back side of the hand (Figure 3).

8.2 Connect the Finger Cuff to the System

Step	Procedure
1	Plug the finger cuff connector into the pressure controller (Figure 3). Refer to the pressure controller instructions for use for more details.
2	If using a heart reference sensor (HRS), attach the finger end of the HRS to the HRS tab on top of the finger cuff (Figure 3). Refer to the HRS instructions for use for more details.
3	Initiate monitoring with the connected Edwards noninvasive monitoring system. To enable advanced parameters, refer to the monitoring system operator's manual for more details.
4	If using double cuff monitoring, repeat all previous steps to apply the second finger cuff.

9.0 Recommendations

The Edwards noninvasive monitoring system will pause blood pressure measurements periodically per the monitoring system configuration settings. At this time, the fingertip should be checked for perfusion and potential tissue damage. The fingertip should also be checked periodically according to hospital protocol.

10.0 MRI Safety

This device is MR-unsafe and poses hazards in the MRI environment. This device contains metallic components, which can experience RF-induced heating in the MRI environment.

11.0 Storage

Store in a cool, dry place.

12.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

14.0 Specification Table

Model	Size Range
AIQCA	43-71 mm

15.0 Warranty

The Acumen IQ finger cuff adult is for single patient use only. The Acumen IQ finger cuff adult is warranted at time of delivery to the end user only. The Acumen IQ finger cuff adult is a non-serviceable part.

16.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): ... 800.822.9837
Outside the U.S. and Canada
(24 hours): ... 949.250.2222
In the UK: ... 0870 606 2040 - Option 4
In Ireland: ... 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Refer to the latest version of the Edwards noninvasive monitoring system operator's manual for more information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

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

Figure 3: Connect the finger cuff to the system.

- Figure 3 : Connecter le manchon de doigt au système.
- Abbildung 3: Fingermanschette am System anschließen.
- Figura 3: Conecte el manguito para el dedo al sistema.
- Figura 3: Collegamento della fascetta per dito al sistema.
- Afbeelding 3: Sluit de vingermanschet aan op het systeem.
- Figur 3: Slut fingermanschetten til systemet.
- Figur 3: Anslut fingermanschetten till systemet.
- Εικόνα 3: Συνδέστε την περιχειρίδα δακτύλου στο σύστημα.
- Figura 3: Conecte a dedeira pletismográfica ao sistema.
- Obrázek 3: Připojte prstovou manžetu k systému.
- 3. ábra: Az ujjmandzsetta csatlakoztatása a rendszerhez
- Rysunek 3: Podłączanie mankietu na palec do systemu.
- Obrázok 3: Pripojenie prstovej manžety k systému.
- Figur 3: Koble fingermansjetten til systemet.
- Kuva 3: Liitä sormimansetti järjestelmään.
- Фигура 3: Свържете маншета за пръст към системата.
- Figura 3: Conectarea manșetei pentru deget la sistem.
- Joonis 3. Ühendage sõrmemansett süsteemiga.
- 3 pav. Piršto manžetės prijungimas prie sistemos.
- 3. attēls. Savienojiet pirksta manšeti ar sistēmu.
- Şekil 3: Parmak manşonunu sisteme bağlayın.
- Рис. 3. Подсоединение пальцевой манжеты к системе
- Slika 3: Spojite manžetnu za prst sa sistemom.
- Slika 3: Povežite naprstak sa sustavom.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
 ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly
 ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları
 ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español
	Fragile, handle with care	Fragile, manipuler avec précaution	Vorsicht! Zerbrechlich!	Frágil; manipular con cuidado
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco
	Manufacturer	Fabricant	Hersteller	Fabricante
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación
LOT	Lot Number	N° du lot	Chargenbezeichnung	Número de lote
	Non-sterile	Non stérile	Unsteril	No estéril
EC REP	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea
QTY	Quantity	Quantité	Menge	Cantidad

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español
	Caution	Avertissement	Vorsicht	Aviso
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM
MD	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario
	Separate collection for electrical and electronic equipment in accordance with EC Directive 2012/19/EU	Tri sélectif des équipements électriques et électroniques conformément à la Directive CE 2012/19/UE	Getrennte Sammlung von Elektro- und Elektronikgeräten gemäß EU-Richtlinie 2012/19/EU	Recogida por separado de equipos eléctricos y electrónicos según la directiva de la CE 2012/19/UE
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)
	Importer	Importateur	Importeur	Importador
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Simbolų seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español
	Model Number	Référence	Modellnummer	Número de modelo
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo

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	Italiano	Nederlands	Dansk	Svenska
	Rappresentante autorizzato nella Comunità Europea/Unione Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen
	Quantità	Hoeveelheid	Mængde	Antal
	Attenzione	Let op	Forsigtig	Var försiktig
	Consultare le istruzioni per l'uso	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen
	Seguire le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Non compatibile con RM	MRI-onveilig	MR-usikker	MR-farlig
	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr	Medicinsk utrustning
	Raccolta separata per attrezzature elettriche ed elettroniche in conformità alla Direttiva CE 2012/19/UE	Gescheiden inzameling voor elektrische en elektronische apparatuur in overeenstemming met EG-richtlijn 2012/19/EU	Separat indsamling af elektrisk og elektronisk udstyr iht. EF-direktivet 2012/19/EU	Separat insamling för elektrisk och elektronisk utrustning i enlighet med EU-direktiv 2012/19/EU
	Conformité Européenne (Marchio CE)	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)
	Importatore	Importeur	Importør	Importör

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

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


	Latviešu	Türkçe	Русский	Srpski
	Pilvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji
	Daudzums	Miktar	Количество	Količina
	Uzmanību!	Dikkat	Предостережение	Oprez
	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu
	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Nedrikt lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo
	Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES Direktīvu 2012/19/ES	2012/19/EU sayılı EC direktifi uyarinca elektrikli ve elektronik ekipmanlar aynı toplanmalıdır	Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU	Odvajeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
 ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly
 ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на
 символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları
 ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski
	Importētājs	İthalatçı	Импортер	Uvoznik
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu
	Modeļa numurs	Model Numarası	Номер модели	Broj modela
	Unikāls ierices identifikators	Benzersiz cihaz tanımlayicisi	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva

	Hrvatski
	Lomljivo, pažljivo rukovati
	Držati suhim
	Proizvođač
	Rok upotrebe
	Nemojte ponovno upotrebljavati

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Hrvatski
	Datum proizvodnje
	Broj serije
	Nesterilno
	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Količina
	Oprez
	Pridržavajte se uputa za upotrebu
 eifu.edwards.com +1 888 570 4016	Pridržavajte se uputa za upotrebu na internetskoj stranici
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Nije sigurno kod pregleda MR-om
	Medicinski uređaj

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
 ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly
 ■ Jelmagyarazat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите
 ■ Legendă de simboluri ■ Simbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembolu Açıklamaları
 ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Hrvatski
	Odvojeno prikupljanje električne i elektroničke opreme u skladu s direktivom EZ-a 2012/19/EU
	Conformité Européenne (oznaka CE)
	Uvoznik
	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Pohranite na hladnom, suhom mjestu
	Broj modela
	Jedinstveni identifikator proizvoda

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Edwards

Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

Made in Dominican Republic
CCT Critical Care Technologies, S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

12/24
10058043002 A / DOC-0255432 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU