Edwards

Directory	
English (en).....	1
Français (fr).....	3
Deutsch (de).....	4
Español (es).....	5
Italiano (it).....	6
Nederlands (nl).....	7
Dansk (da).....	8
Svenska (sv).....	9
Ελληνικά (el).....	10
Português (pt).....	11
Česky (cs).....	12
Magyar (hu).....	13
Polski (pl).....	14
Slovensky (sk).....	15
Norsk (no).....	16
Suomi (fi).....	17
Български (bg).....	18
Română (ro).....	19
Eesti (et).....	20
Lietuvių (lt).....	21
Latviešu (lv).....	22
Türkçe (tr).....	23
Русский (ru).....	24
Srpski (sr).....	25
Hrvatski (hr).....	26

Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike.....	27
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	28

**English**

**Fogarty Corkscrew Catheter**

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

**CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.**

**For Single Use Only**

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

**1.0 Description**

The device is a flexible catheter consisting of a polyvinylchloride outer catheter body surrounding a movable, radiopaque stainless steel spring body. A second radiopaque, spiral-shaped stainless steel cable covered by a latex membrane is attached to the distal end of the catheter. A handle to modulate the spiral cable diameter is attached to the inner spring body at the proximal end of the catheter. The spiral cable design of the Fogarty Corkscrew catheter increases the surface area for entrapping fibrous material.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

This device aids in clot management and is used for restoration of circulatory blood flow to distal arterial system and downstream tissues.

The device is to be used by medical professionals who are trained in the safe use and clinical usage of endovascular catheter technologies as part of their respective institutional guidelines.

**2.0 Intended Use/Purpose**

The Fogarty Corkscrew catheter is intended to be used for the removal of emboli and thrombi from either native vessels or synthetic grafts in the arterial system.

**3.0 Indications**

The Fogarty Corkscrew Catheter is indicated for use in adult patients with emboli/thrombi in the peripheral arterial vasculature in either native vessels or synthetic grafts.

**4.0 Contraindications**

- The Fogarty Corkscrew catheter must not be used outside the arterial system.
- The Fogarty Corkscrew catheter is not recommended for endarterectomy procedures.
- The Fogarty Corkscrew catheter is not intended for use as a vessel dilator.
- The Fogarty Corkscrew catheter must not be used in the venous system.

**5.0 Warnings**

- **Exposure to the atmosphere, handling during insertion, and plaque and other deposits within the blood vessel may cause latex membrane degradation.**
- **To minimize the risk of vessel or membrane damage, do not exceed the maximum recommended pull force (see Catheter Specifications).**

**6.0 Instructions**

**6.1 Technique**

Step	Procedure
1	Inspect the catheter.
2	With the membrane extended, with the control button in the locked position, place the catheter in the artery or graft so that the membrane is past the clot to be removed.

Step	Procedure
3	Contract the membrane by pushing down the control button and pulling it backwards. This causes the membrane to expand into a helical configuration so that it grips or engages the clot.
4	Withdraw the entire catheter from the vessel while the membrane maintains slight contact with luminal wall, removing the occluding material as the catheter is withdrawn. <b>Note: A sensitive touch in controlling the contracted membrane during catheter withdrawal is important so that the membrane diameter can be adjusted to varying luminal diameters. Contraction of the membrane is associated with a feeling of resistance.</b> <b>Note: Because of the complexity and variation in the procedure, the choice of surgical technique, appropriately modified in accordance with the previously described warnings, precautions, and technique, is left to the discretion of the individual surgeon.</b>

**7.0 MRI Information**

This product has not been tested for MRI compatibility. This device is not intended to be used in an MRI environment and should not be used in this environment as it contains metal.

**8.0 Complications**

Potential complications include, but are not limited to, vascular damage, including perforation, vessel rupture, dissection and intimal and/or membrane disruption, bleeding/hemorrhage, local hematomas, distal embolization, arterial embolization, thrombosis, air embolization, local or systemic infection, adverse reaction to device materials, arterial spasm, device breakage requiring intervention to remove, aneurysms and arteriovenous fistula formation.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Edwards, Edwards Lifesciences, the stylized E logo, Corkscrew, Fogarty and Fogarty Corkscrew are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

## 9.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use.

## 10.0 Storage

Store in a cool, dry place.

## 11.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

## 12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada  
(24 hours): ..... 800.822.9837  
Outside the U.S. and Canada  
(24 hours): ..... 949.250.2222  
In the UK: ..... 0870 606 2040 - Option 4  
In Ireland: ..... 01 8211012 - Option 4

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**

## 13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

**WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.**

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

## 14.0 Specification Table

Model	140806	140808	1408010
French Size	4 F	5 F	6 F
(mm)	1.33	1.67	2.0
Usable Length (cm)	80	80	80
Maximum Diameter of Extended Membrane (mm)	2.7	3.1	3.3
Diameter of Contracted Membrane (mm)	6	8	10
Maximum Recommended Pull Force on Contracted Membrane (lb)	0.9 (0.4 kg)	1.3 (0.6 kg)	1.8 (0.8 kg)

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

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Afmetingen	Størrelse	Storlek	Μέγεθος	Tamanho
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Op een koele en droge plaats bewaren.	Skal opbevares køligt og tørt.	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο.	Guardar num local fresco e seco.
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστεριώνετε	Não voltar a esterilizar
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogénico
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marcação CE)
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό	Manter seco
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo

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




	Suomi	Български	Română	Eesti	Lietuvių
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks	Vienkartinio naudojimo
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Katso käyttöohjeet	Направте справка с инструкциите за употреба	Consultați instrucțiunile de utilizare	Palun lugege kasutusjuhiseid!	Žr. naudojimo instrukcijas
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Katso käyttöohjeet verkkosivustolta	Направте справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Viimeinen käyttöpäivä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Naudoti iki nurodytos datos
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Sisältää luonnokumilateksia tai sen jäämiä	Съдържание или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural	Sisaldab naturaalsed kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/ Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Koko	Размер	Mărime	Suurus	Dydis
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Säilytettävä kuivassa ja viileässä.	Да се съхранява на хладно, сухо място.	Stocați într-un loc rece și uscat.	Säilitada jahedas ja kuivas kohas.	Laikykite vésioje, sausoje vietoje.
	Älä steriloitu uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepirogeenne	Nepirogeniškas
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcăj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направте справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarazás ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
	Pidä kuivana	Да се пази сухо	A se păstra uscat	Hoida kuivalt	Laikykite sausoje vietoje
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusis priemonės identifikatorius

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjljiva dužina	Upotrebjljiva dužina
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
	Daudzums	Miktar	Количество	Količina	Količina
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Skatīt lietošanas instrukciju	Kullanım talimatlarına başvurun	См. инструкции по применению	Pogledajte uputstvo za upotrebu	Pogledajte upute za upotrebu
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidadom	Sterilizirano etilen-oksidadom
	Skatīt lietošanas instrukciju timekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na internetskoj stranici
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj uniji
	Izmērs	Boyut	Размер	Veličina	Veličina
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
<b>Rx only</b>	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Uzglabāt vēsā, sausā vietā.	Serin, kuru yerde saklayın.	Хранить в прохладном, сухом месте.	Čuvati na hladnom i suvom mestu.	Pohranite na hladnom, suhom mjestu.
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanim talimatlarını başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda

**Note:** Not all symbols may be included in the labeling of this product. ■ **Remarque:** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykietcie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus:** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminių etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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10/25  
10062731001 A / DOC-0256860 A  
© Copyright 2025, Edwards Lifesciences LLC  
All rights reserved.

  
**Edwards Lifesciences LLC**  
One Edwards Way  
Irvine, CA 92614 USA  
Made in USA

**Edwards Lifesciences B.V.**  
Verlengde Poolseweg 16  
4818 CL Breda, Netherlands

Telephone 949.250.2500  
800.424.3278  
FAX 949.250.2525

Web IFU

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**Edwards Lifesciences GmbH**