Edwards

Directory	
English (en).....	1
Français (fr).....	3
Deutsch (de).....	4
Español (es).....	5
Italiano (it).....	6
Nederlands (nl).....	7
Dansk (da).....	8
Svenska (sv).....	9
Ελληνικά (el).....	10
Português (pt).....	11
Česky (cs).....	12
Magyar (hu).....	13
Polski (pl).....	14
Slovensky (sk).....	15
Norsk (no).....	16
Suomi (fi).....	17
Български (bg).....	18
Română (ro).....	19
Eesti (et).....	20
Lietuvių (lt).....	21
Latviešu (lv).....	22
Türkçe (tr).....	23
Русский (ru).....	24
Srpski (sr).....	25
Hrvatski (hr).....	26
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Simbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	27

English

Fogarty Venous Thrombectomy Catheter

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

Fogarty venous thrombectomy catheters have a long, soft tip which permits advancement past the venous valves without undue trauma. The flexibility required for the venous procedure is supplied by the spring-wound body. The spring is covered with a braiding for additional strength. The gate valve is designed for single-handed operation and to minimize the possibility of leakage. The gate valve has an arrow to indicate the "open" and "closed" positions. There are three available catheter sizes (the 6 French (2.0 mm), 8 French (2.7 mm), and 8/10 French (2.7/3.3 mm)). Each may be quickly identified by the color-coded body. The size and filling capacity are printed on each catheter.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

This device aids in clot management and restoration of circulatory blood flow to veins of the non-central circulatory system and temporary occlusion of blood vessels.

The device is intended to be used by medical professionals who are trained in the safe use and clinical usage of endovascular catheter technologies as part of their respective institutional guidelines.

2.0 Intended Use/Purpose

Fogarty venous thrombectomy catheters are intended for use in the performance of venous thrombectomy and for the temporary occlusion of blood vessels.

3.0 Indications

The Fogarty Venous Thrombectomy Catheters are indicated for use in the peripheral venous vasculature of adult patients for the

removal of emboli and thrombi, and may be used for temporary vessel occlusion.

4.0 Contraindications

Venous thrombectomy catheters are not designed to withstand the additional strain encountered during endarterectomy procedures.

5.0 Warnings

- **Air should not be used for inflation in instances where balloon rupture could produce a dangerous air embolus.**
- **The maximum recommended volume should not be exceeded, as overinflation increases the possibility of balloon rupture. (See specification table).**
- **Undue manipulation or overdistension of the balloon can result in significant damage to the vessel lining, increasing the potential for recurrent thrombosis.**
- **Experience, along with the realization of the limitations of the instrument, have been the most significant factors in reducing the incidence of complications.**

6.0 Instructions

6.1 Balloon Inflation

The balloon should be inflated with a sterile, blood-compatible fluid which may include radiopaque solutions if they are highly dilute and nonparticulate.

Carbon dioxide should be substituted in cases where gas is used as the inflation fluid. However, it should be recognized that CO₂ will readily diffuse through the balloon wall, and corrections for gas loss must be made.

Use of a highly viscous or highly particulate contrast medium is not recommended for balloon inflation as the inflation passage may become occluded. Prior to use, the catheter, with the balloon inflated to its maximum recommended volume, should be carefully inspected by the surgeon.

The smallest syringe which will hold the stated maximum fluid capacity should be used, as the smaller diameter facilitates control of the withdrawal force because greater displacement of the plunger is required to make a given change in a balloon diameter.

The amount of fluid in the syringe should be checked before each inflation. If it is more than the maximum capacity, the syringe should be removed from the catheter and the barrel refilled to the proper volume, taking into account the amount of fluid remaining in the lumen of the catheter.

Inflation of the balloon is associated with a feeling of resistance. When no resistance is encountered, it should be assumed that the balloon has ruptured. Inflation should be discontinued at once, and the catheter withdrawn.

7.0 MRI Information

This product has not been tested for MRI compatibility.

8.0 Complications

As with all catheterization procedures, complications may occur. Incidents such as distal embolization, sepsis/infection, bleeding, blood loss, local hematoma, tissue/vascular trauma, soft tissue trauma/ injury, air emboli, adverse reaction to device materials, ischemia and arteriovenous fistula have been reported with the use of thrombectomy catheters.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

9.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

- Inside the U.S. and Canada (24 hours): 800.822.9837
- Outside the U.S. and Canada (24 hours): 949.250.2222
- In the UK: 0870 606 2040 - Option 4
- In Ireland: 01 821 1012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Edwards, Edwards Lifesciences, the stylized E logo, and Fogarty are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE

this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences GmbH
Edwards Lifesciences LLC
Edwards Lifesciences B.V.