Edwards


Directory

English (en).....	2
Français (fr).....	4
Deutsch (de).....	7
Español (es).....	10
Italiano (it).....	13
Nederlands (nl).....	16
Dansk (da).....	19
Svenska (sv).....	21
Ελληνικά (el).....	23
Português (pt).....	26
Česky (cs).....	29
Magyar (hu).....	31
Polski (pl).....	34
Slovensky (sk).....	37
Norsk (no).....	40
Suomi (fi).....	42
Български (bg).....	45
Română (ro).....	48
Eesti (et).....	51
Lietuvių (lt).....	53
Latviešu (lv).....	55
Türkçe (tr).....	57
Русский (ru).....	59
Srpski (sr).....	62
Hrvatski (hr).....	65
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	67

Fogarty Arterial Embolectomy Catheter

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

FG1002

Fogarty Arterial Embolectomy Catheter

1. Stylet (applicable models only)
2. Inflation Hub
3. Catheter Shaft
4. Latex Balloon

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For single use only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Fogarty embolectomy catheters provide a means of clearing emboli and thrombi from the arterial circulatory system through inflation of the balloon and engagement with the arterial wall beyond the vascular obstruction followed by gentle withdrawal of the catheter thereby removing the obstruction from its position.

To remove fibrous or adherent material, alternative devices such as the Fogarty adherent clot and graft thrombectomy catheter are recommended.

The device is intended to be used by medical professionals (surgeons, interventional clinicians such as radiologists), who are trained in the safe use and clinical usage of endovascular catheter technologies as part of their respective institutional guidelines.

This device aids in clot management and facilitates the removal of emboli and thrombi in the arterial vasculature of the non-central circulatory system.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

2.0 Intended Use/Purpose

The Fogarty arterial embolectomy catheters are intended for the removal of fresh, soft emboli and thrombi from vessels in the arteries of the noncentral circulatory system.

3.0 Indications

The Fogarty Arterial Embolectomy Catheters are indicated for use in adult patients for the removal of fresh, soft emboli and thrombi from vessels in the peripheral arterial system.

4.0 Contraindications

- The arterial embolectomy catheter should not be used outside the arterial system.
- The arterial embolectomy catheter is not recommended for the removal of fibrous, adherent, or calcified material (e.g. chronic clot, atherosclerotic plaque). The catheter is not designed to withstand the additional pull force needed to remove these materials.
- The arterial embolectomy catheter is not intended for use as a vessel dilator.
- The arterial embolectomy catheter should not be used in the venous system.

5.0 Warnings

- Balloon rupture and catheter separation as a result of excessive pull force applied to remove adherent material are the most frequent causes of reported failures. The possibility of balloon rupture must be taken into account when considering the risks involved in any embolectomy procedure.
- To minimize the risk of vessel damage, balloon rupture, or tip detachment, do not exceed the maximum recommended inflation volume and pull force for each size catheter (see **Specification Table**).

- Use of a highly viscous or particulate contrast medium is not recommended for balloon inflation because the inflation lumen may become occluded.
- To avoid air embolus in case of balloon rupture, air should not be used for balloon inflation. Carbon dioxide is the only recommended gas. See Instructions.

6.0 Instructions

6.1 Prepare

Using sterile technique, remove the cap or cap and stiffening stylet from the catheter hub.

6.2 Purge

Inflate the balloon with sterile fluid or gas to the maximum recommended volume. (It is recommended that the 2F (0.67 mm) catheter be inflated with carbon dioxide gas.) Pull a vacuum on the syringe. Repeat until all air is removed.

Note: For all inflations use the smallest syringe capable of holding the stated maximum fluid capacity.

6.3 Inspect

The catheter should be inspected with the balloon inflated during purging. A balloon that does not inflate, leaks, or inflates in a grossly asymmetric (eccentric) manner should not be used.

CAUTION: The amount of fluid in the syringe should be checked before each inflation. Do not exceed the recommended maximum inflation volume. See the specification table.

6.4 Place

With the balloon deflated, insert the catheter into the vessel and beyond the obstructing material.

6.5 Inflate

With the catheter suitably positioned, inflate the balloon with sterile fluid or gas. Balloon inflation should be stopped when the balloon can be felt to engage the arterial wall.

Note: The 2F (0.67 mm) catheter should be inflated with carbon dioxide gas. Gas seepage through the intact balloon will require frequent adjustment of the inflation volume.

CAUTION: To minimize lateral wall pressures and shear forces on the inner surface of the artery, use the smallest inflated balloon diameter that will remove the obstructing material.

CAUTION: Air should not be used for inflation in instances where balloon rupture could produce a dangerous air embolus.

6.6 Withdraw

Remove the occlusive material by gently withdrawing the catheter. During withdrawal, it is important to adjust the balloon diameter to the varying arterial diameters by controlling the inflation volume.

CAUTION: Do not exceed the Maximum Recommended Pull Force. (See Specification Table).

CAUTION: Inflation of the balloon is associated with a feeling of resistance to the fluid or gas injection. When no resistance is encountered it should be assumed that the balloon has ruptured. Discontinue inflation and remove the catheter at once.

7.0 MRI Information

This product has not been tested for MRI compatibility.

8.0 Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection/sepsis, adverse reaction to device materials, blood loss, local hematomas, hemorrhage, intimal disruption, arterial dissection, perforation and vessel rupture, arterial thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm, arteriovenous fistula formation, and balloon rupture with fragmentation, tip separation and distal embolization.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

9.0 How Supplied

Contents sterile and nonpyrogenic if package is undamaged or unopened. Do not use if package is opened or damaged. Do not resterilize. Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): . . 800.822.9837
Outside the U.S. and Canada
(24 hours): . . 949.250.2222
In the UK: . . . 0870 606 2040 - Option 4
In Ireland: . . . 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

STERILE	EO
---------	----

14.0 Specification Table

Catheter French Size (mm)	2 F 0.67	3 F 1.00	4 F 1.33	5 F 1.67	6 F 2.0	7 F 2.3
Maximum Liquid Capacity (ml)	N/A	0.2	0.75	1.5	2.0	2.5
Maximum Gas Capacity (ml)	0.2	0.6	1.7	3.0	4.5	5.0
Maximum Recommended Pull Force on Inflated Balloon (lb)	0.5	0.7	1.5	2.0	2.5	3.5
(kg)	0.23	0.3	0.68	0.9	1.13	1.59
Diameter of Inflated Balloon (mm)	4	5	9	11	13	14
Maximum French Size of Deflated Balloon (mm)	3.4 F 1.13	4.3 F 1.43	5.0 F 1.67	6.0 F 2.0	7.0 F 2.3	8.0 F 2.7

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
	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Balloncapaciteit	Ballonkapacitet	Ballongkapacitet	Χωρητικότητα μπαλονιού	Capacidade do balão
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig	Indeholder eller har spor af naturlig gummi-latex	Innehåller eller har spår av naturgummilatem	Περιέχει λάτεξ φυσικού καουτσούκ	Contém ou está presente látex de borracha natural
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης στον ιστότοπο	Consultar as instruções de utilização no site
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Afmetingen	Størrelse	Storlek	Μέγεθος	Tamanho
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrónicas
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogénico
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marcação CE)
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barriersystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

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




	Suomi	Български	Română	Eesti	Lietuvių
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Pallon tilavuus	Вместимост на балона	Capacitatea balonului	Ballooni mahutavus	Balionėlio talpa
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Sisältää luonnokumilateksia tai sen jäämiä	Съдържание или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural	Sisaldab naturaalsed kummilateksid	Turi savo sudėtyje ar yra natūralios gumos latakso
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Koko	Размер	Mărime	Suurus	Dydis
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjljiva dužina	Upotrebjljiva dužina
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Nelietot atkārtoti	Yeniden kullanmayın	Не использовать повторно	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Daudzums	Miktar	Количество	Količina	Količina
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilenoksidom	Sterilizirano etilenoksidom
	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na internet-skoj stranici

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Izmērs	Boyut	Размер	Veličina	Veličina
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilizirati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одиная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykietcie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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10/25
10062754001 A / DOC-0256884 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH