Edwards

Directory	
English (en).....	1
Français (fr).....	5
Deutsch (de).....	8
Español (es).....	11
Italiano (it).....	14
Nederlands (nl).....	17
Dansk (da).....	20
Svenska (sv).....	23
Ελληνικά (el).....	26
Português (pt).....	29
Česky (cs).....	32
Magyar (hu).....	35
Polski (pl).....	38
Slovensky (sk).....	41
Norsk (no).....	44
Suomi (fi).....	47
Български (bg).....	50
Română (ro).....	53
Eesti (et).....	56
Lietuvių (lt).....	59
Latviešu (lv).....	62
Türkçe (tr).....	65
Русский (ru).....	68
Srpski (sr).....	71
Hrvatski (hr).....	74

Symbol Legend ■ Légende des symboles	
■ Zeichenerklärung ■ Significado de los	
símbolos ■ Legenda dei simboli ■ Lijst	
met symbolen ■ Symbolforklaring	
■ Symbolförklaring ■ Υπόμνημα	
συμβόλων ■ Legenda de símbolos	
■ Legenda se symboly ■ Jelmagyarázat	
■ Legenda symboli ■ Vysvetlivky	
k symbolom ■ Symbolforklaring	
■ Merkkien selitykset ■ Легенда на	
символите ■ Legendă de simboluri	
■ Sümbolite seletus ■ Simbolių	
paaiškinimas ■ Simbolu skaidrojums	
■ Sembol Açıklamaları ■ Условные	
обозначения ■ Legenda sa simbolima	
■ Legenda simbola.....	77

English

Fogarty Thru-Lumen Embolectomy Catheter

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The device is a double lumen catheter with a latex balloon at the distal end. One lumen is used for balloon inflation and is accessed through the gate valve at the proximal end. The thru-lumen is used for infusion

of fluids, sampling of blood or guidewire passage through the distal tip. A removable stainless steel stylet(s) is provided.

To remove fibrous or adherent material, alternative devices such as the Fogarty corkscrew catheter and Fogarty graft thrombectomy catheter are recommended.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when

Edwards, Edwards Lifesciences, the stylized E logo, and Fogarty are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

used in accordance with the established Instructions For Use.

This device aids in clot management, temporary occlusion, infusion, and aspiration.

The device has not yet been tested in pediatric population or in pregnant or lactating women.

Intended users include medical professionals who are trained in the safe use and clinical usage of endovascular catheter technologies as part of their respective institutional guidelines.

2.0 Intended Use/Purpose

The intended purpose is for removal of fresh, soft emboli from vessels in the arteries of the non-central circulatory system, temporary occlusion of blood vessels, infusion of fluids, and blood sampling.

3.0 Indications

The Fogarty Thru-Lumen Embolectomy Catheter is indicated for adults requiring the removal of fresh, soft emboli and thrombi from vessels in the arteries of the peripheral vasculature. Additionally, it may also be used for temporary occlusion of blood vessels, infusion of fluids, and blood sampling in general surgical procedures.

4.0 Contraindications

- The thru-lumen embolectomy catheter should not be used outside the arterial system or for endarterectomy procedures.
- The thru-lumen embolectomy catheter is not recommended for the removal of fibrous, adherent, or calcified material (e.g., chronic clot, atherosclerotic plaque). The catheter is not designed to withstand the additional pull force needed to remove these materials.
- The thru-lumen embolectomy catheter is not intended for use as a vessel dilator.
- The thru-lumen embolectomy catheter should not be used in the venous system.

5.0 Warnings

- **Balloon rupture and catheter separation as a result of excessive pull force applied to remove adherent material are the most frequent causes of reported failures. The possibility of balloon rupture must be taken into account when considering the risks involved in any embolectomy procedure.**
- **To minimize the risk of vessel damage, balloon rupture, or tip detachment, do not exceed the maximum recommended inflation volume and pull force for each size catheter (see Specification Table).**

- **Use of a highly viscous or particulate contrast medium is not recommended for balloon inflation because the inflation lumen may become occluded.**
- **To avoid air embolus in case of balloon rupture, air should not be used for balloon inflation. Carbon dioxide is the only recommended gas. See Instructions.**

6.0 Package Contents

Each Catheter is packaged with an off-the-shelf syringe in the tray for balloon inflation.

7.0 Instructions

Precaution: To minimize lateral wall pressures and shear forces on the intimal lining of the vessel wall, the smallest diameter balloon catheter that will remove the embolus should be used.

7.1 Prepare

Remove all stylets from the catheter.

7.2 Purge

Inflate the balloon with sterile fluid or gas to the maximum recommended volume. Pull a vacuum on the syringe. Repeat until all air is removed.

Note: For all inflations, use the smallest syringe capable of holding the stated maximum fluid capacity. Attach another syringe filled with sterile, heparinized saline or other sterile solution to the thru-lumen port and purge air from the thru-lumen.

7.3 Inspect

The catheter should be inspected with the balloon inflated during purging. A balloon that does not inflate, leaks, or inflates in a grossly asymmetric (eccentric) manner should not be used.

CAUTION: The amount of fluid in the syringe should be checked before each inflation. Do not exceed the recommended maximum inflation volume. See Specification Table.

7.4 Place

With the balloon deflated, insert the catheter into the vessel and beyond the obstructing material either directly by cutdown or by percutaneous technique through a suitable introducer. Heparinized saline or other solutions may be flushed through the catheter thru-lumen as needed during the procedure. Do not advance the catheter against resistance.

7.5 Inflate

With the catheter suitably positioned, inflate the balloon with sterile fluid or gas. Balloon inflation should be stopped when the balloon can be felt to engage the vessel wall.

Note: CO₂ is the preferred medium in cases where gas is used for inflation.

It should be noted that CO₂ diffuses through the balloon wall and corrections for gas loss must be made. Use of liquid inflation media will result in slower inflation/deflation rates.

CAUTION: To minimize lateral wall pressures and shear forces on the inner surface of the artery, use the smallest inflated balloon diameter that will remove the obstructing material.

7.6 Withdraw

Remove the occlusive material by gently withdrawing the catheter. During withdrawal, it is important to adjust the balloon diameter to the varying arterial diameters by controlling the inflation volume.

CAUTION: Do not exceed the Maximum Recommended Pull Force. See Specification Table.

CAUTION: Inflation of the balloon is associated with a feeling of resistance to the fluid or gas injection. When no resistance is encountered it should be assumed that the balloon has ruptured. Discontinue inflation and remove the catheter at once.

8.0 MRI Information

This product has not been tested for MRI compatibility.

9.0 Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection/sepsis, adverse reaction to device materials, bleeding, local hematomas, hemorrhage, intimal disruption, thrombosis, arterial dissection, perforation and vessel rupture, arterial thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm, arteriovenous fistula formation, and balloon rupture with fragmentation, tip separation and distal embolization.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

10.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use.

11.0 Storage

Store in a cool, dry place.

12.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

13.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
 (24 hours): 800.822.9837
 Outside the U.S. and Canada
 (24 hours): 949.250.2222
 In the UK: 0870 606 2040 - Option 4
 In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

14.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to

support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

15.0 Specification Table

Catheter French Size (mm)	3 1.00	4 1.33	5.5 1.83	6 2.0	7 2.3
Maximum Liquid Capacity (ml)	0.15	0.5	0.9	1.4	1.6
Maximum Gas Capacity (ml)	0.4	1.2	1.7	2.3	2.8
Maximum Recommended Pull Force on Inflated Balloon (lb)	0.7	1.5	2.0	2.5	3.5
Maximum Recommended Pull Force on Inflated Balloon (kg)	0.32	0.68	0.91	1.13	1.59
Diameter of Inflated Balloon (mm)	5	9	11	13	14
Maximum French Size of Deflated Balloon (mm)	4.5	5	6	7	8
Maximum Size of Deflated Balloon (mm)	1.50	1.67	2.0	2.3	2.7
Compatible Guidewire Diameter (in)	0.018	0.025	0.035	0.035	0.038
Compatible Guidewire Diameter (mm)	0.46	0.64	0.89	0.89	0.97
* Compatible Introducer Size (F)	4.5	5	6	7	8
* Compatible Introducer Size (mm)	1.50	1.67	2.0	2.3	2.7

* Ease of introducing the catheter through hemostasis valve may vary among manufacturers.

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


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marcação CE)
	Vervaldatum	Dato for sidste anvendelse	Sista förbrukningsdag	Ημερομηνία λήξης	Data de vencimento
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Balloncapaciteit	Ballonkapacitet	Ballongkapacitet	Χωρητικότητα μπαλονιού	Capacidade do balão
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/ Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Spuit	Sprøjte	Spruta	Σύριγγα	Seringa
	Afmetingen	Størrelse	Storlek	Μέγεθος	Tamanho
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogênico
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig	Indeholder eller har spor af naturlig gummlatex	Innehåller eller har spår av naturgummilatex	Περιέχει λάτεξ φυσικού καουτσούκ	Contém ou está presente látex de borracha natural
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης στον ιστότοπο	Consultar as instruções de utilização no site

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






	Suomi	Български	Română	Eesti	Lietuvių
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Ruisku	Спринцовка	Seringă	Süstal	Švirkštas
SZ	Koko	Размер	Mărime	Suurus	Dydis
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
LOT	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Sisältää luonnokumilateksia tai sen jäämiä	Съдържание или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural	Sisaldab naturaalsed kummlateksit	Turi savo sudėtyje ar yra natūralios gumos latekso
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
UDI	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба в уебсайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege kasutusjuhiseid veebisaidilt	Žr. naudojimo instrukcijas interneto svetainėje

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjliva dužina	Upotrebjliva dužina
GW	Ieteicamais vadītājstīgas izmērs	Kılavuz Tel için Önerilen Boyut	Рекомендованный размер проводника	Preporučena veličina žice vodice	Preporučena veličina žice vodilice

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Minimālais ievadītāja izmērs	Minimum introdüser boyutu	Минимальный размер интродьюсера	Minimalna veličina uvodnika	Minimalna veličina uvodnog instrumenta
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Nelietot atkārtoti	Yeniden kullanmayın	Не использовать повторно	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Daudzums	Miktar	Количество	Količina	Količina
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksидом	Sterilizirano etilen-oksидом
	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Pilnvarotais pārstāvis Eiropas Kopienā/ Eiropas Savienībā	Avrupa Birliği'nde/ Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Šīrce	Şırınga	Шприц	Špric	Šprica

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmērs	Boyut	Размер	Veličina	Veličina
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одиная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na internet stranici	Pogledajte upute za upotrebu na internetskoj stranici

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítku tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykietcie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забележка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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10/25
10062756001 A / DOC-0256888 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH