Directory

English (en)	1	Português (pt)	21	Română (ro)	38
Français (fr)	3	Česky (cs)	23	Eesti (et)	40
Deutsch (de)	6	Magyar (hu)	25	Lietuvių (lt)	42
Español (es)	8	Polski (pl)	27	Latviešu (lv)	44
Italiano (it)	11	Slovensky (sk)	30	Türkçe (tr)	46
Nederlands (nl)	13	Norsk (no)	32	Русский (ru)	48
Dansk (da)	15	Suomi (fi)	34	Srpski (sr)	51
Svenska (sv)	17	Български (bg)	36	Hrvatski (hr)	53
Ελληνικά (el)	19				

Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised
■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike55

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων
■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Лєгєндє нє символитє ■ Legendă de simboluri
■ Sūbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola56

English

VolumeView™ Set

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

For single use only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 VolumeView™ Set Components

- Edwards VolumeView™ femoral arterial catheter with insertion kit
- Edwards VolumeView™ venous injectate kit: TruWave™ disposable pressure transducer (DPT) and VolumeView™ thermistor manifold.

2.0 Additional Components Required

- Cooled normal saline filled luer lock syringe(s)
- Central venous catheter (CVC)
- EV1000™ clinical platform or HemoSphere Alta™ Advanced Monitoring Platform
- Edwards VolumeView™ sensor or Acumen IQ™ sensor

Refer to the EV1000™ Clinical Platform or HemoSphere Alta™ Advanced Monitoring Platform operator's manual for sensor compatibility and cable connections.

These are general instructions for setting up the disposable components of the VolumeView™ system. Kit configurations and procedures may vary according to institutional preferences.

3.0 Intended Use/Purpose

The VolumeView™ set, when used in conjunction with a compatible hemodynamic platform, is designed for measurement of continuous and intermittent transpulmonary thermodilution (TPTD) parameters. The VolumeView™ set enables measurement of intermittent cardiac output (iCO) and a range of volumetric parameters including extravascular lung water (EVLW) and global end-diastolic volume (GEDV).

4.0 Description

The VolumeView™ set from Edwards Lifesciences in combination with the applicable Edwards hemodynamic monitoring platform enables the display of valuable volumetric parameters.

The benefits of hemodynamic monitoring may lead to improved data driven clinical decision making for medically necessary intervention and/or clinical re-evaluation.

Device risks include thrombosis, sepsis/infection, bleeding, peripheral ischemia, inappropriate/unintended treatment, hematoma, tissue damage, vessel perforation, embolism, adverse reaction to device materials, and/or delay in treatment.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

This device is used by medical professionals who have been trained in the safe use of invasive hemodynamic technologies in accordance with their institutional guidelines.

5.0 MRI Information

MR Conditional

Non-clinical testing has demonstrated that the VolumeView™ catheter is MR Conditional.

It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla (T) or 3.0 Tesla (T)
- Highest spatial gradient field of 720 gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning.
- Normal mode operation, as defined in IEC 60601-2-33, of the MR system.

In non-clinical testing and analysis, the VolumeView™ catheter was determined to produce a maximum temperature rise of 2.0 °C above background for a SAR of 2.0 W/kg for 15 minutes of scanning in a 1.5 T whole body cylindrical bore MR system. The projected rise above background was 2.4 °C for a SAR of 2.0 W/kg in a 3.0 T whole body cylindrical bore MR system. These calculations overestimate the true *in vivo* rise, since the cooling effects of blood are not considered.

The maximum image artifact extends as far as 10 mm from the catheter for spin echo images and gradient echo images.

The catheter has not been evaluated in MR systems other than 1.5 T or 3.0 T.

Non-clinical testing has demonstrated that the VolumeView™ venous injectate kit is MR Conditional according to the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 gauss/cm or less

Note: The VolumeView™ injectate kit has minor metallic components that do not come into contact with patients during the intended use. Furthermore, the metallic components are insulated and isolated from the patient during their intended use and as such, it was only necessary to evaluate this product for magnetic field interactions. MRI-related heating and artifacts were not characterized for this product because it is not a medical implant.

6.0 Indications

The VolumeView™ set when used with a compatible hemodynamic monitor, is indicated for use primarily for critical care patients in which the balance between cardiac function, fluid status and vascular resistance needs continuous or intermittent assessment. Analysis of the thermodilution curve in terms of mean transit time and the shape is used to determine intravascular and extravascular fluid volumes. When connected to an Edwards oximetry catheter, the monitor measures oximetry in adults and pediatrics. The VolumeView™ set may be used in all settings in which critical care is provided.

The VolumeView™ catheter is indicated for femoral artery insertion.

7.0 Contraindications

The VolumeView™ femoral arterial catheter is contraindicated for venous insertion. Relative contraindications include: placement of a central venous catheter or of a femoral arterial catheter in patients who have severe coagulation abnormalities, burned or infected skin at insertion point, or severe atherosclerosis.

8.0 Warnings

- **Please read this manual before use. Retain for future reference.**
- **Use sterile technique when connecting the fluid path components of the system.**
- **Do not use if expiration date is passed.**
- **This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.**

9.0 Cautions

- **Femoral catheter usage may be associated with thrombotic, infectious or hemorrhagic complications and distal limb ischemia.**
- **Assess limb distal to femoral arterial catheter before and after placement for potential compromised limb perfusion per institutional policies.**
- **You may receive erroneous results if pressure transducer is not zeroed.**

- **Erroneous measurements can be caused by an incorrect catheter position or by electromagnetic disturbance (such as an electric blanket).**
- **Repeat the measurement if the results are implausible.**
- **Do not withdraw a PTFE-coated guidewire through a metal cannula needle as this may damage the guidewire coating.**
- **Some models may contain the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0. Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.**
- **Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.**

10.0 Instructions

Refer to individual components' instructions for use: VolumeView™ sensor or Acumen IQ™ sensor and central venous catheter.

Note: A central venous catheter with a minimum burst pressure of 130 PSI (896 kPa) is recommended.

10.1 Insert VolumeView™ Femoral Arterial Catheter

Step	Procedure
1	Recommendation: Set up and prime VolumeView™ sensor or Acumen IQ™ sensor before insertion of VolumeView™ femoral arterial catheter.
2	Ensure that all connections are secure.
3	The VolumeView™ femoral arterial catheter can be introduced using standard percutaneous Seldinger technique. - Prep and drape puncture site per institutional policy. - Administer local anesthetic per clinician discretion. - Locate femoral artery through palpation or ultrasound guidance. - Introduce 18ga (1.2 mm) thin-wall needle into artery. Confirm placement by observing bright red pulsatile blood flow. - Insert the guidewire through the 18ga (1.2 mm) thin-wall needle. Gentle manipulation may be needed to insert the guidewire. The guidewire should never be forced. Remove the 18ga (1.2 mm) thin-wall needle, leaving the guidewire in place. - A skin nick may be made with the scalpel. - The insertion site may be enlarged by threading the dilator over the guidewire. - Holding the guidewire in place, remove the dilator and thread the VolumeView™ catheter over the guidewire. - Remove the guidewire and confirm placement by observing bright red pulsatile blood flow or presence of arterial waveform on the bedside monitor once connected to the VolumeView™ or Acumen IQ™ sensor.
4	Connect catheter to pre-primed VolumeView™ sensor, Acumen IQ™ sensor, or arterial pressure transducer, then aspirate and flush system confirming patency and that the system is free from air.
5	Once in position, after the guidewire has been removed, secure catheter in place with suture or suitable anchoring device.
6	Confirm appropriate arterial waveform and pressures on patient monitor(s).
7	Connect VolumeView™ femoral arterial catheter's thermistor to the VolumeView™ temperature cable or Patient CCO Cable-HemoSphere Alta™ Advanced Monitoring Platform.
8	Confirm appropriate temperature is displayed on the EV1000™ monitor or HemoSphere Alta™ Advanced Monitor.
9	Level and zero the VolumeView™ or Acumen IQ™ sensor as per the Edwards EV1000™ monitor or HemoSphere Alta™ Advanced Monitor instructions.

10.2 Maintenance and Use of the VolumeView™ Femoral Arterial Catheter *in situ*

Catheter insertion site maintenance to mitigate infection should be per institutional protocol. VolumeView™ sensor or Acumen IQ™ sensor maintenance is required to avoid catheter occlusion and/or avoid inaccurate pressure measurement or cardiac output calculations. Proper VolumeView™ sensor or Acumen IQ™ sensor maintenance includes assuring: proper leveling of sensor, pressure bag inflated to 300 mmHg, adequate flush volume, assessment of waveform quality, and periodic assessment of frequency response.

10.3 VolumeView™ Venous Injectate Kit Setup

Note: Attaching the VolumeView™ venous injectate kit to the central venous catheter is required to obtain VolumeView™ system parameters.

Step	Procedure
1	Insert the central venous catheter according to manufacturer's instructions and institutional policy.

Step	Procedure
2	Prime, level, and zero the TruWave™ transducer: - Flush normal saline through the TruWave™ transducer tubing to remove all air. (Heparin anticoagulation therapy to be administered per institutional protocol.) - Connect the TruWave™ transducer cable to the Edwards EV1000™ monitor or connect HemoSphere Alta™ monitor Pressure Cable to the HemoSphere Alta™ Advanced Monitor as well as the bedside physiologic monitor. - Adjust the level of the sensor port vent (the fluid-air interface) to the level of the right atrium. This is at the phlebostatic axis, determined by the intersection of the midaxillary line and the fourth intercostal space. - Remove the non-vented cap and open the vent port to the atmosphere. - Zero both the Edwards EV1000™ monitor or HemoSphere Alta™ Advanced Monitor as well as the bedside physiologic monitor to atmospheric air. - Close the vent port to the atmosphere and replace the non-vented cap. - The system is ready to begin monitoring pressure.
3	(Optional) If central venous pressure monitoring is desired on patient's bedside physiologic monitor, attach pressure cable between bedside patient monitor and the Edwards EV1000™ monitor or HemoSphere Alta™ Advanced Monitor. CAUTION: Do not connect the set to the catheter during the priming procedure; air can be infused into the patient. CAUTION: To guarantee a constant flush rate of 3 mL/h in the pressure monitoring system, maintain a pressure of 300 mmHg in the pressure cuff.

10.4 VolumeView™ Thermistor Manifold

Note: It is recommended to monitor the balance between oxygen delivery and consumption by monitoring continuous ScvO₂ with an Edwards central venous oximetry catheter.

Note: All components illustrated in Figure 1 on page 59 are required for the system to function properly.

10.5 VolumeView™ Thermistor Manifold Setup

Step	Procedure
1	Per institution guidelines, clean surface of swabbable luer site.
2	Using a luer lock syringe flush normal saline through connections (3 and 4 – Figure 1 on page 59) assuring all air is removed before connecting thermistor manifold to CVC catheter.
3	Connect the swabbable luer site of the VolumeView™ thermistor manifold (3 – Figure 1 on page 59) to the TruWave™ transducer DPT tubing and flush both the TruWave™ transducer and the thermistor manifold with normal saline to prime thermistor manifold removing all air.
4	Connect to CVC lumen (2 - Figure 1 on page 59). CAUTION: It is recommended not to add stopcock(s) in between the CVC lumen and the manifold to avoid damaging/breaking the manifold.
5	Securely attach injectate temperature probe into thermistor manifold (1 – Figure 1 on page 59). Attach the other end of the cable to the Edwards EV1000™ monitor or HemoSphere Alta™ Advanced Monitor via Patient CCO cable HemoSphere Alta™ Advanced Monitoring Platform.
6	Ensure all connections are secure.

10.6 Thermodilution Measurement

Step	Procedure
1	Confirm the cable connections to the Edwards compatible monitor: - Femoral: Blood temperature cable - Venous: Injectate temperature cable - Venous: TruWave™ CV pressure transducer cable or connect HemoSphere Pressure Cable or HemoSphere Alta™ monitor Pressure Cable to the HemoSphere Alta™ Advanced Monitor Note: Use aseptic technique.
2	Select Thermodilution on the EV1000™ monitor or TPTD on the HemoSphere Alta™ Advanced Monitor to start the TPTD procedure (see the manufacturer's manual).
3	Select appropriate syringe volume on the Edwards EV1000™ monitor or HemoSphere Alta™ Advanced Monitor.
4	Select cooled normal saline filled luer lock syringe. CAUTION: Warm injectate may result in inaccurate measurements. Assure timely usage of cooled normal saline filled syringes.
5	Securely attach the luer lock syringe to the swabbable luer site (4 – Figure 1 on page 59).
6	Inject (when prompted) with a smooth, rapid and steady rate while holding and securing the manifold at the base of the swabbable luer (intersection of 3 and 4 in Figure 1 on page 59) to avoid damaging/ breaking the manifold. CAUTION: If significant resistance is encountered during injection, discontinue injection and assess catheter patency and integrity. CAUTION: When the injection has been completed remove the syringe from the swabbable luer site to avoid possible bleed-back into syringe.

11.0 Maintenance and Cleaning

Follow institutional guidelines for care and maintenance of femoral arterial catheters, central venous catheters, and disposable pressure transducers.

Ensure that all lines and connections remain tightly fitted.

12.0 How Supplied

The VolumeView™ set sterilization process utilizes ethylene oxide. Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize. Visually inspect for breaches of the packaging integrity prior to use.

13.0 Storage

Store in a cool, dry place.

14.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

15.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

16.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): . . . 800.822.9837
Outside the U.S. and Canada (24 hours): . . . 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

17.0 Specifications

Model Number	TPTD4F1 TPTD4F5	TPTD5F1 TPTD5F5
Catheter Placement Position, i.e. Femoral, etc.	Femoral	
Outer Diameter (French size and mm)	4 F (1.33 mm)	5 F (1.67 mm)
Usable Length	16 cm	20 cm
Compatible Guidewire Dimensions	0.025" (0.64 mm) 60 cm NiTi	0.032" (0.81 mm) 60 cm NiTi
Material Composition	PU, white, radiopaque	
(Temperature) Measuring Range	15 °C – 45 °C (with EV1000™ clinical platform or HemoSphere Alta™ Advanced Monitor)	
Vessel Dilator Size	4 F (1.4 mm) 10 cm length	5 F (1.7 mm) 15 cm length
Distance Between Length Markers	approximately every 1 cm from 10 cm mark to 16 cm mark	approximately every 1 cm from 10 cm mark to 20 cm mark

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


Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring							
	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op	Forsigtig
	Single use	Ne pas réutiliser	Zum einmaligen Gebrauch	Un solo uso	Monouso	Voor eenmalig gebruik	Til engangsbrug
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid	Mængde
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer	Partinummer
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum	Dato for sidste anvendelse
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/ Den Europæiske Union
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant	Producent
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter	Udvendig diameter
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte	Anvendelig længde
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad	Anbefalet guidewire størrelse
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum	Fremstillingsdato
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterrilizzare	Niet opnieuw steriliseren	Må ikke resteriliseres
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen	Ikke-pyrogen
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata	MRI-veilig onder bepaalde voorwaarden	MR-betinget
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr

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
Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola						
	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Vienkartinio naudojimo	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
	Kiekis	Daudzums	Miktar	Количество	Količina	Količina
	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Naudoti iki nurodytos datos	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе / Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašćeni predstavnik u Europskoj zajednici / Europskoj uniji
	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Išorinis skersmuo	Ārējais diametrs	Dış çap	Внешний диаметр	Spoljašnji prečnik	Vanjski promjer
	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebljiva dužina	Upotrebljiva dužina
	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs	Kılavuz Tel İçin Önerilen Boyut	Рекомендуемый размер проводника	Preporučena veličina žice vodice	Preporučena veličina žice vodilice
	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Žr. naudojimo instrukcijos arba elektronines naudojimo instrukcijas	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению или электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Nepirogeniškas	Nepirogēna	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Sąlyginis MR	Izmantojams MR vidē, ievērojot noteiktus nosacījumus	MR Koşullu	Условно безопасно при проведении МРТ	Uslovno bezbedno za MR	Uvjetno siguran kod pregleda MR-om
	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена, и см. инструкции по применению	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Medicinos priemonė	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski proizvod

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Edwards Lifesciences LLC

One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.

Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

07/25
10063276001 A / D0C-0259128 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

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Edwards Lifesciences GmbH