Edwards Lifesciences

Directory

English (en)	1	Magyar (hu)	45	Русский (ru)	89
Français (fr)	5	Polski (pl)	49	Srpski (sr)	93
Deutsch (de)	9	Slovensky (sk)	53	Hrvatski (hr)	97
Español (es)	13	Norsk (no)	57		
Italiano (it)	17	Suomi (fi)	61		
Nederlands (nl)	21	Български (bg)	65		
Dansk (da)	25	Română (ro)	69		
Svenska (sv)	29	Eesti (et)	73		
Ελληνικά (el)	33	Lietuvių (lt)	77		
Português (pt)	37	Latviešu (lv)	81		
Česky (cs)	41	Türkçe (tr)	85		
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Φιγυρι ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike					101
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola					115

English

Model 93600 CO-Set+ Closed Injectate Delivery System for Cold Injectate

Model 93610 CO-Set+ Closed Injectate Delivery System for Room Temperature Injectate

Model 93520 CO-Set+ Cooling Container

Model 93505 Flow-Through Housing

Model 93650 CO-Set+ Syringe


Model 93522 CO-Set+ In-Line Injectate Temperature Probe

Model 93524 CO-Set+ In-Line Injectate Temperature Probe

Model 93528 CO-Set+ In-Line Injectate Temperature Probe

Model 9850A Injectate Temperature Probe

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Model 93600 CO-Set+ Closed Injectate Delivery System, Cold Injectate

- | | |
|---|--|
| 1. Swan-Ganz Thermodilution Catheter | 10. Snap Clamp |
| 2. Catheter Connecting Cable | 11. Cooling Container |
| 3. Patient Monitor | 12. Tubing Connector |
| 4. Sterile Injectate Solution (User Supplied) | 13. Three Way Stopcock |
| 5. Non-Vented IV Spike | 14. Injectate Temperature Probe |
| 6. Tubing Coil | 15. Flow-Through Check Valve |
| 7. Lid | 16. 10 ml syringe with jacket and contamination shield |
| 8. Thermistor Connector | 17. Cooling Container Bracket |
| 9. Injectate Probe Connector | |

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For single use only

1.0 Description

CO-Set+ is a closed injectate delivery system. It is an accessory to be used with Swan-Ganz catheters that intermittently measure cardiac output (CO). It allows for delivery of injectate for the measurement of Cardiac Output by thermodilution (TD). The unit is designed for use with a patient monitor and a Swan-Ganz thermodilution catheter. These devices are intended to be used by licensed physicians, surgeons, including cardiac, and vascular surgeons.

CO-Set+ closed injectate delivery system for cold injection (Figure 1 on page 101) and for room temperature injectate (Figure 2 on page 103) is connected between a suitable supply of sterile injectate (D5W) and a Swan-Ganz thermodilution catheter.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The measurement of cardiac outputs may improve data driven clinical decision making for medically necessary intervention and/or clinical re-evaluation.

2.0 Intended Use/Purpose

The CO-Set+ closed injectate delivery system is intended to be a pathway to deliver thermodilution injectate (used to make cardiac output determinations) into a thermodilution catheter.

3.0 Indications for Use

The CO-Set+ closed injectate delivery systems used for cold injectate (Model 93600) or room temperature injectate (Model 93610), provide simple and convenient delivery of cold or room temperature injectate for the measurement of cardiac output by thermodilution. The devices are intended for use in adult critically ill patients.

4.0 Contraindications

There are no known contraindications.

5.0 CO-Set+ closed injectate delivery system for cold injectate (Model 93600)

After the system is primed and the cooling coil is immersed in an ice bath, the injectate will reach operating temperature within 5 minutes.

The covered ice bath will remain usable for up to 6 hours; for best results, the ice bath should be checked periodically to ensure that ice is still present.

6.0 Flow-Through Housing (Model 93505)

When used in conjunction with the Injectate Temperature Probe, the Flow-Through Housing (Model 93505) allows direct temperature measurement of the thermal indicator at the proximal lumen infusion site. Measurement of the injectate temperature at the infusion site helps reduce errors frequently encountered in standard injection techniques (for example, injectate warming caused by handling the syringe). The temperature of the injectate is fed directly by the injectate temperature probe into a patient monitor for curve analysis and determination of cardiac output. An appropriate computation constant must be entered into the computer to determine the correct cardiac output. A 3-way stopcock should be placed proximal to the flow-through housing (Figure 11 on page 112) for standard catheter maintenance and fluid withdrawal.

7.0 CO-Set+ Syringe (Model 93650)

The CO-Set+ syringe is a 10 ml syringe for use with the Model 93600 (Cold) or Model 93610 (Room Temperature) closed injectate delivery systems. This individually packaged and sterilized syringe is intended to replace the originally supplied syringe if the need arises.

Additionally, the CO-Set+ syringe may be used with the Model 93505 Flow-Through Housing in place of the user supplied 10 ml syringe described in its Instruction for Use.

8.0 CO-Set+ In-Line Injectate Temperature Probe (Models 93522, 93524, 93528)

The models 93522, 93524, and 93528 CO-Set+ in-line injectate temperature Probes are used in conjunction with injectate delivery systems: model 93600 (Cold) CO-Set+ or model 93610 (Room Temperature) CO-Set+ or model 93505 Flow-Through Housing. The CO-Set+ in-line injectate temperature probe is a three element network consisting of two resistors and one thermistor configured to produce a linear resistance versus temperature curve. The nominal measurement accuracy over the range of 0° to 25 °C is ± 0.3 °C. The CO-Set+ in-line injectate temperature probe is stationed at the point where the syringe injects the thermal indicator into the catheter's proximal lumen hub. The probe is able to precisely measure the injectate temperature as it is being introduced into the catheter.

9.0 Injectate Temperature Probe (Model 9850A)

The Model 9850A Injectate Temperature Probe is a dual thermistor composite that produces a linearized response to temperature change. The nominal measurement accuracy over the range of 0° to 25 °C is ± 0.3 °C. The probe's thermistors are electrically insulated in a tubular stainless steel probe specifically designed for liquid immersion. The probe should not be immersed in liquid past the point where the stainless steel tube attaches to the spring-cord assembly.

The injectate temperature probe attaches to the properly marked branch of the Edwards connecting cable for the Edwards Vigilance monitors. The model 9850A Injectate Temperature Probe is an in-line bath probe used for bolus cardiac output. The probe is 3 1/2 in (8.9 cm) in length and has a 10 ft (3 m) long cable.

10.0 Instructions for Use

Equipment Needed for Thermodilution

- Swan-Ganz thermodilution catheter
- Flow-Through Injectate Temperature Probe
- Sterile D5W in IV container (bag recommended) for use as injectate (500 ml minimum)
- Patient Monitor
- Connecting cable for patient monitor
- Empty flush container for collecting discarded injectate
- Continuous flush device

Equipment Needed for cold injectate:

- CO-Set+ closed injectate delivery system for cold injectate (Model 93600)
- Cooling Container (Model 93520 or equivalent)
- Crushed ice
- Iced water
- Cooling Container IV Pole Mounting Bracket (Model 93521) (optional)

Equipment Needed for room temperature injectate:

- CO-Set+ closed injectate delivery system for room injectate (Model 93610)

10.1 Connecting Injectate Container to Closed System (Figure 4)

Use aseptic technique.

Step	Procedure
1	Suspend injectate container (IV bag or bottle) at the desired height; do not place it higher than three feet (90 cm) above the catheter insertion site.
2	Remove CO-Set+ closed injectate delivery system for cold injection from its package.
3	Close the snap clamp.
4	Remove the protective cap from the IV spike and insert the nonvented spike into the injectate container. Note: If a non-vented glass bottle is used, an airway needle or vented extension set is required to vent the bottle.
5	For cold injectate (Model 93600) only: Pull apart the TUBING COIL (not cooling coil) to obtain the desired tubing length between the injectate container and cooling container. WARNING: Do not pull apart the Cooling Coil. Injectate may not cool to appropriate temperature if Cooling Coil is uncoiled. If additional tubing length is needed, pull apart the tubing coil.

10.2 Connecting Injectate Container to Closed System (Figure 5)

Use Aseptic Technique

Step	Procedure
1	Suspend injectate container (IV bag or bottle) at the desired height; do not place it higher than 3 feet (90 cm) above the catheter insertion site.
2	Remove a CO-Set+ system from its package.
3	Close the snap clamp.
4	Remove the protective cap from the IV spike and insert the non-vented spike into the injectate container.

Note: If a non-vented glass bottle is used, an airway needle or vented extension set is required to vent the bottle.

10.3 Flushing Air from Fluid Path

Use aseptic technique.

Step	Procedure
1	Connect the 10 ml syringe to the flow-through/check valve (Figure 6 on page 108). CAUTION: Do not connect the syringe to the catheter side of the flow-through/check valve. Damage to the flow-through/check valve can occur if the flow-through/check valve is flushed or pressurized from the wrong side. Disabling back-pressurization of the flow-through/check valve can also occur if the 3-way stopcock between the flow-through/check valve and the catheter is not opened to the catheter before injecting. If the flow-through/check valve becomes damaged, then it is possible that blood could be aspirated back into the CO-Set+ system syringe. Should this occur, use of the CO-Set+ system (syringe and flow-through/check valve) should be discontinued and the system replaced.
2	Ensure that all the components are securely connected.
3	Open the snap clamp to allow solution to flow from the injectate container. To collect the discarded solution, hold the flow-through/check valve over an empty container.
4	Slowly pull the syringe plunger to prime the system; then push the plunger in. Repeat 5 or 6 times, or until the system is completely free of air. The system may also be filled by squeezing the IV bag or plastic bottle. Remove the syringe from the flow-through/check valve when filling the system by this method.
5	Return the syringe plunger to the fully depressed position.
6	Close the snap clamp.

10.4 Preparing Cooling Container (Model 93520) (See Fig. 7)

Step	Procedure
1	Pack the container with crushed ice up to the inside ridge in the internal ribs.
2	Seat the cooling coil on the first ridge. For maximum cooling, position the coil with the outlet at the bottom as shown.

Step	Procedure
3	Place the inlet and outlet tubing in the slots provided at the top of the container. Be sure that the white plastic tubing connector is inside the container as shown.
4	Completely cover the cooling coil with crushed ice.
5	Add iced water until it is visible above the ice level.
6	Place any unused length of tubing coil into the ice container to provide additional cold solution.
7	Lock the lid on the cooling container.
8	Let sit for 5 minutes.
9	Flush system once to ensure cold fluid is used for injection (Figure 8 on page 110).

Note: Use the spigot to drain excess water before refilling with more ice. There is no need to disassemble the system.

Note: An optional bracket (Model 93521) is available to mount the container on an IV pole with a diameter between 1/2 and 1 1/2 inches (1.27 and 3.81 cm).

10.5 For Models 93600 and 93610: Connecting System to Catheter and Patient Monitor

Use aseptic technique.

Step	Procedure
1	Verify the absence of air in the fluid path. Attach the flow-through/check valve to the catheter injectate lumen stopcock (Figure 9 on page 110).
2	Insert the injectate temperature probe into the flow-through/check valve until it snaps into place (See Figure 10 on page 111). Secure the system to minimize movement with respect to the patient.
3	Connect the injectate temperature probe cable to the "Injectate Probe" lead of the monitor cable.

10.6 For Models 93505: Connecting System Housing to Catheter and Patient Monitor

Use aseptic technique.

Step	Procedure
1	Verify that the fluid path is completely free of air. Attach the flow-through housing to the catheter injectate lumen stopcock. (Figure 11 on page 112)
2	Insert the reusable injectate temperature probe into the flow-through housing (Figure 12 on page 113). Secure the system to minimize movement with respect to the patient.
3	Connect the injectate temperature probe cable to the "Injectate Probe" lead of the computer cable.

10.7 Prepare Patient Monitor for Operation

Step	Procedure
1	Select injectate volume and catheter size on monitor to determine computation constant or enter the computation constant corresponding to the injectate volume and catheter model being used (see package insert accompanying the catheter for computation constant).
2	Prepare the monitor for operation according to the directions in the appropriate monitor operations manual. CO-Set+ closed injectate delivery system is now ready for use.

10.8 Cardiac Output Measurement

Step	Procedure
1	Open the snap clamp.
2	Turn the stopcock at the catheter injectate hub to close the IV flush and to open the fluid path between the syringe and catheter.
3	Withdraw desired amount (5 to 10 ml) of injectate into syringe.
4	Press start on patient monitor.
5	Immediately perform injections at a rate of at least 2.5 ml/sec. Note: The syringe plunger and the finger grips on the syringe jacket have been designed to help accommodate a wide range of user hand sizes for efficient and comfortable injection of fluid. To ensure uniform fluid injections, care should be taken to avoid pushing in the syringe plunger at an angle to the syringe barrel during injections.
6	Perform cardiac output determination according to instructions in the appropriate patient monitor operations manual. Note: To improve cardiac output reproducibility, space injections approximately one minute apart and increase injectate volume to 10 ml if not already doing so.
7	At the end of the measurements, return the stopcock to its original (flush) position (assuring that the syringe is empty), and close the snap clamp.

10.9 Injectate Temperature

Lower or higher temperatures will affect cardiac output accuracy when using the computation constants provided.

Ensure that the injectate temperature falls within the following ranges:

Model 93600:

Injectate Volume	Temperature Range
10 ml	6 °C - 12 °C
5 ml	8 °C - 16 °C

Model 93610:

Ensure that the injectate temperature falls within the range of 18 °C to 25 °C for 5 or 10 ml injectate.

Model 93505:

Ensure that the injectate temperature falls within 6 °C to 16 °C for cold injectate and 18 °C and 25 °C for room temperature injectate.

Models 93522, 93524, 93528, and 93650:

Please refer to package insert provided with Models 93600, 93610, or 93505 for detailed instructions for use.

10.10 Removing the Injectate Temperature Probe (See Fig. 13)

Grasp the flow-through/check valve, and carefully withdraw (unsnap) the Injectate Temperature Probe. Do not squeeze or pinch the area around the probe connector as this will make probe removal from the flow-through/check valve difficult.

10.11 Model 93505 only: To Remove Flow-Through Temperature Probe (Figure 12)

Grasp the flow-through housing, gently press both tabs, and withdraw the temperature probe.

11.0 Warnings

With the exception of the temperature probes, these devices are designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

12.0 Precautions

Device risks include infection, peripheral air embolism, adverse reaction to device materials, burns or shock, arrhythmia, inappropriate/unintended treatment, and/or delay in monitoring and/or treatment.

Care should be taken to avoid air entrapment during setup and use. The flow-through/check valve and syringe should be secured to minimize movement with respect to the patient.

Models 93522, 93524, 93528, and 9850A:

Do not immerse any electrical connectors of the cable or accessory equipment in liquids.

The recommended Use By Date imprinted on the flag label on the cable represents the recommended date the cable should be replaced or returned to Edwards Lifesciences for evaluation.

Please contact Edwards Technical Support or your local Edwards representative for further assistance.

13.0 How Supplied

Models 93600, 93610, 93650, and 93505:

Contents sterile and nonpyrogenic if package is undamaged or unopened. Do not use if package is opened or damaged. Do not resterilize. Visually inspect for breaches of packaging integrity prior to use.

The CO-Set+ closed injectate delivery system is recommended for single patient use only. It is recommended that CO-Set+ closed injectate delivery system be changed (replaced) for each patient according to hospital procedures for changing tubing sets containing dextrose solutions.

Models 93522, 93524, 93528, and 9850A are reusable and shipped non-sterile. Clean the probe by wiping it with a cleaning solution such as 70% isopropyl alcohol or 10% Bleach Solution.

14.0 Storage

Store in a cool, dry place.

15.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

16.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: This device is designed, intended, and distributed for single use only.

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

17.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Computation Constants for CO-Set+ devices - see package insert accompanying the catheter for computation constant.

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
Figure 2: Model 93610 CO-Set+ Closed Injectate Delivery System (Room Temperature Injectate)

- Figure 2 : Modèle 93610 – Système clos pour administration de solution d'injection CO-Set+ (solution d'injection à température ambiante)
- Abbildung 2: Modell 93610 CO-Set+ geschlossenes Injektatverabreichungssystem (für Injektat mit Raumtemperatur)
- Figura 2: Modelo 93610 del sistema cerrado para inyección CO-Set+ (para inyectables a temperatura ambiente)
- Figura 2: Sistema chiuso per iniettato CO-Set+ (iniettato a temperatura ambiente) modello 93610
- Afbeelding 2: Model 93610 CO-Set+ gesloten injectaatgifteesysteem (voor injectaat op kamertemperatuur)
- Figur 2: Model 93610 CO-Set+ lukket injektatforsyningsssystem (til rumtempereret injektat)
- Figur 2: Modell 93610 CO-Set+ slutet system för administrering av injektat (rumstempererat injektat)
- Εικόνα 2: Κλειστό σύστημα χορήγησης εγχέομένου διαλύματος CO-Set+ (εγχέομενο διάλυμα σε θερμοκρασία δωματίου), μοντέλο 93610
- Figura 2: Sistema de administração de injetado fechado CO-Set+, modelo 93610 (injetado à temperatura ambiente)
- Obrázek 2: Uzavřený zaváděcí systém injektátu CO-Set+ (injektát s pokojovou teplotou), model 93610
- 2. ábra: 93610-es típusú CO-Set+ zárt injektálórendszer (szoba-hőmérsékletű injektátumhoz)
- Rysunek 2: Model 93610: zamknięty system CO-Set+ do podawania roztworu iniekcijnego (roztwór iniekcyjny o temperaturze pokojowej)
- Obrázok 2: Uzavretý aplikačný systém injektátu CO-Set+ (pre injektát s izbovou teplotou), model 93610
- Figur 2: Modell 93610 CO-Set+ lukket injektatleveringsssystem (romtemperert injektat)
- Kuva 2: Mallin 93610 suljettu injektionesteen CO-Set+ -antojärjestelmä (huoneenlämpöinen injektioneste)
- Фигура 2: CO-Set+ затворена система за доставяне на инжектат, модел 93610 (за инжектат със стайна температура)
- Figura 2: Sistemul de administrare a soluției injectate închis CO-Set+, modelul 93610, (soluție injectată la temperatura camerei)
- Joonis 2. Mudel 93610. Süstelahuse suletud manustussüsteem CO-Set+ (toatemperatuuril süstelahus)
- 2 pav. Modelio 93610 „CO-Set+“ uždara švirkščiamojo skysčio leidimo sistema (kambario temperatūros švirkščiamajam skysčiui)
- 2. attēls. Modelis 93610: slēgta injektāta padeves sistēma CO-Set+ (istabas temperatūras injektāts)
- Şekil 2: 93610 Modeli CO-Set+ Kapalı Enjektat İletim Sistemi (Oda Sıcaklığında Enjektat)
- Рис. 2. Система закрытой подачи раствора CO-Set+ модели 93610 (для раствора комнатной температуры)
- Slika 2: Model 93610 CO-Set+ zatvorenog sistema za dopremanje injektata (injektat na sobnoj temperaturi)
- Slika 2: Zatvoreni sustav za uvođenje injektata CO-Set+ model 93610 (injektat na sobnoj temperaturi)

M5033-35



Figure 3 ■ Figure 3 ■ Abbildung 3 ■ Figura 3 ■ Figura 3 ■ Afbeelding 3 ■ Figur 3 ■ Figur 3 ■ Εικόνα 3 ■ Figura 3 ■ Obrázek 3 ■ 3. ábra ■ Rysunek 3
 ■ Obrázok 3 ■ Figur 3 ■ Kuva 3 ■ Φύλλα 3 ■ Figura 3 ■ Joonis 3 ■ 3 pav. ■ 3. attēls. ■ Şekil 3 ■ Рис. 3 ■ Slika 3 ■ Slika 3

M5033-6

1. Sterile Injectate Solution (User supplied) ■ Solution d'injection stérile (fournie par l'utilisateur) ■ Sterile Injektatlösung (vom Benutzer bereitgestellt) ■ Solución de inyectable estéril (suministrada por el usuario) ■ Soluzione di iniettato sterile (fornita dall'utente) ■ Steriele injectaatoplossing (door de gebruiker geleverd) ■ Steril injektatoplsning (brugerleveret) ■ Steril injektatlösning (tillhandhålls av användaren) ■ Στείο εγχυόμενο διάλυμα (παρεχόμενο από τον χρήστη) ■ Solução de injetado esterilizado (fornecida pelo utilizador) ■ Sterilní roztok injektátu (dodaný uživatelem) ■ Steril injektátumoldat (a felhasználó által biztosított) ■ Jałowy roztwór iniekcynny (dostarczany przez użytkownika) ■ Sterilný roztok injektátu (dodaný používateľom) ■ Steril injektatlösning (fra brukeren) ■ Steriili injektionesteliuos (käyttäjän hankkima) ■ Стерильный раствор на инъектат (осигуряван от потребителя) ■ Solutje injectabilā sterilā (asigurată de utilizator) ■ Steriilne süstelahus (kasutaja toodud) ■ Sterilus švirksčiამის ტირპას (პატიკიამს ნაუდოთო) ■ Sterils injektāta šķīdums (nodrošina lietotājs) ■ Steril Enjektat Çözeltilisi (Kullanıcı tarafından sağlanır) ■ Стерильный вводимый раствор (приобретается пользователем) ■ Sterilni rastvor injektata (nabavlja korisnik) ■ Sterilna otopina injektata (nabavlja korisnik)

2. Non-Vented IV Spike ■ Perforateur pour poche de perfusion sans évent ■ Nicht belüfteter Infusionsdorn ■ Punta i.v. sin ventilación ■ Perforatore per endovenosa senza sfiato ■ Infuuspin zonder ontluuchting ■ Ikke-ventileret IV-spids ■ Ikke-ventilerad droppspets ■ Μη αεριζόμενο ρύγχος διάτρησης ασκοῦ ενδοφλέβιου διαλύματος ■ Perforador IV sem ventilação ■ Infuzni hrot bez odvzdušňovacího otvoru ■ Nem szellőző infúziós beszűrőtű ■ Kolec do przekłuwania pojemników do wlewów dożylnych, bez odpowietrznika ■ IV bodec bez odvzdušňovacích otvorov ■ Ikke-ventilert IV-spike ■ Ilma-aukoton pistokärki ■ Спайк за интравенозни разтвори без обезвъздушител ■ Trocar de perfuzor fără orificiu de aerisire ■ Ōhuavata tilgutipiik ■ Neventilioojamas IV smaigas ■ Noslēgts intravenozās sistēmas smaiks uzgalis ■ Havalandırmasız Serum Katatıcı ■ Невентилируемый шип для емкостей с раствором для внутривенных вливаний ■ Šiljak za infuziju bez otvora ■ IV šiljak bez odzračivanja

3. Tubing Coil ■ Serpentin de tubulure ■ Schlauchwickel ■ Bobina del tubo ■ Matassa di tubo ■ Slangspoel ■ Slangespiral ■ Slangspole ■ Σπειρα σωλήνωσης ■ Rolo do tubo ■ Cívka hadičky ■ Feltekereselt csővezeték ■ Zwój przewodu ■ Hadičková cievka ■ Slangkeveil ■ Letkukela ■ Намотка за тръбички ■ Bobinā de tubulaturā ■ Voolikuspiraal ■ Vamzdelio ritē ■ Caurulītes spirāle ■ Hortum Bobini ■ Трубочатая спираль ■ Kalem infuzione linije ■ Namotaj cijevi

4. Cooling Coil ■ Serpentin de refroidissement ■ Kühlspule ■ Bobina de refrigeración ■ Serpentina di raffreddamento ■ Koelspoel ■ Kølespiral ■ Kylspole ■ Στοιχείο ψύξης ■ Bobina de refrigeração ■ Chladicí cívka ■ Hűtőtekercs ■ Wężownica chłodząca ■ Chładiaca cievka ■ Kjølespole ■ Jäähdytyskela ■ Охлаждаща намотка ■ Bobinā de răcire ■ Jahutusspiraal ■ Šaldymo ritē ■ Dzesešanas spirāle ■ Soğutma Bobini ■ Охлаждающий змеевик ■ Kalem za hlađenje ■ Zavojnica za hlađenje

5. Snap Clamp ■ Clamp ■ Schnappklemme ■ Pinza a presión ■ Clamp stringitubo ■ Afsluitklem ■ Sikringsklemme ■ Snäppklämma ■ Κοιμωτός σφιγκτήρας ■ Grampo de retenção ■ Stiskací svorka ■ Pillanatszorító ■ Zacisk ■ Zacvakávacia svorka ■ Smekklekke ■ Puristin ■ Регулираща кламна ■ Clapetă de oprire ■ Plöcksklamber ■ Spaustukas ■ Spaiļu skava ■ Basmali Klempe ■ Фиксирующий зажим ■ Klema ■ Hvataljka

Figure 4: Model 93600 Injectate, IV Spike, Snap Clamp, and Tubing and Cooling Coils

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

Figure 5: Model 93610 View of Injectate Container, IV Spike, Snap Clamp, and Tubing

- Figure 5 : Modèle 93610 Vue du conteneur de solution d'injection, du perforateur pour poche de perfusion, du clamp et de la tubulure
- Abbildung 5: Ansicht des Injektatbehälters, des Infusionsdorns, der Schnappklemme und des Schlauches von Modell 93610
- Figura 5: Modelo 93610: Vista de recipiente de inyectable, punta i.v., pinza a presión y tubos
- Figura 5: Vista di contenitore dell'iniettato, perforatore per endovenosa, clamp stringitubo e tubo, modello 93610
- Afbeelding 5: Model 93610 overzicht van injectaatcontainer, infuuspin, afsluitklem en slangen
- Figur 5: Model 93610 Visning af injektatbeholder, IV-spids, sikringsklemme samt slange
- Figur 5: Modell 93610 Vy av injektatbehållare, droppspets, snäppklämma och slang
- Εικόνα 5: Όψη περιέκτη εγχέομενου διαλύματος, ρύγχους διάτρησης ασκού ενδοφλέβιου διαλύματος, κομπωτού σφιγκτήρα και σωλήνωσης, μοντέλο 93610
- Figura 5: Modelo 93610: Vista do recipiente de injetado, perfurador IV, grampo de retenção e tubagem
- Obrázek 5: Model 93610: pohled na nádobu s injektátem, infuzní hrot, stiskací svorka a hadičku
- 5. ábra: 93610-es típus rajza: Injektátumtartály, infúziós beszúróút, pillanatszorító és csövezeték
- Rysunek 5: Model 93610: pojemnik z roztworem iniekcynnym, kolec do przekłuwania pojemników do wlewów dożylnych, zacisk i przewód
- Obrázok 5: Model 93610: obrázok nádoby na injektát, IV bodca, zcvakávacej svorky a hadičky
- Figur 5: Modell 93610 – visning av injektatbeholder, IV-spike, smekklekke og slange
- Kuva 5: Mallin 93610 injektionestesäiliö, nesteensiirtoletkuston pistokärki, puristin ja letku
- Фигура 5: Модел 93610 – изглед на контейнер за инжектат, спайк за интравенозни разтвори, регулираща клампа и тръбичка
- Figura 5: Ilustrație a modelului 93610 cu recipientul cu soluție injectabilă, trocarul perfuzorului, clapeta de oprire și tubul
- Joonis 5. Mudel 93610. Süstelahuse konteiner, tilgutipiik, plöskklamber ja voolik
- 5 pav. Modelis 93610, švirkščiamajo skysčio talpyklė, IV smaigas, spaustukas ir vamzdeliai
- 5. attēls. Modeļa 93610 injektāta konteiners, IV smaile, piespiedējspāile un caurulītes
- Şekil 5: 93610 Modeli Enjektat Kabı, Serum Kapatıcı, Basmalı Klemp ve Hortumun Görünümü
- Рис. 5. Модель 93610. Изображение емкости с вводимым раствором, шипа для емкости для внутривенного вливания, фиксирующего зажима и трубки
- Slika 5: Model 93610 Prikaz posude injektata, šiljka za infuziju, kleme i cevi
- Slika 5: Model 93610, prikaz spremnika injektata, IV šiljka, hvataljke i cijevi


1. Flow-Through/Check Valve ■ Valve unidirectionnelle / clapet anti-retour ■ Durchfluss-/Rückschlagventil ■ Válvula de retención/paso de flujo ■ Valvola di non ritorno/a flusso passante ■ Doorstroom-/controleklep ■ Gennemløbs-/kontraventil ■ Genomströmnings-/backventil ■ Βαλβίδα συνεχούς ροής/ ανεπίστροφη βαλβίδα ■ Válvula de fluxo/retenção ■ Průtokový/ regulační ventil ■ Átfolyó-/ellenőrzőszelep ■ Zawór przepływowy/zwrotny ■ Prietokový/ spätný ventil ■ Gjennomstrømnings-/tilbakeslagsventil ■ Läpivirtaus-/takaiskuventtiili ■ Поточна/контролна клапа ■ Supapă de pe traseul cu curgere continuă/supapă de control ■ Lävivoolu-/kontrollklapp ■ Srauto tékmés ir (arba) kontrolinis vožtuvas ■ Caurteces vārsts/ pretvārsts ■ Akış Doğrultulu/Çek Valf ■ Проточный/возвратный клапан ■ Ventil za protok/proveru ■ Protočni/jednosmjerni ventil

2. 10 ml syringe with jacket and contamination shield ■ Seringue de 10 ml avec manchon et gaine anticontamination ■ 10- ml-Spritze mit Hülle und Kontaminationsschutz ■ Jeringa de 10 ml con cubierta y funda protectora de contaminación ■ Siringa da 10 ml con camicia e dispositivo anticontaminazione ■ Spuit van 10 ml met omhulsel en verontreinigingsscherm ■ 10 ml-sprøjte med overtræk og kontamineringsbeskytter ■ 10 ml-spruta med mantel och kontamineringsskydd ■ Σύριγγα 10 ml με περίβλημα και προστατευτική διάταξη κατά της μόλυνσης ■ Seringa de 10 ml com revestimento e proteção anticontaminação ■ Injekční stříkačka o objemu 10 ml s pláštěm a ochranou proti kontaminaci ■ 10 mles fecskendő tûsapkával és kontamináció elleni védőburokkal ■ Strzykawka o pojemności 10 ml z płaszczem i osłoną chroniącą przed zanieczyszczeniami ■ 10 ml striekačka s obalom a krytom proti kontaminácii ■ 10 ml sprøyte med kappe og kontamineringsskjold ■ 10 ml:n ruisku, jossa on vaippa ja kontaminaatio suojaus ■ Спринцовка от 10 ml с обивка и щит срещу замърсяване ■ Seringă de 10 ml cu manșon și protecție la contaminare ■ 10 ml süstal koos kesta ja saastetõkkega ■ 10 ml švirkštas su korpusu ir užteršimo skydu ■ 10 ml šļirce ar arvalku un aizsargu pret piesārņojumu ■ Kılıflı ve kontaminasyon korumalı 10 ml şırınga ■ Шприц объемом 10 мл с чехлом и патроном для защиты от загрязнений ■ Špric od 10 ml sa navlakom i štitnikom od kontaminacije ■ Šprica od 10 ml s omotačem i zaštitom od kontaminacije

Figure 6: Model 93600 Syringe and Flow-Through/Check Valve

- Figure 6 : Modèle 93600 – Seringue et valve unidirectionnelle / clapet anti-retour
- Abbildung 6: Modell 93600 Spritze und Durchfluss-/Rückschlagventil
- Figura 6: Modelo 93600 Jeringa y válvula de retención/paso de flujo
- Figura 6: Siringa e valvola di non ritorno/a flusso passante modello 93600
- Afbeelding 6: Model 93600 spuit en doorstroom-/controleklep
- Figur 6: Model 93600 sprøjte og gennemløbs-/kontraventil
- Figur 6: Modell 93600 Spruta och genomströmnings-/backventil
- Εικόνα 6: Σύριγγα και βαλβίδα συνεχούς ροής/ανεπίστροφη βαλβίδα, μοντέλο 93600
- Figura 6: Seringa, modelo 93600, e válvula de fluxo/retenção
- Obrázek 6: Injekční stříkačka a průtokový/regulační ventil, model 93600
- 6. ábra: 93600-as típusú fecskendő és átfolyó-/ellenőrzőszelep
- Rysunek 6: Model 93600: strzykawka i zawór przepływowy/zwrotny
- Obrázok 6: Striekačka a prietokový/spätný ventil, model 93600
- Figur 6: Modell 93600 – sprøyte og gjennomstrømnings-/tilbakeslagsventil
- Kuva 6: Mallin 93600 ruisku ja läpivirtaus-/takaiskuventtiili
- Фигура 6: Спринцовка и поточна/контролна клапа, модел 93600
- Figura 6: Seringa și supapă de pe traseul cu curgere continuă/supapă de control, modelul 93600
- Joonis 6. Mudel 93600. Süstal ja läbivoolu-/kontrollklapp
- 6 pav. Modelis 93600, švirkštas ir srauto tékmés ir (arba) kontrolinis vožtuvas
- 6. attēls. Modelis 93600: šļirce un caurteces vārsts/pretvārsts
- Şekil 6: 93600 Modeli Şırınga ve Akış Doğrultulu/Çek Valf
- Рис. 6. Шприц и проточный/возвратный клапан модели 93600
- Slika 6: Model 93600 Špric i ventil za protok/proveru
- Slika 6: Model 93600, šprica i protočni/jednosmjerni ventil

1. Lid ■ Couvercle ■ Deckel ■ Tapa ■ Coperchio ■ Deksel ■ Låg ■ Lock ■ Πώμα ■ Tampa ■ Viko ■ Fedél ■ Pokrywa ■ Veko ■ Lokk ■ Kansi ■ Kanak ■ Capac ■ Kaas ■ Dangtis ■ Vāks ■ Kapak ■ Крышка ■ Poklopac ■ Poklopac
2. Cooling Coil ■ Serpentin de refroidissement ■ Kühlspeule ■ Bobina de refrigeración ■ Serpentina di raffreddamento ■ Koelspoel ■ Kølespiral ■ Kylspole ■ Στοιχείο ψύξης ■ Bobina de refrigeração ■ Chladicí cívka ■ Hűtőtekercs ■ Wężownica chłodząca ■ Chładiaca cievka ■ Kjølespole ■ Jäähdytyskela ■ Охлаждаща намотка ■ Bobinã de răcire ■ Jahutusspiraal ■ Šaldymo ritė ■ Dzesēšanas spirāle ■ Soğutma Bobini ■ Охлаждающий змеевик ■ Kalem za hlađenje ■ Zavojnica za hlađenje
3. Mounting Tab ■ Patte de fixation ■ Befestigungslasche ■ Pestaña de montaje ■ Linguetta di fissaggio ■ Bevestigingslipje ■ Monteringsstap ■ Monteringsflik ■ Πτερύγιο στερέωσης ■ Encaixe para montagem ■ Montážní otvor ■ Rögzítőfül ■ Element montażowy ■ Upevňovací ploška ■ Monteringsflik ■ Kiinnityskieleke ■ Монтажен палец ■ Ureche de montare ■ Kinnitussakk ■ Tvirtinimo fiksatorius ■ Piestiprināšanas izcilnis ■ Montaj Tırnağı ■ Крепление ■ Jezičak za montiranje ■ Jezičak za postavljanje
4. Thermal Tubing ■ Tubulure thermique ■ Thermoschlauch ■ Tubo térmico ■ Tubo termico ■ Thermische slang ■ Termoslange ■ Termisk slang ■ Θερμική σωλήνωση ■ Tubagem térmica ■ Tepelná hadička ■ Hőszigetelő csővezeték ■ Przewód termiczny ■ Termoizolacyjna hadička ■ Varmeslange ■ Lämpöletku ■ Термална тръбичка ■ Tubulatură termică ■ Termotoru ■ Šildymo vamzdeliai ■ Termālā caurulīte ■ Termal Hortum ■ Теплообменная трубка ■ Termalna cev ■ Toplinska cijev
5. Model 93520 Cooling Container ■ Conteneur de refroidissement, modèle 93520 ■ Kühlbehälter Modell 93520 ■ Modelo 93520 del recipiente de refrigeración ■ Contenitore di raffreddamento modello 93520 ■ Koelcontainer model 93520 ■ Model 93520 Kølebeholder ■ Kylbehållare modell 93520 ■ Περιέκτης ψύξης, μοντέλο 93520 ■ Recipiente de refrigeração, modelo 93520 ■ Chladicí nádoba, model 93520 ■ 93520 típuszámú hűtőtartály ■ Pojemnik chłodzący, model 93520 ■ Chładiaca nádoba, model 93520 ■ Modell 93520 – kjølebeholder ■ Mallin 93520 jäähdytysäiliö ■ Охлаждащ контейнер, модел 93520 ■ Recipient de răcire model 93520 ■ Mudel 93520. Jahutuskonteiner ■ Modelis 93520, šaldymo talpyklė ■ Dzesēšanas tvertnes modelis 93520 ■ 93520 Modeli Soğutma Kabı ■ Охлаждающий контейнер модели 93520 ■ Posuda za hlađenje, model 93520 ■ Spremnik za hlađenje model 93520
6. Outlet Slot ■ Fente de sortie ■ Auslassschlitz ■ Ranura de salida ■ Alloggiamento uscita ■ Afvoersleuf ■ Afløbsåbning ■ Skåra för utlopp ■ Υποδοχή εξόδου ■ Ranhura para tomada ■ Otvor pro vývod ■ Kimeneti nyílás ■ Otwór wylotowy ■ Výstupný otvor ■ Utløpspor ■ Poistoletkun aukko ■ Жлеб за извод ■ Fantă pentru tubul de ieșire ■ Väljastuspilu ■ Išleidimo anga ■ Izplūdes caurulītes grope ■ Çıkış Yuvası ■ Паз для выходной трубки ■ Otvor za odvodnu cev ■ Izlazni prorez
7. Spigot ■ Robinet ■ Hahn ■ Espita ■ Rubinetto ■ Tapkraantje ■ Taphane ■ Kran ■ Kánoula ■ Torneira ■ Kohoutek ■ Csap ■ Kranik ■ Vypúšťací kohútik ■ Tapp ■ Hana ■ Кранче ■ Robinet cu cep ■ Kraan ■ Kaištis ■ Kräns ■ Musluk Tapası ■ Қран ■ Slavina ■ Slavina

Figure 7: Model 93600 Cooling Container and Coil ■ Figure 7 : Modèle 93600 Conteneur et serpentin de refroidissement ■ Abbildung 7: Modell 93600 Kühlbehälter und Spule ■ Figura 7: Modelo 93600 Recipiente y bobina de refrigeración ■ Figura 7: Contenitore e serpentina di raffreddamento, modello 93600 ■ Afbeelding 7: Model 93600 koelcontainer en spoel ■ Figur 7: Model 93600 Kølebeholder og spiral ■ Figur 7: Modell 93600 kylbehållare och spole ■ Εικόνα 7: Περιέκτης ψύξης και στοιχείο ψύξης, μοντέλο 93600 ■ Figura 7: Recipiente de refrigeração, modelo 93600 e bobina ■ Obrázek 7: Model 93600: chladicí nádoba a cívka ■ 7. ábra: 93600-as típusú hűtőtartály és -tekercs ■ Rysunek 7: Model 93600: pojemnik i wężownica chłodząca ■ Obrázok 7: Chładiaca nádoba a cievka, model 93600 ■ Figur 7: Modell 93600 – kjølebeholder og -spole ■ Kuva 7: Mallin 93600 jäähdytysäiliö ja -kela ■ Фигура 7: Охлаждащ контейнер и намотка, модел 93600 ■ Figura 7: Recipient și bobină de răcire model 93600 ■ Joonis 7. Mudel 93600. Jahutuskonteiner ja spiraal ■ 7 pav. Modelis 93600, šaldomoji talpyklė ir ritė ■ 7. attēls. Modeļa 93600 dzesēšanas tvertne un spirāle ■ Şekil 7: 93600 Modeli Soğutma Kabı ve Bobini ■ Рис. 7. Модель 93600. Охлаждающий контейнер и змеевик ■ Slika 7: Model 93600, spremnik i zavojnica za hlađenje



Figure 8: Model 93520 ■ Figure 8 : Modèlle 93520 ■ Abbildung 8: Modell 93520 ■ Figura 8: Modelo 93520 ■ Figura 8: Modello 93520 ■ Afbeelding 8: Model 93520 ■ Figur 8: Model 93520 ■ Figur 8: Modell 93520 ■ Εικόνα 8: Μοντέλο 93520 ■ Figura 8: Modelo 93520 ■ Obrázek 8: Model 93520 ■ 8. ábra: 93520-as típus ■ Rysunek 8: Model 93520 ■ Obrázok 8: Model 93520 ■ Figur 8: Modell 93520 ■ Kuva 8: Malli 93520 ■ Фигура 8: Модел 93520 ■ Figura 8: Modelul 93520 ■ Joonis 8. Mudel 93520 ■ 8 pav. Modelis 93520 ■ 8. attēls. Modelis 93520 ■ Şekil 8: 93520 Modeli ■ Рис. 8. Модель 93520 ■ Slika 8: Model 93520 ■ Slika 8: Model 93520

1. Three-Way Stopcock ■ Robinet d'arrêt à 3 voies ■ Dreiwegehahn ■ Llave de paso de tres vías ■ Rubinetto di arresto a 3 vie ■ 3-wegkraantje ■ Trevejsstophane ■ Trevägskran ■ Τριόδη στρόφιγγα ■ Torneira de passagem de 3 vias ■ Trojcestný uzavírací kohout ■ Háromutas zárócsap ■ Trójdrożny zawór odcinający ■ Trojcestný ventil ■ Treveis stoppekran ■ Kolmitiehana ■ Τριπτόνο σπирателно кранче ■ Robinet de închidere cu trei căi ■ Kolmesuunaline sulgurkraan ■ Trikrýptis čiáupas ■ Trisvirzienu noslēgkrāns ■ Üç Yönlü Musluk ■ Трехходовой запорный кран ■ Trokraka slavinica v Trosmjerni zaporni ventil

2. Flow-Through/Check Valve ■ Valve unidirectionnelle / clapet anti-retour ■ Durchfluss-/Rückschlagventil ■ Válvula de retención/paso de flujo ■ Valvola di non ritorno/a flusso passante ■ Doorstroom-/controleklep ■ Gennemløbs-/kontraventil ■ Genomströmnings-/backventil ■ Βαλβίδα συνεχούς ροής/ ανεπίστροφη βαλβίδα ■ Válvula de fluxo/retenção ■ Prútokový/regulační ventil ■ Átfolyó-/ellenőrzőszelép ■ Zawór przepływowy/zwrotny ■ Prietokový/ spätný ventil ■ Gjennomstrømnings-/tilbakeslagsventil ■ Läpivirtaus-/takaiskuventtiili ■ Поточна/контролна клапа ■ Supapă de pe traseul cu curgere continuă/supapă de control ■ Läbivoolu-/kontrollklapp ■ Srauto tėkmės ir (arba) kontrolinis vožtuvas ■ Caurteces vārsts/ pretvārsts ■ Akış Doğrultulu/Çek Valf ■ Проточный/возвратный клапан ■ Ventil za protok/proveru ■ Protočni/jednosmjerni ventil


3. Catheter Proximal Injectate Hub ■ Embase d'injection proximale du cathéter ■ Proximaler Katheterinjektatanschluss ■ Conector del inyectable proximal al catéter ■ Raccordo di iniezione prossimale del catetere ■ Proximale injectaathub katheter ■ Kateterets proksimale injektatmuffe ■ Kateterns proximale injektatfättning ■ Εγγύς ομφαλός εγγεόμενου διαλύματος καθετήρα ■ Extremidade do injetado proximal do cateter ■ Proximální hrdlo katéturu pro injektát ■ Katéter proximális injektátumcsatlakozója ■ Proksymalne złącze iniekcijne cewnika ■ Proximálne hrdlo katérového injektátu ■ Proksimalt injektatnav for kateter ■ Katetrin proksimaalinen injektionestekanta ■ Проксимален адаптер на катетър за инжектат ■ Racord proximal de injectare al cateterului ■ Kateetri proksimaalne süstelahuse jaotur ■ Kateterio proksimalinio švirksčiavimo skyčio spindžio mova ■ Katetra proksimalā injektāta pieslēgvieta ■ Kateter Proksimal Enjektat Göbeği ■ Проксимальный соединитель катетера для подачи вводимого раствора ■ Proksimalni priključak katetera za injektat ■ Proksimalno čvorište za injektat katetera

4. 10 ml syringe with jacket and contamination shield ■ Seringue de 10 ml avec manchon et gaine anticontamination ■ 10- ml-Spritze mit Hülle und Kontaminationsschutz ■ Jeringa de 10 ml con cubierta y funda protectora de contaminación ■ Siringa da 10 ml con camicia e dispositivo anticontaminazione ■ Spuit van 10 ml met omhulsel en verontreinigingsscherm ■ 10 ml-sprøjte med overtræk og kontamineringsbeskytter ■ 10 ml-spruta med mantel och kontamineringskydd ■ Σύριγγα 10 ml με περιβλημα και προστατευτική διάταξη κατά της μόλυνσης ■ Seringa de 10 ml com revestimento e proteção anticontaminação ■ Injekční stříkačka o objemu 10 ml s pláštěm a ochranou proti kontaminaci ■ 10 mles fesckendõ tûsapkával és kontamináció elleni védõburokkal ■ Strzykawka o pojemności 10 ml z płaszczem i osłoną chroniącą przed zanieczyszczeniami ■ 10 ml striekačka s obalom a krytom proti kontaminácii ■ 10 ml sprøyte med kappe og kontamineringskjold ■ 10 ml:n ruisku, jossa on vaippa ja kontaminaatio suojustus ■ Спринцовка от 10 ml с обивка и щит срещу замърсяване ■ Seringã de 10 ml cu manşon și protecție la contaminare ■ 10 ml švirksčtas su korpusu ir užteršimo skydu ■ 10 ml šļirce ar arvalku un aizsargu pret piesārņojumu ■ Kilifli ve kontaminasyon korumalı 10 ml şırınga ■ Шприц объемом 10 мл с чехлом и патроном для защиты от загрязнений ■ Špric od 10 ml sa navlakom i štítnikom od kontaminacije ■ Šprica od 10 ml s omotačem i zaštitom od kontaminacije

Figure 9: Model 93600 Flow-Through/Check Valve Connected to Catheter Lumen


- Figure 9 : Modèlle 93600 – Valve unidirectionnelle / clapet anti-retour raccordé(e) à la lumière du cathéter
- Abbildung 9: Modell 93600 Durchfluss-/Rückschlagventil angeschlossen an das Katheterlumen
- Figura 9: Modelo 93600 Válvula de retención/paso de flujo conectada a la luz del catéter
- Figura 9: Valvola di non ritorno/a flusso passante collegata al lume del catetere, modello 93600
- Afbeelding 9: Model 93600 doorstroom-/controleklep aangesloten op katheterlumen
- Figur 9: Model 93600 Gennemløbs-/kontraventil tilsluttet til kateterlumen
- Figur 9: Modell 93600 Genomströmnings-/backventil ansluten till kateterns lumen
- Εικόνα 9: Βαλβίδα συνεχούς ροής/ανεπίστροφη βαλβίδα συνδεδεμένη στον αυλό καθετήρα, μοντέλο 93600
- Figura 9: Válvula de fluxo/retenção, modelo 93600, ligada ao lúmen do cateter
- Obrázek 9: Prútokový/regulační ventil, model 93600, připojený k lúminu katéturu

- 9. ábra: 93600-as típus, a katéter lumenéhez csatlakoztatott átfolyó-/ellenőrzőszelep
- Rysunek 9: Model 93600: zawór przepływowy/zwrotny podłączony do kanału cewnika
- Obrázok 9: Prietokový/spätný ventil pripojený k lúmenu katétra, model 93600
- Figur 9: Modell 93600 – gjennomstrømnings-/tilbakeslagsventil koblet til kateterets lumen
- Kuva 9: Mallin 93600 läpivirtaus-/takaiskuventtiili, joka on yhdistetty katetrin luumeniin
- Фигура 9: Поточна/контролна клапа, свързана с лумен на катетър, модел 93600
- Figura 9: Supapa de pe traseul cu curgere continuă/Supapa de control conectată la lumenul cateterului, modelul 93600
- Joonis 9. Mudel 93600. Kateetri valendikuga ühendatud läbivoolu-/kontrollklapp
- 9 pav. Modelio 93600 srauto tėkmės ir (arba) kontrolinis vožtuvas, prijungtas prie kateterio spindžio
- 9. attēls. Modelis 93600: caurteces vārsts/pretvārsts, kas savienots ar katetra lūmenu
- Şekil 9: 93600 Modeli Kateter Lümenine Bağlı Akış Doğrultulu/Çek Valf
- Рис. 9. Проточный/возвратный клапан модели 93600, подсоединенный к просвету катетера
- Slika 9: Model 93600 Ventil za protok/proveru povezan sa lumenom katetera
- Slika 9: Model 93600, protočni/jednosmjerni ventil spojen na lumen katetera

1. Injectate Temperature Probe ■ Sonde de température d'injection ■ Injektat-Temperatursonde ■ Sonda de temperatura del inyectable ■ Sonda di temperatura dell'iniettato
 ■ Injectaattemperatuursonde ■ Injektattemperaturprobe ■ Injektattemperatursond ■ Αισθητήρας θερμοκρασίας εγχέομενου διαλύματος ■ Sonda de temperatura do injetado ■ Sonda k měření teploty injektátu ■ Az injektátum hőmérsékletének mérésére szolgáló szonda ■ Sonda temperatury roztworu iniekcynego ■ Sonda na meranie teploty injektátu
 ■ Injektattemperaturprobe ■ Injektionsteen lämpötilasondi ■ Термосонда за инжектат ■ Sondă de temperatură a soluției injectate ■ Süstelahuse temperatuurisond ■ Švirkščiamoji skysčio temperūros zondas ■ Injektata temperatūras zonde ■ Enjektat Sıcaklık Probu ■ Зонд для измерения температуры вводимого раствора ■ Temperatura sonda injektata
 ■ Sonda za temperaturu injektata


- Figure 10: Temperature Probe Connection ■ Figure 10 : Raccordement de la sonde de température ■ Abbildung 10: Anschluss der Temperatursonde ■ Figura 10: Conexión de la sonda de temperatura ■ Figura 10: Collegamento della sonda di temperatura ■ Afbeelding 10: Verbinding temperatuursonde ■ Figur 10: Tilslutning af temperaturprobe ■ Figur 10: Anslutning för temperatursond ■ Εικόνα 10: Σύνδεση αισθητήρα θερμοκρασίας ■ Figura 10: Ligação da sonda de temperatura ■ Obrázek 10: Připojení sondy k měření teploty ■ 10. ábra: Hőmérsékletmérő szonda csatlakozója ■ Rysunek 10: Złącze sondy temperatury ■ Obrázok 10: Pripojenie sondy na meranie teploty ■ Figur 10: Temperaturprobekobling ■ Kuva 10: Lämpötilasondin liittäminen ■ Фигура 10: Свързване на термосонда ■ Figura 10: Conectarea sondei de temperatură ■ Joonis 10. Temperatuurisondi ühendamine ■ 10 pav. Temperatūros zondo jungtis ■ 10. attēls. Temperatūras zondes savienojums ■ Şekil 10: Sıcaklık Probu Bağlantısı ■ Рис. 10. Подсоединение зонда для измерения температуры ■ Slika 10: Povezivanje temperature sonde ■ Slika 10: Priključak sonde za temperaturu injektata

1. 10 ml syringe with jacket and contamination shield ■ Seringe de 10 ml avec manchon et gaine anticontamination ■ 10- ml-Spritze mit Hülle und Kontaminationsschutz ■ Jeringa de 10 ml con cubierta y funda protectora de contaminación ■ Siringa da 10 ml con camicia e dispositivo anticontaminazione ■ Spuit van 10 ml met omhulsel en verontreinigingscherm ■ 10 ml-sprøjte med overtræk og kontamineringsbeskytter ■ 10 ml-spruta med mantel och kontamineringskydd ■ Σύριγγα 10 ml με περίβλημα και προστατευτική διάταξη κατά της μόλυνσης ■ Seringa de 10 ml com revestimento e proteção anticontaminação ■ Injekční stříkačka o objemu 10 ml s pláštěm a ochranou proti kontaminaci ■ 10 mles fecskendő túsapkával és kontamináció elleni védőburokkal ■ Strykawka o pojemności 10 ml z płaszczem i osłoną chroniącą przed zanieczyszczeniami ■ 10 ml striekačka s obalom a krytom proti kontaminácii ■ 10 ml sprøyte med kappe og kontaminerings skjold ■ 10 ml:n ruisku, jossa on vaippa ja kontaminaatio suojaus ■ Спринцовка от 10 ml с обвивка и щит срещу замърсяване ■ Seringă de 10 ml cu manșon și protecție la contaminare ■ 10 ml süstal koos kesta ja saastetõkkega ■ 10 ml švirkštas su korpusu ir užtersimo skydu ■ 10 ml širice ar arvalku un aizsargu pret piesārņojumu ■ Kılıflı ve kontaminasyon korumalı 10 ml şiringa ■ Шприц объемом 10 мл с чехлом и патроном для защиты от загрязнений ■ Špic od 10 ml sa navlakom i štitnikom od kontaminacije ■ Šprica od 10 ml s omotačem i zaštitom od kontaminacije
2. Flow-Through Housing ■ Dispositif intermédiaire ■ Durchflussgehäuse ■ Carcasa de paso de flujo ■ Alloggiamento a flusso passante ■ Doorstroombehuizing ■ Gennemløbskabinet ■ Genomströmningshölje ■ Περίβλημα διάταξης συνεχούς ροής ■ Fluxómetro ■ Průtokový kryt ■ Átfolyóburkolat ■ Złącze przepływowe ■ Prietokové puzdro ■ Gjennomstrømningshus ■ Läpivirtauskotelo ■ Резервоар за поток ■ Teacă pentru pasajul soluției injectate ■ Kılıflı ve kontaminasyon korumalı 10 ml şiringa ■ Шприц объемом 10 мл с чехлом и патроном для проточный переходник ■ Kučičte za protok ■ Protočno kučičte
3. Proximal Injectate Hub ■ Embout proximal d'injection ■ Proximaler Injektatanschluss ■ Conector del inyectable proximal ■ Raccordo di iniezione prossimale ■ Proximale injectaat hub ■ Proksimal injektatmuffe ■ Proximal injektatfatning ■ Εγγύς ομφαλός εγγεόμενου διαλύματος ■ Extremidade do injetado proximal ■ Proximální hrdlo pro injektát ■ Proximális injek-tátumcsatlakozó ■ Proksymalne złącze iniekcyjne ■ Proximálne hrdlo na injektát ■ Proksimalt injektatnav ■ Proksimaalinen injektionestekanta ■ Проксимален адаптер за инжектат ■ Racord proximal de injectare ■ Proksimaalne süstelahuse jaotur ■ Proksimalinio švirkščiavimo skysčio spindžio mova ■ Proksimālā injektāta pieslēgvieta ■ Proksimal Enjektat Göbeği ■ Проксимальный соединитель для введения раствора ■ Proksimalno čvorište za injektat ■ Proksimalno čvorište za injektat
4. 3-Way Stopcock ■ Robinet d'arrêt à 3 voies ■ Dreivegehahn ■ Llave de paso de tres vías ■ Rubinetto di arresto a 3 vie ■ 3-wegkraantje ■ Trevejsstoppane ■ Trevägskran ■ Τρίοδη στρόφιγγα ■ Torneira de passagem de 3 vias ■ Trojcestný uzavírací kohout ■ Háromutas zárócsap ■ Trójdrożny zawór odcinający ■ 3-cestný ventil ■ Treveis stoppekrän ■ Kolmitiehana ■ 3-пътно спирателно кранче ■ Robinet de închidere cu 3 căi ■ 3-suunaline sulgurkraan ■ Trikyrtis čiaupas ■ Trisvirzienu noslēgkrāns ■ 3 Yönlü Musluk ■ Трехходовой запорный кран ■ Trokraka slavnica ■ Tromjerni zaporni ventil


Figure 11: Model 93505: Flow-Through Housing Connected to Catheter Injectate Lumen

- Figure 11 : Modèle 93505 : dispositif intermédiaire connecté à la lumière d'injection du cathéter
- Abbildung 11: Modell 93505: Durchflussgehäuse angeschlossen an das Katheter-Injektatlumen
- Figura 11: Modelo 93505: Carcasa de paso de flujo conectada a la luz del catéter del inyectable
- Figura 11: Modello 93505: alloggiamento a flusso passante collegato al lume di iniezione del catetere
- Afbeelding 11: Model 93505: doorstroombehuizing aangesloten op katheterinjectaatlumen
- Figur 11: Model 93505: Gennemløbskabinet tilsluttet kateterinjektatlumen
- Figur 11: Modell 93505: Genomströmningshölje anslutet till kateters injektatlumen
- Εικόνα 11: Μοντέλο 93505: περίβλημα διάταξης συνεχούς ροής συνδεδεμένο στον αυλό εγγεόμενου διαλύματος του καθετήρα
- Figura 11: Modelo 93505: fluxómetro ligado ao lúmen de injetado do cateter
- Obrázek 11: Model 93505: průtokový kryt připojený k lúminu injektátu na katétru
- 11. ábra: 93505-ös típus: katéter injektátumlumenéhez csatlakoztatott átfolyóburkolat
- Rysunek 11: Model 93505: złącze przepływowe podłączone do kanału iniekcijnego w cewniku
- Obrázok 11: Model 93505: prietokové puzdro pripojené k lúmenu pre katérový injektát
- Figur 11: Modell 93505: Gjennomstrømningshus koblet til kateterets injektatlumen
- Kuva 11: Malli 93505: läpivirtauskotelo, joka on yhdistetty katetrin injektionestelumlumeniin
- Фигура 11: Модел 93505: Резервоар за поток, свързан с лумена на катетър за инжектат
- Figura 11: Modelul 93505: teacă pentru pasajul soluției injectate, conectată la lumenul de injectare al cateterului
- Joonis 11. Mudel 93505. Kateetri süstelahuse valendikuga ühendatud läbivoolumkorpus
- 11 pav. Modelis 93505: srauto tékmės korpusas, prijungtas prie kateterio švirkščiamojo skysčio spindžio
- 11. attēls. Modelis 93505: caurteces korpusas, kas pieslēgts katetra injektāta lūmenam
- Şekil 11: 93505 Modeli: Kateter Enjektat Lümenine Bağlı Akış Doğrultulu Muhafaza
- Рис. 11. Модель 93505: проточный переходник, подсоединенный к просвету катетера для подачи вводимого раствора
- Slika 11: Model 93505: kučičte za protok povezano sa lumenom katetera za injektat
- Slika 11: Model 93505, protočno kučičte spojeno na lumen za injektat katetera

1. 10 ml syringe with jacket and contamination shield ■ Seringe de 10 ml avec manchon et gaine anticontamination ■ 10- ml-Spritze mit Hülle und Kontaminationsschutz ■ Jeringa de 10 ml con cubierta y funda protectora de contaminación ■ Siringa da 10 ml con camicia e dispositivo anticontaminazione ■ Spuit van 10 ml met omhulsel en verontreinigingsscherm ■ 10 ml-sprøjte med overtræk og kontamineringsbeskytter ■ 10 ml-spruta med mantel och kontamineringskydd ■ Σύριγγα 10 ml με περίβλημα και προστατευτική διάταξη κατά της μόλυνσης ■ Seringa de 10 ml com revestimento e proteção anticontaminação ■ Injekční stříkačka o objemu 10 ml s pláštěm a ochranou proti kontaminaci ■ 10 mles feckendõ tüsapkával és kontamináció elleni védõburokkal ■ Strzykawka o pojemności 10 ml z płaszczem i osłoną chroniącą przed zanieczyszczeniami ■ 10 ml striekačka s obalom a krytom proti kontaminácii ■ 10 ml sprøyte med kappe og kontamineringskjold ■ 10 ml:n ruisku, jossa on vaippa ja kontaminaatio suojaus ■ Спринцовка от 10 ml с обвивка и щит срещу замърсяване ■ Seringă de 10 ml cu manșon și protecție la contaminare ■ 10 ml süstal koos kesta ja saastetõkkega ■ 10 ml švirkštas su korpusu ir užteršimo skydu ■ 10 ml šjirce ar arvalku un aizsargu pret piesārņojumu ■ Kılıflı ve kontaminasyon korumalı 10 ml şiringa ■ Шприц объемом 10 ml с чехлом и патроном для защиты от загрязнений ■ Špic od 10 ml sa navlakom i štitnikom od kontaminacije ■ Šprica od 10 ml s omotačem i zaštitom od kontaminacije
2. Flow-Through Housing ■ Dispositif intermédiaire ■ Durchflussgehäuse ■ Carcasa de paso de flujo ■ Alloggiamento a flusso passante ■ Doorstroombehuizing ■ Gennemløbskabinet ■ Genomströmningshölje ■ Περίβλημα διάταξης συνεχούς ροής ■ Fluxómetro ■ Průtokový kryt ■ Átfolyóburkolat ■ Złącze przepływowe ■ Prietokové puzdro ■ Gjennomstrømningsshus ■ Läpivirtauskotelo ■ Резервоар за поток ■ Teacă pentru pasajul soluției injectate ■ Lävivoolumkorpus ■ Srauto tēkmēs korpusas ■ Caurteces korpus ■ Akış Doğrultulu Muhafaza ■ Проточный переходник ■ Kućište za protok ■ Protočno kućište
3. Proximal Injectate Hub ■ Embout proximal d'injection ■ Proximaler Injektatanschluss ■ Conector del inyectable proximal ■ Raccordo di iniezione prossimale ■ Proximale injectaathub ■ Proksimal injektatmuffe ■ Proximal injektatfatning ■ Εγγύς ομφαλός εγγεόμενου διαλύματος ■ Extremidade do injetado proximal ■ Proximální hrdlo pro injektát ■ Proximális injecktátumcsatlakozó ■ Proksymalne złącze iniekcyjne ■ Proximálne hrdlo na injektát ■ Proksimalt injektatnav ■ Proksimaalinen injektionestekanta ■ Проксимален адаптер за инжектат ■ Racord proximal de injectare ■ Proksimaalne süstelahuse jaotur ■ Proksimalinio švirkščiavimo skyčio spindžio mova ■ Proksimālā injektāta pieslēgvieta ■ Proksimal Enjektat Göbeği ■ Проксимальный соединитель для введения раствора ■ Proksimalno čvorište za injektat ■ Proksimalno čvorište za injektat
4. 3-Way Stopcock ■ Robinet d'arrêt à 3 voies ■ Dreivegehahn ■ Llave de paso de tres vías ■ Rubinetto di arresto a 3 vie ■ 3-wegkraantje ■ Trevejsstophane ■ Trevägskran ■ Τριόδη στρόφιγγα ■ Torneira de passagem de 3 vias ■ Trojcestný uzavírací kohout ■ Háromutas zárócsap ■ Trójdrożny zawór odcinający ■ 3-cestný ventil ■ Treveis stoppekrän ■ Kolmitiehana ■ 3-пътно спирателно кранче ■ Robinet de închidere cu 3 căi ■ 3-suunaline sulgurkraan ■ Trikryptis čiaupas ■ Trisvirzienu noslēgkrāns ■ 3 Yönlü Musluk ■ Трехходовой запорный кран ■ Trokraka slavnica ■ Tromsjerni zaporni ventil
5. Temperature Probe ■ Sonde de température ■ Temperatursonde ■ Sonda de temperatura ■ Sonda di temperatura ■ Temperatuursonde ■ Temperaturprobe ■ Temperatursond ■ Αισθητήρας θερμοκρασίας ■ Sonda de temperatura ■ Teplotní sonda ■ Hőmérsékletmérő sonda ■ Sonda temperatury ■ Sonda na meranie teploty ■ Temperaturprobe ■ Lämpötilasondi ■ Термосонда ■ Sondă de temperatură ■ Temperatuurisond ■ Temperatūros zondas ■ Temperatūras zonde ■ Sıcaklık Probu ■ Зонд для измерения температуры ■ Temperaturna sonda ■ Sonda za temperaturu

Figure 12: Temperature Probe and Flow-Through Housing ■ Figure 12 : Sonde de température et dispositif intermédiaire ■ Abbildung 12: Temperatursonde und Durchflussgehäuse ■ Figura 12: Sonda de temperatura y carcasa de paso de flujo ■ Figura 12: Sonda di temperatura e alloggiamento a flusso passante ■ Afbeelding 12: Temperatuursonde en doorstroombehuizing ■ Figur 12: Temperaturprobe og gennemløbskabinet ■ Figur 12: Temperatursond och genomströmningshölje ■ Εικόνα 12: Αισθητήρας θερμοκρασίας και περίβλημα διάταξης συνεχούς ροής ■ Figura 12: Sonda de temperatura e fluxómetro ■ Obrázek 12: Teplotní sonda a průtokový kryt ■ 12. ábra: Hőmérsékletmérő sonda és átfolyóburkolat ■ Rysunek 12: Sonda temperatury i złącze przepływowe ■ Obrázok 12: Sonda na meranie teploty a prietokové puzdro ■ Figur 12: Temperaturprobe og gjennomstrømningsshus ■ Kuva 12: Lämpötilasondi ja läpivirtauskotelo ■ Фигура 12: Термосонда и резервоар за поток ■ Figura 12: Sonda de temperatură și teacă pentru pasajul soluției injectate ■ Joonis 12. Temperatuurisond ja lävivoolumkorpus ■ 12 pav. Temperatūros zondas ir srauto tēkmēs korpusas ■ 12. attēls. Temperatūras zonde un caurteces korpus ■ Şekil 12: Sıcaklık Probu ve Akış Doğrultulu Muhafaza ■ Рис. 12. Зонд для измерения температуры и проточный переходник ■ Slika 12: Temperaturna sonda i kućište za protok ■ Slika 12: Sonda za temperaturu i protočno kućište


1. Pull (to unsnap) ■ Tirer (pour déclipser) ■ Herausziehen ■ Tira (desencaje) ■ Tirare (per sganciare) ■ Trekken (om los te klikken) ■ Træk (for at frakoble) ■ Dra (för att snäppa loss) ■ Τραβήξτε (για να ξεκουμπώσετε) ■ Puxe (para soltar) ■ Vytáhněte (odpojte) ■ Húzza meg (hogy kipattanjon) ■ Pociągnąć (aby odłączyć od złącza zatraskowego) ■ Potiahnut' (vycvaknúť) ■ Trekk (for å ta den ut av lås) ■ Vedā (napsauta irti) ■ Издърпване (за освобождаване) ■ Trageti (pentru a desprinde) ■ Tõmmake (eemaldamiseks) ■ Patraukite (kad atsikabintų) ■ Vilklt (lai atbrīvotu) ■ (Bağlantıyı kesmek için) çekin ■ Потяните (чтобы отсоединить) ■ Povucite (da biste otkáčili) ■ Povucite (za otkvačivanje)

Figure 13: Removing the Temperature Probe

- Figure 13 : Retrait de la sonde de température
- Abbildung 13: Entfernen der Temperatursonde
- Figura 13: Extracción de la sonda de temperatura
- Figura 13: Rimozione della sonda di temperatura
- Afbeelding 13: De temperatuursonde verwijderen
- Figur 13: Udtagning af temperaturproben
- Figur 13: Ta bort temperatursonden
- Εικόνα 13: Αφαίρεση του αισθητήρα θερμοκρασίας
- Figura 13: Remover a sonda de temperatura
- Obrázek 13: Odstranění sondy k měření teploty
- 13. ábra: A hőmérsékletmérő szonda eltávolítása
- Rysunek 13: Odłączanie sondy temperatury
- Obrázok 13: Odstránenie sondy na meranie teploty
- Figur 13: Fjerne temperaturproben
- Kuva 13: Lämpötilasondin irrottaminen
- Фигура 13: Отстраняване на термосондата
- Figura 13: Îndepărtarea sondei de temperatură
- Joonis 13. Temperatuurisondi eemaldamine
- 13 pav. Temperatūros zondo ištraukimas
- 13. attēls. Temperatūras zondes noņemšana
- Şekil 13: Sıcaklık Probenun Çıkarılması
- Рис. 13. Отсоединение зонда для измерения температуры
- Slika 13: Uklanjanje temperature sonde
- Slika 13: Uklanjanje sonde za temperaturu

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español	Italiano
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Quantity	Quantité	Menge	Cantidad	Quantità
	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	Single Use	Ne pas réutiliser	Zum einmaligen Gebrauch	Un solo uso	Monouso
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Do not re sterilize	Ne pas restériliser	Nicht re sterilisieren	No volver a esterilizar	Non risterilizzare
	Store in a cool, dry place.	Conserver dans un endroit frais et sec.	Das Produkt an einem kühlen, trockenen Ort lagern.	Guárdese en un lugar fresco y seco.	Conservare in un luogo fresco e asciutto.
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Consultare le istruzioni per l'uso

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring
 ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom
 ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas
 ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español	Italiano
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web
	Importer	Importateur	Importeur	Importador	Importatore
	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union	Representante autorizado en la Comunidad Europea/Unión Europea	Rappresentante autorizzato nella Comunità Europea/Unione Europea
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarriersystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring
 ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat / Legenda symboli ■ Vysvetlivky k symbolom
 ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolu paašikinimas
 ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

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
	Suomi	Български	Română	Eesti	Lietuvių
	Mallinnumero	Номер на модела	Număr model	Mudeli number	Modelio numeris
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks.	Vienkartinio naudojimo
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge reesteriliseerige	Kartotinai nesterilizuokite
	Säilytettävä kuivassa ja viileässä.	Да се съхранява на хладно, сухо място.	Stocați într-un loc rece și uscat.	Säilitada jahedas ja kuivas kohas.	Laiykite vėsioje, sausoje vietoje.
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
 eifu.edwards.com + 1 888 570 4016	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring
 ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom
 ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas
 ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
 eifu.edwards.com + 1 888 570 4016	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiata veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană / Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Viimeinen käyttöpäivä	Срок на годност	Data expirării	Kõlblik kuni	Naudoti iki nurodytos datos
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepirogeenne	Nepirogeniškas
	Epästeriili	Нестерилно	Nesteril	Mittesteriilne	Nesterilus
	Sisältää luonnokumilateksia tai sen jäämiä	Съдържание или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural	Sisaldab naturaalsed kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusis priemonės identifikatorius

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring
 ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom
 ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas
 ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Daudzums	Miktar	Количество	Količina	Količina
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Vienreizējai lietošanai	Tek Kullanımlık	Только для одноразового использования	Za jednokratnu upotrebu	Za jednokratnu upotrebu
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Uzglabāt vēsā, sausā vietā.	Serin, kuru yerde saklayın.	Хранить в прохладном, сухом месте.	Čuvati na hladnom i suvom mestu.	Pohranite na hladnom, suhom mjestu.
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na web-stranici
	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
 eifu.edwards.com + 1 888 570 4016	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu na web-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Pilvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašćeni predstavnik u Europskoj zajednici / Europskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Nesterils	Steril Değildir	Нестерильно	Nesterilno	Nesterilno
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilenoksidom	Sterilizirano etilenoksidom
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykietcie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei vältämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminių etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** Na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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Made in Dominican Republic
CCT Critical Care Technologies, S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

02/26
10056131002 A / DOC-0561902 A
© Copyright 2026, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH