Edwards


Directory

English (en).....	2
Français (fr).....	5
Deutsch (de).....	8
Español (es).....	11
Italiano (it).....	14
Nederlands (nl).....	17
Dansk (da).....	20
Svenska (sv).....	23
Ελληνικά (el).....	26
Português (pt).....	29
Česky (cs).....	32
Magyar (hu).....	35
Polski (pl).....	38
Slovensky (sk).....	41
Norsk (no).....	44
Suomi (fi).....	47
Български (bg).....	50
Română (ro).....	53
Eesti (et).....	56
Lietuvių (lt).....	59
Latviešu (lv).....	62
Türkçe (tr).....	65
Русский (ru).....	68
Srpski (sr).....	71
Hrvatski (hr).....	74
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Pav. ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike.....	77
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförlaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola	79

English

Acumen IQ Plus finger cuff

AIQCA2

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Acumen IQ Plus finger cuff, when used with an appropriate monitoring system, provides continuous, noninvasive hemodynamic monitoring and the applicable derived parameters. The Acumen IQ Plus finger cuff utilizes the volume-clamp method to measure blood pressure with an inflatable bladder wrapped around the middle phalanx of the index, middle or fourth/ring finger. This device is supplied non-sterile. Potential risks include inappropriate/unintended treatment, adverse reaction to device materials, patient or clinician burns or electrical shock, peripheral ischemia, and/or minor tissue damage.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

Device is used by clinicians who have been trained in safe use of noninvasive hemodynamic technologies in accordance with their institutional guidelines.

2.0 Intended Use/Purpose

The intended purpose is to noninvasively measure blood pressure and use the information to drive

hemodynamic parameters when connected to the Smart Pressure Controller (PC1Q) and a compatible hemodynamic platform.

3.0 Indications

The Acumen IQ Plus finger cuff is indicated for patients over 18 years of age to continuously noninvasively measure blood pressure and associated hemodynamic parameters when used with a compatible hemodynamic platform.

4.0 Contraindications

In some patients with extreme contraction of the smooth muscle in the arteries and arterioles in the lower arm and hand, such as may be present in patients with Raynaud's disease, blood pressure measurement can become impossible.

5.0 Single Patient Use

The Acumen IQ Plus finger cuff is designed for single patient use. Upon starting measurement, the finger cuff can be used and re-applied for up to 72 hours on one patient. After 8 hours of continuous monitoring on a single finger, the finger cuff should be re-applied to another finger.

Do not attempt to clean and reuse the finger cuff on more than a single patient, such action could lead to illness or an adverse event as the device may not function as originally intended.

6.0 Warnings

- **Improper placement or alignment of the finger cuff can lead to inaccurate monitoring.**
- **Do not apply the finger cuff on injured skin as this may cause further injury.**
- **To reduce the risk of skin irritation and tissue damage, do not monitor longer than 8 hours continuously on a**

Edwards, Edwards Lifesciences, the stylized E logo, Acumen, Acumen IQ, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

single finger. To continue to monitor beyond 8 hours, move the cuff in use to another finger.

- **Measurement on one finger in contradiction with the instructions for use may affect patient comfort and/or lead to minor injuries.**
- **Do not apply finger cuff on a hand or finger when external constriction (that may prevent circulation to the hand or finger) is present.**
- **Do not apply finger cuff on a hand or finger when a second blood pressure measurement device is actively monitoring on the same arm (or hand or finger).**
- **Do not use the finger cuff with magnetic resonance imaging.**

7.0 Cautions

- **Do not use a damaged finger cuff. This may result in inaccurate measurements or may damage the monitoring system.**
- **Never bend a finger cuff to a flat shape as it will damage the finger cuff and affect measurement accuracy.**
- **Excessive ambient lighting may interfere with the finger cuff measurement. Avoid using the finger cuff under close, direct lighting.**
- **Always disconnect the finger cuff when it is not applied to a finger, to prevent damage by accidental overinflation.**
- **The effectiveness of the finger cuffs have not been established in pre-eclamptic patients.**

8.0 Instructions for Use

Refer to Figure 1 on page 77 through Figure 3 on page 78 for figures corresponding to the steps below.

8.1 Apply the Finger Cuff

Step	Procedure
1	Gently open the alignment tab and the finger cuff and place the finger cuff on the middle phalanx of the index, middle or fourth/ring finger. Ensure the finger cuff is centered between the second and third knuckles, and the center line on the inside of the finger cuff (Figure 1). Note: Do not apply the finger cuff on the thumb, small finger, or previously fractured fingers.
2	Align the distal end of the finger with the center line of the alignment tab on the finger cuff (Figure 1).
3	Allow the finger cuff to close around the finger (Figure 2). Ensure that the finger cuff remains aligned and does not rotate.

Step	Procedure
4	Use hook and loop to secure around the finger (Figure 2-A). Secure the distal end of the alignment tab and secure to the top of the finger cuff (Figure 2-B).
5	Lead the finger cuff cable between two fingers to the back side of the hand (Figure 3).

8.2 Connect the Finger Cuff to the System

Step	Procedure
1	Plug the finger cuff connector into the pressure controller (Figure 3). Refer to the pressure controller instructions for use for more details.
2	Initiate monitoring with the connected noninvasive monitoring system. To enable advanced parameters, refer to the monitoring system operator's manual for more details.

9.0 Recommendations

The noninvasive monitoring system will pause blood pressure measurements periodically per the monitoring system configuration settings. At this time, the fingertip should be checked for perfusion and potential tissue damage. The fingertip should also be checked periodically according to hospital protocol.

10.0 MRI Safety

This device is MR-unsafe and poses hazards in the MRI environment. This device contains metallic components, which can experience RF-induced heating in the MRI environment.

11.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

13.0 Storage

Store in a cool, dry place.

14.0 Specification Table

Model	Size Range
AIQCA2	43-71 mm

15.0 Warranty

The Acumen IQ Plus finger cuff is for single patient use only. The Acumen IQ Plus finger cuff is warrantied at time of delivery to the end user only. The Acumen IQ Plus finger cuff is a non-serviceable part.

16.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): . . 800.822.9837

Outside the U.S. and Canada

(24 hours): . . 949.250.2222

In the UK: . . . 0870 606 2040 - Option 4

In Ireland: . . . 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.


Refer to the latest version of the noninvasive monitoring system operator's manual for more information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

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


Figure 3: Connect the finger cuff to the system.

- **Figure 3 : Connecter le manchon de doigt au système.**
- **Abbildung 3: Fingermanschette am System anschließen.**
- **Figura 3: Conecte el manguito para el dedo al sistema.**
- **Figura 3: Collegamento della fascetta per dito al sistema.**
- **Afbeelding 3: Sluit de vingermanschet aan op het systeem.**
 - **Figur 3: Slut fingermanschetten til systemet.**
 - **Figur 3: Anslut fingermanschetten till systemet.**
- **Εικόνα 3: Συνδέστε την περιχειρίδα δακτύλου στο σύστημα.**
- **Figura 3: Ligue o sensor digital ao sistema.**
- **Obrázek 3: Připojte prstovou manžetu k systému.**
- **3. ábra: Az ujjmandzsetta csatlakoztatása a rendszerhez**
- **Rysunek 3: Podłączanie mankietu na palec do systemu.**
- **Obrázok 3: Pripojenie prstovej manžety k systému.**
 - **Figur 3: Koble fingermansjetten til systemet.**
 - **Kuva 3: Liitä sormimansetti järjestelmään.**
- **Фигура 3: Свържете маншета за пръст към системата.**
- **Figura 3: Conectarea manșetei pentru deget la sistem.**
 - **Joonis 3. Ühendage sõrmemansett süsteemiga.**
- **3 pav. Piršto manžetės prijungimas prie sistemos.**
- **3. attēls. Savienojiet pirksta manšeti ar sistēmu.**
- **Şekil 3: Parmak manşonunu sisteme bağlayın.**
- **Рис. 3. Подсоединение пальцевой манжеты к системе**
 - **Slika 3: Spojite manžetnu za prst sa sistemom.**
 - **Slika 3: Povežite naprstak sa sustavom.**

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español	Italiano
	Fragile, handle with care	Fragile, manipuler avec précaution	Vorsicht! Zerbrechlich!	Frágil; manipular con cuidado	Fragile, maneggiare con cura
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile
	Quantity	Quantité	Menge	Cantidad	Quantità
	Caution	Avertissement	Vorsicht	Aviso	Attenzione

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español	Italiano
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Consultare le istruzioni per l'uso
 eifu.edwards.com + 1 888 570 4016	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM	Non compatibile con RM
	Separate collection for electrical and electronic equipment in accordance with EC Directive 2012/19/EU	Tri sélectif des équipements électriques et électroniques conformément à la Directive CE 2012/19/UE	Getrennte Sammlung von Elektro- und Elektronikgeräten gemäß EU-Richtlinie 2012/19/EU	Recogida por separado de equipos eléctricos y electrónicos según la directiva de la CE 2012/19/UE	Raccolta separata per attrezzature elettriche ed elettroniche in conformità alla Direttiva CE 2012/19/UE
	Importer	Importateur	Importeur	Importador	Importatore
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

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

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Niet-steriel	Ikke-steril	Icke-steril	Μη αποστειρωμένο	Não esterilizado
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen	Ακολουθείτε τις οδηγίες χρήσης	Seguir as instruções de utilização
	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen	Ακολουθείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Seguir as instruções de utilização no site
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	MRI-onveilig	MR-usikker	MR-farlig	Μη ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)	Utilização não segura em ambiente de RM
	Gescheiden inzameling voor elektrische en elektronische apparatuur in overeenstemming met EG-richtlijn 2012/19/EU	Separat indsamling af elektrisk og elektronisk udstyr iht. EF-direktivet 2012/19/EU	Separat insamling för elektrisk och elektronisk utrustning i enlighet med EU-direktiv 2012/19/EU	Ξεχωριστή συλλογή για ηλεκτρικό και ηλεκτρονικό εξοπλισμό σύμφωνα με την Οδηγία του Ευρωπαϊκού Κοινοβουλίου 2012/19/ΕΕ	Recolha em separado para o equipamento elétrico e eletrônico, de acordo com a Diretiva CE 2012/19/ UE.
	Importeur	Importør	Importör	Εισαγωγέας	Importador

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola



	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (σήμανση CE)	Conformité Européenne (Marcação CE)
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/ Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/União Europeia
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico

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
	Suomi	Български	Română	Eesti	Lietuvių
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis
	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiате veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Ei sovi magneettikuvaukseen	Небезопасно при MR	Incompatibil RM	Ohtlik magnetresonantsto mograafias	MR nesaugus
	Erottele sähköiset ja elektroniset laitteet EY-direktiivin 2012/19/EU mukaisesti	Разделно събиране на електрическо и електронно оборудване според Директивата на ЕС 2012/19/ЕС	Colectarea separată a echipamentelor electrice și electronice, în conformitate cu directiva CE 2012/19/UE	Elektri- ja elektroonikaseadmet e eraldi kogumine kooskõlas EÜ direktiiviga 2012/19/EL	Atskiras elektros ir elektroninės įrangos atliekų surinkimas pagal EB direktyvą 2012/19/ES
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Älä käyttää, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Mallinnumero	Номер на модела	Număr model	Mudeli number	Modelio numeris
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Nesterils	Steril Değildir	Нестерильно	Nesterilno	Nesterilno
	Daudzums	Miktar	Количество	Količina	Količina
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu
	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Nedrīkst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR	Nije sigurno kod pregleda MR-om

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES Direktīvu 2012/19/ES	2012/19/EU sayılı EC direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır	Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU	Odvojeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU	Odvojeno prikupljanje električne i elektroničke opreme u skladu s direktivom EZ-a 2012/19/EU
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Pilnvarotais pārstāvis Eiropas Kopienā/ Eiropas Savienībā	Avrupa Birliği'nde/ Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto produktu nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykiecie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забележка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almayış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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Made in Dominican Republic
CCT Critical Care Technologies, S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

09/25
10061115001 A/DOC-0568424 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH