Edwards

KONECT RESILIA Aortic Valved Conduit, Model 11060A Accessory Care and Sterilization Instructions for Tray Model TRAY1190

Directory	
English (en).....	1
Français (fr).....	4
Deutsch (de).....	6
Español (es).....	8
Italiano (it).....	10
Nederlands (nl).....	12
Dansk (da).....	14
Svenska (sv).....	16
Ελληνικά (el).....	18
Português (pt).....	20
Česky (cs).....	22
Magyar (hu).....	24
Polski (pl).....	26
Slovensky (sk).....	28
Norsk (no).....	30
Suomi (fi).....	32
Български (bg).....	34
Română (ro).....	36
Eesti (et).....	38
Lietuvių (lt).....	40
Latviešu (lv).....	42
Türkçe (tr).....	44
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları.....	46

English

KONECT RESILIA Aortic Valved Conduit, Model 11060A Accessory Care and Sterilization Instructions for Tray Model TRAY1190

Rx only

1.0 Tray Model TRAY1190

The Instructions For Use supplied with KONECT RESILIA aortic valved conduit (AVC), Model 11060A provide information on the use of this tray.

2.0 Device Description

TRAY1190 is an accessory to KONECT RESILIA aortic valved conduit (AVC), Model 11060A. It does not have direct patient contact. The product is comprised of a tray base, lid, inserts and latch and latch pin. The material for tray base, lid and inserts is made out of Radel polyphenylsulfone. The lid is secured onto the tray by means of two stainless steel latches and latch pins. This device is reusable, supplied non-sterile, and must be cleaned and sterilized prior to its use. Tray holds the sizers during sterilization, storage, and transfer to the surgical site.

3.0 Intended Use and Indications for Use

The Tray Model TRAY1190 is intended to hold the six Sizers Model 1190 during sterilization and storage and allows for multiple sterilization cycles (re-use) by steam (moist heat).

TRAY1190 is designed to contain, store and hold Sizer Model 1190 during product sterilization and transfer.

4.0 Contraindications

There are no known contraindications with the use of TRAY1190.

5.0 How Supplied

The tray is provided non-sterile and must be cleaned and sterilized before each use. It cannot be sterilized in its original packaging.

6.0 Intended User

Primary users are staff responsible for preparation of KONECT RESILIA aortic valved conduit devices prior to implant; operating room (OR) nurses and the cardiac surgeons who perform valve sizing and implantation during valve and/or aortic root repair/replacement procedures.

7.0 Cleaning Instructions

All trays and lids must be cleaned and sterilized prior to each use in accordance with appropriate hospital procedures and the following recommended instructions and parameters.

Accessories should be retired from use when signs of cracking, crazing, or other deterioration is evident. Contact your Edwards Lifesciences sales representative to obtain appropriate replacements.

CAUTION: Examine tray for signs of wear, such as dullness, cracking or crazing. Replace tray if any deterioration is observed.

7.1 Instructions for Manual Cleaning:

1. Rinse the tray under cold, running tap water to remove gross soil.
2. Prepare a neutral pH cleaning solution (such as Enzol or equivalent) per the manufacturer's recommendation.
3. Fully immerse the tray in the cleaning solution and soak for a minimum of 5 minutes.
4. After the 5-minute soak time, use a soft-bristled brush and scrub the tray for a minimum of 2 minutes to remove any visible soil.
5. Remove the tray from the cleaning solution and rinse with cold, running tap water, then in reverse osmosis/deionized water for a minimum of 1 minute.
6. If soil is visible, repeat cleaning process.

7.2 Instructions for Automated Cleaning:

1. Rinse the tray under cold, running tap water to remove gross soil.
2. Prepare a neutral pH cleaning solution (such as Enzol or equivalent) per the manufacturer's recommendation.
3. Fully immerse the tray in the cleaning solution and soak for a minimum of 1 minute. Avoid contact between the sizers.
4. Load the tray with instruments as indicated by the tray markings.
5. Place the loaded tray into the automated washer with the lid detached.
6. Set motor speed to HIGH, if applicable, and select the cycle and cycle parameters per the following recommendations:

PHASE	RECIRCULATION TIME (MINUTES)	TEMPERATURE	DETERGENT TYPE AND CONCENTRATION
Pre-wash	2:00	Cold tap water	N/A
*Enzyme wash	2:00	Hot tap water	Neutral pH enzymatic detergent (Enzol or equivalent)
Wash	2:00	50-60 °C (122-140 °F)	Alkaline detergent (Neodisher Mediclean Forte or equivalent)
Rinse	1:00	60-70 °C (140-158 °F)	N/A
Drying	7:00	90-115 °C (194-239 °F)	N/A

*The enzyme wash phase is optional and may vary based on the washer model.

8.0 Disinfection

8.1 Instructions for Thermal Disinfection:

1. Load tray with sizers into washer/disinfector with the lid detached.
2. Set the motor speed to HIGH, if applicable, and select the following recommended parameters:

PHASE	RECIRCULATION TIME (MINUTES)	WATER TEMPERATURE	DETERGENT TYPE AND CONCENTRATION
Thermal rinse	5:00	90 °C (194 °F)	N/A

8.2 Instructions for High Level Disinfection:

1. Equilibrate a bath of high-level disinfectant (Cidex OPA or equivalent) to a minimum of 20 °C (68 °F), or follow manufacturer's recommendations.
2. Fully immerse the tray in the disinfectant solution, and ensure all air bubbles are removed from the surface of the tray by wiping with a sterile lint-free cloth.
3. Allow the tray to soak for a minimum of 12 minutes.
4. Thoroughly rinse the tray by fully immersing in sterile purified water or equivalent, agitating, and soaking for a minimum of 1 minute.
5. Dry the tray using a sterile lint-free cloth.

9.0 Sterilization Instructions

CAUTION: Do not sterilize any of the accessories in their original packaging. Each institution should use procedures that include biological indicators to determine the effectiveness of the sterilization procedure.

CAUTION: Do not stack trays during sterilization.

The user is responsible for the qualification of any deviations from the recommended method of cleaning or sterilization.

Sterilizer Type	Configuration	Preconditioning	Minimum Temperature	Exposure Time	Minimum Dry Time
Gravity Displacement	Double Wrapped	N/A	132-137 °C (270-279 °F)	15-18 minutes	40 minutes
Pre-vacuum	Double Wrapped	4 pulses	132-137 °C (270-279 °F)	3-18 minutes	30 minutes

Examine all instruments prior to use. If particles are visible, repeat the cleaning and sterilization process.

10.0 Device Disposal

Devices that have signs of wear, cracking or crazing may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

11.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed.

The following table contains the Basic UDI-DI:

Product	TRAY
Models	TRAY1190
Basic UDI-DI	0690103D002OBT000UR

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at http://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

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


	ISO Reg. No. ¹	English	Français	Deutsch	Español	Italiano
	2493	Catalogue Number	Référence catalogue	Artikelnummer	Número de catálogo	Numero di catalogo
	N/A	Quantity	Quantité	Menge	Cantidad	Quantità
	2492	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
Rx only	N/A	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
 eifu.edwards.com + 1 888 570 4016	1641	Consult instructions for use on the website	Consulter le mode d'emploi sur le site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web
	0434A	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	N/A	Unique Device Identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	2606	Do not use if package is damaged	Ne pas utiliser si le conditionnement est endommagé	Bei beschädigter Verpackung nicht verwenden	No lo utilice si el envase está dañado	Non utilizzare se la confezione è danneggiata
	N/A	Authorized representative in the European Community	Représentant autorisé dans la Communauté européenne	Bevollmächtigter in der Europäischen Gemeinschaft	Representante autorizado en la Comunidad Europea	Rappresentante autorizzato nella Comunità Europea
	3082	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	2497	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	2609	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile
	N/A	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	N/A	Importer	Importateur	Importeur	Importador	Importatore


	ISO Reg. No. ¹	Nederlands	Dansk	Svenska	Ελληνικά	Português
	2493	Catalogusnummer	Katalognummer	Katalognummer	Αριθμός καταλόγου	Número de catálogo
	N/A	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	2492	Lotnummer	Partinummer	Lotnummer	Αριθμός παρτίδας	Número de lote

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları

	ISO Reg. No. ¹	Nederlands	Dansk	Svenska	Ελληνικά	Português
Rx only	N/A	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	1641	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης στον ιστότοπο	Consultar as instruções de utilização no site
	0434A	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
UDI	N/A	Unieke hulpmiddelidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	2606	Niet gebruiken als de verpakking is beschadigd	Må ikke anvendes, hvis emballagen er beskadiget	Använd inte om förpackningen är skadad	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά	Não utilizar se a embalagem estiver danificada
EC REP	N/A	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap	Autoriseret repræsentant i Det Europæiske Fællesskab	Auktoriserad representant inom Europeiska gemenskapen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia
	3082	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	2497	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	2609	Niet-steriel	Ikke-steril	Icke-steril	Μη αποστειρωμένο	Não esterilizado
MD	N/A	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	N/A	Importeur	Importør	Importör	Εισαγωγέας	Importador

	ISO Reg. No. ¹	Česky	Magyar	Polski	Slovensky	Norsk
REF	2493	Katalogové číslo	Katalógusszám	Numer katalogowy	Katalógové číslo	Katalognummer
#	N/A	Množství	Mennyiség	Ilość	Množstvo	Antall
LOT	2492	Číslo šarže	Tételszám	Numer serii	Číslo šarže	Lotnummer
Rx only	N/A	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	1641	Postupujte podle návodu k použití, který naleznete na webových stránkách	Olvassa el a használati utasítást a weboldalon	Należy zapoznać się z instrukcją użycia na stronie internetowej	Prečítajte si návod na použitie na webovej stránke	Se bruksanvisningen på nettstedet


Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen
 ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli
 ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus
 ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları

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
Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

05/24
10046356001 B / DOC-0172701 B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Web IFU

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Edwards Lifesciences GmbH