Edwards

Edwards PASCAL Transcatheter Valve Repair System

Directory	
English (EN).....	1
Français (FR).....	5
Deutsch (DE).....	10
Español (ES).....	15
Italiano (IT).....	19
Nederlands (NL).....	24
Dansk (DA).....	28
Svenska (SV).....	32
Ελληνικά (EL).....	36
Português (PT).....	41
Česky (CS).....	46
Magyar (HU).....	50
Polski (PL).....	54
Slovensky (SK).....	59
Norsk (NO).....	63
Suomi (FI).....	67
Български (BG).....	72
Română (RO).....	76
Eesti (ET).....	81
Lietuvių (LT).....	85
Latviešu (LV).....	89
Türkçe (TR).....	93
Русский (RU).....	97
Hrvatski (HR).....	102
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrazky ■ Ábrák ■ Rysunki ■ Obrazky ■ Figur ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike.....	106
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de simbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola.....	112

English

Edwards PASCAL Transcatheter Valve Repair System

Instructions for Use

The Edwards PASCAL transcatheter valve repair system includes the following configurations (herein referred to as the PASCAL system):

Table 1: Model 10000

Model Number	Device	Compatible Devices	
		Model Numbers	Device
10000IS	PASCAL Implant System	10000GS 10000ST	Guide Sheath Stabilizer
10000ISM	PASCAL Ace Implant System	20000ST 10000T	Stabilizer Rail System Table

Table 2: Model 15000

Model Number	Device	Compatible Devices	
		Model Numbers	Device
15000IS	PASCAL Implant System	10000GS 20000ST	Guide Sheath Stabilizer Rail System
15000ISM	PASCAL Ace Implant System	10000T	Table

Implant System

The Implant System consists of the Steerable Catheter (outermost layer), the Implant Catheter (innermost layer), and the Implant (hereinafter refers to implants from Model 10000IS, Model 10000ISM, Model 15000IS and Model 15000ISM). The Implant System percutaneously delivers the Implant to the valve via a femoral vein access using a transvenous approach.

Implant (Figures 1-3)

The Implant is deployed and secured to the leaflets of the valve, acting as a filler in the regurgitant orifice. The primary components of the Implant are the Spacer, Paddles, and Clasps constructed from Nitinol and covered in polyethylene terephthalate. The 10000IS

and 15000IS Implants also comprise a titanium nut and bolt, PEEK bushing, and a silicone seal. The 10000ISM and 15000ISM Implants also comprise a titanium nut, bolt, distal and proximal plate, and a silicone seal and are smaller size implants.

The Implant has four main paddle positions: Elongated, Closed, Leaflet-Capture-Ready, and Leaflet-Captured.

Steerable Catheter (Figure 4)

The Steerable Catheter has a rotational control knob that actuates the flexion mechanism to navigate and position the Implant to the target location. A radiopaque marker band located on the distal portion of the catheter indicates the end of the flex section.

Implant Catheter (Figure 4)

The Implant is provided attached to the Implant Catheter by sutures and a threaded shaft. The Implant Catheter controls the deployment of the Implant. The three primary controls are the Sliders, the threaded Actuation Knob, and the Release Knob. The Sliders control the Implant Clasps (retracting the Sliders raises the Clasps and advancing the Sliders lowers the Clasps). The threaded Actuation Knob controls the Implant Paddles (retracting the Actuation Knob closes the Paddles and advancing the Actuation Knob opens the Paddles). The Release Knob controls the release of the Implant from the Implant Catheter. The Implant Catheter is provided assembled within the Steerable Catheter.

Guide Sheath (Figure 5)

The Guide Sheath set includes a steerable Guide Sheath and Introducer. The Guide Sheath provides atrial access. It has a hydrophilic coating and a rotational control knob which actuates the flexion mechanism to position the Guide Sheath at the target location. The Introducer is compatible with a 0.035 inch (0.89 mm) guidewire.

Stabilizer (Figures 6 and 7)

The Stabilizer is indicated to aid with positioning and stabilization of the PASCAL system during implantation procedures. The Stabilizer can be attached to the system as needed any time during the procedure. The use of the Stabilizer is optional.

Table (Figure 8)

The Table is used outside of the sterile field to provide a stable platform for the Implant System, Guide Sheath, and Stabilizer. The Table is height-adjustable. The use of the Table is optional.

Loader (Figure 9)

The Peel Away Loader is used to introduce the Implant and delivery catheters through the Guide Sheath seals. The Loader is included in the Implant System and/or Guide Sheath packaging for user convenience.

1.0 Indications

The Edwards PASCAL transcatheter valve repair system is indicated for the percutaneous reconstruction of an insufficient mitral and/or tricuspid valve through tissue approximation.

2.0 Contraindications

The PASCAL system is contraindicated in patients with:

- Patient in whom a TEE is contraindicated or screening TEE is unsuccessful
- Echocardiographic evidence of intracardiac mass, thrombus, or vegetation
- Presence of an occluded or thrombosed IVC filter that would interfere with the delivery catheter, or ipsilateral deep vein thrombosis is present
- Known hypersensitivity to nitinol (nickel or titanium) or contraindication to procedural medications which cannot be adequately managed medically
- History of bleeding diathesis or coagulopathy or patient who refuses blood transfusions

Additionally, the PASCAL system is contraindicated in mitral patients with contraindication to transseptal catheterization.

3.0 Warnings

3.1 Anatomic Considerations

For optimal results, the following anatomic patient characteristics should be considered. The safety and effectiveness of the PASCAL system outside of these conditions has not been established. Use outside these conditions may interfere with placement of the Implant or native valve leaflet insertion.

Applicable to mitral and tricuspid patients:

- Evidence of moderate to severe calcification in the grasping area
- Evidence of severe calcification in the annulus or subvalvular apparatus
- Presence of significant cleft or perforation in the grasping area
- Leaflet mobility length < 8 mm

Applicable to mitral patients only:

- Flail width > 15 mm and/or flail gap > 10 mm
- Transseptal puncture height < 3.5 cm
- LA diameter ≤ 35 mm
- Presence of two or more significant jets
- Presence of one significant jet in the commissural area
- Mitral valve area (MVA) < 4.6 cm²
- LVEDD > 8.0 cm

Applicable to tricuspid patients only:

- Presence of primary non-degenerative tricuspid disease

Edwards, Edwards Lifesciences, the stylized E logo, CLASP, PASCAL, and PASCAL Ace are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

3.2 Device Handling

- The devices are designed, intended, and distributed for single use only. There are no data to support the sterility, non-pyrogenicity, and functionality of the devices after reprocessing.
- Devices should be handled using standard sterile technique to prevent infection.
- Do not expose any of the devices to any solutions, chemicals, etc., except for the sterile physiological and/or heparinized saline solution. Irreparable damage to the device, which may not be apparent under visual inspection, may result.
- Do not use any of the devices in the presence of combustible or flammable gases, anesthetics, or cleaners/disinfectants.
- Do not use the devices if the expiration date has elapsed.
- Do not use if the packaging seal is broken or if the packaging is damaged for sterile devices.
- Do not use if any of the devices were dropped, damaged or mishandled in any way.
- Standard flushing and de-airing technique should be used during preparation and throughout procedure to prevent air embolism.

3.3 Clinical Warnings

- As with any implanted medical device, there is a potential for an adverse immunological response.
- Serious adverse events, sometimes leading to surgical intervention and/or death, may be associated with the use of this system ("Potential Adverse Events"). A full explanation of the benefits and risks should be given to each prospective patient before use.
- Careful and continuous medical follow-up is advised so that implant-related complications can be diagnosed and properly managed.
- Anticoagulation therapy must be determined by the physician per institutional guidelines.

4.0 Precautions

4.1 Precautions Prior to Use

- Patient selection should be performed by a multi-disciplinary heart team specializing in the treatment of mitral and/or tricuspid regurgitation to assess patient risk and anatomical suitability.

4.2 Precautions After Use

- Long-term durability has not been established for the implant. Regular medical follow-up is advised to evaluate implant performance.
- Short-term anticoagulation therapy may be necessary after valve repair with the PASCAL device. Prescribe anticoagulation and other medical therapy per institutional guidelines.

5.0 Potential Adverse Events

Complications associated with standard cardiac catheterization, the use of anesthesia and use of the PASCAL system could lead to the following outcomes: conversion to open surgery, emergent or non-emergent reoperation, explant, permanent disability, or death. Physicians are encouraged to report suspected device related events to Edwards or the assigned hospital authorities.

The following anticipated adverse events have been identified as possible complications of the PASCAL procedure:

- Abnormal lab values
- Allergic reaction to anesthetic, contrast, heparin, Nitinol
- Anemia or decreased Hgb, may require transfusion
- Aneurysm or pseudoaneurysm
- Angina or chest pain
- Anaphylactic shock
- Arrhythmias – atrial (i.e. AF, SVT)
- Arrhythmias – ventricular (i.e. VT, VF)
- Arterio-venous fistula
- Atrial septal injury requiring intervention
- Bleeding
- Cardiac arrest
- Cardiac failure
- Cardiac injury, including perforation
- Cardiac tamponade/pericardial effusion
- Cardiogenic shock
- Chordal entanglement or rupture that may require intervention
- Coagulopathy, coagulation disorder, bleeding diathesis
- Conduction system injury which may require permanent pacemaker
- Deep vein thrombosis (DVT)
- Deterioration of native valve (e.g. leaflet tearing, retraction, thickening)
- Dislodgement of previously deployed implant
- Dyspnea
- Edema
- Electrolyte imbalance
- Emboli/embolization including air, particulate, calcific material, or thrombus
- Endocarditis
- Esophageal irritation
- Esophageal perforation or stricture
- Exercise intolerance or weakness
- Failure to retrieve any PASCAL system components
- Fever
- Gastrointestinal bleeding or infarct
- Heart failure
- Hematoma
- Hemodynamic compromise
- Hemolysis
- Hemorrhage requiring transfusion or intervention
- Hypertension
- Hypotension
- Implant deterioration (wear, tear, fracture, or other)
- Implant embolization
- Implant malposition or failure to deliver to intended site
- Implant migration
- Implant thrombosis
- Infection
- Inflammation
- LVOT obstruction

- Mesenteric ischemia
- Multi-system organ failure
- Myocardial infarction
- Nausea and/or vomiting
- Nerve injury
- Neurological symptoms, including dyskinesia, without diagnosis of TIA or stroke
- Non-neurological thromboembolic events
- Pain
- Papillary muscle damage
- Paralysis
- PASCAL system component(s) embolization
- Peripheral ischemia
- Pleural effusion
- Pulmonary edema
- Pulmonary embolism
- Reaction to anti-platelet or anticoagulation agents
- Renal failure
- Renal insufficiency
- Respiratory compromise, respiratory failure, atelectasis, pneumonia- may require prolonged ventilation
- Retroperitoneal bleed
- Septal damage or perforation
- Septicemia, sepsis
- Skin burn, injury or tissue changes due to exposure to ionizing radiation
- Single leaflet device attachment (SLDA)
- Stroke
- Syncope
- Transient ischemic attack (TIA)
- Urinary tract infection and/or bleeding
- Valve injury
- Valve stenosis
- Valvular regurgitation
- Vascular injury or trauma, including dissection or occlusion
- Vessel spasm
- Ventricular wall damage or perforation
- Wound dehiscence, delayed or incomplete healing
- Worsening of heart failure
- Worsening regurgitation / valvular insufficiency

6.0 How Supplied

6.1 Packaging

The Guide Sheath, Implant System, and Stabilizer are individually packaged and ethylene oxide sterilized. The Table is packaged and provided non-sterile.

6.2 Storage

The PASCAL system should be stored in a cool, dry place.

7.0 Directions for Use

7.1 Physician Training

The implanting physician shall be experienced in transcatheter techniques and trained on the PASCAL system and implant procedure. The final decision for PASCAL device implantation should be made by physicians specializing in treatment of mitral and/or tricuspid regurgitation in specialized centers who can determine a reasonable chance of significant clinical improvement should be expected based on stage of disease and comorbidity.

7.2 Equipment and Materials

- Standard cardiac catheterization lab equipment
- Fluoroscopy system
- Transesophageal echocardiography (TEE) capabilities (2D and 3D)
- Pigtail catheter for contrast injection (with compatible sheath)
- Venous puncture kit
- Transseptal needle, sheath, and guidewire
- Exchange length 0.035 inch (0.89 mm) guidewire
- Basins
- 50-60 cc syringes with Luer fitting
- Heparinized saline
- Hemostat
- Surgical towels (e.g. size 43 x 69 cm)
- Optional: Step-up dilators
- Optional: Continuous physiological saline drip (Rolling IV pole, IV tubing with thumbwheel occluders, 1-liter bags of heparinized sterile saline solution)
- Optional: Pressure monitoring device

7.3 Device Preparation

7.3.1 Table

Step	Procedure
1	Remove the Table(s) from packaging and inspect for damage.
2	Assemble the Table(s) as seen in Figure 8.

7.3.2 Stabilizer

Step	Procedure
1	Remove Stabilizer components from packaging and inspect for damage.
2	Assemble the Stabilizer as necessary as seen in Figure 6.

7.3.3 Guide Sheath

Step	Procedure
1	Remove the Guide Sheath, Loader, and Introducer from packaging and inspect for gross damage.
2	While keeping distal tip raised, flush and de-air Guide Sheath with heparinized saline.
3	While keeping the distal tip raised, insert Introducer into Guide Sheath. Flush the Introducer and wipe Guide Sheath with heparinized saline prior to use.

7.3.4 Implant System

Step	Procedure
1	Remove Implant System and Loader from packaging and inspect for gross damage. Verify both Slider Stopcocks are in open position. WARNING: If Slider Stopcocks are not in open position, use of the device may result in infection.
2	Advance the Actuation Knob (rotate the Actuation Knob counterclockwise or press the Actuation Button to push forward the Actuation Knob) until it is flush against the Clasp Positioning Tool.
3	Remove the Slider Pin and remove suture slack. Lock Slider Stopcocks and secure Slider Pin. Remove Clasp Positioning Tool.
4	Fully retract and advance Sliders to confirm proper Clasp motion and close the Implant (rotate the Actuation Knob clockwise or press the Actuation Button to retract the Actuation Knob).
5	Advance the Steerable Catheter. Ensure Sliders are fully retracted and Actuation Knob fully retracted. Orient Implant Catheter Handle vertically so the Release Knob is against the table.
6	Flush heparinized saline through the Implant Catheter.
7	Once Saline is seen exiting from the distal end of the Implant Catheter, lower the Implant Catheter Handle and raise the distal end of the Implant Catheter while continuing to flush with heparinized saline.
8	Fully retract Steerable Catheter. Advance the Sliders and the Actuation Knob to set Implant in Elongated position.
9	Remove Loader cap and guide the Loader cap onto the Implant System.
10	Insert the Implant through the proximal end of the Loader until it exits the distal end. Connect the Loader and Loader cap.
11	While keeping Loader and distal tip raised, flush heparinized saline through the Steerable Catheter.
12	Gradually retract Implant Catheter into Steerable Catheter and Implant into the Loader while continuing to flush through Steerable Catheter until the distal end of the Implant is fully in the Loader.

7.4 Implant Procedure

Delivery of the implant should be performed under general anesthesia with hemodynamic monitoring in an operating room, hybrid operating room or cath lab with fluoroscopic and echocardiographic imaging capabilities.

Note: Prior to implant procedure, refer to Anatomic Considerations (Section 3.1) as use outside of stated conditions may interfere with placement of the Implant or native valve leaflet insertion.

CAUTION: During the procedure, heparin should be administered so that the ACT is maintained at ≥ 250 sec.

CAUTION: Excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

7.4.1 Patient Preparation

Step	Procedure
1	Prior to sterile draping the patient, assemble and position the Table(s) between the legs of the patient, adjusting height of the Table(s) as needed. Use towels as support between Table(s) and patient's legs. WARNING: The Table is provided non-sterile; introduction of the Table into the sterile field may result in infection.
2	After sterile draping, assemble and attach Stabilizer as needed any time during procedure.

7.4.2 Femoral Vein Access and Sheath Introduction

Step	Procedure
1	Access the common femoral vein using conventional percutaneous puncture methods.
2	<u>For PASCAL Mitral Procedures:</u> Access the left atrium via transvenous, transseptal techniques using conventional percutaneous methods and place guidewire in left atrium. Dilate the vessel, as needed. CAUTION: (for mitral procedures only) Inappropriate puncture may result in cardiac structure damage requiring surgical repair or other intervention. <u>For PASCAL Tricuspid Procedures:</u> Access the right atrium using conventional percutaneous methods and place guidewire in right atrium. Dilate the vessel, as needed.

Step	Procedure
3	For PASCAL Mitral Procedures: Insert Guide Sheath with Introducer over guidewire until Guide Sheath tip is securely across the septum, using flex mechanism as needed. <u>For PASCAL Tricuspid Procedures:</u> Insert Guide Sheath with Introducer over guidewire until Guide Sheath tip is inside the right atrium. CAUTION: Excessive manipulation may result in dislodgement or disturbance of a previously implanted device, cardiac structure damage requiring surgical repair or other intervention.
4	Remove the Introducer and guidewire. Do not aspirate and flush the Guide Sheath until the Implant System is inserted. CAUTION: Aspiration or connection of a continuous physiological saline flush to the Guide Sheath prior to Implant System insertion may result in air embolism.

7.4.3 Implant Delivery

Step	Procedure
1	Insert the Implant System with the Loader into the Guide Sheath.
2	Advance Implant System until the Implant exits the Loader. Retract and peel away the Loader.
3	Aspirate and flush Guide Sheath with heparinized saline. Utilizing the specified syringe, aspirate a minimum of 45 cc. CAUTION: Failure to fully aspirate Guide Sheath may result in air embolism.
4	If desired, connect the continuous physiological saline drip to the Implant Catheter. CAUTION: Connection of the continuous physiological saline drip to the Implant System prior to aspiration may result in air embolism.
5	Advance Implant System until the Implant exits the distal end of the Guide Sheath.
6	Retract the Actuation Knob to get the Implant in Closed position. Retract the Sliders.
7	Adjust Guide Sheath as needed.
8	At the discretion of the treating physician, if pressure monitoring is used to continuously assess atrial pressure during procedure, please follow the pressure monitor manufacturer's instruction of use. Connect a fluid filled pressure monitoring device to the steerable catheter. Aspirate and then calibrate at the patient's heart level before obtaining measurement. Note: Pressure monitoring should be used in conjunction with echo. Pressure should be reconciled against echo and Doppler readings. When assessing atrial pressure, ensure that the distal tip of the Implant Catheter is fully exposed from the Steerable Catheter.
9	Advance Implant System as needed. Manipulate Steerable Catheter and Guide Sheath (flex-unflex, torque in opposing directions, advance-retract) as needed until Implant is centered in the target coaptation zone with the appropriate trajectory. CAUTION: Excessive manipulation may result in dislodgement or disturbance of a previously implanted device, cardiac structure damage requiring surgical repair, or other intervention. Note: The radiopaque marker band on the Steerable Catheter indicates the end of the flex section and can be visualized on fluoroscopy.
10	Advance the Actuation Knob to get the Implant into Leaflet-Capture-Ready position. Note: For a tricuspid procedure, after placing implant in Leaflet-Capture-Ready Position, pull the Slider Pin and move one slider to identify which Implant Clasp it controls via imaging and secure Slider Pin.
11	Torque the Implant Catheter, as needed, to orient the Paddles.
12	Advance the Implant through the valve until Paddles are below the free edge of the leaflets.
13	Verify location and orientation of the Implant and adjust position slightly as needed. CAUTION: Excessive manipulation of the Implant below the leaflets may cause the Implant to be entangled in the chords; chordal entanglement may result in cardiac injury, worsening regurgitation, difficulty or inability to remove the Implant requiring additional intervention.
14	Under imaging guidance, retract the Implant until leaflets are positioned between Paddles and Clasps.
15	Advance Slider(s) so the leaflet(s) are secured between the Clasps and Paddles. This can be performed for both leaflets simultaneously (Slider Pin engaged to move both Clasps) or each leaflet individually (Slider Pin disengaged to move individual Clasp).
16	Verify leaflet insertion with imaging. If leaflets are not secured between Clasps and Paddles, retract the Sliders to release the leaflets and reattempt.
17	Once leaflets are secured between the Clasps and Paddles, close the Implant.
18	Advance Implant Catheter slightly to release tension on leaflets.

Step	Procedure
19	<p>Assess regurgitation, and reposition as needed. Once the Implant position is confirmed, ensure Implant is closed.</p> <p>If repositioning within the ventricle is needed, retract the Sliders and advance the Actuation Knob to set the Implant in Leaflet-Capture-Ready position. Adjust Clasps and Implant orientation as needed.</p> <p>If repositioning within the atrium is needed, retract the Sliders and advance the Actuation Knob to elongate the Implant slowly under fluoroscopic guidance while ensuring that the Actuation Wire does not bend, and retract the Implant back into the atrium.</p> <p>CAUTION: Failure to elongate the Implant when retracting into the atrium during repositioning may result in leaflet damage or chordal entanglement.</p> <p>CAUTION: Failure to release leaflets from Clasps and Paddles prior to repositioning may result in leaflet damage.</p>
20	<p>To release the Implant from the catheter:</p> <ol style="list-style-type: none"> Ensure that the distal tip of the Implant Catheter is fully exposed from the Steerable Catheter. Cut outside suture on the proximal end of each Slider. Open both Slider Stopcocks to unlock sutures. Pull the Slider Pin to fully remove sutures. Close both Slider Stopcocks after removal of sutures. Remove the Release Pin. Rotate counterclockwise and retract the Release Knob until the Implant is released, as confirmed via imaging. <p>Note: Prior to Implant release, if needed, it is possible to retrieve the Implant System back into the Guide Sheath for removal. To retrieve:</p> <ol style="list-style-type: none"> Retract the Sliders. Elongate the Implant slowly under fluoroscopic guidance while ensuring that the Actuation Wire does not bend. Then retract Implant into the atrium. Set the Implant in Closed position. Unflex the Steerable Catheter and retract the Implant System until the Implant is adjacent to the tip of the Guide Sheath. Advance the Sliders. Set the Implant in Elongated position. Retract the Sliders to open the Clasps to approximately 45° on each side. Retract entire Implant System through the Guide Sheath. <p>CAUTION: Failure to cut suture in prescribed location may result in inability to release implant or introduction of fiber that may lead to micro-embolism.</p> <p>CAUTION: Failure to follow prescribed release steps may result in difficulty or inability to release Implant, requiring additional intervention.</p> <p>CAUTION: Releasing the Implant prior to confirmation that leaflets are securely captured between Paddles and Clasps may result in Implant movement or dislodgement leading to a single leaflet device attachment (SLDA) or other potential adverse events requiring additional intervention.</p> <p>WARNING: Re-use of the devices (including Implant System and Guide Sheath) after retrieval may cause embolism of foreign material or infection. Device may malfunction if re-use is attempted.</p> <p>Note: If an additional implant is placed at the decision of the treating physician, caution should be taken to avoid dislodgement of the previously placed implant. Crossing the valve in a low profile implant configuration can minimize interaction with the previously placed implant.</p> <p>CAUTION: Excessive manipulation may result in dislodgement or disturbance of a previously implanted device, cardiac structure damage requiring surgical repair or other intervention.</p>

7.4.4 Device Removal and Closure

Step	Procedure
1	<p>Retract Implant Catheter completely into Steerable Catheter. Gradually unflex and remove Implant System. Gradually unflex and remove Guide Sheath.</p> <p>CAUTION: Failure to unflex devices prior to removal may result in vessel damage.</p>
2	<p>Perform standard percutaneous closure of access site.</p>

8.0 Magnetic Resonance (MR) Safety

Non-clinical testing has demonstrated that the PASCAL implant is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic fields of 1.5 T and 3.0 T
- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode).

Under the scan conditions defined above, the Implant is expected to produce a maximum temperature rise of less than 4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device in a worst case multiple implant configuration extends up to 15 mm from the implant when imaged in the worst case gradient echo pulse sequence in a 3.0 T MRI system.

9.0 Recovered Implant and Device Disposal

Edwards Lifesciences is interested in obtaining recovered clinical specimens of the Implant for analysis. A written report summarizing our findings will be provided upon completion of our evaluation. Please contact Edwards for return of the recovered Implant.

If you do decide to return any of the devices, please follow the following instructions:

- Unopened Package with Sterile Barrier Intact:**

If the pouches have not been opened, return the device in its original packaging.

- Package Opened but Not Implanted:**

If a pouch is opened, the device is no longer sterile. Please return the device in its original packaging.

- Explanted Implant:**

The explanted implant should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to Edwards.

9.1 Disposal

Used devices may be handled and disposed of in the same manner as hospital waste and bio-hazardous materials in accordance with local regulations as there are no special risks related to the disposal of these devices.

10.0 Summary of Clinical Experience

10.1 CLASP Study

Clinical data in this section includes information obtained from the CLASP clinical study which studied the PASCAL transcatheter valve repair system in the mitral valve.

A multi-center, multi-national, prospective, single-arm study (CLASP) was conducted to assess the safety, performance, and clinical outcomes of the PASCAL system. All enrolled study patients were assessed for clinical follow-up at 30 days, 6 months, 1 year, and will continue annually for 5 years post-implant procedure.

The primary safety endpoint of the CLASP study was a composite of major adverse events (MAEs) at 30 days. The MAEs include: cardiovascular mortality, stroke, myocardial infarction, new need for renal replacement therapy, severe bleeding, and re-intervention for study device-related complications.

The primary performance endpoints of the study include device success, procedural success, and clinical success. The secondary endpoints of the study include clinical, safety, and functional outcomes at 30-day, 6 month, 1 year, and annual follow-up time points.

Device success is defined as device deployment as intended and successful delivery system retrieval as intended at the time of the patient's exit from the cardiac catheterization laboratory. Analysis of device success was performed per device.

Procedural success is defined as device success with MR severity $\leq 2+$ at discharge (Echo Core Lab-evaluated) and without the need for a surgical or percutaneous intervention prior to hospital discharge. Procedural success was analyzed per patient.

Clinical success is defined as procedural success with evidence of MR reduction MR $\leq 2+$ and without MAEs at 30 days (analyzed per patient).

An independent core lab assessed all echocardiographic data. An independent clinical events committee (CEC) adjudicated safety events and a data safety monitoring board (DSMB) independently reviewed aggregate safety data and evaluated trends of adverse events and their effect on trial conduct and device risk assessment.

10.1.1 CLASP Study Results

The mean age of the patients treated was 75.4 years and 55.0% were male. All patients had NYHA Class II, III or IV heart failure. The mean Logistic EuroSCORE I, EuroSCORE II and STS Mortality Score were 14.4%, 5.8% and 4.7%, respectively. At baseline, 50.9% had moderate-severe MR, and 48.1% had severe MR.

10.1.2 Performance

The performance endpoint included three components of success: device, procedural, and clinical. Device success was achieved in 92.0% of devices attempted. Procedural success was achieved in 93.5% of patients. Clinical success was achieved in 86.0% of patients.

10.1.3 Safety

At 30 days, the composite MAE rate was 8.0%. Nine patients experienced 11 MAEs prior to 30-day follow-up. The CEC Adjudicated MAE at 30 Days by Counts are shown below.

Major Adverse Event (MAE)	Summary Statistics
Cardiovascular Mortality	0.9% (1/112)
Stroke	0.9% (1/112)
Myocardial Infarction	0.0% (0/112)
New Need for Renal Replacement Therapy	0.0% (0/112)
Severe Bleeding	7.1% (8/112)
Re-Intervention for Study Device Related Complications	0.9% (1/112)
Composite MAE Rate	8.0% (9/112)

Note: Categorical measures - % (n/Total no)

In the overall population, MR grade reduction (i.e. MR $\leq 2+$) was observed in 95.3% of patients at discharge, 96.1% at 30 days, 98.8% at 6 months, and 100% at 1 year.

10.1.4 Study Conclusion

The data collected in the CLASP study supports the safety and performance of the PASCAL system in patients with mitral regurgitation. The number of patients with follow-up of one year or more is limited and long-term follow-up data is collected by means of a Post Marketing Clinical Follow-up study.

10.2 CLASP TR Study

Clinical data in this section includes information obtained from the tricuspid clinical study (CLASP TR) of the PASCAL transcatheter valve repair system.

A prospective, single arm, multi-center study (CLASP TR) was conducted to assess the safety, performance, and clinical outcomes of the PASCAL system. All enrolled study patients were assessed for clinical follow-up at 30 days, 6 months, 1 year, and will continue annually for 5 years post-implant procedure.

The primary safety endpoint of the CLASP TR study was a composite of major adverse events (MAEs) at 30 days. The MAEs include: cardiovascular mortality, myocardial infarction, stroke, renal complications requiring unplanned dialysis or renal replacement therapy, severe bleeding, unplanned or emergency re-intervention (either percutaneous or surgical) related to the device, and major access site and vascular complications requiring intervention.

The primary performance endpoints of the study include device success, procedural success, and clinical success. The secondary endpoints of the study include clinical, safety, and functional outcomes at 30-day, 6-month, 1-year, and annual follow-up time points.

Device success is defined as device deployed as intended and successful delivery system retrieval as intended at the time of the patient's exit from the cardiac catheterization laboratory. Analysis of device success was performed per device.

Procedural success is defined as device success with evidence of a reduction in TR grade by at least one grade (scale: non/trace, mild, moderate, severe, massive, torrential) at the end of the

procedure, and without the need for a surgical or percutaneous intervention prior to hospital discharge. Procedural success was analyzed per patient.

Clinical success is defined as procedural success without MAEs at 30 days (analyzed per patient).

An independent core lab assessed all echocardiographic data. An independent clinical events committee (CEC) adjudicated safety events and a data safety monitoring board (DSMB) independently reviewed aggregate safety data and evaluated trends of adverse events and their effect on trial conduct and device risk assessment.

10.2.1 CLASP TR Study Results

The mean age of the patients treated was 76.3 years and 53% were female. All patients had NYHA Class II or III heart failure. The mean EuroSCORE II and STS Mortality Score were 5.3% and 7.3%, respectively. At baseline, 81% of patients presented with severe or greater TR.

10.2.2 Performance

The performance endpoint included three components of success: device, procedural, and clinical. Device success was achieved in 82.2% of devices attempted. Of the implanted patients, procedural success was achieved in 95.5% of patients and clinical success was achieved in 86.4% of patients.

10.2.3 Safety

At 30 days, the composite MAE rate was 5.9%. Two patients experienced 3 MAEs prior to 30-day follow-up. The CEC-adjudicated MAEs at 30 Days by Counts are shown below.

Major Adverse Event (MAE)	Summary Statistics
Cardiovascular Mortality	0.0% (0/34)
Myocardial Infarction (MI)	0.0% (0/34)
Stroke	0.0% (0/34)
Renal Complications requiring unplanned dialysis or renal replacement therapy	0.0% (0/34)
Severe Bleeding	5.9% (2/34)
Unplanned or emergency re-intervention (either percutaneous or surgical) related to the device	0.0% (0/34)
Major access site and vascular complications requiring intervention	0.0% (0/34)
Composite MAE Rate	5.9% (2/34)

Note: Categorical measures - % (n/Total no)

In the implanted population, a reduction in TR grade by at least one grade was observed in 88.9% and 82.6% of patients at discharge and 30 days, respectively. A reduction in TR grade by at least two grades was observed in 65.2% of patients at 30 days.

10.2.4 Study Conclusion

The data collected, to date, in the CLASP TR study supports the safety and performance of the PASCAL system in patients with tricuspid regurgitation. The number of patients with follow-up of one year or more is limited and long-term follow-up data is collected by means of a Post Marketing Clinical Follow-up study.

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






	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk
	Catalogue Number	Référence catalogue	Artikelnummer	Número de catálogo	Numero di catalogo	Catalogusnummer	Katalognummer
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer	Partinummer
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid	Mængde
	Contents	Contenu	Inhalt	Contenido	Contenuti	Inhoud	Indhold
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte	Anvendelig længde
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken	Må ikke genanvendes
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op	Forsigtig
	Warning	Mise en garde	Warnung	Advertencia	Avvertenza	Waarschuwing	Advarsel
	Note	Remarque	Hinweis	Nota	Nota	Opmerking	Bemærk
	Consult instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen
	Do not use if package is damaged	Ne pas utiliser si le conditionnement est endommagé	Bei beschädigter Verpackung nicht verwenden	No lo utilice si el envase está dañado	Non utilizzare se la confezione è danneggiata	Niet gebruiken als de verpakking is beschadigd	Må ikke anvendes, hvis emballagen er beskadiget
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt
	Keep away from sunlight	Tenir à l'abri de la lumière du soleil	Vor Sonneneinstrahlung schützen	Mantener alejado de la luz del sol	Tenere lontano dalla luce solare	Uit het zonlicht houden	Beskyttes mod sollys
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantener asciutto	Droog houden	Opbevares tørt
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling
	Sterilized using liquid chemical	Stérilisé à l'aide d'un produit chimique liquide	Mit flüssigen Chemikalien sterilisiert	Esterilizado con sustancia química líquida	Sterilizzato con sostanze chimiche liquide	Gesteriliseerd met vloeibare chemische stof	Steriliseret med flydende kemikalie
	Do not re-sterilize	Ne pas restériliser	Nicht reesterilisieren	No volver a esterilizar	Non risterilizzare	Niet opnieuw steriliseren	Må ikke reesteriseres

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola

	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile	Niet-steriel	Ikke-steril
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen	Ikke-pyrogen
	Non-DEHP	Sans DEHP	Nicht DEHP-haltig	Sin-DEHP	Non contiene DEHP	Niet-DEHP	Ikke-DEHP
	Use-by date	Date d'expiration	Verfallsdatum	Fecha de caducidad	Data di scadenza	Vervaldatum	Sidste anvendelsesdato
	Serial Number	Numéro de série	Serienummer	Número de serie	Numero di serie	Serienummer	Serienummer
	Authorized representative in the European Community	Représentant autorisé dans la Communauté européenne	Bevollmächtigter in der Europäischen Gemeinschaft	Representante autorizado en la Comunidad Europea	Rappresentante autorizzato nella Comunità Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap	Autoriseret repræsentant i Det Europæiske Fællesskab
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant	Producent
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum	Fremstillingsdato
	For use with size 44 mm Edwards transcatheter heart valve	À utiliser avec une valve cardiaque transcathéter Edwards de 44 mm	Zur Verwendung mit einer 44-mm-Transkatheter-Herzklappe von Edwards	Para su uso con una válvula cardíaca transcatéter Edwards de 44 mm	Da utilizzare con valvole cardiache transcatetere Edwards da 44 mm	Voor gebruik met de Edwards transkatheterhartklep van 44 mm	Til brug med størrelse 44 mm Edwards transkateterklap
	For use with size 48 mm Edwards transcatheter heart valve	À utiliser avec une valve cardiaque transcathéter Edwards de 48 mm	Zur Verwendung mit einer 48-mm-Transkatheter-Herzklappe von Edwards	Para su uso con una válvula cardíaca transcatéter Edwards de 48 mm	Da utilizzare con valvole cardiache transcatetere Edwards da 48 mm	Voor gebruik met de Edwards transkatheterhartklep van 48 mm	Til brug med størrelse 48 mm Edwards transkateterklap
	Temperature limit	Limite de température	Temperaturgrenze	Límite de temperatura	Limiti di temperatura	Temperatuurlimiet	Temperaturgrænse
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter	Udvendig diameter
	Inner diameter	Diamètre interne	Innendurchmesser	Diámetro interior	Diametro interno	Binnendiameter	Indre diameter
	Recommended guidewire length	Longueur recommandée du fil-guide	Empfohlene Führungsdrahtlänge	Longitud recomendada de la guía	Lunghezza consigliata per il filo guida	Aanbevolen voerdraadlengte	Anbefalet guidewirelængde
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad	Anbefalet guidewirestørrelse
	Guidewire compatibility	Compatibilité du fil-guide	Führungsdraht-Kompatibilität	Compatibilidad con la guía	Compatibilità del filo guida	Compatibiliteit voerdraad	Guidewirekompatibilitet
	Size	Taille	Größe	Tamaño	Misura	Afmetingen	Størrelse
	Catheter shaft size	Taille du corps du cathéter	Katheterschaftgröße	Tamaño del eje del catéter	Misura dello stelo del catetere	Afmetingen katheterschacht	Kateterrørsstørrelse
	Balloon diameter	Diamètre du ballonnet	Ballondurchmesser	Diámetro del balón	Diametro del palloncino	Ballondiameter	Ballondiameter
	Balloon working length	Longueur utile du ballonnet	Arbeitslänge des Ballons	Longitud útil del balón	Lunghezza utile del palloncino	Werklengte ballon	Ballonens arbejdslængde

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola

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

	Norsk	Suomi	Български	Română	Eesti	Lietuvių	Latviešu
	Holdes vekk fra sollys	Säilytä auringonvalolta suojattuna	Да се пази от слънчева светлина	A se feri de lumina solară	Hoida eemal päikesevalgusest	Saugoti nuo saulės šviesos	Sargāt no saules gaismas
	Oppbevares tørt	Pidä kuivana	Да се пази сухо	A se păstra uscat	Hoida kuivalt	Laikykite sausoje vietoje	Uzglabāt sausu
	Medisinsk utstyr	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos prietaisas	Mediciniska ierice
	Se bruksanvisningen på nettsiden	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje	Skatīt lietošanas instrukciju tīmekļa vietnē
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Sterilisert med etylenoksid	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksīdu
	Sterilisert med stråling	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant	Sterilizēts apstarojot
	Sterilisert med flytende kjemikalie	Steriloitu nestemäisellä kemikaalilla	Стерилизирано с течни химикали	Sterilizat cu substanțe chimice lichide	Steriliseeritud vedelat kemikaali kasutades	Sterilizuota skystomis cheminėmis medžiagomis	Sterilizēts, izmantojot šķidru ķīmisko vielu
	Må ikke resteriliseres	Älä steriloitu uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti
	Usteril	Epästeriili	Нестерилно	Nesteril	Mittesteriilne	Nesterilus	Nesterils
	Ikke-pyrogen	Pyrogeeniton	Непирогенно	Non pyrogen	Mittepürogeenne	Nepirogeniškas	Nepirogēns
	Ikke-DEHP	Ei sisällä DEHP:tä	Не съдържа DEHP	Nu conține DEHP	Mitte-DEHP	Ne DEHP	DEHP nesaturōšs
	Utløpsdato	Viimeinen käyttöpäivämäärä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Tinkamumo terminas	Derīguma termiņš
	Serienummer	Sarjanumero	Сериен номер	Număr de serie	Seerianumber	Serijos numeris	Sērijas numurs
	Autorisert representant i Det europeiske fellesskap	Valtuutettu edustaja Euroopan yhteisössä	Упълномощен представител в Европейската общност	Reprezentant autorizat în Comunitatea Europeană	Volitatud esindaja Euroopa Ühenduses	Igaliotasis atstovas Europos Bendrijoje	Pilnvarotais pārstāvis Eiropas Kopienā
	Produsent	Valmistaja	Производител	Producător	Tootja	Gamintojas	Ražotājs
	Produksjonsdato	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data	Izgatavošanas datums
	Til bruk med Edwards transkateter hjerteklaff, størrelse 44 mm	Käytetään 44 mm:n kokoisten katetrilla asennettavien Edwards - sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа на Edwards с размер 44 mm	A se utiliza cu valva cardiacă transcateter Edwards de mărimea 44 mm	Kasutamiseks ettevõtte Edwards 44 mm transkateetriga südameklapiga	Skirta naudoti su 44 mm dydžio transkateeteriniu širdies vožtuvu „Edwards“	Lietošanai ar 44 mm izmēra Edwards transkatetra sirds vārstuli
	Til bruk med Edwards transkateter hjerteklaff, størrelse 48 mm	Käytetään 48 mm:n kokoisten katetrilla asennettavien Edwards - sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа на Edwards с размер 48 mm	A se utiliza cu valva cardiacă transcateter Edwards de mărimea 48 mm	Kasutamiseks ettevõtte Edwards 48 mm transkateetriga südameklapiga	Skirta naudoti su 48 mm dydžio transkateeteriniu širdies vožtuvu „Edwards“	Lietošanai ar 48 mm izmēra Edwards transkatetra sirds vārstuli
	Temperaturrense	Lämpötilaraja	Температурно ограничение	Limită de temperatură	Temperatuuripiirang	Temperatūros apribojimas	Temperatūras ierobežojums

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola

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Edwards

EC REP

Edwards Lifesciences Services GmbH

Edisonstr. 6
85716 Unterschleissheim
Germany

03/2021
10038791002 A
© Copyright 2021, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC

One Edwards Way
Irvine, CA 92614 USA

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU