Edwards

Edwards PASCAL Precision Transcatheter Valve Repair System

Directory	
English (en).....	1
Français (fr).....	6
Deutsch (de).....	11
Español (es).....	16
Italiano (it).....	21
Nederlands (nl).....	26
Dansk (da).....	32
Svenska (sv).....	37
Ελληνικά (el).....	42
Português (pt).....	47
Česky (cs).....	52
Magyar (hu).....	57
Polski (pl).....	62
Slovensky (sk).....	68
Norsk (no).....	73
Suomi (fi).....	78
Български (bg).....	83
Română (ro).....	88
Eesti (et).....	93
Lietuvių (lt).....	98
Latviešu (lv).....	103
Türkçe (tr).....	108
Русский (ru).....	113
Hrvatski (hr).....	119
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунок ■ Slike.....	124
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Simbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola.....	130

English

Edwards PASCAL Precision Transcatheter Valve Repair System

Instructions for Use

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

The Edwards PASCAL Precision transcatheter valve repair system includes the following model numbers (herein referred to as the PASCAL Precision system):

Model Number	Device
20000IS	PASCAL Precision system – implant system
20000ISM	PASCAL Precision system – PASCAL Ace implant system
20000GS	PASCAL Precision system – guide sheath

The PASCAL Precision system is compatible with the following:

Model Number	Device
10000T	PASCAL system – table
20000ST	PASCAL system – stabilizer rail system

• Implant System (Figure 4)

The implant system consists of the steerable catheter (outermost layer), the implant catheter (innermost layer), and the implant (hereinafter refers to the PASCAL and PASCAL Ace implants). The implant system percutaneously delivers the implant to the valve via a femoral vein access using a transvenous, transseptal [mitral] and transvenous [tricuspid] approach.

• Implant (Figures 1-3)

The implant is deployed and secured to the leaflets of the valve, acting as a filler in the regurgitant orifice. The primary components of the implant are the spacer, paddles, and clasps constructed from Nitinol and covered in polyethylene terephthalate. The PASCAL Ace implant is a smaller size to provide physicians options. Suggested considerations for selecting PASCAL Ace implant include smaller landing zones and dense chordal regions. Refer to the PASCAL Precision system Physician Training Materials for additional information on the differences in size of spacer and other components for PASCAL and PASCAL Ace implants, including additional considerations for implant selection.

The implant has four main paddle positions: elongated, closed, leaflet-capture-ready, and leaflet-captured.

Note: PASCAL Ace implant is a naming convention that refers to an additional implant size with the same indication for use as the PASCAL implant.

• Implant Catheter (Figure 4)

The implant is provided attached to the implant catheter by sutures and a threaded shaft. The implant catheter controls the deployment of the implant. The four primary controls are the clasp sliders, the paddle knob, the implant release knob, and the suture locks. The clasp sliders control the clasps (retracting the clasp sliders raises the clasps and advancing the clasp sliders lowers the clasps). The paddle knob controls the paddles (rotating the paddle knob clockwise closes the paddles and rotating the paddle knob counterclockwise opens the paddles). The implant release knob controls the release of the implant from the implant catheter. The suture locks control release of the sutures from the clasps. The implant catheter is provided assembled within the steerable catheter.

• Steerable Catheter (Figure 4)

The steerable catheter has a rotational control knob (flex knob) that actuates the flexion mechanism to navigate and position the implant to the target location. A radiopaque marker band located on the distal portion of the catheter indicates the end of the flex section.

• Guide Sheath (Figure 5)

The guide sheath is used to provide atrial access. It has a hydrophilic coating and a rotational control knob (flex knob) which actuates the flexion mechanism to position the guide sheath at the target location.

• Introducer (Figure 5)

The introducer is used to facilitate tracking the guide sheath to the desired location. The introducer is compatible with a 0.035 inch (0.89 mm) guidewire. The introducer is included in the guide sheath packaging.

• Loader (Figure 8)

The loader is used to introduce the implant and delivery catheters through the guide sheath seals. The loader is included in the implant system and guide sheath packaging for user convenience.

• Stabilizer Rail System (Figure 6)

The stabilizer rail system is intended to aid with positioning and stabilization of the PASCAL Precision system during implantation procedures. The stabilizer can be attached to the guide sheath and implant system as needed any time during the procedure. The use of the stabilizer rail system is optional.

• Table (Figure 7)

The table is intended to provide a stable platform for the implant system, guide sheath, and stabilizer rail system of the PASCAL Precision system. The table is used outside of the sterile field. The table is height-adjustable. The use of the table is optional.

1.0 Indications for Use

1.1 Intended Use

The PASCAL Precision system is intended to repair an insufficient mitral and/or tricuspid valve via percutaneous reconstruction through tissue approximation. The PASCAL Precision system percutaneously delivers the implant to the valve via a femoral vein access using a transvenous, transseptal (mitral) and transvenous (tricuspid) approach.

1.2 Target Patient Population

The PASCAL Precision system is intended for:

- Adult patients with clinically significant, symptomatic mitral regurgitation (moderate-to-severe or severe MR)
- Adult patients with clinically significant, symptomatic tricuspid regurgitation (severe or greater TR) despite medical therapy

Patient selection should be performed by a multidisciplinary expert heart team specializing in the treatment of mitral and tricuspid regurgitation. The heart team should weigh the benefits and risks of all possible interventions prior to treating patients with the PASCAL Precision system. Patients are candidates for the PASCAL Precision system if they are not deemed suitable for cardiac surgery, including minimally invasive cardiac surgery, and are anatomically suitable for treatment with the PASCAL Precision system.

2.0 Contraindications

The PASCAL Precision system is contraindicated in mitral and tricuspid patients with:

- Patient in whom a TEE is contraindicated or screening TEE is unsuccessful
- Echocardiographic evidence of intracardiac mass, thrombus, or vegetation
- Presence of an occluded or thrombosed IVC filter that would interfere with the delivery catheter, or ipsilateral deep vein thrombosis is present
- Known hypersensitivity to nitinol (nickel or titanium) or contraindication to procedural medications which cannot be adequately managed medically
- History of bleeding diathesis or coagulopathy or patient who refuses blood transfusions

Additionally, the PASCAL Precision system is contraindicated in mitral patients with contraindication to transseptal catheterization.

3.0 Warnings

3.1 Anatomic Considerations

For optimal results, the following anatomic patient characteristics should be considered. Valve anatomy which might limit proper PASCAL Precision system access, use and/or deployment, or

sufficient reduction in mitral or tricuspid regurgitation should be considered by a multi-disciplinary heart team. Safety and effectiveness have not been established for patients with anatomic characteristics including, but not limited to, the following:

Applicable to mitral and tricuspid patients:

- Evidence of moderate to severe calcification in the grasping area
- Evidence of severe calcification in the annulus or subvalvular apparatus
- Presence of significant cleft or perforation in the grasping area
- Leaflet mobility length < 8 mm

Applicable to mitral patients only:

- Flail width > 15 mm and/or flail gap > 10 mm
- Transseptal puncture height < 3.5 cm
- LA diameter ≤ 35 mm
- Presence of two or more significant jets
- Presence of one significant jet in the commissural area
- Mitral valve area (MVA) < 4.0 cm²
- LVEDD > 8.0 cm

Applicable to tricuspid patients only:

- Presence of primary non-degenerative tricuspid disease

3.2 Device Handling

- The devices are designed, intended, and distributed for single use only. There are no data to support the sterility, non-pyrogenicity, and functionality of the devices after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.
- Devices should be handled using standard sterile technique to prevent infection.
- Do not expose any of the devices to any solutions, chemicals, etc., except for the sterile physiological and/or heparinized saline solution. Irreparable damage to the device, which may not be apparent under visual inspection, may result.
- Do not use any of the devices in the presence of combustible or flammable gases, anesthetics, or cleaners/disinfectants.
- Do not use the devices if the expiration date has elapsed.
- Do not use if the packaging seal is broken or if the packaging is damaged for sterile devices.
- Do not use if any of the devices were dropped, damaged, or mishandled in any way.
- Standard flushing and de-airing technique should be used during preparation and throughout procedure to prevent air embolism.

3.3 Clinical Warnings

- As with any implanted medical device, there is a potential for an adverse immunological response.
- Serious adverse events, sometimes leading to surgical intervention and/or death, may be associated with the use of this system ("Potential Adverse Events"). A full explanation of the benefits and risks should be given to each prospective patient before use.
- Careful and continuous medical follow-up is advised so that implant-related complications can be diagnosed and properly managed.
- Anticoagulation therapy must be determined by the physician per institutional guidelines.
- The PASCAL Precision system has not been evaluated in pregnant or pediatric patients.

4.0 Precautions

4.1 Precautions Prior to Use

- Patient selection should be performed by a multi-disciplinary heart team specializing in the treatment of mitral and/or tricuspid regurgitation to assess patient risk and anatomical suitability after considering the various device treatment options.

4.2 Precautions After Use

- Long-term durability has not been established for the implant. Regular medical follow-up is advised to evaluate implant performance.
- Short-term anticoagulation therapy may be necessary after valve repair with the PASCAL Precision system. Prescribe anticoagulation and other medical therapy per institutional guidelines.

5.0 Potential Adverse Events

Complications associated with standard cardiac catheterization, the use of anesthesia, and use of the PASCAL Precision system could lead to the following outcomes: conversion to open surgery, emergent or non-emergent reoperation, explant, permanent disability, or death. Physicians are encouraged to report suspected device related events to Edwards or the assigned hospital authorities.

The following anticipated adverse events have been identified as possible complications of the PASCAL Precision procedure:

- Abnormal lab values
- Allergic reaction to anesthetic, contrast, heparin, Nitinol
- Anemia or decreased Hgb, may require transfusion
- Aneurysm or pseudoaneurysm
- Angina or chest pain
- Anaphylactic shock
- Arrhythmias – atrial (i.e. AF, SVT)
- Arrhythmias – ventricular (i.e. VT, VF)
- Arterio-venous fistula
- Atrial septal injury requiring intervention
- Bleeding
- Cardiac arrest
- Cardiac failure
- Cardiac injury, including perforation
- Cardiac tamponade/pericardial effusion
- Cardiogenic shock
- Chordal entanglement or rupture that may require intervention
- Coagulopathy, coagulation disorder, bleeding diathesis
- Conduction system injury which may require permanent pacemaker
- Deep vein thrombosis (DVT)
- Deterioration of native valve (e.g. leaflet tearing, retraction, thickening)
- Dislodgement of previously deployed implant
- Dyspnea
- Edema

- Electrolyte imbalance
- Emboli/embolization including air, particulate, calcific material, or thrombus
- Endocarditis
- Esophageal irritation
- Esophageal perforation or stricture
- Exercise intolerance or weakness
- Failure to retrieve any PASCAL Precision system components
- Fever
- Gastrointestinal bleeding or infarct
- Heart failure
- Hematoma
- Hemodynamic compromise
- Hemolysis
- Hemorrhage requiring transfusion or intervention
- Hypertension
- Hypotension
- Implant deterioration (wear, tear, fracture, or other)
- Implant embolization
- Implant malposition or failure to deliver to intended site
- Implant migration
- Implant thrombosis
- Infection
- Inflammation
- LVOT obstruction
- Mesenteric ischemia
- Multi-system organ failure
- Myocardial infarction
- Nausea and/or vomiting
- Nerve injury
- Neurological symptoms, including dyskinesia, without diagnosis of TIA or stroke
- Non-neurological thromboembolic events
- Pain
- Papillary muscle damage
- Paralysis
- PASCAL Precision system component(s) embolization
- Peripheral ischemia
- Pleural effusion
- Pulmonary edema
- Pulmonary embolism
- Reaction to anti-platelet or anticoagulation agents
- Renal failure
- Renal insufficiency
- Respiratory compromise, respiratory failure, atelectasis, pneumonia - may require prolonged ventilation
- Retroperitoneal bleed
- Septal damage or perforation
- Septicemia, sepsis
- Skin burn, injury, or tissue changes due to exposure to ionizing radiation
- Single leaflet device attachment (SLDA)
- Stroke
- Syncope
- Transient ischemic attack (TIA)
- Urinary tract infection and/or bleeding
- Valve injury
- Valve stenosis
- Valvular regurgitation
- Vascular injury or trauma, including dissection or occlusion
- Vessel spasm
- Ventricular wall damage or perforation
- Wound dehiscence, delayed or incomplete healing
- Worsening of heart failure
- Worsening regurgitation / valvular insufficiency

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

6.0 How Supplied

6.1 Packaging

The implant system, guide sheath, and stabilizer rail system are individually packaged and ethylene oxide sterilized. The table is packaged and provided non-sterile.

6.2 Storage

The PASCAL Precision system should be stored in a cool, dry place.

7.0 Directions for Use

7.1 Physician Training

The PASCAL Precision system is intended for use by interventional cardiologists and cardiac surgeons, with other support staff that are trained to assist in mitral and/or tricuspid heart valve repair.

The implanting physician shall be experienced in transcatheter techniques and trained on the PASCAL Precision system and procedure. All physicians conducting a procedure with the PASCAL Precision system must be trained per Edwards training requirements summarized below:

- PASCAL Precision system Physician Training Manual Didactic Session: device design, procedural imaging, procedural steps and challenging situations
- Hands-on Benchtop Model: practical exercise of procedural steps
- Physiologic Simulation Model: practical exercise of procedural steps with procedural imaging

The final decision for PASCAL Precision device implantation should be made by physicians specializing in treatment of mitral and/or tricuspid regurgitation in specialized centers who can

determine a reasonable chance of significant clinical improvement should be expected based on stage of disease and comorbidity.

7.2 Equipment and Materials

- Standard cardiac catheterization lab equipment
- Fluoroscopy system
- Transesophageal echocardiography (TEE) capabilities (2D and 3D)
- Venous puncture kit
- Transseptal needle, sheath, and guidewire (for mitral procedures only)
- Exchange length 0.035 inch (0.89 mm) guidewire
- Basins
- 50-60 cc syringes with Luer fitting
- Heparinized saline
- Hemostat
- Surgical towels (e.g. size 43 x 69 cm)
- Optional: pigtail catheter for contrast injection (with compatible sheath)
- Optional: step-up dilators
- Optional: continuous physiological saline drip (rolling IV pole, IV tubing with thumbwheel occluders, 1-liter bags of heparinized sterile saline solution)
- Optional: pressure monitoring device

7.3 Device Preparation

7.3.1 Table

Step	Procedure
1	Remove the table from packaging and inspect for damage.
2	Assemble the table as seen in Figure 7.

7.3.2 Stabilizer Rail System

Step	Procedure
1	Remove stabilizer rail system components from packaging and inspect for damage.

7.3.3 Guide Sheath

Step	Procedure
1	Remove the guide sheath, loader, and introducer from packaging and inspect for gross damage.
2	While keeping distal tip raised, flush and de-air guide sheath with heparinized saline.
3	While keeping the distal tip raised, insert introducer into guide sheath. Flush the introducer and wipe guide sheath with heparinized saline prior to use.

7.3.4 Implant System – System Check, Clasp Check, and Resetting

Step	Procedure
1	Remove implant system and loader from packaging and inspect for gross damage. CAUTION: If vented cap is not present on the implant catheter flush port, use of the device may result in infection.
2	Fully elongate implant. Fully retract and advance clasp sliders to confirm proper clasp motion.
3	If clasps do not move properly, follow steps below to reset. If clasps do move properly, continue to next section "Implant System – Flushing and Preparation."
4	Ensure implant is fully closed. Loosen and remove suture locks from suture lock base. Note: Ensure free end of suture is not pulled into handle while loosening suture lock.
5	Fully retract clasp sliders and place clasp setting tool flush with suture locks, suture lock bases, and implant release knob.
6	Pull free end of suture on one suture lock base to remove suture slack. Release tension on free end of suture, replace and tighten the suture lock. Repeat for second suture lock.
7	Remove clasp setting tool. Fully elongate implant. Fully advance and retract clasp sliders to confirm proper clasp motion.

7.3.5 Implant System – Flushing and Preparation

Step	Procedure
1	Close implant.
2	Ensure clasp sliders are fully retracted and implant is fully closed.
3	Remove vented cap from implant catheter flush port. Raise distal end of the implant catheter and flush with heparinized saline.
4	Attach flush port cap to implant catheter flush port.
5	Attach implant release cover to implant catheter handle.
6	Fully retract implant catheter. Advance the clasp sliders and set implant in elongated position.
7	Remove loader cap and guide the loader cap onto the implant system.
8	Insert the implant through the proximal end of the loader until it exits the distal end. Connect the loader and loader cap.
9	Advance implant catheter fully so implant exits loader.
10	While keeping loader and distal tip raised, flush heparinized saline through the steerable catheter.
11	Gradually retract implant catheter into steerable catheter and implant into the loader while continuing to flush through steerable catheter until the distal end of the implant is fully in the loader.

7.4 Implant Procedure

All implant procedure steps apply to both PASCAL mitral and PASCAL tricuspid procedures unless otherwise specified.

Delivery of the implant should be performed under general anesthesia with hemodynamic monitoring in an operating room, hybrid operating room or cath lab with fluoroscopic and echocardiographic imaging capabilities.

Note: Prior to implant procedure, refer to Anatomic Considerations (Section 3.1) as use outside of stated conditions may interfere with placement of the implant or native valve leaflet insertion.

CAUTION: During the procedure, heparin should be administered so that the ACT is maintained at ≥ 250 sec.

CAUTION: Excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

7.4.1 Patient Preparation

Step	Procedure
1	Prior to sterile draping the patient, assemble and position the table between the legs of the patient, adjusting height of the table as needed. Use towels as support between table and patient's legs. CAUTION: The table is provided non-sterile; introduction of the table into the sterile field may result in infection.
2	After sterile draping, assemble and attach stabilizer rail system as needed any time during procedure.

7.4.2 Femoral Vein Access and Sheath Introduction

Step	Procedure
1	Access the common femoral vein using conventional percutaneous puncture methods.
2	<u>For PASCAL Mitral Procedures:</u> Access the left atrium via transvenous, transseptal techniques using conventional percutaneous methods and place guidewire in left atrium. Dilate the vessel, as needed. CAUTION: (for mitral procedures only) Inappropriate puncture may result in cardiac structure damage requiring surgical repair or other intervention. <u>For PASCAL Tricuspid Procedures:</u> Access the right atrium using conventional percutaneous methods and place guidewire in right atrium. Dilate the vessel, as needed.
3	<u>For PASCAL Mitral Procedures:</u> Insert guide sheath with introducer over guidewire until guide sheath tip is securely across the septum, using flex mechanism as needed. <u>For PASCAL Tricuspid Procedures:</u> Insert guide sheath with introducer over guidewire until guide sheath tip is inside the right atrium. CAUTION: Excessive manipulation may result in dislodgement or disturbance of a previously implanted device, cardiac structure damage requiring surgical repair or other intervention.
4	Remove the introducer and guidewire. Do not aspirate and flush the guide sheath until the implant system is inserted. CAUTION: Aspiration or connection of a continuous physiological saline flush to the guide sheath prior to implant system insertion may result in air embolism.

7.4.3 Implant Navigation and Placement

Step	Procedure
1	Insert the implant system with the loader into the guide sheath.
2	Advance implant system until the implant exits the loader. Retract and peel away the loader.
3	Aspirate and flush guide sheath with heparinized saline. Utilizing the specified syringe, aspirate a minimum of 45 cc. CAUTION: Failure to fully aspirate guide sheath or aspiration without the presence of the flush port cap on the implant catheter flush port may result in air embolism.
4	If desired, connect the continuous physiological saline drip to the implant catheter. CAUTION: Connection of the continuous physiological saline drip to the implant system prior to aspiration may result in air embolism.
5	Advance implant system until the implant exits the distal end of the guide sheath.
6	Set the implant in closed position. Retract the clasp sliders.
7	Adjust guide sheath as needed.
8	At the discretion of the treating physician, if pressure monitoring is used to continuously assess atrial pressure during procedure, please follow the pressure monitor manufacturer's instruction of use. Connect a fluid filled pressure monitoring device to the steerable catheter. Aspirate and then calibrate at the patient's heart level before obtaining measurement. Refer to the PASCAL Precision system Physician Training Material for additional guidance on atrial pressure monitoring, including limitations. Note: Pressure monitoring should be used in conjunction with echo. Pressure should be reconciled against echo and Doppler readings. When assessing atrial pressure, ensure that the distal tip of the implant catheter is fully exposed from the steerable catheter.

Step	Procedure
9	Advance implant system as needed. Manipulate steerable catheter and guide sheath (flex-unflex, torque in opposing directions, advance-retract) as needed until implant is centered in the target coaptation zone with the appropriate trajectory. CAUTION: Excessive manipulation may result in dislodgement or disturbance of a previously implanted device, cardiac structure damage requiring surgical repair, or other intervention. Note: The radiopaque marker band on the steerable catheter indicates the end of the flex section and can be visualized on fluoroscopy.
10	Rotate the paddle knob to set the implant into leaflet-capture-ready position.
11	Torque the implant catheter, as needed, to orient the paddles.
12	Move one clasp slider to identify which clasp it controls via imaging. Once identified, ensure sliders are fully retracted.
13	Advance the implant through the valve until paddles are below the free edge of the leaflets.
14	Verify location and orientation of the implant and adjust position slightly as needed. CAUTION: Excessive manipulation of the implant below the leaflets may cause the implant to be entangled in the chords; chordal entanglement may result in cardiac injury, worsening regurgitation, difficulty or inability to remove the implant requiring additional intervention.
15	Under imaging guidance, retract the implant until leaflets are positioned between paddles and clasps.
16	Advance clasp slider(s) so the leaflet(s) are secured between the clasps and paddles. This can be performed for both leaflets simultaneously (clasp lock engaged to move both clasps) or each leaflet individually (clasp lock disengaged to move individual clasp).
17	Verify leaflet insertion with imaging. If leaflet(s) are not secured between clasps and paddles, retract the clasp slider(s) to release the leaflet(s) and reattempt.
18	Once leaflets are secured between the clasps and paddles, close the implant.
19	Advance implant catheter slightly to release tension on leaflets.
20	Assess regurgitation, and reposition as needed. Once the implant position is confirmed, ensure implant is closed. If repositioning within the ventricle is needed, retract the clasp sliders and set the implant in leaflet-capture-ready position. Adjust clasps and implant orientation as needed. If repositioning within the atrium is needed, retract the clasp sliders and elongate the implant slowly under fluoroscopic guidance while ensuring that the actuation wire does not bend, and retract the implant back into the atrium. CAUTION: Failure to elongate the implant when retracting into the atrium during repositioning may result in leaflet damage or chordal entanglement. CAUTION: Failure to release leaflets from clasps and paddles prior to repositioning may result in leaflet damage.

7.4.4 Implant Retrieval (if needed)

Prior to implant release, if needed, it is possible to retrieve the implant system back into the guide sheath for removal. Follow the steps below to retrieve the implant.

Refer to the PASCAL Precision system Physician Training Materials for additional considerations on implant retrieval maneuvers.

CAUTION: Excessive manipulation may result in dislodgement or disturbance of a previously implanted device, cardiac structure damage requiring surgical repair or other intervention.

Step	Procedure
1	Retract the clasp sliders.
2	Elongate the implant slowly under fluoroscopic guidance while ensuring that the actuation wire does not bend. Then retract implant into the atrium. Set the implant in closed position.
3	Unflex the steerable catheter and retract the implant system until the implant is adjacent to the tip of the guide sheath.
4	Advance the clasp sliders.
5	Set the implant in elongated position.
6	Retract the clasp sliders to open the clasps to approximately 45° on each side.
7	Retract entire implant system through the guide sheath.

7.4.5 Implant Release

To release the implant follow the steps below:

CAUTION: Failure to follow prescribed release steps may result in difficulty or inability to release implant, requiring additional intervention.

CAUTION: Releasing the implant prior to confirmation that leaflets are securely captured between paddles and clasps may result in implant movement or dislodgement leading to a single leaflet device attachment (SLDA) or other potential adverse events requiring additional intervention.

CAUTION: Re-use of the devices (including implant system and guide sheath) after retrieval may cause embolism of foreign material or infection. Device may malfunction if re-use is attempted.

Note: If an additional implant [PASCAL or PASCAL Ace] is placed at the decision of the treating physician, caution should be taken to avoid dislodgement of the previously placed implant. Crossing the valve in a low profile implant configuration can minimize interaction with the previously placed implant.

CAUTION: Excessive manipulation may result in dislodgement or disturbance of a previously implanted device, cardiac structure damage requiring surgical repair or other intervention.

Step	Procedure
1	Ensure that the distal tip of the implant catheter is fully exposed from the steerable catheter.
2	Unscrew and remove implant release cover from implant catheter handle.
3	Unthread and remove one suture lock from the suture lock base.
4	Pull suture lock away from handle to fully remove suture.
5	Repeat steps for other suture lock.
6	Rotate counterclockwise and retract the implant release knob until the implant is released, as confirmed via imaging.
7	Replace suture locks, as needed.

7.4.6 Device Removal and Closure

Step	Procedure
1	Retract implant catheter completely into steerable catheter. Gradually unflex and remove implant system. Gradually unflex and remove guide sheath. CAUTION: Failure to unflex devices prior to removal may result in vessel damage.
2	Perform standard percutaneous closure of access site.

8.0 Magnetic Resonance (MR) Safety

Non-clinical testing has demonstrated that the PASCAL and PASCAL Ace implants are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic fields of 1.5 T and 3.0 T
- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode).

Under the scan conditions defined above, the implant is expected to produce a maximum temperature rise of less than 4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device in a worst case multiple implant configuration extends up to 15 mm from the implant when imaged in the worst case gradient echo pulse sequence in a 3.0 T MRI system.

9.0 Recovered Implant and Device Disposal

Edwards Lifesciences is interested in obtaining recovered clinical specimens of the implant for analysis. A written report summarizing our findings will be provided upon completion of our evaluation. Please contact Edwards for return of the recovered implant.

If you do decide to return any of the devices, please follow the following instructions:

- **Unopened Package with Sterile Barrier Intact:**
If the pouches have not been opened, return the device in its original packaging.
- **Package Opened but Not Implanted:**
If a pouch is opened, the device is no longer sterile. Please return the device in its original packaging.
- **Explanted Implant:**
The explanted implant should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to Edwards.

9.1 Disposal

Use universal precautions for biohazards and sharps to avoid user injury. Used devices (includes all devices that come in contact with patients) should be handled and disposed of in accordance with institutional guidelines on biohazardous materials and hospital waste to avoid possible cross-contamination.

10.0 Summary of Clinical Experience

A high-level summary of all clinical experience is provided below. A detailed summary of all clinical experience can be found in the SSCP.

10.1 CLASP Study

Results from the CLASP study, a multi-center, multi-national, prospective, single-arm study, show high rates of device, procedural, and clinical success, a reduction in MR grade through 2 years, an improvement in quality of life, exercise capacity, and functional status through 1 year, and an acceptable MAE rate.

10.2 CLASP IID/IIF (IID Roll-In Cohort) Study

Results from the DMR roll-in cohort of the pivotal, randomized, controlled CLASP IID/IIF study, show an improvement in MR grade and clinical outcomes through 30 days along with an MAE rate within the anticipated range.

10.3 MiCLASP PMCF Study

Results from the MiCLASP PMCF study, a multi-center, single-arm, prospective post market clinical follow up study, show a reduction in MR grade and an improvement in functional status, exercise capacity, and quality of life through 1 year. Evaluation of the MAE rate at 1 year indicates an acceptable safety profile.

10.4 CLASP TR Study

Results from the CLASP TR study, a prospective, single arm, multi-center study, shows a high rate of device, procedural, and clinical success, a reduction in TR grade through 6 months, and an improvement in clinical outcomes through 1 year. The 1-year MAE rate indicates an acceptable safety profile.

10.5 Additional Studies

The following PASCAL implant studies (MR and TR) including TriCLASP PMCF study (Tricuspid), PASCAL implant post-market Registry, CLASP ITR (Tricuspid), and CLASP IID/IIF (Mitral) Registry and Randomized Cohorts are currently ongoing and have yet to reach their primary endpoints; therefore, results are not presented here.

11.0 Summary of Safety and Clinical Performance (SSCP)

The SSCP (Summary of Safety and Clinical Performance) has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which the CE certification has been granted. The SSCP contains a relevant summary of the same information.

The Notified Body has taken notice of and agreed with the benefit-risk rationales for the short- and long-term safety and effectiveness of the PASCAL Precision system.

Conformity of the PASCAL Precision system of the Performance Requirements (GSPR) for safety (MDR GSPR 1), performance (MDR GSPR 1), acceptability of side-effects (MDR GSPR 8), usability (MDR GSPR 5), device lifetime (MDR GSPR 6), acceptable benefit-risk profile (MDR GSPR 8) has been established for the labelled indications.

Refer to <https://meddeviceinfo.edwards.com/> for a SSCP for this medical device.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

12.0 Clinical Benefits

The clinical benefits of the PASCAL Precision system for the treatment of MR include the following:

- Effective and stable reduction in mitral regurgitation.
- Enables minimally invasive percutaneous treatment option for mitral regurgitation.
- Improvement in functional status, exercise capacity, and quality of life.

The clinical benefits of the PASCAL Precision system for the treatment of TR include the following:

- Effective and stable reduction in tricuspid regurgitation.
- Enables minimally invasive percutaneous treatment option for tricuspid regurgitation.
- Improvement in functional status, exercise capacity, and quality of life.

There are no clinical benefits specific to the table, as the table is an optional non-patient contacting accessory. Benefits of the table are functional in nature and related to the intended use of the accessory to support the PASCAL Precision system.

There are no clinical benefits specific to the stabilizer rail system, as the stabilizer rail system is an optional non-patient contacting accessory used only during the implant procedure. Benefits of the stabilizer rail system are functional in nature and related to the intended use of the accessory to support the PASCAL Precision system.

13.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed. The Basic UDI-DI for the PASCAL Precision system can be used to locate the SSCP.

The following table contains the Basic UDI-DIs for the PASCAL Precision system and compatible devices:

Product	Model	Basic UDI-DI
PASCAL Precision system – implant system	20000IS	0690103S004PAS000BC
PASCAL Precision system – PASCAL Ace implant system	20000ISM	0690103S004PAS000BC
PASCAL Precision system – guide sheath	20000GS	0690103S004PAS000BC
PASCAL system – stabilizer rail system	20000ST	0690103D004PAC000S6
PASCAL system – table	10000T	0690103D004PAC000S6

14.0 Expected Lifetime of the Device

The PASCAL and PASCAL Ace implants are subject to rigorous pre-clinical durability testing per the testing requirements and were successfully tested to a minimum of 5 years of simulated wear. The actual lifetime performance in humans depends on multiple biological factors and varies widely from patient to patient. Specific activities or conditions that could shorten or lengthen the device lifetime have not been established.

15.0 Patient Information

A patient implant card is provided with each implant system. After implantation, please complete all requested information and provide the implant card to the patient. This implant card allows patients to inform healthcare providers what type of implant they have when they seek care.

16.0 Performance Characteristics

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The design of the PASCAL Precision delivery catheters is based on user feedback. The claim for "Precision" will be confirmed during post-market monitoring.

17.0 Qualitative and Quantitative Information related to the PASCAL and PASCAL Ace Implants

The PASCAL implant comprises a titanium nut and bolt, PEEK bushing, and a silicone seal. The PASCAL Ace implant comprises a titanium nut, bolt, distal and proximal plate, and a silicone seal.

The following table shows the qualitative and quantitative information on the materials and substances:

Substance	CAS	Model Mass Range (mg)
Titanium	7440-32-6	254 - 324
Nickel	7440-02-0	235 - 258
Polyethylene terephthalate	25038-59-9	82.5 - 98.1
Polyethylene	9002-88-4	14.3 - 42.6
Polyether ether ketone	29658-26-2	0 - 23.2
Aluminium	7429-90-5	3.39 - 9.13
Vanadium	7440-62-2	2.16 - 6.32
Silicon dioxide	7631-86-9	5.08 - 5.40
Polydimethylsiloxane	63148-62-9	3.86 - 4.16
Perfluoropolyether	69991-67-9	2.44 - 2.52
Polytetrafluoroethylene	9002-84-0	1.15 - 1.22
Iron	7439-89-6	0 - 0.567
Titanium dioxide	13463-67-7	0.180 - 0.541
Oxygen	7782-44-7	0 - 0.355
Carbon	7440-44-0	0 - 0.285
Cobalt	7440-48-4	0 - 0.226
Antimony trioxide	1309-64-4	0.0847 - 0.118
Niobium	2023505	0 - 0.113
Nitrogen	7727-37-9	0 - 0.0918
Chromium	7440-47-3	0 - 0.0452
Copper	7440-50-8	0 - 0.0452
Hydrogen	1333-74-0	0 - 0.0384
4,4'-Difluorobenzophenone	345-92-6	0 - 0.00141
Diphenyl sulfone	127-63-9	0 - 0.00114
Erucamide	112-84-5	0.000516 - 0.00102
4-Dodecylbenzenesulfonic acid	121-65-3	0.0000906 - 0.00857
Decamethylcyclopentasiloxane; D5	541-02-6	0 - 0.000698
Dodecamethylcyclohexasiloxane; D6	540-97-6	0 - 0.000698
Octamethylcyclotetrasiloxane; D4	556-67-2	0 - 0.000651

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








	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk
	Catalogue Number	Référence catalogue	Artikelnummer	Número de catálogo	Numero di catalogo	Catalogusnummer	Katalognummer
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer	Partinummer
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello	Modelnummer	Modelnummer
	Contents	Contenu	Inhalt	Contenido	Contenuti	Inhoud	Indhold
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte	Anvendelig længde
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken	Må ikke genanvendes
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op	Forsigtig
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt
	Keep away from sunlight	Tenir à l'abri de la lumière du soleil	Vor Sonneneinstrahlung schützen	Mantener alejado de la luz del sol	Tenere lontano dalla luce solare	Uit het zonlicht houden	Beskyttes mod sollys
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto	Droog houden	Opbevares tørt
	Importer	Importateur	Importeur	Importador	Importatore	Importeur	Importør
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo	Unieke instrumentidentificatie	Unik udstyrsidentifikation
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid	Mængde
	Date of implantation	Date de l'implantation	Implantationsdatum	Fecha de implantación	Data di impianto	Datum van implantatie	Implantationsdato
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarriersystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem
	Consult instructions for use on the website	Consulter le mode d'emploi sur le site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola


	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling
	Sterilized using liquid chemical	Stérilisé à l'aide d'un produit chimique liquide	Mit flüssigen Chemikalien sterilisiert	Esterilizado con sustancia química líquida	Sterilizzato con sostanze chimiche liquide	Gesteriliseerd met vloeibare chemische stof	Steriliseret med flydende kemikalie
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare	Niet opnieuw steriliseren	Må ikke resteriliseres
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile	Niet-steriel	Ikke-steril
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen	Ikke-pyrogen
	Non-DEHP	Sans DEHP	Nicht DEHP-haltig	Sin-DEHP	Non contiene DEHP	Niet-DEHP	Ikke-DEHP
	Use-by date	Date d'expiration	Verfallsdatum	Fecha de caducidad	Data di scadenza	Vervaldatum	Sidste anvendelsesdato
	Serial Number	Numéro de série	Serienummer	Número de serie	Numero di serie	Serienummer	Serienummer
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant	Producent
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum	Fremstillingsdato
	Patient information website	Site web d'information du patient	Website mit Patienteninformationen	Sitio web de información para el paciente	Sito Web con le informazioni per i pazienti	Patiënteninformatiewebsite	Websted med patientinformation
	Patient Name	Nom du patient	Patientenname	Nombre del paciente	Nome del paziente	Patiëntnaam	Patientens navn
	For use with size 44 mm Edwards transcatheter heart valve	À utiliser avec une valve cardiaque transcathéter Edwards de 44 mm	Zur Verwendung mit einer 44-mm-Transkatheter-Herzklappe von Edwards	Para su uso con una válvula cardíaca transcatéter Edwards de 44 mm	Da utilizzare con valvole cardiache transcatetere Edwards da 44 mm	Voor gebruik met de Edwards transcatheter-hartklep van 44 mm	Til brug med størrelse 44 mm Edwards transkateterklap
	For use with size 48 mm Edwards transcatheter heart valve	À utiliser avec une valve cardiaque transcathéter Edwards de 48 mm	Zur Verwendung mit einer 48-mm-Transkatheter-Herzklappe von Edwards	Para su uso con una válvula cardíaca transcatéter Edwards de 48 mm	Da utilizzare con valvole cardiache transcatetere Edwards da 48 mm	Voor gebruik met de Edwards transcatheter-hartklep van 48 mm	Til brug med størrelse 48 mm Edwards transkateterklap
	Temperature limit	Limite de température	Temperaturgrenze	Limite de temperatura	Limiti di temperatura	Temperatuurlimiet	Temperaturgrænse
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter	Udvendig diameter
	Inner diameter	Diamètre interne	Innendurchmesser	Diámetro interior	Diametro interno	Binnendiameter	Indre diameter
	Recommended guidewire length	Longueur recommandée du fil-guide	Empfohlene Führungsdrahtlänge	Longitud recomendada de la guía	Lunghezza consigliata per il filo guida	Aanbevolen voerdraadlengte	Anbefalet guidewirelængde

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola

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
	Norsk	Suomi	Български	Română	Eesti	Lietuvių	Latviešu
	Medisinsk utstyr	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiineseade	Medicinos priemonė	Mediciniska ierice
	Enkelt, sterilt barrieresystem	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema	Viena sterilitātes aizsargslāpa sistēma
	Se bruksanvisningen på nettstedet	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба в уебсайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege kasutusjuhiseid veebisaidilt	Žr. naudojimo instrukcijas interneto svetainėje	Skatīt lietošanas instrukciju tīmekļa vietnē
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Sterilisert med etylenoksid	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksīdu
	Sterilisert med stråling	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiirgamist kasutades	Sterilizuota švitinant	Sterilizēts apstarojot
	Sterilisert med flytende kjemikalie	Steriloitu nestemäisellä kemikaalilla	Стерилизирано с течни химикали	Sterilizat cu substanțe chimice lichide	Steriliseeritud vedelat kemikaali kasutades	Sterilizuota skystomis cheminėmis medžiagomis	Sterilizēts, izmantojot šķidru ķīmisko vielu
	Må ikke resteriliseres	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti
	Usteril	Epästeriili	Нестерилно	Nesteril	Mittesteriilne	Nesterilus	Nesterils
	Ikke-pyrogen	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas	Nepirogēna
	Ikke-DEHP	Ei sisällä DEHP:tä	Не съдържа DEHP	Nu conține DEHP	Mitte-DEHP	Ne DEHP	DEHP nesaturošs
	Utløpsdato	Viimeinen käyttöpäivämäärä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Tinkamumo terminas	Derīguma termiņš
	Serienummer	Sarjanumero	Сериен номер	Număr de serie	Seerianumber	Serijos numeris	Sērijas numurs
	Autorisert representant i Det europeiske fellesskap / EU	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană / Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjunga	Pilnvarotais pārstāvis Eiropas Kopienā / Eiropas Savienībā
	Produsent	Valmistaja	Производител	Producător	Tootja	Gamintojas	Ražotājs
	Produksjonsdato	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data	Izgatavošanas datums
	Nettsted for pasientinformasjon	Potilastietojen verkkosivusto	Уеб сайт с информация за пациента	Site web cu informații pentru pacient	Patsiendiandmete veebileht	Informacijos pacientui interneto svetainė	Pacienta informācijas tīmekļa vietne
	Navn på pasient	Potilaan nimi	Име на пациента	Numele pacientului	Patsiendi nimi	Paciento vardas, pavardė	Pacienta vārds un uzvārds
	Til bruk med Edwards transkateter hjerteklaff, størrelse 44 mm	Käytetään 44 mm:n kokoisten katetrilla asennettavien Edwards - sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа на Edwards с размер 44 mm	A se utiliza cu valvă cardiacă transcateter Edwards de mărimea 44 mm	Kasutamiseks ettevõtte Edwards 44 mm transkateetriga südameklapiga	Skirta naudoti su 44 mm dydžio transkateriniu širdies vožtuvu „Edwards“	Lietošanai ar 44 mm izmēra Edwards transkatetra sirds vārstuli
	Til bruk med Edwards transkateter hjerteklaff, størrelse 48 mm	Käytetään 48 mm:n kokoisten katetrilla asennettavien Edwards - sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа на Edwards с размер 48 mm	A se utiliza cu valvă cardiacă transcateter Edwards de mărimea 48 mm	Kasutamiseks ettevõtte Edwards 48 mm transkateetriga südameklapiga	Skirta naudoti su 48 mm dydžio transkateriniu širdies vožtuvu „Edwards“	Lietošanai ar 48 mm izmēra Edwards transkatetra sirds vārstuli

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola

	Norsk	Suomi	Български	Română	Eesti	Lietuvių	Latviešu
	Temperaturgrense	Lämpötilaraja	Температурно ограничение	Limită de temperatură	Temperatuuriipiirang	Temperatūros apribojimas	Temperatūras ierobežojums
	Ytre diameter	Ulkhalkaisija	Външен диаметър	Diametru exterior	Välisläbimõõt	Išorinis skersmuo	Ārējais diametrs
	Innvendig diameter	Sisäläpimitta	Вътрешен диаметър	Diametru interior	Sisemine läbimõõt	Vidinis skersmuo	Iekšējais diametrs
	Anbefalt ledevaierlengde	Suosittelu ohjainlangan pituus	Препоръчителна дължина на теления водач	Lungimea recomandată a firului de ghidaj	Soovituslik juhtetraadi pikkus	Rekomenduojamas kreipiamosios vielos ilgis	Ieteicamais vadītājstīgas garums
GW	Anbefalt ledevaierstørrelse	Ohjainlangan suositeltu koko	Препоръчителен размер телен водач	Mărimea recomandată a firului de ghidaj	Juhtetraadi soovituslik suurus	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs
GWC	Ledevaierkompatibilitet	Ohjainlangan yhteensopivuus	Съвместимост на теления водач	Compatibilitatea firului de ghidaj	Juhtetraadi ühilduvus	Kreipiamosios vielos suderinamumas	Vadītājstīgas saderība
SZ	Størrelse	Koko	Размер	Mărime	Suurus	Dydis	Izmērs
	Kateterskaftets størrelse	Katetrin varren koko	Размер на шифта на катетъра	Mărime a tijei cateterului	Kateetri varre suurus	Kateterio vamzdelio dydis	Katetra ass izmērs
	Ballongdiameter	Pallon halkaisija	Диаметър на балона	Diametru balon	Ballooni läbimõõt	Balionėlio skersmuo	Balona diametrs
	Ballongens arbelsslengde	Pallon käyttöpituu	Работна дължина на балона	Lungimea de lucru a balonului	Ballooni tööpikkus	Balionėlio darbinis ilgis	Balona darba garums
	[Kun for implantat] Implantatenheten har blitt fastslått å være MR-sikker under spesifiserte forhold når den brukes under forholdene som er angitt i bruksanvisningen.	[Vain implantti] Implanttilaitteen katsottiin olevan ehdollisesti turvallinen magneettikuvaussessa, kun sitä käytetään ohjeissa lueteltujen ehtojen mukaisesti.	[Само за импланта] За имплантираното изделие е определено, че е безопасно при MR при определени условия, когато се използва при условията, посочени в инструкциите.	[Numai pentru implanturi] S-a stabilit că dispozitivul implantabil este condiționat RM atunci când este utilizat în condițiile enumerate în instrucțiunile.	[Ainult implantaat] Implantaatseade on määratud ohutuks magnetresonantsomograafias teatud tingimustel, kui seda kasutatakse juhendis nimetatud tingimustel.	[Tik implantas] Buvo nustatyta, kad implantu prietaisas yra santykinai saugus MR aplinkoje, kai naudojamas pagal šioje naudojimo instrukcijoje nurodytas sąlygas.	[Tikai implantam] Tika noteikts, ka implantēšanas ierīce ir izmantojama MR vidē, ievērojot noteiktus nosacījumus, ja tā tiek izmantota saskaņā ar instrukcijā minētajiem nosacījumiem.
	MR-usikker	Ei sovi magneettikuvausseen	Небезопасно при MR	Incompatibil RM	Ohtlik magnetresonants-tomograafias	MR nesaugus	Nedrīkst lietot MR vidē
	Navn og adresse til implanterende helseinstitusjon/ helsepersonell	Implantoivan terveydenhuoltolaitoksen / terveyspalvelujen tarjoajan nimi ja osoite	Име и адрес на здравното заведение/ доставчика, извършващи имплантирането	Numele și adresa instituției/furnizorului de servicii medicale pentru implanturi	Implanteerimist teostava tervishoiuasutuse/ teenuse osutaja nimi ja aadress	Implantuojančios sveikatos priežiūros įstaigos / teikėjo pavadinimas ir adresas	Veselības aprūpes iestādes, kas veic implantēšanu / implantēšanas nodrošinātāja nosaukums un adrese

	Türkçe	Русский	Hrvatski
REF	Katalog Numarası	Номер по каталогу	Kataloški broj
LOT	Lot Numarası	Номер партии	Broj serije
#	Model Numarası	Номер модели	Broj modela
	İçindekiler	Содержимое	Sadržaj
	Kullanılabilir uzunluk	Рабочая длина	Upotrebjliva dužina
	Yeniden kullanmayın	Не использовать повторно	Nemojte ponovno upotrebljavati

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символах ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola

	Türkçe	Русский	Hrvatski
	Dikkat	Предостережение	Oprez
	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Pohranite na hladnom, suhom mjestu
	Güneş ışığından uzak tutun	Беречь от солнечных лучей	Držite dalje od sunčeve svjetlosti
	Kuru halde tutun	Беречь от влаги	Držati suhim
	İthalatçı	Импортёр	Uvoznik
UDI	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator proizvoda
QTY	Miktar	Количество	Količina
31	İmplantasyon tarihi	Дата имплантации	Datum ugradnje
MD	Tıbbi cihaz	Медицинское устройство	Medicinski uređaj
	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sustav jedne sterilne barijere
	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte upute za upotrebu na internetskoj stranici
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
STERILE EO	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizirano etilen-oksidadom
STERILE R	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizirano zračenjem
STERILE LC	Sıvı kimyasal kullanılarak sterilize edilmiştir	Стерилизовано с использованием жидкого химического вещества	Sterilizirano tekućom kemikalijom
	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovno sterilizirati
	Steril değildir	Нестерильно	Nesterilno
	Nonpirogeniktir	Апирогенно	Nepirogeno

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda simbola

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12/2022
10047576002 A
© Copyright 2022, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences Services GmbH