Edwards

Edwards EVOQUE Tricuspid Valve Replacement System

Directory	
English (en).....	1
Français (fr).....	12
Deutsch (de).....	23
Español (es).....	34
Italiano (it).....	45
Svenska (sv).....	56
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Figurer.....	66
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Symbolförklaring.....	68

English

Edwards EVOQUE Tricuspid Valve Replacement System

Instructions for Use

Caution: Implantation of the transcatheter tricuspid valve should only be performed by physicians who have received training for the Edwards EVOQUE tricuspid valve replacement system.

- Do not attempt to use the Edwards EVOQUE valve (herein referred to as the EVOQUE valve), delivery system, or associated accessories before completely reading and understanding the information contained in this booklet.
- Failure to follow these instructions, warnings, and precautions may lead to device damage or patient injury. Use of the EVOQUE system should be restricted to those physicians trained to perform invasive endovascular procedures and to those physicians trained in the proper use of the system.
- Consult with authorized Edwards personnel for appropriate selection of EVOQUE valve size.
- The EVOQUE valve, delivery system, dilator kit, loading system, and stabilizer are supplied STERILE. The base and plate are supplied NON-STERILE.
- All devices are supplied for single use only. After use, dispose in accordance with hospital administrative and/or government policy. Do not re-sterilize.

1.0 Introduction

The Edwards EVOQUE tricuspid valve replacement system (herein referred to as the EVOQUE system) is designed to replace the native tricuspid valve in patients with symptomatic, severe tricuspid valve regurgitation without the need for conventional open-heart surgery. The EVOQUE system consists of four (4) elements and includes three (3) optional accessories, as outlined below:

Product Name	44 mm	48 mm	52 mm
Edwards EVOQUE Valve (EVOQUE Valve)	9850EV44	9850EV48	9850EV52
Edwards EVOQUE Tricuspid Delivery System	9850TDS		
Edwards EVOQUE Dilator Kit	9850DK		
Edwards EVOQUE Loading System	9850LS		

Optional Accessories	Model Number
Edwards EVOQUE Stabilizer	9850SB
Edwards EVOQUE Stabilizer Base	9850BA
Edwards EVOQUE Stabilizer Plate	9850PT

Note: The devices that comprise the EVOQUE system are intended to be used together and are not intended to be used as standalone devices. The information in this IFU is applicable to all the EVOQUE system devices.

2.0 Device Description

Edwards EVOQUE Valve (Figure 1)

The EVOQUE valve is intended to be implanted in the tricuspid position as an alternative to surgery. The EVOQUE valve consists of a trileaflet bovine pericardial tissue valve, nitinol frame, and fabric skirt, and is packaged and terminally sterilized in glutaraldehyde.

Valve size recommendations are based on native valve annulus size, as measured by computed tomography (CT). Patient anatomical factors and imaging modalities should be considered during valve size selection.

Device Diameter (Recommended Valve Size)	Systole		Diastole	
	Recommended Treatable Perimeter - Derived Diameter Range (mm)	Maximum Treatable Annulus Length (mm)	Recommended Treatable Perimeter - Derived Diameter Range (mm)	Maximum Treatable Annulus Length (mm)
44	36.5 - 43	45.5	39.6 - 45.5	50
48	40 - 47	49.5	43.2 - 49.5	54
52	45 - 51	53.5	46.8 - 53.5	58

Edwards EVOQUE Tricuspid Delivery System (Figure 2)

The delivery system has an outer diameter of 28 F and is intended to deliver the EVOQUE valve in the crimped position via the transfemoral venous approach. The delivery system handle contains a primary flex knob, secondary flex knob, and depth knob to facilitate EVOQUE valve alignment and positioning in the native valve, and a capsule knob and release knob to control the expansion and release of the EVOQUE valve.

- **Edwards EVOQUE Dilator Kit (Figure 3)**

The 24 F, 28 F, and 33 F diameter hydrophilic-coated dilators are intended to dilate the access site, facilitating delivery system insertion. All dilators accommodate a 0.035 inch (0.89 mm) guidewire and are tapered to minimize access site trauma.

- **Edwards EVOQUE Loading System (Figure 4)**

The loading system, which consists of multiple components, is intended to facilitate loading and attachment of the EVOQUE valve onto the delivery system. The loading system assists in crimping the EVOQUE valve to the appropriate diameter, which allows the outer capsule to advance over the EVOQUE valve.

- **Edwards EVOQUE System Optional Accessories: Stabilizer, Base, and Plate (Figure 5)**

The stabilizer, base, and plate are intended to secure the delivery system at an angle appropriate for the transfemoral venous approach and to enable fine adjustments of the position of the delivery system during the implantation procedure. The base is height-adjustable to accommodate patient lower extremities and is intended to provide a stable platform for the stabilizer. The plate is intended to provide a stable, flat surface for the base on the operating table.

3.0 Facilities and Training Requirement

Facilities that intend to perform an implantation procedure utilizing the EVOQUE system must have access to cine fluoroscopy and transesophageal echocardiography (TEE) throughout the procedure. In addition, the implanting physicians must have prompt access to facilities with the necessary equipment, instruments, supplies, and personnel to perform emergency tricuspid valve surgery, if required.

A comprehensive training program is provided by Edwards Lifesciences and must be completed by the implanting physicians before use of the EVOQUE system. The implanting physicians should have advanced technical knowledge and experience in related catheter-based procedures.

4.0 Indications for Use

4.1 Intended Use

The EVOQUE tricuspid valve replacement system is intended for the reduction of tricuspid regurgitation (TR) for patients who remain symptomatic on medical therapy and in whom no other surgical or transcatheter treatment option exists per heart team decision. The EVOQUE system percutaneously delivers the Implant to the tricuspid valve via femoral vein access using a transvenous approach.

4.2 Target Patient Population

- Patients with symptomatic, severe TR who do not have any of the following:
 - Severe aortic, mitral and/or pulmonic valve stenosis and/or regurgitation
 - LVEF < 25%
 - Severe right ventricular dysfunction as measured by echocardiography
 - Pacemaker dependency¹ on trans-tricuspid lead without alternative pacing option
 - Need for emergent or urgent surgery for any reason or any planned cardiac surgery within the next 12 months
 - Resting systolic blood pressure < 90 or > 160 mmHg after repeated measurements
 - Refractory heart failure requiring advanced intervention (i.e., left ventricular assist device, transplantation) (ACC/AHA/ESC/EACTS Stage D heart failure)
 - An occluded or thrombosed IVC filter that would interfere with the delivery catheter
- Patient selection should be performed by a multi-disciplinary heart team specializing in the treatment of tricuspid regurgitation to assess patient risk and anatomical suitability.

¹ Pacemaker dependency may increase the risk of serious injury or death from sudden pacemaker failure

5.0 Contraindications

The EVOQUE valve is contraindicated for patients with any of the following conditions:

- Active endocarditis within 3 months or infection requiring antibiotic therapy (oral or intravenous) within 2 weeks of the scheduled implant.
- Untreatable hypersensitivity or contraindication to any of the following: all antiplatelets, all anticoagulants, nitinol alloys (nickel and titanium), bovine tissue, glutaraldehyde, contrast media, or transesophageal echocardiography.
- Tricuspid valve anatomy that precludes proper device deployment and functionality based on CT and echocardiographic evaluation.

Note: Patient must be able to tolerate at least one antiplatelet medication AND one anticoagulant medication.

6.0 Warnings

- Do not use the EVOQUE valve outside of the labeled indication.
- The EVOQUE valve, delivery system, loading system, dilator kit, and stabilizer are designed, intended, and distributed for STERILE single use only. The base and plate are provided nonsterile for single use only. Do not resterilize or reuse any of the devices. There are no data to support the sterility, nonpyrogenicity, or functionality of the devices after reprocessing.
- Ensure proper sterile techniques are utilized during the preparation, transfer, and use of the devices.
- Do not use the valve if the tamper evident seal is broken, the storage solution does not completely cover the valve, the temperature indicator has been activated, the valve is damaged, or the expiration date has elapsed. The EVOQUE valve must remain hydrated at all times. The valve cannot be exposed to solutions, antibiotics, or chemicals other than its shipping storage solution and sterile physiologic saline solution. This will help prevent leaflet damage that may impact valve functionality. Keep the EVOQUE valve hydrated with normal saline until ready for implantation.
- Ensure the correct valve size is selected. Implantation of the improper size (i.e., undersizing or oversizing) may lead to paravalvular leak (PVL), migration, embolization, and/or annular damage.
- Patients with previously-implanted devices (e.g., IVC filter) should be carefully assessed prior to insertion of the delivery system to avoid potential damage to vasculature or a previously-implanted device.
- Patients with pre-existing cardiac leads should be carefully assessed prior to implantation to avoid potential adverse interaction between devices.
- Care should be taken when implanting cardiac leads after EVOQUE valve implantation to avoid potential adverse interaction between the devices.
- Patients implanted with the EVOQUE valve should maintain anticoagulant/antiplatelet therapy as determined by their physicians in accordance with current guidelines, to minimize the risk of valve thrombosis or thromboembolic events.
- There are no data to support device safety and performance if the patient has:
 - Echocardiographic evidence of severe right ventricular dysfunction;
 - Pulmonary arterial systolic pressure (PASP) > 70 mmHg by echo Doppler;
 - A trans-tricuspid pacemaker or defibrillator lead that has been implanted in the RV within the last 3 months;
 - Dependency on a trans-tricuspid pacemaker without alternative pacing options.

7.0 Precautions

7.1 Precautions Prior to Use

- The patient's eligibility depends on the anatomic conditions based on CT scan.
- It is advised that a multi-disciplinary heart team be of the opinion that EVOQUE valve implantation is preferable to alternative percutaneous device solutions, including minimally-invasive open heart surgery.
- It is advised that a multi-disciplinary heart team takes into consideration the severity of disease and the chances of reversibility of right heart failure based on a complete hemodynamic assessment.
- The patient should have a life expectancy of at least 12 months.

7.2 Precautions

- The EVOQUE valve is to be used only with the 9850TDS delivery system and 9850LS loading system.
- The procedure should be conducted under appropriate imaging modalities, such as transesophageal echocardiography (TEE), fluoroscopy, and/or intracardiac echocardiography (ICE).
- Glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Safety Data Sheet available from Edwards Lifesciences.
- Conduction disturbances may occur before, during, or following implantation of the EVOQUE valve, which may require continuous ECG monitoring before hospital discharge. If a patient has confirmed or suspected conduction disturbances, consider patient monitoring and/or electrophysiology evaluation.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Long-term durability has not been established for the EVOQUE valve. Clinical data are reflective of short-term follow-up. Regular medical follow-up is advised to evaluate EVOQUE valve performance.
- Implantation of the EVOQUE valve should be postponed in patients with (1) a history of myocardial infarction within one month (30 days) of planned intervention, (2) pulmonary emboli within 3 months (90 days) of planned intervention, (3) cerebrovascular accident (stroke or transient ischemic attack [TIA]) within 3 months (90 days) of planned intervention, (4) active upper gastrointestinal bleeding within 3 months (90 days) prior to procedure requiring transfusion.

8.0 Potential Adverse Events

Potential adverse events related to standard cardiac catheterization, use of anesthesia, the EVOQUE valve, and the implantation procedure include:

- Abnormal lab values
- Allergic reaction to anesthesia, contrast media, anti-coagulation medication, or device materials
- Anaphylactic shock

- Anemia or decreased hemoglobin (Hgb), may require transfusion
- Aneurysm or pseudoaneurysm
- Angina or chest pain
- Arrhythmia – atrial (i.e., atrial fibrillation, supraventricular tachycardia)
- Arrhythmias – ventricular (i.e., ventricular tachycardia, ventricular fibrillation)
- Arterio-venous fistula
- Bleeding
- Cardiac arrest
- Cardiac (heart) failure
- Cardiac injury, including perforation
- Cardiac tamponade / pericardial effusion
- Cardiogenic shock
- Chordal entanglement or rupture that may require intervention
- Coagulopathy, coagulation disorder, bleeding diathesis
- Conduction system injury, which may require implantation of a pacemaker (temporary or permanent)
- Conversion to open heart surgery
- Coronary artery occlusion
- Damage to or interference with function of pacemaker or implantable cardioverter defibrillator (ICD)
- Death
- Edema
- Electrolyte imbalance
- Embolization including air, particulate, calcific material, or thrombus
- Emergent cardiac surgery
- Endocarditis
- Esophageal irritation
- Esophageal perforation or stricture
- EVOQUE system component(s) embolization
- Failure to retrieve any EVOQUE system components
- Fever
- Gastrointestinal bleeding
- Hematoma
- Hemodynamic compromise
- Hemolysis / hemolytic anemia
- Hemorrhage requiring transfusion/surgery
- Hypertension
- Hypotension
- Inflammation
- Injury to the tricuspid apparatus including chordal damage, rupture, papillary muscle damage
- Local and systemic infection
- Mesenteric ischemia or bowel infarction
- Multi-system organ failure
- Myocardial infarction
- Nausea and/or vomiting
- Nerve injury
- Neurological symptoms, including dyskinesia, without diagnosis of TIA or stroke
- Non-emergent reoperation
- Pain
- Pannus formation
- Paralysis
- Percutaneous valve intervention
- Peripheral ischemia
- Permanent disability
- Pleural effusion
- Pneumonia
- Pulmonary edema
- Pulmonary embolism
- Reaction to anti-platelet or anticoagulation agents
- Rehospitalization
- Renal failure
- Respiratory failure, atelectasis – may require prolonged intubation
- Retroperitoneal bleed
- Right ventricular outflow tract (RVOT) obstruction
- Septicemia, sepsis
- Skin burn, injury, or tissue changes due to exposure to ionizing radiation
- Stroke
- Structural deterioration (wear, fracture, calcification, leaflet tear, leaflet thickening, stenosis of implanted device, or new leaflet motion disorder)
- Thromboembolism
- Transient ischemic attack (TIA)
- Valve dislodgement/embolization
- Valve endocarditis
- Valve explant
- Valve leaflet entrapment
- Valve malposition
- Valve migration
- Valve paravalvular leak (PVL)
- Valve regurgitation (new or worsening tricuspid, aortic, mitral, pulmonary)
- Valve thrombosis
- Vascular injury or trauma, including dissection or occlusion
- Vessel spasm
- Wound dehiscence, delayed or incomplete healing

For a patient/user/third party in the European Economic area: if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

9.0 Additional Equipment

The implantation procedure requires additional equipment that is not supplied with the EVOQUE system. The additional equipment is provided below.

9.1 Equipment for EVOQUE Valve Loading

Note: Volumes reflect adequate amount for preparing 1 implant.

- 3500 ml (minimum) sterile physiological saline solution at ambient temperature (~ 23 °C)

- 500 ml (minimum) heparinized saline solution (2 units/ml) at ambient temperature (~ 23 °C)
- 4 sterile bowls (≥ 500 ml, ≥ 7 cm depth, plastic)
- 1 large sterile bowl (≥ 2 L, ≥ 10 cm depth, plastic)
- 1 scalpel, no. 11 scalpel blade
- 1 luer locking syringe (≥ 20 cc)
- Gauze pads
- Blunt tip forceps
- Blunt tip scissors
- Sterile towels

9.2 Equipment for Access, Procedure, and Monitoring

- Standard cardiac catheterization lab equipment
- Femoral vessel introducer sheath
- Fluoroscopy (fixed, mobile, or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal echocardiography capabilities
- Steerable introducer sheath
- Exchange length 0.035 inch (0.89 mm) max guidewire
- Extra small curve 0.035 inch (0.89 mm) max guidewire
- Right coronary arterial catheter and guidewire
- Sterile table for EVOQUE valve and device preparation

9.3 Standby Equipment

- Arterial bypass cannula (~18 F)
- Cardiopulmonary bypass machine
- Compliant balloon (> 20 mm diameter, 9 cc contrast volume)
- Diluted radiopaque contrast medium (15:85 medium to saline dilution)
- High pressure contrast injector
- Intra-aortic balloon pump and appropriately sized balloon
- Intracardiac echocardiography (ICE) equipment
- Pigtail angiographic catheter
- Venous return bypass cannula (~18 F)
- Transthoracic echocardiographic (TTE) equipment
- Vascular access lubricant
- Temporary pacing equipment

10.0 INSTRUCTION FOR USE

10.1 Prior to Use Inspection

Before using the EVOQUE system, visually examine each element and accessory for evidence of gross damage (e.g., a cracked jar or lid, leakage, broken or missing seals) that may have compromised the packaging sterility (if applicable) or functionality of the components.

WARNING: Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging and/or sterile barriers and any components have been opened or damaged, or the expiration date has elapsed, as sterility and/or function may be compromised.

WARNING: Do not mishandle the EVOQUE valve or use device/container if found to be damaged, leaking, or without adequate sterilant (not fully submerged in glutaraldehyde, or missing intact seals). The EVOQUE valve must not be used for implantation, as sterility may be compromised.

WARNING: Do not use the EVOQUE valve if the expiration date has elapsed, as either sterility or valve function may be compromised.

WARNING: Do not use the EVOQUE valve if the temperature indicator has been activated, as valve function may be compromised.

WARNING: Do not use the EVOQUE valve if the tamper evident seal is broken, as sterility may be compromised.

10.2 Base & Plate Set Up

If desired, the plate is placed on the operating table beneath the patient's leg to support the base during the procedure.

After the patient has been positioned on the operating table, the base is placed over the patient's leg, on top of the plate, at the desired distance from the mid-sternum.

Establish sterile barrier.

The base adapter is then placed in-line with the access site and attached to the front of the base using a clamp. The base adapter and clamp are both supplied with the stabilizer.

10.3 Device Preparation

All device preparations will be performed by authorized Edwards personnel.

WARNING: Do not mishandle the EVOQUE valve. EVOQUE valve leaflets mishandled or damaged during any part of the loading process will require replacement of the EVOQUE valve.

WARNING: The EVOQUE valve should not remain fully crimped for a period of time longer than 120 minutes, as it may impair valve functionality.

CAUTION: Do not place the jar or the pouch of the delivery system into the sterile field. The outside of the jar and the pouch are not sterile, and the contents within the jar and pouch should be handled using standard aseptic techniques to prevent contamination.

CAUTION: To reduce the risk of contamination, do not open the EVOQUE valve jar until implantation is certain.

CAUTION: Ensure that the entire suture is removed when removing the serial number tag from the EVOQUE valve, as it may lead to emboli.

CAUTION: Do not allow the EVOQUE valve to come in contact with any sharp instrument, as it may impair valve function.

WARNING: Adequate rinsing with normal saline must be performed before implantation to reduce glutaraldehyde concentration, as it may result in glutaraldehyde toxicity.

CAUTION: Avoid contact of the leaflet tissue or the rinse solution with towels, linens, or other sources of lint or particulate matter that could be transferred to the leaflet tissue, as it may lead to emboli.

10.4 EVOQUE Valve Implant

10.4.1 Guidewire Placement

Prepare femoral venous access using standard interventional techniques.

WARNING: Do not use excessive force and/or manipulation during advancement and positioning of the guidewire, as it may result in perforation/dissection of arteries, veins, and/or other cardiac structures. This may also result in cardiac arrhythmias and conduction disturbances.

Step	Procedure
1	Advance steerable sheath into the right atrium at the exit of the IVC.
2	Insert a guidewire through the steerable sheath.
3	Advance the guidewire across the tricuspid valve. Note: Other interventional devices and techniques (e.g., guide catheters) may be used to assist the guidewire in crossing the tricuspid valve.
4	Establish proper guidewire path and confirm no entanglement with cardiac structures.

10.4.2 EVOQUE Valve Delivery

WARNING: Avoid excessive movement of the delivery system while executing the procedure in order to protect vasculature or cardiac structures. Avoid excessive rotation of the delivery system to maintain functionality of the delivery system.

Note: Flush the delivery system with heparinized saline as needed throughout the procedure.

Step	Procedure
1	Ensure the hydrophilic coating on dilators and delivery system is activated prior to use. Dilate access site. An Edwards EVOQUE dilator should be used as needed.
2	Insert the delivery system over the guidewire.
3	Advance the delivery system until the distal end of the tapered tip is positioned at the junction between IVC and right atrium.

Step	Procedure
4	Using fluoroscopy, ensure the delivery system is oriented in the proper orientation. WARNING: The primary flex on the delivery system flexes in the direction of the flush ports; care should be taken to ensure that the delivery system is orientated correctly at this point.
5	Retract the sheath.
6	Flex and orient the delivery system towards the tricuspid valve.
7	Advance the delivery system to cross the tricuspid valve. Note: Delivery system flex, delivery system rotation, and guidewire position may be adjusted during valve crossing to optimize crossing position.
8	Using echocardiography and fluoroscopic guidance, verify that the delivery system has crossed through the tricuspid valve into the right ventricle.
9	Dock the stabilizer onto the base adapter and secure it to the base, if applicable.
10	Attach the delivery system and sheath to the stabilizer, if applicable.
11	Adjust the delivery system as needed to ensure hemodynamic stability.
12	Using pre-operative MSCT data (if available), position the C-arm at the optimal viewing projection.
13	Position the delivery system to be coaxial to the tricuspid annulus while minimizing contact with the native anatomy.
14	Using echocardiography and fluoroscopic guidance, confirm that the EVOQUE valve is positioned at the appropriate depth and coaxiality with relation to the native valve. CAUTION: Once the capsule is retracted to expose the EVOQUE valve anchors, the valve is unable to be retrieved or recaptured into the delivery system. WARNING: Maintain central position of the delivery system within the native valve during deployment to ensure proper positioning of the valve.
15	Retract the outer capsule until the anchors are exposed.
16	Adjust the position of the EVOQUE valve so that the anchors are positioned within native leaflets as dictated by patient anatomy.
17	Retract the outer capsule until the desired EVOQUE valve diameter is achieved.
18	Capture the leaflets.
19	Confirm positioning of the EVOQUE valve by using echo imaging to assess engagement of the leaflets. Adjust the EVOQUE valve positioning as needed.
20	Using echo imaging, observe the motion of the native leaflets and adjust the position of the EVOQUE valve as required to fully engage native tricuspid valve leaflets.
21	Once full engagement has been confirmed, ensure that the EVOQUE valve is perpendicular to the tricuspid annular plane.
22	Retract the tapered tip until it is positioned within the EVOQUE valve.
23	Rotate the release knob to withdraw the inner capsule until the EVOQUE valve is released from the delivery system. CAUTION: Care should be taken during final release of the EVOQUE valve using the release knob, as premature release may impact performance of the EVOQUE valve.
24	Using echo and fluoroscopic imaging, assess final position and functionality of the EVOQUE valve.

10.4.3 Delivery System Removal

WARNING: Take care to maintain position of delivery system central within the EVOQUE valve during delivery system removal, as failure to do so may impact valve functionality or lead to dislodgement of the valve.

Note: Delivery system can be removed from the stabilizer at any point during removal, if applicable.

Step	Procedure
1	Fully retract the tapered tip.
2	Unflex and retract delivery system as needed until the tapered tip is above the locking tabs of the EVOQUE valve. Adjust the guidewire as needed to maintain a central position relative to the EVOQUE valve. Ensure that the locking ring is free from the EVOQUE valve.
3	Rotate the release knob so that the inner capsule is in contact with the tapered tip.
4	Unflex and retract the delivery system as needed.
5	Rotate the capsule retraction knob until the outer capsule is in contact with the inner capsule.
6	Ensure the delivery system is fully unflexed and remove the delivery system from the access site. Note: A sheath may be used to seal the femoral vein following system removal.
7	Perform femoral closure as applicable using standard interventional techniques.
8	Perform a ventriculogram if needed to assess final position of the EVOQUE valve.

11.0 How Supplied

11.1 Sterilization and Packaging

The EVOQUE valve is provided sterile by terminal liquid sterilization and is non-pyrogenic. It is packaged and sterilized in a glutaraldehyde solution inside a jar to which a seal has been applied. The outer surface of the jar is not sterile and must not be placed in the sterile field. The EVOQUE valve is supplied with a temperature indicator and should not be used if the indicator has been activated.

The delivery system, dilator kit, and loading system are supplied sterilized by ethylene oxide and provided non-pyrogenic. The components are secured on a card and packaged in a pouch and shelf box.

The stabilizer is provided ethylene oxide sterilized. The components are secured on a card and packaged in a pouch and shelf box.

The base and plate are provided non-sterile. The components are packaged in individual shipper boxes.

11.2 Storage

The EVOQUE valve should be stored between 10 °C and 25 °C (50 °F and 77 °F). Stock inspection and rotation at regular intervals are recommended so that the EVOQUE valve with an earlier expiration date is used first.

The delivery system, dilator kit, loading system, stabilizer, base, and plate should be stored in a cool, dry place that is contamination free.

12.0 Magnetic Resonance Imaging (MRI) Safety Information

Non-clinical testing has demonstrated that the Edwards EVOQUE valve is MR Conditional. A patient with the valve can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 and 3 T only
- Maximum spatial gradient magnetic field of 3000 gauss/cm (30.0 T/m) or less
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg
- Normal operating mode of the MR system for both gradients and SAR

Under the scan conditions defined above, the EVOQUE valve is expected to produce a maximum temperature rise of 4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the EVOQUE valve extends approximately 0.8 cm from the device when imaged with a gradient echo or spin echo pulse sequence and a 3 T MRI system.

13.0 Expected Lifetime of the Device

The EVOQUE valve has been subjected to rigorous pre-clinical durability testing per the valve testing requirements. The valves were successfully tested to 5 years of simulated wear. Actual lifetime performance in humans depends on many factors and varies from patient to patient.

14.0 Patient Information

A patient implant card is provided with each EVOQUE valve. After implantation, please complete all requested information and provide the implant card to the patient. The serial number is found on the EVOQUE valve package. This implant card allows patients to inform healthcare providers what type of implant they have when they seek care.

15.0 Recovered Implant and Device Disposal

Edwards Lifesciences is interested in obtaining recovered clinical specimens of the EVOQUE valve for analysis. A written report summarizing our findings will be provided upon completion of our evaluation. Please contact Edwards for return of the recovered valve.

If you do decide to return any of the devices, please follow the following instructions:

- **Unopened Package with Sterile Barrier Intact:** If the pouches have not been opened, return the device in its original packaging.
- **Package Opened but Not Implanted:** If a pouch is opened, the device is no longer sterile. Please return the device in its original packaging.
- **Explanted Implant:** The explanted implant should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit for return to Edwards.

15.1 Disposal

Use universal precautions for biohazards and sharps to avoid user injury. Used devices (includes all devices that come in contact with patients) should be handled and disposed of in accordance with institutional guidelines on biohazardous materials and hospital waste to avoid possible cross-contamination.

16.0 Qualitative and Quantitative Information Related to the EVOQUE Valve

For the EVOQUE valve, the following table shows the qualitative and quantitative information on the materials and substances:

EVOQUE Valve

Substance	CAS	Model Mass Range (mg)
Nickel	7440-02-0	1250 - 1470
Titanium	7440-32-6	981 - 1171
Polyethylene terephthalate	25038-59-9	529 - 583
Polyethylene	9002-88-4	285
Polyurethane	9009-54-5	280 - 281
Collagens, bovine, polymers with glutaraldehyde	2370819-60-4	113
Titanium dioxide	13463-67-7	1.22 - 3.26
Cobalt	7440-48-4	0 - 1.29
Iron	7439-89-6	0.00 - 1.29
Carbon	7440-44-0	0 - 1.03
Oxygen	7782-44-7	0 - 1.03
Siloxanes and silicones, di-methyl, polymers with methyl silsesquioxanes, polyethylene-polypropylene glycol, and polyethylene-polypropylene glycol monobutyl ether	68554-68-7	0.762
Antimony trioxide	1309-64-4	0.632 - 0.697
Niobium	7440-03-1	0 - 0.643
Silicon dioxide	7631-86-9	0.598
Polybutadiene diacrylate	9003-17-2	0.452
Polybutylate	24936-97-8	0.333
Chromium	7440-47-3	0 - 0.257
Copper	7440-50-8	0 - 0.257
Hydrogen	1333-74-0	0 - 0.129
Nitrogen	7727-37-9	0 - 0.129
4-Dodecylbenzenesulfonic acid	121-65-3	0.0511 - 0.0531
D&C Green No. 6	128-80-3	0.0480 - 0.0502
Erucamide	112-84-5	0.0208 - 0.0249
Glycerol	56-81-5	0.00924 - 0.0168
Bis(2-dimethylaminoethyl) ether	3033-62-3	0.00550 - 0.00790
Triethylenediamine	280-57-9	0.00580 - 0.00611

EVOQUE Delivery System:

This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless-steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

17.0 Summary of Clinical Experience

Clinical data obtained from the TRISCEND II pivotal trial and TRISCEND study, two clinical trials of the EVOQUE system, are presented in this section.

The clinical benefits of the EVOQUE system for the treatment of tricuspid regurgitation (TR) include the following:

- Effectively reduces tricuspid regurgitation.
- Offers a safe treatment option for tricuspid regurgitation.
- Enables a minimally invasive tricuspid valve replacement procedure.
- Enables clinical improvements in quality of life, functional status, and exercise capacity.

Treatment with the EVOQUE system may lead to elimination or a significant reduction of TR, which can result in clinically meaningful improvements in quality of life, cardiac function, and exercise capacity. The clinical improvement depends on the seriousness and reversibility of the tricuspid regurgitation-associated disease. Follow-up data will be collected and shared with physicians by means of regular updates.

Clinical experience with the EVOQUE system in the TRISCEND II pivotal trial, supported by the TRISCEND study, confirms sufficient safety, efficacy, and performance of the EVOQUE system. Both studies utilized an independent core lab to assess all echocardiographic data. An independent CEC adjudicated safety events, and a DSMB independently reviewed aggregate safety data and evaluated trends of adverse events and their effect on trial conduct and device risk assessment. Long-term follow-up data will be collected as part of the 5-year follow-up for patients enrolled in both clinical trials to confirm safety and effectiveness of the EVOQUE system.

17.1 TRISCEND II Pivotal Trial Results

The TRISCEND II pivotal trial is a prospective, global, multi-center randomized controlled clinical trial conducted to evaluate the safety and effectiveness of the EVOQUE system with optimal medical therapy (OMT) (EVOQUE+OMT) compared to OMT alone in the treatment of patients with at least severe TR. All enrolled study patients are followed at 30 days, 6 months, 1 year, and annually for 5 years. Data described below represent the primary analysis cohort, which includes the first 150 randomized patients, who were randomized 2:1 to EVOQUE+OMT (n=96) or OMT alone (n=54).

17.1.1 Key Demographics and Baseline Characteristics

Enrolled patients had a mean age of 79 years, and 80% were female. The mean EuroSCORE II was 5.3% for EVOQUE+OMT patients and 5.4% for OMT alone patients. Atrial fibrillation was reported for >96% of both treatment groups, and pulmonary hypertension was reported for >70% of both treatment groups. Pre-existing pacemakers or ICDs were present in both treatment groups (36.5% for EVOQUE+OMT, 42.6% for OMT alone).

17.1.2 Safety

The TRISCEND II pivotal trial met the primary safety endpoint. The composite major adverse event (MAE) rate at 30 days for EVOQUE+OMT patients was 27.4%, with a one sided 97.5% upper confidence limit of 36.9%, which is less than the pre-specified performance goal of 70%, based on surgical predicates.

Twenty-six patients in the EVOQUE+OMT group had 35 MAEs prior to the 30-day follow-up visit. The CEC-adjudicated MAEs at 30 days are shown below.

Table 1: Safety Endpoint at 30 Days (Primary Analysis Cohort) - mITT (Safety) Population

Major Adverse Event	Per Patient (%) N=95
Cardiovascular mortality	3.2%
Myocardial infarction	1.1%
Stroke	0.0%
New need for renal replacement therapy	1.1%
Severe bleeding ¹	10.5%
Non-elective tricuspid valve re-intervention, percutaneous or surgical	0.0%
Major access site and vascular complications	3.2%

Major Adverse Event	Per Patient (%) N=95
Major cardiac structural complications due to access-related issues	2.1%
Device-related pulmonary embolism	1.1%
Arrhythmia and conduction disorder requiring permanent pacing	14.7%
Composite MAE	27.4%
¹ Fatal, life-threatening, extensive, or major bleeding, as defined by MVARC	


Beyond the 30-day periprocedural period, the overall rates of site-reported serious adverse events were comparable between the two treatment groups (39.1% of EVOQUE+OMT patients compared to 32.7% of OMT alone patients). Immediately following implantation of the EVOQUE valve, there were no incidences of valve embolization or migration, and no patients required treatment for PVL. As determined by the echocardiographic core lab, over 85% of patients in the EVOQUE+OMT group of the Primary Analysis Cohort had none/trace PVL at discharge, 30 days, and 6 months. Mild PVL was reported to be 7.9% at discharge, 13.1% at 30 days, and 8.8% at 6 months. Moderate PVL was reported to be 4.5% at discharge, 1.2% at 30 days, and 2.5% at 6 months. No patients were reported to have severe or greater PVL at any time point. Currently, no evident harms to patients have been identified for patients with mild valvular and paravalvular leakage; however, potential long-term sequelae are a subject of investigation in post-approval activities.

17.1.3 Primary Effectiveness Endpoints

17.1.3.1 TR Grade Reduction

The TRISCEND II pivotal trial met the first co-primary effectiveness endpoint of TR grade reduction to none/trace, mild, or moderate at 6 months. Treatment with EVOQUE+OMT was demonstrated to be superior to OMT alone in the reduction of severe or greater TR at baseline to moderate or less TR at 6 months. The proportion of patients with severe or greater TR reduced to none/trace, mild, or moderate TR at 6 months was 98.8% (80/81) in the EVOQUE+OMT group and 21.6% (8/37) in the OMT alone group, a difference of 77.1% between the two groups ($p < 0.001$).


Figure 6: TR Severity Grade by Visit – mITT (Effectiveness) Population

17.1.3.2 Hierarchical Composite Endpoint

The TRISCEND II pivotal trial met the second co-primary effectiveness composite endpoint consisting of clinically meaningful quality of life, functional status, and exercise capacity outcomes. Treatment of patients with at least severe TR with EVOQUE+OMT was shown to be superior to OMT alone ($p < 0.001$) in a hierarchical analysis of patient-reported (Kansas City Cardiomyopathy Questionnaire [KCCQ]), clinician-reported (New York Heart Association [NYHA] functional classification), and performance (Six-Minute Walk Distance [6MWD]) outcomes. The unmatched win ratio analysis demonstrated that the EVOQUE system is 4.6 times more likely to improve symptomatic and functional outcomes compared to the medical therapy alone (the lower bound of the one-sided 97.5% confidence interval is 2.6).


Figure 7: Win Ratio – mITT (Effectiveness) Population

17.1.3.2.1 Quality of Life Improvement

Compared to baseline, patients in the EVOQUE+OMT group experienced a mean improvement in their KCCQ overall summary scores of 17.0 points at 30 days and 21.5 points at 6 months. Patients in the OMT alone group experienced a mean improvement in their KCCQ overall summary scores of 1.2 points at 30 days and 3.7 points at 6 months post-implant compared to baseline.


Figure 8: KCCQ by Visit – mITT (Effectiveness) Population

17.1.3.2.2 NYHA Functional Classification Improvement

The proportion of patients in NYHA class I/II in the EVOQUE+OMT group increased from 20.8% (20/96) at baseline to 86.7% (72/83) at 30 days and 90.0% (72/80) at 6 months. The proportion of patients in NYHA class I/II in the OMT alone group was 29.6% at baseline, 30.4% at 30 days, and 34.1% at 6 months.

Figure 9: NYHA Functional Class – mITT (Effectiveness) Population

17.1.3.2.3 Six-Minute Walk Distance Improvement

In paired analysis of 6MWD, the median distance for the EVOQUE+OMT group improved between baseline and 30 days (+14 m; n=81) and 6 months (+18 m; n=73). The median distance for the OMT alone group worsened at 30 days (-30 m; n=39) and 6 months (-16 m; n=34).

17.1.4 TRISCEND II Pivotal Trial Conclusion

Data collected in the TRISCEND II pivotal trial demonstrate the EVOQUE system is a safe and effective treatment option for patients with at least severe tricuspid regurgitation. The trial met the primary safety and effectiveness endpoints, demonstrating (1) an acceptable 30-day safety profile for the EVOQUE system compared to the surgical predicate, (2) superiority of EVOQUE+OMT to OMT alone in TR grade reduction at 6 months, and (3) superiority of EVOQUE+OMT to OMT alone in a hierarchical assessment of clinically meaningful symptomatic and functional improvement (KCCQ, NYHA, and 6MWD) at 6 months.

17.2 TRISCEND Study Results

Results from the single-arm TRISCEND study are an important component of the totality of data for the EVOQUE system, as they provide informative supporting data on the safety, effectiveness, and potential long-term benefits of the EVOQUE system.

The TRISCEND study is a single-arm, prospective, multicenter trial for patients with at least moderate tricuspid regurgitation. As of the data extract, 176 of the 186 total patients who received the EVOQUE valve in the TRISCEND study had completed follow-up through 30 days, with 118 patients completing their 1-year follow-up visit and 31 patients completing 2 years of follow-up. While there are limitations to drawing conclusions from a single-arm study, data from the TRISCEND study are useful for the overall assessment of the EVOQUE system, especially when considering the potential benefits on long-term clinical outcomes.

The TRISCEND study shows that the EVOQUE system is safe and provides durable TR reduction with symptomatic and functional improvement for patients. Prior to treatment with the EVOQUE system, 88.2% of patients had severe or greater TR (i.e., TR grade ≥ 3 , assessed by TTE). At discharge, 96.1% of patients had none/trace or mild TR (i.e., TR grade ≤ 1). The reduction in TR to none/trace or mild was sustained at 30 days (97.0%), 6 months (99.2%), 1 year (97.9%), and 2 years (100.0%). Less than 4% of patients had moderate TR at discharge through 2 years, and no patients were reported to have severe TR after EVOQUE valve implantation. As determined by the echocardiographic core lab, over 82% of patients had none/trace PVL, 11% - 18% had mild PVL, 0% - 4% had moderate PVL, and no patients had severe or greater PVL through 2 years follow-up. Most importantly, treated patients have experienced clinically meaningful improvement in quality of life, functional status, and exercise capacity that is sustained through 1 year, and available data suggest a trend for sustained improvement through 2 years. The percentage of patients in NYHA functional class I/II improved from less than 30% at baseline to over 80% during follow-up (81.5% at 30 days, 90.6% at 6 months, 93.0% at 1 year, and 87.0% at 2 years).

Table 2: Summary of TRISCEND Study Clinical Results

Outcome Measure	Time Point			
	30 Days	6 Months	1 Year	2 Years
Patients Treated	n=176	n=142	n=118	n=31
Kaplan-Meier estimate of composite MAEs**	17.3%	21.3%	25.1%	25.1%
Proportion of patients with mild or none/trace TR	97.0% (162/167)	99.2% (125/126)	97.9% (94/96)	100.0% (17/17)
Mean improvement in KCCQ score from baseline (points)	20.5**	26.0**	24.7**	26.5**
Improvement by ≥ 1 NYHA functional class	69.0% (116/168)	79.5% (101/127)	85.0% (85/100)	78.3% (18/23)
Mean improvement in 6MWD from baseline (m)	39.8**	49.5**	51.0**	58.6 [†]
Kaplan-Meier estimate of all-cause mortality [‡]	2.6%	7.7%	10.3%	13.2%
Kaplan-Meier estimate of heart failure hospitalization [‡]	3.7%	8.5%	14.7%	16.9%
* Composite MAE rate for TRISCEND study does not include event category of "arrhythmia and conduction disorder requiring permanent pacing."				
** p<0.001 compared to baseline				
[†] p=0.036				
[‡] CEC-adjudicated rate				

Data collected in the TRISCEND study support the safety and performance of the EVOQUE system in patients with tricuspid regurgitation. Patients treated with the EVOQUE system experienced significant reductions in TR with clinically meaningful and statistically significant improvement in quality of life, functional status, and exercise capacity. Improvements were sustained through 1 year, and available data suggest a trend for sustained improvement through 2 years.

17.3 Clinical Conclusion

In conclusion, the totality of data demonstrates that the benefits of the EVOQUE system outweigh the risks of its use and provides a reasonable assurance of safety and effectiveness. Data from the primary analysis cohort of the TRISCEND II pivotal trial, supported by the results of the TRISCEND study, establish the EVOQUE system as a safe and effective treatment option for patients with at least severe TR.

18.0 Summary of Safety and Clinical Performance (SSCP)

Refer to <https://meddeviceinfo.edwards.com/> for a SSCP for this medical device.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device. This is the location of the SSCP after the launch of European Database on Medical Devices/Eudamed.

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information. The Notified Body has taken notice of and agreed with the benefit-risk rationales for the short- and long-term safety and effectiveness of the EVOQUE system.

Conformity with the entire EVOQUE system of the Performance Requirements (GSPR) for safety (MDR GSPR1), performance (MDR GSPR1), acceptability of side-effects (MDR GSPR8), usability (MDR GSPR5), device lifetime (MDR GSPR6), and acceptable benefit-risk profile (MDR GSPR8) has been established for the labelled indications.

19.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed. The Basic UDI-DI for the valves, delivery system, and accessories can be used to locate the SSCP. The following table contains the Basic UDI-DIs for the EVOQUE system:

Product	Model			Basic UDI-DI
Edwards EVOQUE Valve	9850EV44	9850EV48	9850EV52	0690103D004EVV000Z3
Edwards EVOQUE Tricuspid Delivery System	9850TDS			0690103D004EVD000V5
Edwards EVOQUE Dilator Kit	9850DK			0690103D004EVK000WN
Edwards EVOQUE Loading System	9850LS			0690103D004EVL000VV
Edwards EVOQUE Stabilizer	9850SB			0690103D004EVA000UG
Edwards EVOQUE Stabilizer Base	9850BA			0690103D004EVA000UG
Edwards EVOQUE Stabilizer Plate	9850PT			0690103D004EVA000UG

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© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences GmbH
Edwards Lifesciences LLC
Edwards Lifesciences B.V.