Edwards

SAPIEN M3 Transcatheter Mitral Valve Replacement System

Directory	
English (en).....	1
Français (fr).....	12
Deutsch (de).....	23
Español (es).....	35
Italiano (it).....	46
Nederlands (nl).....	58
Dansk (da).....	69
Svenska (sv).....	80
Ελληνικά (el).....	91
Português (pt).....	103
Česky (cs).....	114
Polski (pl).....	125
Norsk (no).....	137
Suomi (fi).....	148
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Rysunki ■ Figurer ■ Kuvat.....	159
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset.	165

English

SAPIEN M3 Transcatheter Mitral Valve Replacement System

Instructions for Use

Implantation of the SAPIEN M3 Dock and SAPIEN M3 Transcatheter Heart Valve should be performed only by physicians who have received Edwards Lifesciences training.

1.0 Device Description

The SAPIEN M3 transcatheter mitral valve replacement system (also known as SAPIEN M3 system) consists of the following devices:

Product Name	Model/REF
SAPIEN M3 Transcatheter Heart Valve	9880TFX29M
SAPIEN M3 Dock Steerable Catheter ^[1]	9880DDS
Edwards Commander M Delivery System ^[2]	9880CM29
SAPIEN M3 Crimper	9880CR
SAPIEN Stabilizer Rail System	9880SRS
^[1] Includes a steerable catheter, dock, and preparation accessories	
^[2] Includes a loader and crimp stopper	

The SAPIEN M3 system is intended to be used with:

Product Name	Model/REF
Edwards 23F Guide Sheath ^[1]	9880GS
Edwards Locking Syringe	96406
^[1] Includes an introducer	

The SAPIEN M3 system is intended to be used with one or the other supporting devices:

Edwards single-use accessory:

Product Name	Model/REF
Edwards Table	9880T

Edwards reusable accessories:

Product Name	Model/REF
Edwards Reusable Platform	10000UP
Edwards Reusable Plate	10000PT
Edwards Reusable Cradle	10000CR

• SAPIEN M3 Transcatheter Heart Valve (Figure 1)

The SAPIEN M3 transcatheter heart valve (also known as SAPIEN M3 valve or valve) is a bioprosthesis comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, three bovine pericardial tissue leaflets, polyethylene terephthalate (PET) fabric inner skirt, and full-frame outer skirt. The valve frame inflow and outflow are covered with expanded polytetrafluoroethylene (ePTFE). A green suture line indicates the valve frame inflow. All bovine pericardium is treated according to the Carpentier-Edwards TheraFix process.

Valve Size	Valve Height
29 mm	22.5 mm

• SAPIEN M3 Dock Steerable Catheter (Figure 2)

The SAPIEN M3 dock steerable catheter (also known as steerable catheter) is used for delivery of the SAPIEN M3 dock (also known as dock) to its intended location and consists of a dock, steerable catheter, dock handle, a hydrophilic coated removable sleeve, and device preparation accessories. The steerable catheter handle contains two (2) flex wheels, tip flex wheel and back flex wheel, that control steerable catheter flexing and aid in gaining access to the left ventricle. Two (2) radiopaque rings are provided to define the tip flex and back flex regions of the steerable catheter. There is a radiopaque tip on the steerable catheter. The dock handle aids in advancing and/or retrieving the dock and can be locked to the steerable catheter. The dock encircles the native mitral leaflets which applies an inward force to the mitral apparatus, pulling the leaflets and chordae to the dock center and approximating the papillary muscles, to provide a landing zone for the SAPIEN M3 valve. The dock is attached to the dock handle by a suture. The dock has a nitinol core that is covered with ePTFE tubing and PET braid. There is a self-expanding PVL guard made of a PET covered nitinol braid which expands at the medial commissure of the native mitral valve. The dock includes one (1) radiopaque marker to visualize proper deployment positioning within the anatomy, and two (2) radiopaque markers to visualize final deployment positioning of the PVL guard. The removable sleeve covers the dock during encircling and is designed to facilitate encircling of the mitral apparatus and is removed after the dock is positioned in the anatomy. The sleeve has a radiopaque marker that aids in visualization. A release assembly at the proximal end of the dock handle allows for suture slack or suture tension. The suture can be cut to release the dock from the dock handle. Additional provided accessories (two white luer caps, one red stopcock, one extension tubing, and one aspiration adapter) are to be used during device preparation.

Model	Shaft Effective Length	Shaft O.D.
9880DDS	113 cm	18 F (6.2 mm)

• Edwards Commander M Delivery System (Figure 3)

The Edwards Commander M delivery system (also known as valve delivery system) is used for delivery of the SAPIEN M3 valve and consists of a balloon catheter and flex catheter. The valve delivery system includes a tapered tip to facilitate crossing of the native valve/dock. The handle contains a flex wheel and a balloon lock. The balloon catheter has radiopaque markers defining the crimp location for the valve. A radiopaque double marker proximal to the balloon indicates the flex catheter and guide sheath position during deployment. A balloon cover and crimp alignment tool are packaged on the valve delivery system. A stylet is included within the guidewire lumen. The crimp alignment tool assists with crimping the valve onto the balloon in the proper location. A peel-away loader and 2-piece crimp stopper are included in the packaging. The inflation parameters for valve deployment are:

Model	Shaft Effective Length	Shaft O.D.	Inflation Volume	Balloon Outer Diameter
9880CM29	112 cm	16 F (5.4 mm)	33 ml (initial)	30 mm
			37 ml (post-dilation)	31 mm

• SAPIEN Stabilizer Rail System (Figure 4)

The SAPIEN stabilizer rail system (also known as stabilizer) holds the guide sheath and the steerable catheter to support and maintain catheter positioning and must be used in conjunction with either the single-use or reusable supporting accessories.

• SAPIEN M3 Crimper (Figure 5)

The SAPIEN M3 crimper (also known as crimper) reduces the diameter of the valve to mount it onto the Edwards Commander M delivery system. The crimper is comprised of a housing, a compression mechanism that is closed with a handle located on the housing, and a crimp stopper.

• Edwards 23F Guide Sheath

For Edwards 23F guide sheath (also known as guide sheath), refer to Edwards 23F Guide Sheath Instructions for Use (herein referred to as Edwards 23F Guide Sheath IFU).

• Edwards Table

For Edwards Table (also known as table), refer to Edwards Table Instructions for Use (herein referred to as Edwards Table IFU).

• Edwards Reusable Accessories

For Edwards reusable accessories (also known as reusable accessories), refer to Reusable Platform (Model 10000UP), Reusable Plate (Model 10000PT), and Reusable Cradle (Model 10000CR) Instructions for Use (herein referred to as Edwards Reusable Accessories IFU).

• Edwards Locking Syringe

For Edwards locking syringe (also known as inflation device), refer to Edwards Locking Syringe Instructions for Use (herein referred to as Edwards Locking Syringe IFU).

Note: For proper volume sizing, the Edwards Commander M delivery system should be used with the inflation device provided by Edwards Lifesciences.

2.0 Intended User and Intended Use / Intended Purpose

2.1 Intended Users

The intended users of the SAPIEN M3 system include two user groups:

One intended user group is responsible for device preparation and includes operating room nurses or cardiovascular technicians who have received Edwards Lifesciences training, or certified Edwards field representatives.

The second intended user group is responsible for procedural device use and includes interventional cardiologists or cardiovascular surgeons.

2.2 Intended Use / Intended Purpose

The SAPIEN M3 transcatheter heart valve is intended to replace the mitral valve. The SAPIEN M3 dock steerable catheter is intended to deliver the SAPIEN M3 dock to its intended location. The SAPIEN M3 dock, which is attached to the SAPIEN M3 dock steerable catheter, is intended to create a landing zone for the implantation of the SAPIEN M3 valve. The Edwards Commander M delivery system is intended to facilitate the placement of the SAPIEN M3 valve. The SAPIEN M3 crimp and crimp stopper is intended to reduce the diameter of the valve thereby mounting the SAPIEN M3 valve onto the Edwards Commander M delivery system. The SAPIEN stabilizer rail system is intended to hold the Edwards 23F guide sheath and the SAPIEN M3 dock steerable catheter to support and maintain catheter positioning.

3.0 Intended Patient Population and Indication for Use

3.1 Intended Patient Population

The SAPIEN M3 system is intended for adult patients with mitral valve regurgitation.

3.2 Indication for Use

The SAPIEN M3 system is indicated for treatment of patients with symptomatic mitral regurgitation (moderate-to-severe or severe MR) who are deemed unsuitable for surgery or transcatheter edge-to-edge therapy by a heart team.

4.0 Contraindications

Use of the SAPIEN M3 system and accessories is contraindicated in patients with:

- Evidence of intracardiac mass, thrombus, vegetation, active infection or endocarditis
- Inability to tolerate intraprocedural heparin
- Inability to tolerate anticoagulation/antiplatelet therapy

5.0 Warnings

- The SAPIEN M3 system devices are designed, intended, and distributed STERILE for single use only. Do not resterilize or reuse the devices. There are no data to support the sterility, non-pyrogenicity, and functionality of the devices after reprocessing. The table and reusable accessories are non-sterile.
- Do not mishandle the SAPIEN M3 system devices or use them if the packaging or any components are not sterile, have been opened or are damaged (e.g., kinked or stretched), or the expiration date has elapsed.
- Patients with hypersensitivities to cobalt, nitinol (nickel or titanium), chromium, molybdenum, manganese, silicon, bovine tissue, and/or polymeric materials may have an allergic reaction/immunological response to these materials.
- Accelerated deterioration of the valve may occur in patients with altered calcium metabolism.
- Exercise caution when implanting a valve in patients with clinically significant coronary artery disease as it may result in myocardial ischemia.
- Prior to delivery, the valve must always remain hydrated and cannot be exposed to solutions other than its shipping storage solution and sterile physiologic rinsing solution. Valve leaflets mishandled or damaged during any part of the procedure will require replacement of the valve.
- Do not use the valve if the tamper-evident seal is broken, the storage solution does not completely cover the valve, the temperature indicator has been activated, the valve is damaged, or the expiration date has elapsed.
- Do not add or apply antibiotics to the storage solution, rinse solutions, or the valve.
- The physician must verify correct orientation of the valve prior to its implantation.
- The procedure should be conducted under 3D echocardiography and fluoroscopic guidance. Some fluoroscopically guided procedures are associated with a risk of radiation injury to the skin. These injuries may be painful, disfiguring, and long-lasting.
- Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.
- Observation of the pacing lead throughout the procedure is essential to avoid the potential risk of pacing lead perforation.
- In the event of device malfunction or device damage during use (e.g., destructive deformation to the catheter, balloon burst, etc.) safely remove the device(s). If unable to safely remove the device(s), conversion to surgery is recommended.
- Prior to valve deployment, 3D echocardiographic and fluoroscopic (short-axis view) verification must be used to confirm that the guidewire passes through the center of the implanted dock and has unrestricted movement. Failure to do so can result in chordal rupture and/or the valve being deployed outside of target location.
- Incorrect positioning of the dock and/or valve may lead to left ventricular outflow tract obstruction, paravalvular leak (PVL), valve migration, or valve embolization.
- Valve recipients must be on appropriate anticoagulation regimen, determined at the physician's discretion based on individual subject needs for a minimum of 6 months. Failure to anticoagulate and bridge appropriately will lead to valve thrombosis. For subjects receiving vitamin K antagonists, target range for INR is 2.5 to 3.5. After 6 months, continued antithrombotic therapy is recommended as tolerated.

6.0 Precautions

- Glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Material Safety Data Sheet available from Edwards Lifesciences.
- Additional precautions for transseptal replacement of a mitral valve include abnormalities in the caval vein precluding safe transvenous femoral access for transseptal approach, presence of atrial septal occluder device, or calcium preventing safe transseptal access.

- Patients with a pre-existing prosthesis should be evaluated for the location, shape, construction, and characteristics of the prosthesis (e.g., low-deployed aortic prosthesis, rigid or small annuloplasty ring, septal occluder, etc.) as it may interfere with SAPIEN M3 system deployment, functionality, or dock/valve durability.
- Patients with mitral annular calcification should be evaluated for the characteristics of the calcium and mitral pathology as it may interfere with the dock trajectory during deployment, result in malposition of the dock/valve, and/or have an increased risk of PVL.
- Patient's sub-valvular anatomy should be evaluated for the characteristics of papillary muscles, chordae, and ventricular wall as it may interfere with or prevent dock deployment.
- Patients with the following characteristics have an increased risk of PVL which may lead to hemolysis and/or intervention:
 - Compromised leaflet integrity (e.g., perforation, endocarditis, Barlow's syndrome, etc.)
 - Flail or prolapse located at the commissures
 - Flail or prolapse located at P3 leaflet in conjunction with a commissural distance > 42 mm
 - Any large non-commissural flail or prolapse
- To maintain proper valve leaflet coaptation, do not overinflate the deployment balloon.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Long-term durability has not been established for the valve. Regular medical follow-up is advised to evaluate valve performance.
- The safety and effectiveness of the SAPIEN M3 system have not been established for patients who have/are:
 - A left ventricular end-diastolic diameter greater than 75 mm
 - A commissural distance greater than 50 mm
 - A left ventricular ejection fraction below 25%
 - Severe RV dysfunction
 - Severe mitral annular calcification
 - History of heart transplant
 - Severe pulmonary hypertension
 - Blood dyscrasias defined as: leukopenia (WBC < 3000 cells/ μ L), acute anemia (Hb < 9 g/dL), thrombocytopenia (platelet count < 50,000 cells/ μ L), or history of bleeding diathesis or coagulopathy
 - Female who is pregnant or lactating
 - Children

7.0 Potential Adverse Events

The following potential risks are associated with the procedure and device usage including potential access complications associated with standard cardiac catheterization, the potential risks of anesthesia, and the use of angiography.

- Death
- Stroke or neurological dysfunction
- Cardiovascular injury – cardiac structure complications
- Cardiovascular injury – vascular complications
- Cardiovascular injury – access related complications
- Heart failure or low cardiac output / Worsening of heart failure
- Renal insufficiency or renal failure
- Cardiogenic shock
- Cardiac arrest
- Pericardial effusion or cardiac tamponade
- Thromboembolism including air, calcific valve material, or thrombus
- Retroperitoneal bleed
- Arrhythmia
- Hypertension or hypotension
- New or worsening valvular regurgitation
- Bleeding / Hematoma / Hemorrhage
- Hemolysis that may require transfusion or intervention
- Device/Valve thrombosis
- Respiratory insufficiency or respiratory failure
- Paravalvular or transvalvular leak
- Device deterioration (wear, fracture, calcification, or other)
- Reoperation / Reintervention
- Device explants
- Pleural effusion
- LVOT obstruction
- Emergency cardiac surgery
- Conversion to cardiac surgery
- Thoracic bleeding
- Valve stenosis
- Myocardial infarction
- Pulmonary edema
- Transient ischemic attack including clusters
- Device migration, malposition or embolization
- Infection including septicemia and endocarditis
- Allergic reaction to anesthesia, contrast media, or device materials
- Deterioration of native valve (leaflet tear/tearing, leaflet retraction, leaflet thickening, or other)
- Structural valve deterioration (wear, fracture, calcification, leaflet tear/tearing from the stent posts, leaflet retraction, suture line disruption of components of a prosthetic valve, thickening, stenosis)

- Nonstructural valve dysfunction
- Atrial septal defect
- Syncope
- Conduction system defect which may require a permanent pacemaker
- Skin burn
- Mechanical failure of delivery system, and/or accessories
- Valve deployment in an unintended location
- Abnormal lab values (including electrolyte imbalance)
- Angina
- Anemia
- Fever
- Inflammation
- Pain or changes at the access site

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

8.0 Directions for Use

8.1 Required Equipment

- Standard cardiac catheterization lab equipment
- Fluoroscopy (fixed, mobile, or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- 3D Transesophageal echocardiography system
- Instrumentation for transseptal access
- Exchange length 0.035" (0.89 mm) extra-stiff guidewires
- Exchange length (≥ 270 cm) pre-shaped stiff TAVR type wire
- ≥ 120 cm pigtail catheter
- Temporary pacemaker and pacing leads
- Sterile heparinized saline (2,000 IU/L) bags (x2)
- Infusion pumps (x2) and infusion pump tubing
- Sterile extension tubing, length > 36 inches (x2)
- Sterile rinsing basins (x6), physiological saline, heparinized saline, diluted radiopaque contrast medium (15 ml contrast: 85 ml heparinized saline)
- SAPIEN M3 Valve
- Edwards 23F Guide Sheath
- SAPIEN M3 Dock Steerable Catheter
- Edwards Commander M Delivery System
- Edwards Locking Syringe
- SAPIEN Stabilizer Rail System
- SAPIEN M3 Crimper
- Edwards Table or Edwards Reusable Accessories

8.2 Device Handling and Preparation

Follow sterile technique during device preparation and implantation.

8.2.1 Prepare the Devices

The SAPIEN M3 system is intended to be used with either the Edwards Table or Edwards Reusable Accessories.

8.2.1.1 Edwards Table

The device is designed, intended, and distributed for single use. For Edwards Table, refer to the Edwards Table IFU.

8.2.1.2 Edwards Reusable Accessories

The devices are designed, intended, and distributed for multiple uses. For reusable accessories, refer to the Edwards Reusable Accessories IFU.

8.2.1.3 Edwards 23F Guide Sheath

For Edwards 23F guide sheath device preparation, refer to the Edwards 23F Guide Sheath IFU.

8.2.1.4 SAPIEN M3 Dock Steerable Catheter

1. Verify expiration date, model number, and visually inspect for breaches in package integrity prior to opening sterile package.
2. Visually inspect all components for damage.
3. Fill the tray reservoir beneath the dock with heparinized saline to submerge the dock and sleeve.
4. Attach provided red stopcock to dock handle flush port.
5. Using a 20 cc or larger syringe, flush dock handle flush port with 20 ml heparinized saline. Attach provided white luer cap to release assembly flush port.
6. Using a 20 cc or larger syringe, flush the dock handle flush port with an additional 20 ml heparinized saline. Close the red stopcock to the system.
7. Using a 50 cc or larger syringe, flush steerable catheter flush port with heparinized saline until heparinized saline comes out of the sleeve lock. Lock sleeve lock and continue flushing until air is removed from sleeve lumen. Close the stopcock to the system.
8. Set up two continuous, heparinized (2,000 IU/L) saline pumps and sterile IV lines infusing at 200 ml/hr each.
9. Connect one, deaired IV line with extension tubing to steerable catheter flush port. Open the stopcock to the system.
10. Connect another, deaired IV line with extension tubing to dock handle flush port. Open the red stopcock to the system. Attach second provided white luer cap to the red stopcock.
11. Align sleeve and dock tip. Lock sleeve lock and leave $\frac{1}{4}$ turn of dock out of steerable catheter after alignment.
12. Attach provided extension tubing and 50 cc or larger syringe to aspiration adapter and de-air.

13. Insert dock into blue cap of aspiration adapter and tighten cap.
14. Ensure steerable catheter tip is secured in tray and submerged. Perform distal aspiration.
15. Loosen blue cap and remove aspiration adapter from dock.
16. Confirm that heparinized saline is flowing out of the steerable catheter tip and the sleeve tip.

WARNING: Failure to ensure the heparinized saline drips are each infusing at a rate of 200 ml/hr may result in thrombus formation in the steerable catheter.

17. Retract the sleeved dock fully into steerable catheter. Confirm dock position inside steerable catheter.
18. Lock the dock lock.
19. Confirm all connections are tight, IV pumps are flowing, and all flexes are fully relaxed.

8.2.1.5 Edwards Commander M Delivery System

1. Verify expiration date, model number, and visually inspect for breaches in package integrity prior to opening sterile package.
2. Visually inspect all components for damage.
3. Remove all components from the package.
4. Carefully remove the balloon cover from the valve delivery system. Remove the crimp alignment tool and set aside.
5. Flush the valve delivery system with heparinized saline.
6. Place the loader cap onto valve delivery system.
7. Attach a high pressure 3-way stopcock to the balloon inflation port. Partially fill a 50 cc or larger syringe with 15-20 ml diluted contrast medium and attach to the 3-way stopcock.
8. De-air the inflation device with diluted contrast medium. Lock the inflation device and attach to the 3-way stopcock.
9. Close the 3-way stopcock to the inflation device and fully de-air the balloon using a 50 cc or larger syringe. Slowly release the plunger and leave zero-pressure in the system.
10. Ensure 33 ml is in the inflation device. Lock the inflation device.
11. Open stopcock to the inflation device. Remove the 50 cc or larger syringe.

CAUTION: Maintain the inflation device provided by Edwards Lifesciences in the locked position until valve deployment to minimize the risk of premature balloon inflation and subsequent improper valve deployment.

8.2.2 SAPIEN M3 Valve Rinsing Procedure

Before opening the valve jar, verify expiration date, model number, and visually inspect the TagAlert temperature indicator and for breaches in package integrity prior to opening sterile package. Carefully examine for evidence of damage (e.g., a cracked jar or lid, leakage, or broken or missing seals).

CAUTION: Valves from containers found to be damaged, leaking, without adequate sterilant, or missing intact seals must not be used for implantation as sterility may be compromised.

1. Set up three (3) sterile bowls with at least 500 ml of sterile physiological saline.
2. Carefully remove the valve/holder assembly from the jar without touching the tissue. Verify the valve serial number with the number on the jar lid and record in the patient information documents. Inspect the valve for any signs of damage to the frame, outer skirt, or tissue. Do not use if damaged.
3. Rinse the valve three times with agitation for a minimum of a minute in each of the three separate bowls filled with physiological saline. Be sure the saline solution completely covers the valve and holder. Ensure the rinse solution in the first and second bowls are not re-used. The valve should be left in the final rinse solution until needed to prevent the tissue from drying.

WARNING: Adequate rinsing with physiological saline must be performed before implantation to reduce glutaraldehyde concentration. Failure to rinse the valve in the three bowls with agitation for a minimum of one minute in each bowl may result in glutaraldehyde toxicity.

CAUTION: Do not allow the valve to come into contact with the bottom or sides of the rinse bowl during agitation or swirling in the rinse solution. Direct contact between the identification tag and valve is also to be avoided during the rinse procedure. No other objects should be placed in the rinse bowls. The valve should be kept hydrated to prevent the tissue from drying.

8.2.3 Mount and Crimp the SAPIEN M3 Valve onto the Edwards Commander M Delivery System

1. Remove the valve from the holder and remove the ID tag.
2. Verify expiration date, model number, and visually inspect for breaches in package integrity of the SAPIEN M3 Crimper prior to opening sterile package.
3. Remove crimper from packaging and inspect for damage.
4. Attach the 2-piece crimp stopper to the base of the crimper.
5. With the crimper in the open position, gently place the valve into the crimper aperture. Rotate the handle until the crimper aperture contacts the valve.
6. Insert the valve delivery system coaxially within the valve ensuring that the inflow of the valve (green suture line) is oriented towards the valve delivery system handle.
7. Position the valve on the balloon such that the outflow edge of the valve is proximal to the distal shoulder of the valve delivery system.
8. Crimp the valve until it reaches the initial stop located on the 2-piece crimp stopper. Hold for 5 seconds. Remove the valve delivery system from crimper.
9. Place the crimp alignment tool onto stylet tip.
10. Adjust the partially crimped valve until the edge of the valve is touching the crimp alignment tool.
11. Remove crimp alignment tool from stylet while keeping stylet in place.
12. Remove the initial stop from 2-piece crimp stopper, leaving the final stop in place.
13. With the crimper in the open position, gently place the partially crimped valve into the center of the crimper aperture.
14. Fully crimp the valve until it reaches the final crimp stopper. Hold for 5 seconds.
15. Fully crimp two (2) more times, holding 5 seconds for each crimp, for a total of three (3) full crimps.
16. Advance the flex catheter to the edge of the valve by retracting the balloon shaft. Lock the balloon lock.
17. Place the loader over the tapered tip, valve, and flex catheter. Align the distal tip of loader to the distal bump of the balloon.

CAUTION: The valve should not remain fully crimped and/or in the loader for over 15 minutes, as leaflet damage may result, and impact valve functionality.

18. Attach the loader cap to the loader, flush the valve delivery system through the flush port to de-air loader.

Remove the stylet and flush the guidewire lumen of the valve delivery system.

CAUTION: Keep the valve hydrated until ready for implantation to prevent damage to the leaflets which may impact valve functionality.

CAUTION: The physician must verify correct orientation of the valve prior to its implantation to prevent the risk of severe patient harm; the green suture (inflow) of the valve should be oriented towards the proximal (handle) end of the valve delivery system.

8.3 SAPIEN M3 Dock and SAPIEN M3 Valve Delivery Procedure

Delivery of the dock and valve should be performed under general anesthesia with hemodynamic monitoring in an operating room, hybrid operating room, or catheterization laboratory with fluoroscopic and 3D echocardiographic imaging capabilities.

CAUTION: Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

CAUTION: Excessive device manipulation may result in cardiac structure damage requiring surgical repair or other intervention.

CAUTION: Ensure that a 15F or larger catheter/device is across the guide sheath seals when aspirating and flushing the guide sheath to reduce the risk of air embolization.

8.3.1 Patient Preparation

1. If using the Edwards Table:

Prior to sterile draping the patient for standard transseptal catheterization via femoral vein access, position the table between the legs of the patient, adjusting height of the table as needed.

WARNING: The table is provided non-sterile; introduction of the table into the sterile field may result in infection.

OR

If using the Edwards Reusable Accessories:

Prior to sterile draping the patient for standard transseptal catheterization via femoral vein access, assemble and position the reusable accessories around the legs of patient, adjusting the height of the platform as needed and maintaining a level surface. Place the cradle on the platform in-line with the intended femoral vein access site. Refer to the Edwards Reusable Accessories IFU.

CAUTION: Always maintain a level surface with the reusable platform as tilting the platform could increase air embolization risk.

WARNING: The reusable accessories are non-sterile; introduction of the reusable accessories into the sterile field may result in infection.

8.3.2 Baseline Parameters

1. Introduce a pacemaker (PM) lead and position appropriately.
2. Set the stimulation parameters to obtain 1:1 capture, and test pacing.
3. Mark the mitral annular plane using a ventriculogram or anatomical landmarks.

8.3.3 Transseptal Access and Edwards 23F Guide Sheath Introduction

1. Insert the guide sheath according to the Edwards 23F Guide Sheath IFU and follow additional instructions for use with the SAPIEN M3 system:
 - a) Access the common femoral vein using conventional percutaneous puncture methods.
 - b) Access the left atrium via transseptal puncture using conventional percutaneous methods and place guidewire in the left atrium.

CAUTION: Inappropriate puncture may result in cardiac structure damage, requiring surgical repair or other intervention.
 - c) Administer heparin to maintain the ACT at ≥ 300 sec.
2. Attach the stabilizer onto the table or reusable accessories (whichever is in use).
 - a) Insert the guide sheath and introducer with the flushport oriented away from the operator. Advance until guide sheath tip is across the septum.
3. Secure the guide sheath to the stabilizer.
 - a) Unlock introducer from guide sheath and slowly retract the introducer and guidewire.
 - b) With the introducer and wire remaining across the guide sheath seals, aspirate and flush the guide sheath. Remove the introducer and wire slowly.
4. Adjust guide sheath (flex, torque) until the tip of the guide sheath is parallel to the mitral valve plane, as confirmed via fluoroscopy and echocardiography.

8.3.4 SAPIEN M3 Dock Steerable Catheter Positioning

1. Confirm that the heparinized saline is flowing out of the steerable catheter tip.

WARNING: Failure to ensure a continuous heparinized saline drip is infusing at a rate of 200 ml/hr per pump may result in thrombus formation in the steerable catheter.
2. Insert steerable catheter into guide sheath just past the guide sheath seals.
3. Aspirate and flush guide sheath.
4. Continue advancing the steerable catheter until the tip is positioned at the tip of guide sheath. Secure the steerable catheter to the stabilizer.
5. Under echocardiographic and fluoroscopic guidance, manipulate steerable catheter and guide sheath until steerable catheter tip is in the left ventricle through the medial commissure.
6. Under echocardiographic guidance, verify ventricular access at the medial commissure.

8.3.5 SAPIEN M3 Dock Deployment

1. Under echocardiographic and fluoroscopic guidance, advance dock handle to deploy $\frac{1}{4}$ turn of the dock in the ventricle.

WARNING: To avoid damage to cardiac structures, ensure that the steerable catheter tip is not pointing directly towards the wall of the left atrium or left ventricle.
2. Under echocardiographic guidance, verify the $\frac{1}{4}$ turn of the dock has not crossed the aortic valve.
3. Manipulate steerable catheter and guide sheath to orient dock trajectory to be parallel with the mitral plane.
4. Advance dock handle to encircle anterior and posterior leaflets with the first functional turn of the dock.

CAUTION: Do not apply back flex during encircling as it can result in cardiac structure damage requiring surgical repair or other intervention.

WARNING: Do not apply excessive force to the dock. Excessive force applied to the dock can be transferred to the anatomy and may result in cardiac structure damage requiring surgical repair or other intervention.

Note: Excessive force is indicated visually by the dock tip not advancing freely (dock loading), the dock diameter enlarging, the dock handle experiencing pushback, or if the dock functional turns are overlapping.

5. Using 2D echocardiography, verify that both mitral leaflets have been captured by the first functional turn of the dock and are moving freely within the dock.
6. Reposition dock as needed.
7. Continue advancing dock handle until the deployment marker band has exited the tip of the steerable catheter. Reengage the dock lock.
If the dock functional turns are overlapping, retract the dock and re-attempt encircling. If unable to achieve an acceptable dock orientation, remove the device.

CAUTION: Do not apply back flex during encircling as it can result in cardiac structure damage requiring surgical repair or other intervention.

WARNING: Do not apply excessive force to the dock. Excessive force applied to the dock can be transferred to the anatomy and may result in cardiac structure damage requiring surgical repair or other intervention.

Note: Excessive force is indicated visually by the dock tip not advancing freely (dock loading), the dock diameter enlarging, the dock handle experiencing pushback, or if the dock functional turns are overlapping.

8. Retract steerable catheter until straddling markers are inside the guide sheath.
9. Apply posterior torque to the steerable catheter. Unflex guide sheath and steerable catheter. Verify that the dock tip and deployment marker band are aligned in the correct position.
10. Pull the sleeve handle to unsleeve the dock. Ensure the sleeve marker band remains on the dock. Relock the sleeve lock.
11. Unlock the dock lock. Hold dock handle in place and retract steerable catheter to be flush with the guide sheath tip. Relock the dock lock.
12. Seat the PVL guard by advancing the sleeve over the dock until the sleeve marker band is between the two (2) seating marker bands.
13. Retract the sleeve handle to completely unsleeve the dock. Relock the sleeve lock.

8.3.6 SAPIEN M3 Dock Release

1. If needed, minimize length of guide sheath in the atrium.
2. Advance dock handle until dock handle tip exits the steerable catheter and the guide sheath.
3. Provide suture slack. Translate suture to the tip of the dock. Repeat as necessary.
4. Maneuver the guide sheath to reduce distance from the atrial turn of the dock.
5. Allow the dock to settle over several cardiac cycles to verify that the dock deployment marker band continues to remain in the correct location.

WARNING: Incorrect position of dock may result in damage to cardiac structures, paravalvular leak or damage to native mitral leaflets.

6. In commissure-to-commissure view on echo, ensure the dock is less than or equal to 12 mm ventricular from the mitral plane.

WARNING: Failure to ensure correct dock depth may result in paravalvular leak and/or cardiac structure damage requiring surgical repair or other intervention.

7. Turn the red stopcock off to the steerable catheter.

WARNING: Failure to turn off the flow to the dock handle flush port could result in an air embolism.

8. Open the red suture cover, cut the suture, and remove suture from the steerable catheter. Re-attach the release assembly.

CAUTION: Do not cut suture, release the dock, or proceed with remaining steps until dock position and appropriate capture of mitral leaflets has been verified under echocardiography and fluoroscopy.

WARNING: Incorrect position of dock may result in paravalvular leak and/or cardiac structure damage requiring surgical repair or other intervention.

9. Unflex steerable catheter. Slowly remove the steerable catheter while leaving guide sheath in place.
10. Prior to removing the steerable catheter, aspirate and flush the guide sheath. Fully remove the steerable catheter.
11. Verify that the dock remains less than or equal to 12 mm ventricular from the mitral plane.

8.3.7 SAPIEN M3 Valve Delivery

1. Orient guide sheath tip towards middle of implanted dock.
2. Insert a pigtail (≥ 120 cm in length) with a preloaded 0.035" stiff guide wire into guide sheath and advance through the middle of the implanted dock. Verify guidewire location using both echocardiography and fluoroscopy.

WARNING: 3D echocardiographic and fluoroscopic (short-axis view) verification must be used to confirm that the guidewire passes through the middle of the dock and has unrestricted movement. Failure to do so can result in the valve being deployed outside of target location, valve embolization, and/or cardiac structure damage requiring surgical repair or other intervention.

3. Remove the pigtail catheter while maintaining guidewire position in the left ventricle.
4. Confirm valve is crimped in the proper orientation prior to insertion.
5. Insert crimped valve and loader fully into the guide sheath.

CAUTION: To prevent possible leaflet damage and possible impact to valve functionality, the valve should not remain in the guide sheath for over 5 minutes.

6. Aspirate and flush guide sheath.
7. Attach 50 cc or larger syringe with > 40 ml heparinized saline. Leave stopcock open to guide sheath.

WARNING: Failure to attach a 50 cc or larger syringe with > 40 ml heparinized saline to the guide sheath and leave the stopcock open to the syringe can introduce air into the left atrium.

8. Advance crimped valve out of loader.
9. Withdraw and peel away loader.
10. Advance the crimped valve until the valve delivery system tip is at the tip of the guide sheath.
11. Verify the valve is correctly positioned using the valve alignment markers.
12. Close the stopcock, and remove the syringe from the guide sheath.
13. Ensure guide sheath is less than 25% flexed.
14. Advance the valve delivery system until half of the valve exits the guide sheath tip.

15. Apply additional flex (up to 60%) and/or torque to the guide sheath to achieve coaxiality with the dock.
16. Fully advance valve out of guide sheath and position into the dock.
Note: There will be a change in force when the valve exits the guide sheath. Advance the valve slowly to ensure that the valve does not interact with the dock. Use care to minimize interaction between the dock and valve during valve positioning.
17. Retract the flex catheter tip to the double markers. Lock the balloon lock. Unflex and retract the guide sheath tip to the double markers.
18. Verify the correct position of the valve with respect to the dock prior to valve deployment.
19. Begin valve deployment:
 - a) Unlock the inflation device provided by Edwards Lifesciences.
 - b) Hold ventilation and begin rapid pacing.
 - c) Deploy the valve using a slow controlled inflation using the entire volume in the inflation device, hold for 3 seconds.
 - d) Deflate the balloon. When the balloon catheter has been completely deflated, turn off the pacemaker and resume ventilation.
20. Prepare for post-dilation:

Add 4 ml heparinized saline with contrast (85:15) to the inflation device or balloon via the 3-way stopcock.
21. Verify the correct position of the balloon with respect to the valve.
22. Begin post-dilation:
 - a) Unlock the inflation device provided by Edwards Lifesciences.
 - b) Hold ventilation and begin rapid pacing (if desired).
 - c) Post-dilate the valve by inflating the balloon with the entire volume in the inflation device, hold for 3 seconds.
 - d) Deflate the balloon. When the balloon catheter has been completely deflated, turn off the pacemaker and resume ventilation (if used).
23. Verify that the flex catheter tip is locked over the double markers and guide sheath is partially unflexed. Retract the valve delivery system into the guide sheath.
24. If the valve delivery system is not immediately removed, aspirate and flush guide sheath with 30 ml of heparinized saline.

8.3.8 SAPIEN M3 System Removal

1. Remove the valve delivery system from the guide sheath.
2. Remove the guide sheath according to the Edwards 23F Guide Sheath IFU.
3. Close the access site per standard of care.

9.0 How Supplied

The SAPIEN M3 valve is supplied sterilized with glutaraldehyde solution. The valve is supplied nonpyrogenic in packaging to which a tamper evident seal has been applied.

The SAPIEN M3 dock steerable catheter, Edwards Commander M delivery system, SAPIEN M3 crimper, and SAPIEN stabilizer rail system are ethylene oxide sterilized.

9.1 Storage

The SAPIEN M3 valve must be stored between 10 °C to 25 °C (50 °F to 77 °F). Each valve is shipped in an enclosure containing a temperature indicator to detect exposure of the valve to extreme temperature.

The SAPIEN M3 dock steerable catheter, SAPIEN stabilizer rail system, SAPIEN M3 Crimper, Edwards Commander M delivery system, and supporting devices should be stored in a cool, dry place.

10.0 Magnetic Resonance (MR) Safety Information

MR Conditional

Non-clinical testing has demonstrated that the SAPIEN M3 dock, with a deployed SAPIEN M3 valve, is MR conditional. A patient can be scanned safely immediately after placement of this valve under the following conditions:

- Static magnetic field of 1.5 Tesla (T) or 3 T.
- Spatial magnetic gradient field of 4000 Gauss/cm (40 T/m) or less.
- Maximum MR system-reported, whole body averaged specific absorption rate (WB-SAR) of 2.0 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the dock and valve are expected to produce a temperature rise of 2 °C or less after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 1.1 cm from the dock and valve when imaged with a gradient echo pulse sequence and a 3.0 T MR system. The lumen of the valve inside the dock was partially to fully obscured in spin and gradient echo images.

11.0 Qualitative and Quantitative Information related to the SAPIEN M3 Valve and SAPIEN M3 Dock

This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

For the valve and dock, the following tables show the qualitative and quantitative information on the materials and substances:

Table 1: SAPIEN M3 Transcatheter Heart Valve

Substance	CAS	Model Mass Range (mg)
Cobalt	7440-48-4	319 - 427

Substance	CAS	Model Mass Range (mg)
Nickel	7440-02-0	361 - 405
Polyethylene terephthalate	25038-59-9	386
Chromium	7440-47-3	208 - 230
Polytetrafluoroethylene	9002-84-0	203
Collagens, bovine, polymers with glutaraldehyde	2370819-60-4	141
Molybdenum	7439-98-7	98.5 - 115
Polyethylene	9002-88-4	31.3
Iron	7439-89-6	0 - 10.9
Titanium	7440-32-6	0 - 10.9
Titanium dioxide	13463-67-7	1.44-2.13
Manganese	7439-96-5	0 - 1.64
Silicon	7440-21-3	0 - 1.64
Antimony trioxide	1309-64-4	0.461
Carbon	7440-44-0	0 - 0.274
Polybutylate	24936-97-8	0.246
Boron	7440-42-8	0 - 0.164
Phosphorus	7723-14-0	0 - 0.164
Sulfur	7704-34-9	0 - 0.109
D&C Green No. 6	128-80-3	0.0356 - 0.0372
Silicon dioxide	7631-86-9	0.00381
Erucamide	112-84-5	0.00370-0.00379
4-Dodecylbenzenesulfonic acid	121-65-3	0.000399

Table 2: SAPIEN M3 Dock

Substance	CAS	Model Mass Range (mg)
Polytetrafluoroethylene	9002-84-0	1010
Nickel	7440-02-0	710 - 742
Titanium	7440-32-6	557 - 593
Polyethylene terephthalate	25038-59-9	79.7
Platinum	7440-06-4	20.3
Polyethylene	9002-88-4	10.7
Iridium	7439-88-5	2.25
Cobalt	7440-48-4	0 - 0.651
Iron	7439-89-6	0 - 0.651
Carbon	7440-44-0	0 - 0.521
Oxygen	7782-44-7	0 - 0.521
Titanium dioxide	13463-67-7	0.195 - 0.466
Niobium	7440-03-1	0 - 0.326
Chromium	7440-47-3	0 - 0.130
Copper	7440-50-8	0 - 0.130
Antimony trioxide	1309-64-4	0.0964
Hydrogen	1333-74-0	0 - 0.0651
Nitrogen	7727-37-9	0 - 0.0651
Erucamide	112-84-5	0.000813 - 0.000830
4-Dodecylbenzenesulfonic acid	121-65-3	0.000875

12.0 Summary of Safety and Clinical Performance (SSCP)

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information.

The Notified Body has taken notice of and agreed with the benefit-risk rationales for the short- and long-term safety and effectiveness of the SAPIEN M3 system.

Conformity of the entire SAPIEN M3 system with the Performance Requirements (GSPR) for safety (MDR GSPR1), performance (MDR GSPR1), acceptability of side-effects (MDR GSPR8), usability (MDR GSPR5), device lifetime (MDR GSPR6), and acceptable benefit-risk profile (MDR GSPR8) has been established for the labelled indications.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

13.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed. The Basic UDI-DI for the SAPIEN M3 Transcatheter Heart Valve, SAPIEN M3 Dock Steerable Catheter, and Edwards Commander M Delivery System can be used to locate the SSCP.

The following tables contain the Basic UDI-DIs:

Product	Model	Basic UDI-DI
SAPIEN M3 Transcatheter Heart Valve	9880TFX29M	0690103D003SM3000TP
SAPIEN M3 Dock Steerable Catheter	9880DDS	0690103D003M3D000L6
Edwards Commander M Delivery System	9880CM29	0690103D003CM3000M7

14.0 Expected Lifetime of the Device

The expected lifetime of the SAPIEN M3 dock and SAPIEN M3 valve is at least 5 years based on rigorous pre-clinical durability testing and successful completion of 5 years of simulated wear testing performed following internationally recognized valve testing standards. Actual lifetime performance depends on multiple biological factors and can vary from patient to patient.

15.0 Patient Information

Patient implant cards are provided with the SAPIEN M3 valve and SAPIEN M3 dock. After implantation, please complete all requested information and provide the implant cards to the patient. The serial numbers are found on the packages. These implant cards allow patients to inform healthcare providers what type of implants they have when they seek care.

16.0 Recovered Valve and Device Disposal

An explanted SAPIEN M3 valve and/or SAPIEN M3 dock should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an explant kit.

The used devices may be disposed of in the same manner that the hospital waste and biohazardous materials are handled. There are no special or unusual risks related to the disposal of the devices. For disposal of the reusable accessories, refer to the Edwards Reusable Accessories IFU.

17.0 Clinical Studies

For clinical study information, refer to the SSCP for clinical outcomes.

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











	English	Français	Deutsch	Español	Italiano
29 mm	For use with size 29 mm Edwards transcatheter heart valve	À utiliser avec une valve cardiaque transcathéter Edwards de 29 mm	Zur Verwendung mit einer 29-mm-Transkatheter-Herzklappe von Edwards	Para su uso con una válvula cardiaca transcatéter Edwards de 29 mm	Da utilizzare con valvole cardiache transcatetere Edwards da 29 mm
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata
	Contents	Contenu	Inhalt	Contenido	Contenuti
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
MD	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Contains biological material of animal origin	Contient des matières biologiques d'origine animale	Enthält biologische Materialien tierischen Ursprungs	Contiene materiales biológicos de origen animal	Contiene materiale biologico di origine animale
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Time & Temperature Sensitive	Sensible au temps et à la température	Zeit- und temperaturempfindlich	Sensible al tiempo y a la temperatura	Sensibile al tempo e alla temperatura
	Contains hazardous substances	Contient des substances dangereuses	Enthält gefährliche Substanzen	Contiene sustancias peligrosas	Contiene sostanze pericolose
SZ	Size	Taille	Größe	Tamaño	Misura
WO	Work Order	Bon de travail	Arbeitsablauf	Orden de trabajo	Ordine di lavoro

	Nederlands	Dansk	Svenska	Ελληνικά	Português
REF	Nabestellingsnummer	Genbestillingsnummer	Artikelnummer	Αριθμός επανάληψης παραγγελίας	Número de referência do produto
#	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Bruikbare lengte	Anvendelig længde	Brukbar längd	Ωφέλιμο μήκος	Comprimento útil
	Niet hergebruiken	Må ikke genanvendes	Får inte återanvändas	Μην επαναχρησιμοποιείτε	Não reutilizar
LOT	Lotnummer	Partinummer	Lotnummer	Αριθμός παρτίδας	Número de lote
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Se bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης	Consultar as instruções de utilização
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Consultar as instruções de utilização no site

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos
■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Buitendiameter	Udvendig diameter	Yttre diameter	Εξωτερική διάμετρος	Diâmetro exterior
	Binnendiameter	Indre diameter	Invändig diameter	Εσωτερική διάμετρος	Diâmetro interno
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό	Manter seco
	Uit het zonlicht houden	Beskyttes mod sollys	Skyddas mot solljus	Φυλάσσετε μακριά από την άμεση ηλιακή ακτινοβολία	Manter afastado da luz solar
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Μοναδικό αναγνωριστικό συσκευής	Identificador único de dispositivo
	Temperatuurlimiet	Temperaturgrænse	Temperaturgräns	Όριο θερμοκρασίας	Limite de temperatura
	Steriel	Steril	Steril	Αποστειρωμένο	Esterilizado
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met vloeibare chemische stof	Steriliseret med flydende kemikalie	Steriliserad med flytande kemikalie	Αποστειρωμένο με χρήση υγρών χημικών	Esterilizado por produto químico líquido
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrieresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Enkelvoudig steriel barrièresysteem met beschermende verpakking aan de binnenzijde	Enkelt sterilt barrieresystem med beskyttende indre emballage	Enkelt sterilt barriärsystem med skyddsförpackning inuti	Σύστημα μονού αποστειρωμένου φραγμού που περιέχει προστατευτική συσκευασία	Sistema de barreira única esterilizada com embalagem de proteção no interior
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Vervaldatum	Sidste anvendelsesdato	Sista förbrukningsdag	Ημερομηνία λήξης	Data de expiração
	Serienummer	Serienummer	Serienummer	Αριθμός σειράς	Número de série
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Compatibiliteit voedraad	Guidewirekompatibilitet	Ledarkompatibilitet	Συμβατότητα οδηγού σύρματος	Compatibilidade do fio-guia
	Nominale druk	Nominelt tryk	Nominellt tryck	Ονομαστική πίεση	Pressão nominal
	Nominale barstdruk	Nominelt sprængningstryk	Beräknat bristningstryck	Ονομαστική πίεση ρήξης	Pressão de rutura nominal
	Aanbevolen voedraadlengte	Anbefalet guidewirelængde	Rekommenderad ledarlängd	Συνιστώμενο μήκος οδηγού σύρματος	Comprimento recomendado do fio-guia
	Minimale afmetingen schacht	Minimal hylsterstørrelse	Minsta hylsstorlek	Ελάχιστο μέγεθος θηκαριού	Tamanho mínimo da bainha
	Afmetingen katheterschacht	Katerrørsstørrelse	Kateterskiftstorlek	Μέγεθος στελέχους καθετήρα	Tamanho do eixo do cateter
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Ballondiameter	Ballondiameter	Ballongdiameter	Διάμετρος μπαλονιού	Diâmetro do balão
	Werklengte ballon	Ballonens arbejdslængde	Ballongens arbetslängd	Ωφέλιμο μήκος μπαλονιού	Comprimento útil do balão
	Voor gebruik met de Edwards transkatheterhartklep van 20 mm	Til brug med størrelse 20 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 20 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 20 mm	Para utilização com a válvula cardíaca transcateter Edwards de 20 mm
	Voor gebruik met de Edwards transkatheterhartklep van 23 mm	Til brug med størrelse 23 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 23 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 23 mm	Para utilização com a válvula cardíaca transcateter Edwards de 23 mm
	Voor gebruik met de Edwards transkatheterhartklep van 26 mm	Til brug med størrelse 26 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 26 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 26 mm	Para utilização com a válvula cardíaca transcateter Edwards de 26 mm
	Voor gebruik met de Edwards transkatheterhartklep van 29 mm	Til brug med størrelse 29 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 29 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 29 mm	Para utilização com a válvula cardíaca transcateter Edwards de 29 mm
	MRI-veilig onder bepaalde voorwaarden	MR-betinget	MR-villkorlig	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις	Utilização condicionada em ambiente de RM
	Inhoud	Indhold	Innehåll	Περιεχόμενο	Conteúdo
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogênico
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Bevat biologisch materiaal van dierlijke oorsprong	Indeholder biologisk materiale af animalsk oprindelse	Innehåller biologiskt material av animaliskt ursprung	Περιέχει βιολογικό υλικό ζωικής προέλευσης	Contém material biológico de origem animal
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Gevoelig voor tijd en temperatuur	Følsom overfor tid og temperatur	Tids- och temperaturkänslig	Ευαίσθητο ως προς τον χρόνο και τη θερμοκρασία	Sensível à passagem do tempo e à temperatura
	Bevat gevaarlijke stoffen	Indeholder farlige stoffer	Innehåller farliga ämnen	Περιέχει επικίνδυνες ουσίες	Contém substâncias perigosas
SZ	Afmetingen	Størrelse	Storlek	Μέγεθος	Tamanho
WO	Werkorder	Arbejdsorde	Arbetsorder	Εντολή εργασίας	Ordem de trabalho

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
Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

2025-03
10061250001 B
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU