Edwards

Edwards 23F Guide Sheath

Directory	
English (en).....	1
Français (fr).....	5
Deutsch (de).....	8
Español (es).....	11
Italiano (it).....	14
Nederlands (nl).....	17
Dansk (da).....	20
Svenska (sv).....	23
Ελληνικά (el).....	26
Português (pt).....	29
Česky (cs).....	32
Polski (pl).....	35
Norsk (no).....	38
Suomi (fi).....	41
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset.....	44

English

Edwards 23F Guide Sheath

Instructions For Use

The product should only be used by physicians trained and experienced in interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.

1.0 Device Description

The Edwards 23F guide sheath (also known as guide sheath and shown in Figure 1) includes an articulating hydrophilic coated guide sheath and a hydrophilic coated introducer. The guide sheath has a radiopaque soft tip and a flex wheel which flexes the guide sheath towards the flush port. The introducer is compatible with a 0.035 inch (0.89 mm) guidewire.


Figure 1: Edwards 23F Guide Sheath

Edwards 23F Guide Sheath	Model: 9880GS
Sheath I.D.	23F (7.6 mm)
Sheath O.D.	29F (9.6 mm)
Sheath Effective Length	77 cm
Introducer I.D.	0.89 mm (0.035 in)
Introducer O.D.	24F (7.9 mm)
Introducer Effective Length	103 cm

2.0 Intended User and Intended Use / Intended Purpose

2.1 Intended Users

The intended users of the guide sheath include two user groups:

One intended user group is responsible for device preparation and includes operating room nurses or cardiovascular technicians who have received Edwards Lifesciences training, or certified Edwards field representatives.

The second intended user group is responsible for procedural device use and includes interventional cardiologists or cardiovascular surgeons.

2.2 Intended Use / Intended Purpose

The Edwards 23F guide sheath is intended to provide access to the venous vasculature.

3.0 Intended Patient Population and Indication for Use

3.1 Intended Patient Population

The Edwards 23F guide sheath is intended for adult patients undergoing a SAPIEN M3 transcatheter mitral valve replacement procedure.

3.2 Indication for Use

The Edwards 23F guide sheath is indicated to provide venous vascular access to cardiac structures enabling the introduction and removal of SAPIEN M3 transcatheter mitral valve replacement devices.

4.0 Contraindications

There are no known contraindications.

5.0 Warnings

- The devices are designed, intended, and distributed for single use only. Do not resterilize or reuse the devices. There are no data to support the sterility, non-pyrogenicity, and functionality of the devices after reprocessing.
- Do not mishandle the device or use it if the packaging or any components are not sterile, have been opened or are damaged (i.e., kinked or stretched, etc.), or the expiration date has elapsed.
- Procedures should be conducted under echocardiographic and fluoroscopic guidance. Some fluoroscopically guided procedures are associated with a risk of radiation injury to the skin. These injuries may be painful, disfiguring, and long-lasting.
- The minimum ID of the guide sheath is 23F (7.6 mm). Characteristics of the device(s) to be inserted into the guide sheath should be evaluated to prevent damage to the interior liner of the guide sheath, damage to the device(s) being inserted, and/or injury to the patient.
- In the event of device malfunction or device damage during use (e.g. destructive deformation to the catheter) safely remove the device(s). If unable to safely remove the device(s), conversion to surgery is recommended.
- Patient injury could occur if the guide sheath is not unflexed prior to removal.

6.0 Precautions

- Abnormalities in the caval vein, the presence of an atrial septal occluder device, or presence of calcium may preclude safe transseptal, transvenous femoral access.
- Use caution in tortuous or calcified vessels that would prevent safe entry of the guide sheath and introducer.
- The safety and effectiveness of the Edwards 23F guide sheath have not been established for patients who have/are:
 - A known hypersensitivity or contraindication to heparin
 - Female who is pregnant or lactating
 - Children

7.0 Potential Adverse Events

The following potential risks are associated with the device usage including potential access complications associated with standard cardiac catheterization, the potential risks of anesthesia, and the use of angiography.

- Death
- Stroke/TIA or nerve injury
- Cardiovascular injury – cardiac structure complications
- Cardiovascular injury – vascular complications
- Cardiovascular injury – access related complications
- Cardiac arrest
- Pericardial effusion or cardiac tamponade
- Thromboembolism including air, calcific valve material, or thrombus
- Arrhythmia
- Bleeding / Hematoma / Hemorrhage that may require transfusion or intervention
- Pleural effusion
- Emergency cardiac surgery

- Myocardial infarction
- Infection including septicemia and endocarditis
- Allergic reaction to anesthesia, contrast media, or device materials
- Deterioration of native valve (leaflet tear/tearing, leaflet retraction, or other)
- Atrial septal defect that may require intervention
- Conduction system defect which may require a permanent pacemaker
- Skin burn
- Vessel spasm
- Catheter entrapment
- Fever
- Inflammation
- Pain or changes at the access site

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

8.0 Directions for Use

8.1 Device Handling and Preparation

1. Verify expiration date, model number, and visually inspect for breaches in package integrity prior to opening sterile package.
2. Visually inspect guide sheath and introducer for damage.
3. While keeping distal tip raised, flush guide sheath with heparinized saline.
4. Hydrate the length of the introducer with heparinized saline.
5. Insert introducer into guide sheath partially.
6. While keeping distal tip raised, flush guide sheath with heparinized saline.
7. Advance introducer and twist to lock to guide sheath.
8. Flush introducer with heparinized saline. Hydrate the length of the guide sheath with heparinized saline.

8.2 Device Use

CAUTION: Excessive device manipulation may result in cardiac structure damage requiring surgical repair or other intervention.

1. Access the common femoral vein using conventional percutaneous puncture methods.
2. Access the left atrium via transseptal puncture using conventional percutaneous methods and place a guidewire in the left atrium.

CAUTION: Inappropriate puncture may result in cardiac structure damage, requiring surgical repair or other intervention.

3. Administer heparin to maintain the ACT at ≥ 300 sec.
4. Insert the guide sheath and introducer with the flushport oriented away from the operator. Advance until guide sheath tip is at the desired location.
5. Unlock introducer from guide sheath and slowly retract the introducer (and guidewire if applicable).
6. With the introducer and wire remaining across the guide sheath seals, aspirate and flush the guide sheath. Remove the introducer (and guidewire if applicable) slowly.
7. Insert the device(s) into the sheath. When inserting / removing device(s), aspirate and flush the guide sheath.

CAUTION: Ensure a 15F or larger catheter/device is across the guide sheath seals when aspirating and flushing the guide sheath to reduce the risk of air embolization.

8. After the completion of the procedure and removal of device(s), fully unflex and retract the guide sheath. Retract guide sheath into the right atrium and assess the residual atrial septal defect.

CAUTION: Patient Injury could occur if the guide sheath is not unflexed prior to removal.

9. Remove the guide sheath without torquing.
10. Close the access site per standard of care.

9.0 How Supplied

The Edwards 23F Guide Sheath is sterilized with ethylene oxide.

10.0 Storage

The Edwards 23F Guide Sheath should be stored in a cool, dry place.

11.0 Device Disposal

Used devices may be handled and disposed of in the same manner as hospital waste and biohazardous materials in accordance with local regulations. There are no special risks related to the disposal of these devices.

12.0 Hazardous Substances

This device does not contain hazardous substances.

13.0 Summary of Safety and Clinical Performance (SSCP)

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information.

The Notified Body has taken notice of and agreed with the benefit-risk rationales for the safety and effectiveness of the Edwards 23F guide sheath.

Conformity of the Edwards 23F guide sheath with the Performance Requirements (GSPR) for safety (MDR GSPR1), performance (MDR GSPR1), acceptability of side-effects (MDR GSPR8), usability (MDR GSPR5), device lifetime (MDR GSPR6), and acceptable benefit-risk profile (MDR GSPR8) has been established for the labelled indications.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

14.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed. The Basic UDI-DI for the Edwards 23F Guide Sheath can be used to locate the SSCP.

The following tables contain the Basic UDI-DIs:

Product	Model	Basic UDI-DI
Edwards 23F Guide Sheath	9880GS	0690103D003EGS000SQ

15.0 Clinical Studies

For clinical study information, refer to the SSCP for clinical outcomes.

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





	English	Français	Deutsch	Español	Italiano
29 mm	For use with size 29 mm Edwards transcatheter heart valve	À utiliser avec une valve cardiaque transcathéter Edwards de 29 mm	Zur Verwendung mit einer 29-mm-Transkatheter-Herzklappe von Edwards	Para su uso con una válvula cardíaca transcatéter Edwards de 29 mm	Da utilizzare con valvole cardiache transcaterere Edwards da 29 mm
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata
	Contents	Contenu	Inhalt	Contenido	Contenuti
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
MD	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Contains biological material of animal origin	Contient des matières biologiques d'origine animale	Enthält biologische Materialien tierischen Ursprungs	Contiene materiales biológicos de origen animal	Contiene materiale biologico di origine animale
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Time & Temperature Sensitive	Sensible au temps et à la température	Zeit- und temperaturempfindlich	Sensible al tiempo y a la temperatura	Sensibile al tempo e alla temperatura
	Contains hazardous substances	Contient des substances dangereuses	Enthält gefährliche Substanzen	Contiene sustancias peligrosas	Contiene sostanze pericolose
SZ	Size	Taille	Größe	Tamaño	Misura
WO	Work Order	Bon de travail	Arbeitsablauf	Orden de trabajo	Ordine di lavoro

	Nederlands	Dansk	Svenska	Ελληνικά	Português
REF	Nabestellingsnummer	Genbestillingsnummer	Artikelnummer	Αριθμός επανάληψης παραγγελίας	Número de referência do produto
#	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Bruikbare lengte	Anvendelig længde	Brukbar längd	Ωφέλιμο μήκος	Comprimento útil
	Niet opnieuw gebruiken	Må ikke genbruges	Får inte återanvändas	Μην επαναχρησιμοποιείτε	Não reutilizar
LOT	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Se bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης	Consultar as instruções de utilização
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Consultar as instruções de utilização no site

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos
■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Buitendiameter	Udvendig diameter	Yttre diameter	Εξωτερική διάμετρος	Diâmetro exterior
	Binnendiameter	Indre diameter	Invändig diameter	Εσωτερική διάμετρος	Diâmetro interno
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό	Manter seco
	Uit het zonlicht houden	Beskyttes mod sollys	Skyddas mot solljus	Φυλάσσετε μακριά από την άμεση ηλιακή ακτινοβολία	Manter afastado da luz solar
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Μοναδικό αναγνωριστικό συσκευής	Identificador único de dispositivo
	Temperatuurlimiet	Temperaturgrænse	Temperaturgräns	Όριο θερμοκρασίας	Limite de temperatura
	Steriel	Steril	Steril	Αποστειρωμένο	Esterilizado
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Gesteriliseerd met vloeibare chemische stof	Steriliseret med flydende kemikalie	Steriliserad med flytande kemikalie	Αποστειρωμένο με χρήση υγρών χημικών	Esterilizado por produto químico líquido
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrieresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Enkelvoudig steriel barrièresysteem met beschermende verpakking aan de binnenzijde	Enkelt sterilt barrieresystem med beskyttende indre emballage	Enkelt sterilt barriärsystem med skyddsförpackning inuti	Σύστημα μονού αποστειρωμένου φραγμού που περιέχει προστατευτική συσκευασία	Sistema de barreira única esterilizada com embalagem de proteção no interior
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Vervaldatum	Sidste anvendelsesdato	Sista förbrukningsdag	Ημερομηνία λήξης	Data de expiração
	Serienummer	Serienummer	Serienummer	Αριθμός σειράς	Número de série
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos
■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Compatibiliteit voedraad	Guidewirekompatibilitet	Ledarkompatibilitet	Συμβατότητα οδηγού σύρματος	Compatibilidade do fio-guia
	Nominale druk	Nominelt tryk	Nominellt tryck	Ονομαστική πίεση	Pressão nominal
	Nominale barstdruk	Nominelt sprængningstryk	Beräknat bristningstryck	Ονομαστική πίεση ρήξης	Pressão de rutura nominal
	Aanbevolen voedraadlengte	Anbefalet guidewirelængde	Rekommenderad ledarlängd	Συνιστώμενο μήκος οδηγού σύρματος	Comprimento recomendado do fio-guia
	Minimale afmetingen schacht	Minimal hylsterstørrelse	Minsta hylsstorlek	Ελάχιστο μέγεθος θηκαριού	Tamanho mínimo da bainha
	Afmetingen katheterschacht	Katerrørsstørrelse	Kateterskafstorlek	Μέγεθος στελέχους καθετήρα	Tamanho do eixo do cateter
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Ballondiameter	Ballondiameter	Ballongdiameter	Διάμετρος μπαλονιού	Diâmetro do balão
	Werklengte ballon	Ballonens arbejdslængde	Ballongens arbetslängd	Ωφέλιμο μήκος μπαλονιού	Comprimento útil do balão
	Voor gebruik met de Edwards transkatheterhartklep van 20 mm	Til brug med størrelse 20 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 20 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 20 mm	Para utilização com a válvula cardíaca transcateter Edwards de 20 mm
	Voor gebruik met de Edwards transkatheterhartklep van 23 mm	Til brug med størrelse 23 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 23 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 23 mm	Para utilização com a válvula cardíaca transcateter Edwards de 23 mm
	Voor gebruik met de Edwards transkatheterhartklep van 26 mm	Til brug med størrelse 26 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 26 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 26 mm	Para utilização com a válvula cardíaca transcateter Edwards de 26 mm
	Voor gebruik met de Edwards transkatheterhartklep van 29 mm	Til brug med størrelse 29 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 29 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 29 mm	Para utilização com a válvula cardíaca transcateter Edwards de 29 mm
	MRI-veilig onder bepaalde voorwaarden	MR-betinget	MR-villkorlig	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις	Utilização condicionada em ambiente de RM
	Inhoud	Indhold	Innehåll	Περιεχόμενο	Conteúdo
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogênico
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Bevat biologisch materiaal van dierlijke oorsprong	Indeholder biologisk materiale af animalsk oprindelse	Innehåller biologiskt material av animaliskt ursprung	Περιέχει βιολογικό υλικό ζωικής προέλευσης	Contém material biológico de origem animal
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Gevoelig voor tijd en temperatuur	Følsom overfor tid og temperatur	Tids- och temperaturkänslig	Ευαίσθητο ως προς τον χρόνο και τη θερμοκρασία	Sensível à passagem do tempo e à temperatura
	Bevat gevaarlijke stoffen	Indeholder farlige stoffer	Innehåller farliga ämnen	Περιέχει επικίνδυνες ουσίες	Contém substâncias perigosas
SZ	Afmetingen	Størrelse	Storlek	Μέγεθος	Tamanho
WO	Werkorder	Arbejdsorde	Arbetsorder	Εντολή εργασίας	Ordem de trabalho

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
Edwards

Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

2025-01
10061251001 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU