Edwards

Edwards Table

Directory	
English (en).....	1
Français (fr).....	3
Deutsch (de).....	5
Español (es).....	7
Italiano (it).....	9
Nederlands (nl).....	11
Dansk (da).....	13
Svenska (sv).....	15
Ελληνικά (el).....	17
Português (pt).....	19
Česky (cs).....	21
Polski (pl).....	23
Norsk (no).....	25
Suomi (fi).....	27
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Rysunki ■ Figurer ■ Kuvat.....	29
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Legenda symboli ■ Symbolforklaring ■ Merkkien selitykset.....	30

English

Edwards Table

Instructions for Use

1.0 Device Description

Edwards Table (Figure 1)

The Edwards Table (also known as table) is a single-use, height-adjustable accessory device to aid the positioning and stabilization of interventional catheter devices and/or accessories. The table is non-sterile, non-patient contacting and used outside the sterile field.

Edwards single-use accessory:

Product Name	Model/REF
Edwards Table	9880T

2.0 Intended User and Intended Use / Intended Purpose

2.1 Intended Users

The Edwards table is intended to be used by physicians and support staff trained to these instructions for use.

2.2 Intended Use / Intended Purpose

The Edwards table is intended to aid the positioning and stabilization of interventional catheter devices and/or accessories.

3.0 Intended Patient Population

The Edwards table is intended for patients undergoing an interventional procedure.

4.0 Contraindications

- There are no known contraindications.

5.0 Warnings

- The Edwards table is non-sterile. Introduction of the device into the sterile field may cause infection.
- Do not mishandle the Edwards table or use it if the packaging has been opened or the product is damaged, or the expiration date has elapsed.

6.0 Precautions

- Before use, visually inspect the device to verify that no damage has occurred during shipping and handling.
- This device has not been used or tested to support any devices other than interventional catheter devices and/or accessories from Edwards Lifesciences.

7.0 Potential Adverse Events

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

8.0 Directions for Use

8.1 Prepare the Device - Edwards Table Assembly

1. Verify expiration date and model number prior to opening package.
2. Visually inspect the table for damage.
3. Attach the extension pieces to the bases as necessary.
4. Attach the bases to the table.

8.2 Device Use and Patient Preparation

1. Prior to sterile draping the patient for standard transseptal catheterization via femoral vein access, position the table between the legs of the patient, adjusting the height of the table as needed.

WARNING: The table is provided non-sterile; introduction of the table into the sterile field may result in infection.

2. For use with interventional catheter devices provided by Edwards, please refer to individual system IFUs provided by Edwards.

9.0 How Supplied

The Edwards table is packaged and provided non-sterile.

10.0 Storage

The Edwards table should be stored in a cool, dry place.

11.0 Device Disposal

After use, the table may be a potential biohazard. The table should be disposed of in the same manner that the hospital waste and biohazardous materials are handled. There are no special or unusual risks related to the disposal of the table.

12.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed system. The following table contains the basic UDI-DI for the Edwards table:

Product	Model	Basic UDI-DI
Edwards Table	9880T	0690103D004PAC000S6

13.0 Hazardous Substances

This device does not contain hazardous substances.

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


	English	Français	Deutsch	Español	Italiano
29 mm	For use with size 29 mm Edwards transcatheter heart valve	À utiliser avec une valve cardiaque transcathéter Edwards de 29 mm	Zur Verwendung mit einer 29-mm-Transkatheter-Herzklappe von Edwards	Para su uso con una válvula cardíaca transcatéter Edwards de 29 mm	Da utilizzare con valvole cardiache transcaterere Edwards da 29 mm
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata
	Contents	Contenu	Inhalt	Contenido	Contenuti
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
MD	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Contains biological material of animal origin	Contient des matières biologiques d'origine animale	Enthält biologische Materialien tierischen Ursprungs	Contiene materiales biológicos de origen animal	Contiene materiale biologico di origine animale
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Time & Temperature Sensitive	Sensible au temps et à la température	Zeit- und temperaturempfindlich	Sensible al tiempo y a la temperatura	Sensibile al tempo e alla temperatura
	Contains hazardous substances	Contient des substances dangereuses	Enthält gefährliche Substanzen	Contiene sustancias peligrosas	Contiene sostanze pericolose
SZ	Size	Taille	Größe	Tamaño	Misura
WO	Work Order	Bon de travail	Arbeitsablauf	Orden de trabajo	Ordine di lavoro

	Nederlands	Dansk	Svenska	Ελληνικά	Português
REF	Nabestellingsnummer	Genbestillingsnummer	Artikelnummer	Αριθμός επανάληψης παραγγελίας	Número de referência do produto
#	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Bruikbare lengte	Anvendelig længde	Brukbar längd	Ωφέλιμο μήκος	Comprimento útil
	Niet hergebruiken	Må ikke genanvendes	Får inte återanvändas	Μην επαναχρησιμοποιείτε	Não reutilizar
LOT	Lotnummer	Partinummer	Lotnummer	Αριθμός παρτίδας	Número de lote
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Se bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης	Consultar as instruções de utilização
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Consultar as instruções de utilização no site

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos
■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Buitendiameter	Udvendig diameter	Yttre diameter	Εξωτερική διάμετρος	Diâmetro exterior
	Binnendiameter	Indre diameter	Invändig diameter	Εσωτερική διάμετρος	Diâmetro interno
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό	Manter seco
	Uit het zonlicht houden	Beskyttes mod sollys	Skyddas mot solljus	Φυλάσσετε μακριά από την άμεση ηλιακή ακτινοβολία	Manter afastado da luz solar
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Μοναδικό αναγνωριστικό συσκευής	Identificador único de dispositivo
	Temperatuurlimiet	Temperaturgrænse	Temperaturgräns	Όριο θερμοκρασίας	Limite de temperatura
	Steriel	Steril	Steril	Αποστειρωμένο	Esterilizado
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met vloeibare chemische stof	Steriliseret med flydende kemikalie	Steriliserad med flytande kemikalie	Αποστειρωμένο με χρήση υγρών χημικών	Esterilizado por produto químico líquido
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barriersystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Enkelvoudig steriel barrièresysteem met beschermende verpakking aan de binnenzijde	Enkelt sterilt barriersystem med beskyttende indre emballage	Enkelt sterilt barriärsystem med skyddsförpackning inuti	Σύστημα μονού αποστειρωμένου φραγμού που περιέχει προστατευτική συσκευασία	Sistema de barreira única esterilizada com embalagem de proteção no interior
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Vervaldatum	Sidste anvendelsesdato	Sista förbrukningsdag	Ημερομηνία λήξης	Data de expiração
	Serienummer	Serienummer	Serienummer	Αριθμός σειράς	Número de série
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli
■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos
■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Compatibiliteit voedraad	Guidewirekompatibilitet	Ledarkompatibilitet	Συμβατότητα οδηγού σύρματος	Compatibilidade do fio-guia
	Nominale druk	Nominelt tryk	Nominellt tryck	Ονομαστική πίεση	Pressão nominal
	Nominale barstdruk	Nominelt sprængningstryk	Beräknat bristningstryck	Ονομαστική πίεση ρήξης	Pressão de rutura nominal
	Aanbevolen voedraadlengte	Anbefalet guidewirelængde	Rekommenderad ledarlängd	Συνιστώμενο μήκος οδηγού σύρματος	Comprimento recomendado do fio-guia
	Minimale afmetingen schacht	Minimal hylsterstørrelse	Minsta hylsstorlek	Ελάχιστο μέγεθος θηκαριού	Tamanho mínimo da bainha
	Afmetingen katheterschacht	Katerrørsstørrelse	Kateterskiftstorlek	Μέγεθος στελέχους καθετήρα	Tamanho do eixo do cateter
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Ballondiameter	Ballondiameter	Ballongdiameter	Διάμετρος μπαλονιού	Diâmetro do balão
	Werklengte ballon	Ballonens arbejdslængde	Ballongens arbetslängd	Ωφέλιμο μήκος μπαλονιού	Comprimento útil do balão
	Voor gebruik met de Edwards transkatheterhartklep van 20 mm	Til brug med størrelse 20 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 20 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 20 mm	Para utilização com a válvula cardíaca transcateter Edwards de 20 mm
	Voor gebruik met de Edwards transkatheterhartklep van 23 mm	Til brug med størrelse 23 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 23 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 23 mm	Para utilização com a válvula cardíaca transcateter Edwards de 23 mm
	Voor gebruik met de Edwards transkatheterhartklep van 26 mm	Til brug med størrelse 26 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 26 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 26 mm	Para utilização com a válvula cardíaca transcateter Edwards de 26 mm
	Voor gebruik met de Edwards transkatheterhartklep van 29 mm	Til brug med størrelse 29 mm Edwards -transkateterherteklap	För användning med Edwards kateterburna hjärtklaff, storlek 29 mm	Για χρήση με διακαθετηριακή καρδιακή βαλβίδα Edwards μεγέθους 29 mm	Para utilização com a válvula cardíaca transcateter Edwards de 29 mm
	MRI-veilig onder bepaalde voorwaarden	MR-betinget	MR-villkorlig	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις	Utilização condicionada em ambiente de RM
	Inhoud	Indhold	Innehåll	Περιεχόμενο	Conteúdo
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogênico
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Bevat biologisch materiaal van dierlijke oorsprong	Indeholder biologisk materiale af animalsk oprindelse	Innehåller biologiskt material av animaliskt ursprung	Περιέχει βιολογικό υλικό ζωικής προέλευσης	Contém material biológico de origem animal
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Legenda simboli ■ Symbolforklaring ■ Merkkien selitykset


	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Gevoelig voor tijd en temperatuur	Følsom overfor tid og temperatur	Tids- och temperaturkänslig	Ευαίσθητο ως προς τον χρόνο και τη θερμοκρασία	Sensível à passagem do tempo e à temperatura
	Bevat gevaarlijke stoffen	Indeholder farlige stoffer	Innehåller farliga ämnen	Περιέχει επικίνδυνες ουσίες	Contém substâncias perigosas
	Afmetingen	Størrelse	Storlek	Μέγεθος	Tamanho
	Werkorder	Arbejdsorde	Arbetsorder	Εντολή εργασίας	Ordem de trabalho

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
2025-01
10061252001 A
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU

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Edwards Lifesciences GmbH