Edwards

Edwards EVOQUE Tricuspid Valve Replacement System

Directory	
English (en).....	1
Français (fr).....	8
Deutsch (de).....	15
Español (es).....	22
Italiano (it).....	29
Nederlands (nl).....	36
Dansk (da).....	43
Svenska (sv).....	50
Ελληνικά (el).....	57
Português (pt).....	64
Česky (cs).....	71
Magyar (hu).....	78
Polski (pl).....	85
Slovenský (sk).....	92
Norsk (no).....	99
Suomi (fi).....	106
Български (bg).....	113
Română (ro).....	120
Eesti (et).....	127
Lietuvių (lt).....	134
Latviešu (lv).....	141
Türkçe (tr).....	148
Русский (ru).....	155
Srpski (sr).....	162
Hrvatski (hr).....	169
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike.....	176
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sūbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	179

English

Edwards EVOQUE Tricuspid Valve Replacement System

Instructions for Use

Caution: Implantation of the transcatheter tricuspid valve should only be performed by physicians who have received training for the Edwards EVOQUE tricuspid valve replacement system.

- Do not attempt to use the Edwards EVOQUE valve (herein referred to as the EVOQUE valve), delivery system, or associated accessories before completely reading and understanding the information contained in this booklet.
- Failure to follow these instructions, warnings, and precautions may lead to device damage or patient injury. Use of the EVOQUE system should be restricted to those physicians trained to perform invasive endovascular procedures and to those physicians trained in the proper use of the system.
- Consult with authorized Edwards personnel for appropriate selection of EVOQUE valve size.
- The EVOQUE valve, delivery system, dilator kit, loading system, and stabilizer are supplied STERILE. The single-use accessories (base and plate) and reusable accessories (rail, platform, and plate) are supplied NON-STERILE.
- For the Edwards reusable platform and plate, refer to the reprocessing instructions in the Edwards Reusable Accessories IFU for wipe down steps. For the Edwards reusable rail assembly, refer to the reprocessing instructions in the Edwards Reusable Rail Assembly IFU. All other devices are supplied for single-use only. After use, dispose in accordance with hospital administrative and/or government policy.

1.0 Introduction

The Edwards EVOQUE tricuspid valve replacement system (herein referred to as the EVOQUE system) is designed to replace the native tricuspid valve in patients with symptomatic, severe tricuspid valve regurgitation without the need for conventional open-heart surgery. The EVOQUE system consists of four (4) elements which are used with either of the two (2) accessories sets, as outlined below:

Product Name	44 mm	48 mm	52 mm	56 mm
Edwards EVOQUE Valve (EVOQUE Valve)	9850EV44	9850EV48	9850EV52	9850EV56
Edwards EVOQUE Loading System or Edwards EVOQUE Loading System and Trimmer		9850LS or 9850LSB		9850LSB
Edwards EVOQUE Tricuspid Delivery System		9850TDS		
Edwards EVOQUE Dilator Kit		9850DK		

Accessories Set 1	
Product Name	Model Number
Edwards EVOQUE Stabilizer	9850SB
Edwards EVOQUE Stabilizer Base	9850BA
Edwards EVOQUE Stabilizer Plate	9850PT

Accessories Set 2	
Product Name	Model Number
Edwards EVOQUE Stabilizer	9850SZ
Edwards Reusable Rail	10500RL
Edwards Reusable Platform	10000UP
Edwards Reusable Plate	10000PT

Note: The EVOQUE system is designed for use with either Accessories set 1 or Accessories set 2. Accessories set 1 is considered an optional and single-use set for use with the EVOQUE system. Components from set 1 and set 2 are not interchangeable.

Note: The devices that comprise the EVOQUE system are intended to be used together and are not intended to be used as standalone devices. The information in this IFU is applicable to all the EVOQUE system devices.

2.0 Device Description

• Edwards EVOQUE Valve (Figure 1)

The EVOQUE valve is intended to be implanted in the tricuspid position as an alternative to surgery. The EVOQUE valve consists of a trileaflet bovine pericardial tissue valve, nitinol frame, and fabric skirt, and is packaged and terminally sterilized in glutaraldehyde.

Valve size recommendations are based on native valve annulus size, as measured by computed tomography (CT). Patient anatomical factors and imaging modalities should be considered during valve size selection.

Device Diameter (Recommended Valve Size)	Systole		Diastole	
	Recommended Treatable Perimeter - Derived Diameter Range (mm)	Maximum Treatable Annulus Length (mm)	Recommended Treatable Perimeter - Derived Diameter Range (mm)	Maximum Treatable Annulus Length (mm)
44	36.5 - 43	45.5	39.6 - 45.5	50
48	40 - 47	49.5	43.2 - 49.5	54
52	45 - 51	53.5	46.8 - 53.5	58
56	48 - 55	57.5	50.4 - 57.5	62

• Edwards EVOQUE Tricuspid Delivery System (Figure 2)

The delivery system has an outer diameter of 28 F and is intended to deliver the EVOQUE valve in the crimped position via the transfemoral venous approach. The delivery system handle contains a primary flex knob, secondary flex knob, and depth knob to facilitate EVOQUE valve alignment and positioning in the native valve, and a capsule knob and release knob to control the expansion and release of the EVOQUE valve.

• Edwards EVOQUE Dilator Kit (Figure 3)

The 24 F, 28 F, and 33 F diameter hydrophilic-coated dilators are intended to dilate the access site, facilitating delivery system insertion. All dilators accommodate a 0.035 inch (0.89 mm) guidewire and are tapered to minimize access site trauma.

• Edwards EVOQUE Loading System (Figure 4) / Loading System and Trimmer (Figure 5)

The loading system, which consists of multiple components, is intended to facilitate loading and attachment of the EVOQUE valve onto the delivery system. The loading system assists in crimping the EVOQUE valve to the appropriate diameter, which allows the outer capsule to advance over the EVOQUE valve.

• Edwards EVOQUE System Accessories Set 1: Stabilizer (Model 9850SB), Base, and Plate (Figure 6)

The single-use stabilizer with single-use base and plate are intended to secure the delivery system at an angle appropriate for the transfemoral venous approach and to enable fine adjustments of the position of the delivery system during the implantation procedure. The base is height-adjustable to accommodate patient lower extremities and is intended to provide a stable base for the stabilizer. The plate is intended to provide a stable, flat surface for the base on the operating table.

• Edwards EVOQUE System Accessories Set 2: Stabilizer (Model 9850SZ) and Reusable Rail, Platform, and Plate (Figure 7)

The reusable accessories compatible with the EVOQUE system include the reusable rail, platform, and plate. The reusable rail is a stable, sterilizable rail that provides a surface to secure compatible stabilizer to the reusable platform. The reusable platform allows for positioning and stabilization of Edwards transcatheter replacement delivery systems through adjustable leg height and a stable platform. The reusable plate is placed below the reusable platform, providing a flat and stable surface during use. The single-use stabilizer is intended to aid with positioning and stabilization of the EVOQUE system during implantation procedures. For the Accessories set utilizing reusable accessories (Figure 7), refer to the Edwards Reusable Accessories IFU for set up of the base platform and plate. If reusable rail is non-sterile or if unpackaging for the first time, refer to the Edwards Reusable Rail Assembly IFU.

3.0 Facilities and Training Requirement

Facilities that intend to perform an implantation procedure utilizing the EVOQUE system must have access to cine fluoroscopy and transesophageal echocardiography (TEE) throughout the procedure. In addition, the implanting physicians must have prompt access to facilities with the necessary equipment, instruments, supplies, and personnel to perform emergency tricuspid valve surgery, if required.

A comprehensive training program is provided by Edwards Lifesciences and must be completed by the implanting physicians before use of the EVOQUE system. The implanting physicians should have advanced technical knowledge and experience in related catheter-based procedures.

4.0 Indications for Use

4.1 Intended Use

The EVOQUE tricuspid valve replacement system is intended for the reduction of tricuspid regurgitation (TR) for patients who remain symptomatic on medical therapy and in whom no other surgical or transcatheter treatment option exists per heart team decision. The EVOQUE system percutaneously delivers the Implant to the tricuspid valve via femoral vein access using a transvenous approach.

4.2 Target Patient Population

- Patients with symptomatic, severe TR who do not have any of the following:
 - Severe aortic, mitral and/or pulmonic valve stenosis and/or regurgitation
 - LVEF < 25%
 - Severe right ventricular dysfunction as measured by echocardiography
 - Pacemaker dependency¹ on trans-tricuspid lead without alternative pacing option
 - Need for emergent or urgent surgery for any reason or any planned cardiac surgery within the next 12 months
 - Resting systolic blood pressure < 90 or > 160 mmHg after repeated measurements
 - Refractory heart failure requiring advanced intervention (i.e., left ventricular assist device, transplantation) (ACC/AHA/ESC/EACTS Stage D heart failure)
 - An occluded or thrombosed IVC filter that would interfere with the delivery catheter
- Patient selection should be performed by a multi-disciplinary heart team specializing in the treatment of tricuspid regurgitation to assess patient risk and anatomical suitability.

¹ Pacemaker dependency may increase the risk of serious injury or death from sudden pacemaker failure

5.0 Contraindications

The EVOQUE valve is contraindicated for patients with any of the following conditions:

- Active endocarditis within 3 months or infection requiring antibiotic therapy (oral or intravenous) within 2 weeks of the scheduled implant.
- Untreatable hypersensitivity or contraindication to any of the following: all antiplatelets, all anticoagulants, nitinol alloys (nickel and titanium), bovine tissue, glutaraldehyde, contrast media, or transesophageal echocardiography.
- Tricuspid valve anatomy that precludes proper device deployment and functionality based on CT and echocardiographic evaluation.

Note: Patient must be able to tolerate at least one antiplatelet medication AND one anticoagulant medication.

6.0 Warnings

- Do not use the EVOQUE valve outside of the labeled indication.
- The EVOQUE valve, delivery system, loading system, dilator kit, and stabilizer are designed, intended, and distributed for STERILE single-use only. The EVOQUE stabilizer base and EVOQUE stabilizer plate are provided nonsterile for single-use only. Do not resterilize or reuse any of the single-use devices. There are no data to support the sterility, nonpyrogenicity, or functionality of the devices after reprocessing.
- Ensure proper sterile techniques are utilized during the preparation, transfer, and use of the devices.

- Do not use the valve if the tamper evident seal is broken, the storage solution does not completely cover the valve, the temperature indicator has been activated, the valve is damaged, or the expiration date has elapsed. The EVOQUE valve must remain hydrated at all times. The valve cannot be exposed to solutions, antibiotics, or chemicals other than its shipping storage solution and sterile physiologic saline solution. This will help prevent leaflet damage that may impact valve functionality. Keep the EVOQUE valve hydrated with normal saline until ready for implantation.
- Ensure the correct valve size is selected. Implantation of the improper size (i.e., undersizing or oversizing) may lead to paravalvular leak (PVL), migration, embolization, and/or annular damage.
- Patients with previously-implanted devices (e.g., IVC filter) should be carefully assessed prior to insertion of the delivery system to avoid potential damage to vasculature or a previously-implanted device.
- Patients with pre-existing cardiac leads should be carefully assessed prior to implantation to avoid potential adverse interaction between devices.
- Care should be taken when implanting cardiac leads after EVOQUE valve implantation to avoid potential adverse interaction between the devices.
- Patients implanted with the EVOQUE valve should maintain anticoagulant/antiplatelet therapy as determined by their physicians in accordance with current guidelines, to minimize the risk of valve thrombosis or thromboembolic events.
- There are no data to support device safety and performance if the patient has:
 - Echocardiographic evidence of severe right ventricular dysfunction;
 - Pulmonary arterial systolic pressure (PASP) > 70 mmHg by echo Doppler;
 - A trans-tricuspid pacemaker or defibrillator lead that has been implanted in the RV within the last 3 months;
 - Dependency on a trans-tricuspid pacemaker without alternative pacing options.

7.0 Precautions

7.1 Precautions Prior to Use

- The patient's eligibility depends on the anatomic conditions based on CT scan.
- It is advised that a multi-disciplinary heart team be of the opinion that EVOQUE valve implantation is preferable to alternative percutaneous device solutions, including minimally-invasive open heart surgery.
- It is advised that a multi-disciplinary heart team takes into consideration the severity of disease and the chances of reversibility of right heart failure based on a complete hemodynamic assessment.
- The patient should have a life expectancy of at least 12 months.

7.2 Precautions

- The EVOQUE valve is to be used only with the EVOQUE delivery system and EVOQUE loading system.
- The procedure should be conducted under appropriate imaging modalities, such as transesophageal echocardiography (TEE), fluoroscopy, and/or intracardiac echocardiography (ICE).
- Glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Safety Data Sheet available from Edwards Lifesciences.
- Conduction disturbances may occur before, during, or following implantation of the EVOQUE valve, which may require continuous ECG monitoring before hospital discharge. If a patient has confirmed or suspected conduction disturbances, consider patient monitoring and/or electrophysiology evaluation.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Long-term durability has not been established for the EVOQUE valve. Clinical data are reflective of short-term follow-up. Regular medical follow-up is advised to evaluate EVOQUE valve performance.
- Implantation of the EVOQUE valve should be postponed in patients with (1) a history of myocardial infarction within one month (30 days) of planned intervention, (2) pulmonary emboli within 3 months (90 days) of planned intervention, (3) cerebrovascular accident (stroke or transient ischemic attack [TIA]) within 3 months (90 days) of planned intervention, (4) active upper gastrointestinal bleeding within 3 months (90 days) prior to procedure requiring transfusion.

8.0 Potential Adverse Events

Potential adverse events related to standard cardiac catheterization, use of anesthesia, the EVOQUE valve, and the implantation procedure include:

- Abnormal lab values
- Allergic reaction to anesthesia, contrast media, anti-coagulation medication, or device materials
- Anaphylactic shock
- Anemia or decreased hemoglobin (Hgb), may require transfusion
- Aneurysm or pseudoaneurysm
- Angina or chest pain
- Arrhythmia – atrial (i.e., atrial fibrillation, supraventricular tachycardia)
- Arrhythmias – ventricular (i.e., ventricular tachycardia, ventricular fibrillation)
- Arterio-venous fistula
- Bleeding
- Cardiac arrest
- Cardiac (heart) failure
- Cardiac injury, including perforation
- Cardiac tamponade / pericardial effusion
- Cardiogenic shock
- Chordal entanglement or rupture that may require intervention
- Coagulopathy, coagulation disorder, bleeding diathesis
- Conduction system injury, which may require implantation of a pacemaker (temporary or permanent)
- Conversion to open heart surgery
- Coronary artery occlusion
- Damage to or interference with function of pacemaker or implantable cardioverter defibrillator (ICD)
- Death
- Edema
- Electrolyte imbalance
- Embolization including air, particulate, calcific material, or thrombus
- Emergent cardiac surgery
- Endocarditis
- Esophageal irritation
- Esophageal perforation or stricture
- EVOQUE system component(s) embolization
- Failure to retrieve any EVOQUE system components
- Fever
- Gastrointestinal bleeding
- Hematoma
- Hemodynamic compromise
- Hemolysis / hemolytic anemia
- Hemorrhage requiring transfusion/surgery
- Hypertension
- Hypotension
- Inflammation
- Injury to the tricuspid apparatus including chordal damage, rupture, papillary muscle damage
- Local and systemic infection
- Mesenteric ischemia or bowel infarction
- Multi-system organ failure
- Myocardial infarction
- Nausea and/or vomiting
- Nerve injury
- Neurological symptoms, including dyskinesia, without diagnosis of TIA or stroke
- Non-emergent reoperation
- Pain
- Pannus formation

- Paralysis
- Percutaneous valve intervention
- Peripheral ischemia
- Permanent disability
- Pleural effusion
- Pneumonia
- Pulmonary edema
- Pulmonary embolism
- Reaction to anti-platelet or anticoagulation agents
- Rehospitalization
- Renal failure
- Respiratory failure, atelectasis – may require prolonged intubation
- Retroperitoneal bleed
- Right ventricular outflow tract (RVOT) obstruction
- Septicemia, sepsis
- Skin burn, injury, or tissue changes due to exposure to ionizing radiation
- Stroke
- Structural deterioration (wear, fracture, calcification, leaflet tear, leaflet thickening, stenosis of implanted device, or new leaflet motion disorder)
- Thromboembolism
- Transient ischemic attack (TIA)
- Valve dislodgement/embolization
- Valve endocarditis
- Valve explant
- Valve leaflet entrapment
- Valve malposition
- Valve migration
- Valve paravalvular leak (PVL)
- Valve regurgitation (new or worsening tricuspid, aortic, mitral, pulmonary)
- Valve thrombosis
- Vascular injury or trauma, including dissection or occlusion
- Vessel spasm
- Wound dehiscence, delayed or incomplete healing

For a patient/user/third party in the European Economic area: if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

9.0 Additional Equipment

The implantation procedure requires additional equipment that is not supplied with the EVOQUE system. The additional equipment is provided below.

9.1 Equipment for EVOQUE Valve Loading

Note: Volumes reflect adequate amount for preparing 1 implant.

- 3500 ml (minimum) sterile physiological saline solution at ambient temperature (~ 23 °C)
- 500 ml (minimum) heparinized saline solution (2 units/ml) at ambient temperature (~ 23 °C)
- 4 sterile bowls (≥ 500 ml, ≥ 7 cm depth, plastic)
- 1 large sterile bowl (≥ 2 L, ≥ 10 cm depth, plastic)
- 1 scalpel, no. 11 scalpel blade
- 1 luer locking syringe (≥ 20 cc)
- Gauze pads
- Blunt tip forceps
- Blunt tip scissors
- Sterile towels

9.2 Equipment for Access, Procedure, and Monitoring

- Standard cardiac catheterization lab equipment
- Femoral vessel introducer sheath
- Fluoroscopy (fixed, mobile, or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal echocardiography capabilities
- Steerable introducer sheath
- Exchange length 0.035 inch (0.89 mm) max guidewire
- Extra small curve 0.035 inch (0.89 mm) max guidewire
- Right coronary arterial catheter and guidewire
- Sterile table for EVOQUE valve and device preparation

9.3 Standby Equipment

- Arterial bypass cannula (~18 F)
- Cardiopulmonary bypass machine
- Compliant balloon (> 20 mm diameter, 9 cc contrast volume)
- Diluted radiopaque contrast medium (15:85 medium to saline dilution)
- High pressure contrast injector
- Intra-aortic balloon pump and appropriately sized balloon
- Intracardiac echocardiography (ICE) equipment
- Pigtail angiographic catheter
- Venous return bypass cannula (~18 F)
- Transthoracic echocardiographic (TTE) equipment
- Vascular access lubricant
- Temporary pacing equipment

10.0 INSTRUCTION FOR USE

10.1 Prior to Use Inspection

Before using the EVOQUE system, visually examine each element and accessory for evidence of gross damage (e.g., a cracked jar or lid, leakage, broken or missing seals) that may have compromised the packaging sterility (if applicable) or functionality of the components.

WARNING: Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging and/or sterile barriers and any components have been opened or damaged, or the expiration date has elapsed, as sterility and/or function may be compromised.

WARNING: Do not mishandle the EVOQUE valve or use device/container if found to be damaged, leaking, or without adequate sterilant (not fully submerged in glutaraldehyde, or missing intact seals). The EVOQUE valve must not be used for implantation, as sterility may be compromised.

WARNING: Do not use the EVOQUE valve if the expiration date has elapsed, as either sterility or valve function may be compromised.

WARNING: Do not use the EVOQUE valve if the temperature indicator has been activated, as valve function may be compromised.

WARNING: Do not use the EVOQUE valve if the tamper evident seal is broken, as sterility may be compromised.

10.2 Accessories Setup

10.2.1 Accessories Set 1 Setup

Step	Procedure
1	Place the plate on the operating table beneath the patient's leg to support the base during the procedure.
2	After the patient has been positioned on the operating table, place single-use base over the patient's leg, on top of the plate, at the desired distance from the mid-sternum.
3	Establish sterile barrier.
4	After sterile draping, place single-use base adapter in-line with the access site and attached to the front of the base using a clamp (Figure 6).

Note: In the Accessories set 1 configuration, the base adapter and clamp are both supplied with the stabilizer.

10.2.2 Accessories Set 2 Setup

For the Accessories set 2 configuration utilizing reusable accessories (Figure 7), refer to the Edwards Reusable Accessories IFU for set up of the platform and plate. If reusable rail is non-sterile or if unpackaging for the first time, refer to Edwards Reusable Rail Assembly IFU. Once reusable rail is sterile, refer to below instructions.

Step	Procedure
1	Remove Edwards reusable rail from sterilization wrap (blue wrap) and inspect for damage.
2	Prior to sterile draping the patient, assemble and position the reusable plate and reusable platform around the legs of patient, adjusting the height and angle of the platform as needed. Refer to the Edwards Reusable Accessories IFU.
3	Establish sterile barrier.
4	After sterile draping, assemble and attach Edwards reusable rail and stabilizer to the reusable platform.
5	Following the procedure, remove the reusable rail from the reusable platform. Perform standard wipe down of reusable rail before sending for reprocessing. Refer to the Edwards Reusable Rail Assembly IFU.

10.3 Device Preparation

All device preparations will be performed by authorized Edwards personnel.

WARNING: Do not mishandle the EVOQUE valve. EVOQUE valve leaflets mishandled or damaged during any part of the loading process will require replacement of the EVOQUE valve.

WARNING: The EVOQUE valve should not remain fully crimped for a period of time longer than 120 minutes, as it may impair valve functionality.

CAUTION: Do not place the jar or the pouch of the delivery system into the sterile field. The outside of the jar and the pouch are not sterile, and the contents within the jar and pouch should be handled using standard aseptic techniques to prevent contamination.

CAUTION: To reduce the risk of contamination, do not open the EVOQUE valve jar until implantation is certain.

CAUTION: Ensure that the entire suture is removed when removing the serial number tag from the EVOQUE valve, as it may lead to emboli.

CAUTION: Do not allow the EVOQUE valve to come in contact with any sharp instrument, as it may impair valve function.

WARNING: Adequate rinsing with normal saline must be performed before implantation to reduce glutaraldehyde concentration, as it may result in glutaraldehyde toxicity.

CAUTION: Avoid contact of the leaflet tissue or the rinse solution with towels, linens, or other sources of lint or particulate matter that could be transferred to the leaflet tissue, as it may lead to emboli.

10.4 EVOQUE Valve Implant

10.4.1 Guidewire Placement

Prepare femoral venous access using standard interventional techniques.

WARNING: Do not use excessive force and/or manipulation during advancement and positioning of the guidewire, as it may result in perforation/dissection of arteries, veins, and/or other cardiac structures. This may also result in cardiac arrhythmias and conduction disturbances.

Step	Procedure
1	Advance steerable sheath into the right atrium at the exit of the IVC.
2	Insert a guidewire through the steerable sheath.
3	Advance the guidewire across the tricuspid valve. Note: Other interventional devices and techniques (e.g., guide catheters) may be used to assist the guidewire in crossing the tricuspid valve.
4	Establish proper guidewire path and confirm no entanglement with cardiac structures.

10.4.2 EVOQUE Valve Delivery

WARNING: Avoid excessive movement of the delivery system while executing the procedure in order to protect vasculature or cardiac structures. Avoid excessive rotation of the delivery system to maintain functionality of the delivery system.

Note: Flush the delivery system with heparinized saline as needed throughout the procedure.

Step	Procedure
1	Ensure the hydrophilic coating on dilators and delivery system is activated prior to use. Dilate access site. An Edwards EVOQUE dilator kit should be used as needed.
2	Insert the delivery system over the guidewire.
3	Advance the delivery system until the distal end of the tapered tip is positioned at the junction between IVC and right atrium.
4	Using fluoroscopy, ensure the delivery system is oriented in the proper orientation. WARNING: The primary flex on the delivery system flexes in the direction of the flush ports; care should be taken to ensure that the delivery system is orientated correctly at this point.
5	Retract the sheath.
6	Flex and orient the delivery system towards the tricuspid valve.
7	Advance the delivery system to cross the tricuspid valve. Note: Delivery system flex, delivery system rotation, and guidewire position may be adjusted during valve crossing to optimize crossing position.
8	Using echocardiography and fluoroscopic guidance, verify that the delivery system has crossed through the tricuspid valve into the right ventricle.
9	If using the Accessories set 1 configuration, dock the stabilizer onto the base adapter and secure it to the base. If using the Accessories set 2 configuration, dock the rail onto the platform and secure it. Next, dock the stabilizer onto the rail and secure it.
10	Attach the delivery system and sheath to the stabilizer.
11	Adjust the delivery system as needed to ensure hemodynamic stability.
12	Using pre-operative MSCT data (if available), position the C-arm at the optimal viewing projection.
13	Position the delivery system to be coaxial to the tricuspid annulus while minimizing contact with the native anatomy.
14	Using echocardiography and fluoroscopic guidance, confirm that the EVOQUE valve is positioned at the appropriate depth and coaxiality with relation to the native valve. CAUTION: Once the capsule is retracted to expose the EVOQUE valve anchors, the valve is unable to be retrieved or recaptured into the delivery system. WARNING: Maintain central position of the delivery system within the native valve during deployment to ensure proper positioning of the valve.
15	Retract the outer capsule until the anchors are exposed.
16	Adjust the position of the EVOQUE valve so that the anchors are positioned within native leaflets as dictated by patient anatomy.
17	Retract the outer and inner capsule until the desired EVOQUE valve diameter is achieved.
18	Engage the leaflets.
19	Confirm positioning of the EVOQUE valve by using echo imaging to assess engagement of the leaflets. Adjust the EVOQUE valve positioning as needed.
20	Using echo imaging, observe the motion of the native leaflets and adjust the position of the EVOQUE valve as required to fully engage native tricuspid valve leaflets.

Step	Procedure
21	Once full engagement has been confirmed, ensure that the EVOQUE valve is perpendicular to the tricuspid annular plane.
22	Retract the tapered tip until it is positioned within the EVOQUE valve.
23	Retract the inner capsule until the EVOQUE valve is released from the delivery system. CAUTION: Care should be taken during final release of the EVOQUE valve using the release knob, as premature release may impact performance of the EVOQUE valve.
24	Using echo and fluoroscopic imaging, assess final position and functionality of the EVOQUE valve.

10.4.3 Delivery System Removal

WARNING: Take care to maintain position of delivery system central within the EVOQUE valve during delivery system removal, as failure to do so may impact valve functionality or lead to dislodgement of the valve.

Note: Delivery system can be removed from the stabilizer at any point during removal, if applicable.

Step	Procedure
1	Fully retract the tapered tip.
2	Unflex and retract delivery system as needed until the tapered tip is above the locking tabs of the EVOQUE valve. Adjust the guidewire as needed to maintain a central position relative to the EVOQUE valve. Ensure that the locking ring is free from the EVOQUE valve.
3	Rotate the release knob so that the inner capsule is in contact with the tapered tip.
4	Unflex and retract the delivery system as needed.
5	Rotate the capsule retraction knob until the outer capsule is in contact with the inner capsule.
6	Ensure the delivery system is fully unflexed and remove the delivery system from the access site. Note: A sheath may be used to seal the femoral vein following system removal.
7	Perform femoral closure as applicable using standard interventional techniques.
8	Perform a ventriculogram if needed to assess final position of the EVOQUE valve.

11.0 How Supplied

11.1 Sterilization and Packaging

The EVOQUE valve is provided sterile by terminal liquid sterilization and is non-pyrogenic. It is packaged and sterilized in a glutaraldehyde solution inside a jar to which a seal has been applied. The outer surface of the jar is not sterile and must not be placed in the sterile field. The EVOQUE valve is supplied with a temperature indicator and should not be used if the indicator has been activated.

The delivery system, dilator kit, and loading system are supplied sterilized by ethylene oxide and provided non-pyrogenic. The components are secured on a card and packaged in a pouch and shelf box.

The stabilizer is provided ethylene oxide sterilized. The components are secured on a card and packaged in a pouch and shelf box.

The single-use base and plate or the reusable rail, platform, and plate are provided non-sterile.

The components are packaged in individual shipper boxes.

11.2 Storage

The EVOQUE valve should be stored between 10 °C and 25 °C (50 °F and 77 °F). Stock inspection and rotation at regular intervals are recommended so that the EVOQUE valve with an earlier expiration date is used first.

The delivery system, dilator kit, loading system, single-use stabilizer/base/plate and reusable platform/plate should be stored in a cool, dry place that is contamination free. Refer to the Edwards Reusable Rail Assembly IFU for storage of the reusable rail.

12.0 Magnetic Resonance Imaging (MRI) Safety Information

Non-clinical testing has demonstrated that the Edwards EVOQUE valve is MR Conditional. A patient with the valve can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 and 3 T only
- Maximum spatial gradient magnetic field of 3000 gauss/cm (30.0 T/m) or less
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg
- Normal operating mode of the MR system for both gradients and SAR

Under the scan conditions defined above, the EVOQUE valve is expected to produce a maximum temperature rise of 4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the EVOQUE valve extends approximately 0.8 cm from the device when imaged with a gradient echo or spin echo pulse sequence and a 3 T MRI system.

13.0 Expected Lifetime of the Device

The EVOQUE valve has been subjected to rigorous pre-clinical durability testing per the valve testing requirements. The valves were successfully tested to 5 years of simulated wear. Actual lifetime performance in humans depends on many factors and varies from patient to patient.

14.0 Patient Information

A patient implant card is provided with each EVOQUE valve. After implantation, please complete all requested information and provide the implant card to the patient. The serial number is found on the EVOQUE valve package. This implant card allows patients to inform healthcare providers what type of implant they have when they seek care.

15.0 Recovered Implant and Device Disposal

Edwards Lifesciences is interested in obtaining recovered clinical specimens of the EVOQUE valve for analysis. A written report summarizing our findings will be provided upon completion of our evaluation. Please contact Edwards for return of the recovered valve.

If you do decide to return any of the devices, please follow the following instructions:

- **Unopened Package with Sterile Barrier Intact:** If the pouches have not been opened, return the device in its original packaging.
- **Package Opened but Not Implanted:** If a pouch is opened, the device is no longer sterile. Please return the device in its original packaging.
- **Explanted Implant:** The explanted implant should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit for return to Edwards.

15.1 Disposal

Use universal precautions for biohazards and sharps to avoid user injury. Used devices (includes all devices that come in contact with patients) should be handled and resterilized or disposed of in accordance with institutional guidelines on biohazardous materials and hospital waste to avoid possible cross-contamination. For disposal of the reusable platform and plate accessories, refer to the Edwards Reusable Accessories IFU.

16.0 Qualitative and Quantitative Information Related to the EVOQUE Valve

For the EVOQUE valve, the following table shows the qualitative and quantitative information on the materials and substances:

EVOQUE Valve

Substance	CAS	Model Mass Range (mg)
Nickel	7440-02-0	1250 - 1470
Titanium	7440-32-6	981 - 1171
Polyethylene terephthalate	25038-59-9	529 - 649

Substance	CAS	Model Mass Range (mg)
Polyethylene	9002-88-4	285
Polyurethane	9009-54-5	0 - 281
Collagens, bovine, polymers with glutaraldehyde	2370819-60-4	113
Titanium dioxide	13463-67-7	1.41 - 3.65
Cobalt	7440-48-4	0 - 1.29
Iron	7439-89-6	0 - 1.29
Carbon	7440-44-0	0 - 1.03
Oxygen	7782-44-7	0 - 1.03
Antimony trioxide	1309-64-4	0.632 - 0.777
Siloxanes and silicones, di-methyl, polymers with methyl silsesquioxanes, polyethylene-polypropylene glycol, and polyethylene-polypropylene glycol monobutyl ether	68554-68-7	0 - 0.762
Niobium	7440-03-1	0 - 0.643
Silicon dioxide	7631-86-9	0.00515 - 0.598
Polybutadiene diacrylate	9003-17-2	0 - 0.452
Polybutilite	24936-97-8	0.333
Chromium	7440-47-3	0 - 0.257
Copper	7440-50-8	0 - 0.257
Hydrogen	1333-74-0	0 - 0.129
Nitrogen	7727-37-9	0 - 0.129
4-Dodecylbenzenesulfonic acid	121-65-3	0.000541 - 0.0532
D&C Green No. 6	128-80-3	0.0481 - 0.0503
Erucamide	112-84-5	0.00508 - 0.0256
Glycerol	56-81-5	0 - 0.0168
Bis(2-dimethylaminoethyl) ether	3033-62-3	0 - 0.00790
Triethylenediamine	280-57-9	0 - 0.00611

EVOQUE Delivery System:

This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless-steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

17.0 Summary of Safety and Clinical Performance (SSCP)

Refer to <https://meddeviceinfo.edwards.com/> for a SSCP for this medical device.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device. This is the location of the SSCP after the launch of European Database on Medical Devices/Eudamed.

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information. The Notified Body has taken notice of and agreed with the benefit-risk rationales for the short- and long-term safety and effectiveness of the EVOQUE system.

Conformity with the entire EVOQUE system of the Performance Requirements (GSPR) for safety (MDR GSPR1), performance (MDR GSPR1), acceptability of side-effects (MDR GSPR8), usability (MDR GSPR5), device lifetime (MDR GSPR6), and acceptable benefit-risk profile (MDR GSPR8) has been established for the labelled indications.

18.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed. The Basic UDI-DI for the valves, delivery system, and accessories can be used to locate the SSCP. The following table contains the Basic UDI-DIs for the EVOQUE system:

Product	Model				Basic UDI-DI
Edwards EVOQUE Valve	9850EV44	9850EV48	9850EV52	9850EV56	0690103D004EUV000Z3
Edwards EVOQUE Loading System or Edwards EVOQUE Loading System and Trimmer	9850LS or 9850LSB			9850LSB	0690103D004EVL000VV
Edwards EVOQUE Tricuspid Delivery System	9850TDS				0690103D004EVD000V5
Edwards EVOQUE Dilator Kit	9850DK				0690103D004EVK000WN

Accessories Set 1		
Edwards EVOQUE Stabilizer	9850SB	0690103D004EVA000UG
Edwards EVOQUE Stabilizer Base	9850BA	0690103D004EVA000UG
Edwards EVOQUE Stabilizer Plate	9850PT	0690103D004EVA000UG

Accessories Set 2		
Edwards EVOQUE Stabilizer	9850SZ	0690103D004EVA000UG
Edwards Reusable Rail	10500RL	0690103D004REU000YA
Edwards Reusable Platform	10000UP	0690103D004REU000YA
Edwards Reusable Plate	10000PT	0690103D004REU000YA

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
	Suomi	Български	Română	Eesti	Lietuvių
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Säilytä auringonvalolta suojattuna	Да се пази от слънчева светлина	A se feri de lumina solară	Hoida eemal päikesevalgusest	Saugoti nuo saulės šviesos
	Pidä kuivana	Да се пази сухо	A se păstra uscat	Hoida kuivalt	Laikykite sausoje vietoje
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
STERILE EO	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
STERILE R	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant
STERILE LC	Steriloitu nestemäisellä kemikaalilla	Стерилизирано с течни химикали	Sterilizat cu substanțe chimice lichide	Steriliseeritud vedelat kemikaali kasutades	Sterilizuota skystomis cheminėmis medžiagomis
	Älä steriloi uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Epästeriili	Нестерилно	Nesteril	Mittesteriilne	Nesterilus
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Ei sisällä DEHP:tä	Не съдържа DEHP	Nu conține DEHP	Mitte-DEHP	Ne DEHP
	Viimeinen käyttöpäivämäärä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Tinkamumo terminas
SN	Sarjanumero	Сериен номер	Număr de serie	Seerianumber	Serijos numeris
EC REP	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană / Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
44 mm	Käytetään 44 mm:n kokoisten katetrilla asennettavien Edwards -sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа на Edwards с размер 44 mm	A se utiliza cu valva cardiacă transcateter Edwards de mărimea 44 mm	Kasutamiseks ettevõtte Edwards 44 mm transkateetriga südameklapiga	Skirta naudoti su 44 mm dydžio transkateteriniu širdies vožtuvu „Edwards“
48 mm	Käytetään 48 mm:n kokoisten katetrilla asennettavien Edwards -sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа на Edwards с размер 48 mm	A se utiliza cu valva cardiacă transcateter Edwards de mărimea 48 mm	Kasutamiseks ettevõtte Edwards 48 mm transkateetriga südameklapiga	Skirta naudoti su 48 mm dydžio transkateteriniu širdies vožtuvu „Edwards“
52 mm	Käytetään 52 mm:n kokoisten katetrilla asennettavien Edwards -sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа на Edwards с размер 52 mm	A se utiliza cu valva cardiacă transcateter Edwards de mărimea 52 mm	Kasutamiseks ettevõtte Edwards 52 mm suuruse kateedriklause südameklapiga	Skirta naudoti su 52 mm dydžio transkateteriniu širdies vožtuvu „Edwards“
56 mm	Käytetään 56 mm:n kokoisten katetrilla asennettavien Edwards -sydänläppien kanssa	За употреба с транскатетърна сърдечна клапа Edwards с размер 56 mm	A se utiliza cu valva cardiacă transcateter Edwards de mărimea 56 mm	Kasutamiseks ettevõtte Edwards 56 mm transkateetriga südameklapiga	Skirta naudoti su 56 mm dydžio transkateteriniu širdies vožtuvu „Edwards“
	Lämpötilaraja	Температурно ограничение	Limită de temperatură	Temperatuuripiirang	Temperatūros apribojimas
MD	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola



	Suomi	Български	Română	Eesti	Lietuvių
	Ulkoalkaisija	Външен диаметър	Diametru exterior	Välisläbimõõt	Išorinis skersmuo
	Sisäläpimitta	Вътрешен диаметър	Diametru interior	Sisemine läbimõõt	Vidinis skersmuo
	Suosittelu ohjainlangan pituus	Препоръчителна дължина на теления водач	Lungimea recomandată a firului de ghidaj	Soovituslik juhtetraadi pikkus	Rekomenduojamas kreipiamosios vielos ilgis
GW	Ohjainlangan suositeltu koko	Препоръчителен размер телен водач	Mărimea recomandată a firului de ghidaj	Juhtetraadi soovituslik suurus	Rekomenduojamas kreipiamosios vielos dydis
GWC	Ohjainlangan yhteensopivuus	Съвместимост на теления водач	Compatibilitatea firului de ghidaj	Juhtetraadi ühilduvus	Kreipiamosios vielos suderinamumas
SZ	Koko	Размер	Mărime	Suurus	Dydis
	Katetrin varren koko	Размер на шифта на катетъра	Mărime a tije cateterului	Kateetri varre suurus	Kateterio vamzdelio dydis
	Pallon halkaisija	Диаметър на балона	Diametru balon	Ballooni läbimõõt	Balionėlio skersmuo
	Pallon käyttöpituu	Работна дължина на балона	Lungimea de lucru a balonului	Ballooni tööpikkus	Balionėlio darbinis ilgis
	Ehdollisesti turvallinen magneettikuvauksessa	Безопасно при MR при определени условия	Condiționat RM	Ohutu magnetresonantstomograafias teatud tingimustel	Sąlyginis MR
	Ei sovi magneettikuvaukseen	Небезопасно при MR	Incompatibil RM	Ohtlik magnetresonantstomograafias	MR nesaugus
UDI	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusi priemonės identifikatorius
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Sisältää vaarallisia aineita	Съдържа опасни вещества	Conține substanțe periculoase	Sisaldab ohtlikke aineid	Sudėtyje yra pavojingų medžiagų
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilus barjero sistema
	Sisältää eläinperäistä biologista materiaalia	Съдържа биологичен материал от животински произход	Conține material biologic de origine animală	Sisaldab loomset päritolu bioloogilist materjali	Sudėtyje yra biologinių gyvūninės kilmės medžiagų

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
#	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
REF	Kataloga numurs	Katalog Numarası	Номер по каталогу	Kataloški broj	Kataloški broj
LOT	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
QTY	Daudzums	Miktar	Количество	Količina	Količina

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символах ■ Legenda de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Saturs	İçindekiler	Содержимое	Sadržaj	Sadržaj
	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjliva dužina	Upotrebjliva dužina
	Nelietot atkārtoti	Yeniden kullanmayın	Не использовать повторно	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na internet stranici	Pogledajte upute za upotrebu na internetskoj stranici
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Nelietot produktu, ja attēlota norāde	Ürünü endikasyon gösterilmişse kullanmayın	Не используйте изделие, если отображается индикация	Nemojte koristiti proizvod ako je prikazana indikacija	Nemojte upotrebljavati proizvod ako je prikazana indikacija
	Lietot produktu, ja attēlota norāde	Ürünü endikasyon gösterilmişse kullanın	Используйте изделие, если отображается индикация	Koristite proizvod ako je prikazana indikacija	Upotrebljavajte proizvod ako je prikazana indikacija
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Sargāt no saules gaismas	Güneş ışığından uzak tutun	Беречь от солнечных лучей	Držite podalje od sunčeve svetlosti	Držite dalje od sunčeve svjetlosti
	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Sterilizēts, izmantojot šķidru ķīmisko vielu	Sıvı kimyasal kullanılarak sterilize edilmiştir	Стерилизовано с использованием жидкого химического вещества	Sterilisano korišćenjem tečne hemikalije	Sterilizirano tekućom kemikalijom
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nesterils	Steril Değildir	Нестерильно	Nesterilno	Nesterilno
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	DEHP nesaturošs	DEHP değildir	Не содержит диэтилгексилфталата (ДЭГФ)	Bez DEHP	Ne sadrži DEHP
	Derīguma termiņš	Şu tarihe kadar kullanın	Срок годности	Rok korišćenja	Rok upotrebe

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символах ■ Legenda de simboluri ■ Sümbölitte seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
SN	Sērijas numurs	Seri Numarası	Серийный номер	Serijski broj	Serijski broj
EC REP	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
44 mm	Lietošanai ar 44 mm izmēra Edwards transkatetra sirds vārstuli	44 mm boyutundaki Edwards transkateter kalp kapağı ile kullanım içindir	Для использования с сердечным клапаном для чрескатетерного введения Edwards размером 44 мм	Za upotrebu sa Edwards transkateterskim srčanim zaliskom veličine 44 mm	Za upotrebu s transkateterskim srčanim zaliskom tvrtke Edwards veličine 44 mm
48 mm	Lietošanai ar 48 mm izmēra Edwards transkatetra sirds vārstuli	48 mm boyutundaki Edwards transkateter kalp kapağı ile kullanım içindir	Для использования с сердечным клапаном для чрескатетерного введения Edwards размером 48 мм	Za upotrebu sa Edwards transkateterskim srčanim zaliskom veličine 48 mm	Za upotrebu s transkateterskim srčanim zaliskom tvrtke Edwards veličine 48 mm
52 mm	Lietošanai ar 52 mm izmēra Edwards transkatetra sirds vārstuli	52 mm boyutundaki Edwards transkateter kalp kapağı ile kullanım içindir	Для использования с сердечным клапаном для чрескатетерного введения Edwards размером 52 мм	Za upotrebu sa Edwards transkateterskim srčanim zaliskom veličine 52 mm	Za upotrebu s transkateterskim srčanim zaliskom društva Edwards veličine 52 mm
56 mm	Lietošanai ar 56 mm izmēra Edwards transkatetra sirds vārstuli	56 mm boyutundaki Edwards transkateter kalp kapağı ile kullanım içindir	Для использования с сердечным клапаном для чрескатетерного введения Edwards размером 56 мм	Za upotrebu sa Edwards transkateterskim srčanim zaliskom veličine 56 mm	Za upotrebu s transkateterskim srčanim zaliskom tvrtke Edwards veličine 56 mm
	Temperatūras ierobežojums	Sıcaklık sınırı	Ограничение по температуре	Ograničenje temperature	Ograničenje temperature
MD	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Ārējais diametrs	Dış çap	Внешний диаметр	Spoljašnji prečnik	Vanjski promjer
	Iekšējais diametrs	İç çap	Внутренний диаметр	Unutrašnji prečnik	Unutrašnji promjer
	Ieteicamais vadītājstīgas garums	Kılavuz tel için önerilen uzunluk	Рекомендуемая длина проводника	Preporučena dužina žice vodice	Preporučena dužina žice vodilice
GW	Ieteicamais vadītājstīgas izmērs	Kılavuz Tel İçin Önerilen Boyut	Рекомендованный размер проводника	Preporučena veličina žice vodice	Preporučena veličina žice vodilice
GWC	Vadītājstīgas saderība	Kılavuz tel uyumluluğu	Совместимо с проводником	Kompatibilnost žice vodice	Kompatibilnost žice vodilice
SZ	Izmērs	Boyut	Размер	Veličina	Veličina
	Katetra ass izmērs	Kateter shaft boyutu	Размер стержня катетера	Veličina osovine katetera	Veličina osovine katetera
	Balona diametrs	Balon çapı	Диаметр баллона	Prečnik balona	Promjer balona
	Balona darba garums	Balon çalışma uzunluğu	Рабочая длина баллона	Radna dužina balona	Radna duljina balona
	Izmantojams MR vidē, ievērojot noteiktus nosacījumus	MR Koşullu	Условно безопасно при проведении МРТ	Uslovno bezbedno za MR	Uvjetno siguran kod pregleda MR-om
	Nedrīkst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR	Nije sigurno kod pregleda MR-om

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Satur bīstamas vielas	Tehlikeli maddeler içerir	Содержит вредные вещества	Sadrži opasne supstance	Sadrži opasne tvari
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Satur dzīvnieku izcelsmes bioloģisko materiālu	Hayvan kökenli biyolojik madde içerir	Содержит материалы животного происхождения.	Sadrži biološki materijal životinjskog porekla	Sadrži materijale životinjskog podrijetla

Note: The labeling of this product may not contain every symbol depicted in this legend. ■ **Remarque :** il est possible que certains symboles de cette légende n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Die Kennzeichnung dieses Produkts enthält möglicherweise nicht jedes Symbol, das in dieser Zeichenerklärung abgebildet ist. ■ **Nota:** Es posible que no todos los símbolos mostrados en esta leyenda aparezcan en el etiquetado de este producto. ■ **Nota:** l'etichetta di questo prodotto potrebbe non contenere tutti i simboli riportati in questa legenda. ■ **Opmerking:** De etikettering van dit product bevat mogelijk niet elk symbool dat in deze legenda wordt weergegeven. ■ **Bemærk:** Produktets mærkning omfatter muligvis ikke alle symboler gengivet i denne forklaring. ■ **Obs!** Märkningen av produkten innehåller eventuellt inte alla symboler som avbildas i denna symbolförklaring. ■ **Σημείωση:** Η σήμανση του προϊόντος ενδέχεται να μην περιέχει κάθε σύμβολο που αναφέρεται στο παρόν υπόμνημα. ■ **Nota:** a rotulagem deste produto pode não conter todos os símbolos ilustrados nesta legenda. ■ **Poznámka:** Označení tohoto výrobku nemusí obsahovat všechny symboly uvedené v této legendě. ■ **Megjegyzés:** A termék címkéje nem feltétlenül tartalmazza az ezen jelmagyarázatban szereplő összes szimbólumot. ■ **Uwaga:** Etykieta tego wyrobu może nie zawierać wszystkich symboli przedstawionych w niniejszej legendzie. ■ **Poznámka:** Na štítku tohto produktu sa nemusia nachádzať všetky symboly uvedené v týchto vysvetlivkách. ■ **Merk:** Merkingen av dette produktet inneholder kanskje ikke alle symbolene som er gjengitt i denne symbolforklaringen. ■ **Huomautus:** Tämän tuotteen merkinnät eivät välttämättä sisällä kaikkia tässä selityksessä olevia symboleita. ■ **Забележка:** Етикетите по този продукт може да не съдържат всеки символ, представен в тази легенда. ■ **Notă:** este posibil ca eticheta acestui produs să nu conțină fiecare simbol descris în această legendă. ■ **Märkus:** Tootel esinev märgistus ei pruugi sisaldada kõiki tingimärkide seletuses välja toodud sümboliteid. ■ **Pastaba.** Šio gaminio etiketėje gali nebūti visų šiamė paaiškinime pavaizduotų simbolių. ■ **Piezīme.** Šī produkta marķējuma var nebūt visu simbolu, kuri iekļauti šajā simbolu skaidrojumā. ■ **Not:** Bu ürünün etiketi bu açıklamada verilen her sembolü içermeyebilir. ■ **Примечание.** Маркировка изделия может содержать не все условные обозначения. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** oznaka ovoga proizvoda možda ne sadrži svaki simbol prikazan u ovoj legendi simbola.

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12/2024
10058872003 C
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU

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Edwards Lifesciences GmbH