Edwards

<b>Directory</b>	
English (en).....	1
Français (fr).....	7
Deutsch (de).....	13
Español (es).....	19
Italiano (it).....	25
Nederlands (nl).....	31
Dansk (da).....	37
Svenska (sv).....	43
Ελληνικά (el).....	49
Português (pt).....	55
Česky (cs).....	61
Magyar (hu).....	67
Polski (pl).....	73
Slovensky (sk).....	79
Norsk (no).....	85
Suomi (fi).....	91
Български (bg).....	97
Română (ro).....	103
Eesti (et).....	109
Lietuvių (lt).....	115
Latviešu (lv).....	121
Türkçe (tr).....	127
Русский (ru).....	133
Srpski (sr).....	140
Hrvatski (hr).....	146
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike.....	152
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	156

**English**

**Rx Only**

**Edwards Reusable Rail Set**

**Accessories Care and Reprocessing Instructions for Edwards Reusable Rail (10500RL), Edwards Reusable Rail Extension (10500EX), and Edwards Rail Tray (TRAY10500RL)**

**Instructions for Use**

**1.0 Introduction**

The Edwards reusable rail (rail), reusable rail extension (rail extension), and rail tray (tray), collectively referred to as the Edwards Reusable Rail Set, are reusable and sterilizable accessories indicated for use with compatible Edwards transcatheter cardiac therapies (Figure 1). The rail and rail extension are optional, non-patient contacting accessories that fasten to a compatible reusable platform to serve as a rail for compatible modular stabilizers that support and maintain catheter positioning during implantation procedures. The rail tray is an optional, non-patient contacting accessory which serves as a container for compatible accessories for storage and/or reprocessing.

The Edwards Reusable Rail Set is intended to be used in a cardiac procedure room (e.g. catheterization laboratory, hybrid room, etc.) by physicians and/or support staff and in sterile processing areas by sterile processing department staff.

**2.0 Intended Use**

The reusable rail is intended to aid the positioning and stabilization of delivery systems during intra-cardiac procedures.

The reusable rail extension is intended to aid the positioning and stabilization of delivery systems during intra-cardiac procedures.

The rail tray is intended to facilitate the sterile processing and storage of compatible intra-cardiac device accessories.

**3.0 Precautions**

Read all indications and precautions in these Instructions for Use (IFU) before using the rail, rail extension, and rail tray. The IFU is available electronically, via [www.eifu.edwards.com](http://www.eifu.edwards.com). You may also request a copy by calling Edwards Headquarters at +1.888.570.4016 or by contacting your local clinical representative.

The rail, rail extension, and rail tray are supplied non-sterile. Prior to use, cleaning and sterilization must be performed according to the instructions provided in this document. Failure to properly clean and / or re-sterilize may result in infection.

Original packaging of the rail and rail extension are not intended for sterilization. Processing using the original packaging could result in compromised sterility.

Do not use the devices if damage is observed.

## 4.0 How Supplied

The rail, rail extension, and rail tray are supplied non-sterile. See the *Figures* appendix for reference.

## 5.0 Limitations on Reprocessing

These reprocessing instructions have been validated as being capable of preparing the rail and rail extension for sterile, procedural use. It is the responsibility of the user/hospital/healthcare provider to ensure that reprocessing is performed using the appropriate equipment and materials specified throughout this document, and that personnel have been adequately trained to these reprocessing instructions to achieve the desired level of decontamination; this requires that equipment and processes be validated and routinely monitored. Any deviation by the user/hospital/healthcare provider from the reprocessing methods outlined in these instructions should be evaluated for effectiveness to avoid potential adverse consequences.

End of life is normally determined by wear and damage due to use (see Section 5.1 for inspection instructions).

### 5.1 Instructions for Device Visual and Functional Inspections

#### 5.1.1 Preparation for Visual and Functional Inspections

Step	Procedure
1	Ensure to remove all accessories from the rail and rail extension.
2	If not disassembled, ensure to remove the rail extension from the rail by pulling the rail extension upwards from the rail (Figure 4).

#### 5.1.2 Visual Inspections

Step	Procedure
1	Inspect the rail, rail extension, and rail tray for any unremovable biohazardous stains.
2	Inspect the rail, rail extension, and rail tray for cracking and/or gouges, corrosion, or other deterioration/damage are evident in a well-lit environment. Normal wear and tear is expected due to repeated use.
3	The defective device should be retired from further use and replaced if any inspections are unacceptable.

#### 5.1.3 Functional Inspections

Step	Procedure
1	Inspect that the knob of the rail fully translates in threads by rotating the knob clockwise and counterclockwise.
2	Inspect that the rail and extension both fit within the specified brackets within the rail tray and the tray lid latches to the tray base.
3	The defective device should be retired from further use and replaced if any inspections are unacceptable.

### 5.2 Device Disposal

Damaged devices may be handled and disposed of in the same manner as hospital waste in accordance with local regulations as there are no special risks related to the disposal of these devices.

## 6.0 Cleaning Instructions

The rail, rail extension, and rail tray must be cleaned prior to each use in accordance with appropriate hospital procedure and the following recommended instructions and parameters. After cleaning, the rail and rail extension should be placed in their corresponding brackets within the rail tray (Figure 2); the rail tray should then be wrapped in two layers of standard sterilization wrap until sterilized.

**Note: Do not use metallic brushes, scrub pads, harsh chemicals (such as chlorine or caustic soda), or other abrasive cleaning aids when cleaning the devices. They can cause permanent device damage.**

The devices should be cleaned using one of the following methods: 1) Manual Cleaning with Enzymatic Detergent (Section 6.1), 2) Manual Cleaning with Alkaline Detergent (Section 6.2), or 3) Automated Cleaning with either Enzymatic or Alkaline Detergent (Section 6.3). The procedure for each of these methods is presented below.

**WARNING: Failure to properly place, clean, and sterilize the rail and rail extension in the rail tray according to these instructions for use may result in infection, user injury, or patient injury.**

#### 6.1 Instructions for Manual Cleaning with Enzymatic Detergent

Step	Procedure
1	Ensure to remove all accessories from the rail and rail extension.
2	Disassemble the rail and rail extension. a) If the rail and rail extension are contained within the rail tray, remove the lid from the rail tray then remove the rail and rail extension. b) If the rail and rail extension are assembled together, pull up on the rail extension to remove from the rail (Figure 4). For the rail, ensure that the knob is completely unthreaded from the proximal clamp assembly and the threads are exposed (Figure 3).
3	Rinse the rail, rail extension, and rail tray under running $\geq 20$ °C utility (potable) water for 1 minute to remove gross soil. If soil is visible, use a soft-bristled brush to scrub visible soil. While rinsing the rail: a) Focus on threads, crevices, sliding surfaces, and engraved surfaces.

Step	Procedure
	<p>b) Ensure that the knob remains fully unthreaded from the proximal clamp assembly threads and internal knob threads are fully exposed (Figure 3).</p> <p>While rinsing the rail extension:</p> <p>a) Focus on crevices and engraved surfaces.</p> <p>While rinsing the rail tray:</p> <p>a) Focus on crevices, latching side clips, handles, brackets, and engraved surfaces.</p>
4	Prepare an enzymatic and pH neutral detergent, such as Valsure® Enzymatic Cleaner or equivalent, per the manufacturer's recommendations using utility (potable) water.
5	<p>Fully immerse the separated rail, rail extension, and rail tray in the prepared detergent bath and allow them to soak for a minimum of 2 minutes or per manufacturer's recommendation, whichever is greater. Soak the rail such that the knob threads are fully exposed (Figure 3).</p> <p>When cleaning multiple devices simultaneously, avoid contact between devices.</p>
6	<p>After the soak time, use a soft-bristled brush to scrub the rail, rail extension, and rail tray below the water line for a minimum of 1 minute each to remove all visible soil.</p> <p>For the rail:</p> <p>a) Focus on internal and external threads, crevices, sliding surfaces, and engraved surfaces.</p> <p>b) While brushing, ensure that the knob remains fully unthreaded from the proximal clamp assembly, and the internal knob threads are fully exposed (Figure 3).</p> <p>For the rail extension:</p> <p>a) Focus on crevices and engraved surfaces.</p> <p>For the rail tray:</p> <p>a) Focus on crevices, latching side clips, handles, brackets, and engraved surfaces.</p>
7	While submerged, wipe all surfaces of the rail, rail extension, and rail tray for 1 minute with a clean cloth (e.g., lint-free cloths) moistened with the prepared detergent solution.
8	<p>Remove the rail, rail extension, and rail tray from the detergent bath and rinse them under running utility (potable) water for a minimum of 5 minutes to remove detergent residuals.</p> <p>a) Ensure that running water passes through threads, crevices, latching side clips, handles, sliding surfaces and engraved surfaces of all the devices.</p>
9	Dry the rail, rail extension, and rail tray using lint-free cloths and/or filtered pressurized air.
10	Inspect the rail, rail extension, and rail tray for visible contamination. Repeat steps 3 through 9 if contamination is observed.
11	<p>If the devices are to be disinfected next, disregard this step and proceed to Section 7.0.</p> <p>After cleaning, place the rail and rail extension in their corresponding brackets within the rail tray (Figure 2).</p> <p>Wrap the rail tray (with the rail and rail extension inside) in two layers of standard sterilization wrap until sterilized.</p>

## 6.2 Instructions for Manual Cleaning with Alkaline Detergent

Step	Procedure
1	Ensure to remove all accessories from the rail and rail extension.
2	<p>Disassemble the rail and rail extension.</p> <p>a) If the rail and rail extension are contained within the rail tray, remove the lid from the rail tray then remove the rail and rail extension.</p> <p>b) If the rail and rail extension are assembled together, pull up on the rail extension to remove from the rail (Figure 4).</p> <p>For the rail, ensure that the knob is completely unthreaded from the proximal clamp assembly and the threads are exposed (Figure 3).</p>
3	<p>Rinse the rail, rail extension, and rail tray under running <math>\geq 20</math> °C utility (potable) water for 1 minute to remove gross soil. If soil is visible, use a soft-bristled brush to scrub visible soil.</p> <p>While rinsing the rail:</p> <p>a) Focus on threads, crevices, sliding surfaces, and engraved surfaces.</p> <p>b) Ensure that the knob remains fully unthreaded from the proximal clamp assembly threads and internal knob threads are fully exposed (Figure 3).</p> <p>While rinsing the rail extension:</p> <p>a) Focus on crevices and engraved surfaces.</p> <p>While rinsing the rail tray:</p> <p>a) Focus on crevices, latching side clips, handles, brackets, and engraved surfaces.</p>
4	Prepare an alkaline detergent, such as Neodisher® Mediclean Forte detergent or equivalent, per the manufacturer's recommendations using utility (potable) water.

<b>Step</b>	<b>Procedure</b>
<b>5</b>	Fully immerse the separated rail, rail extension, and rail tray in the prepared detergent bath and allow them to soak for a minimum of 10 minutes or per manufacturer's recommendation, whichever is greater. Soak the rail such that the knob threads are fully exposed (Figure 3). When cleaning multiple devices simultaneously, avoid contact between devices.
<b>6</b>	After the soak time, use a soft-bristled brush to scrub the rail, rail extension, and rail tray below the water line for a minimum of 1 minute each to remove all visible soil. For the rail: a) Focus on internal and external threads, crevices, sliding surfaces, and engraved surfaces. b) While brushing, ensure that the knob remains fully unthreaded from the proximal clamp assembly, and the internal knob threads are fully exposed (Figure 3). For the rail extension: a) Focus on crevices and engraved surfaces. For the rail tray: a) Focus on crevices, latching side clips, handles, brackets, and engraved surfaces.
<b>7</b>	While submerged, wipe all surfaces of the rail, rail extension, and rail tray for 1 minute with a clean cloth (e.g., lint-free cloths) moistened with the prepared detergent solution.
<b>8</b>	Remove the rail, rail extension, and rail tray from the detergent bath and rinse them under running utility (potable) water for a minimum of 5 minutes to remove detergent residuals. a) Ensure that running water passes through threads, crevices, latching side clips, handles, sliding surfaces and engraved surfaces of all the devices.
<b>9</b>	Dry the rail, rail extension, and rail tray using lint-free cloths and/or filtered pressurized air.
<b>10</b>	Inspect the rail, rail extension, and rail tray for visible contamination. Repeat steps 3 through 9 if contamination is observed.
<b>11</b>	If the devices are to be disinfected next, disregard this step and proceed to Section 7.0. After cleaning, place the rail and rail extension in their corresponding brackets within the rail tray (Figure 2). Wrap the rail tray (with the rail and rail extension inside) in two layers of standard sterilization wrap until sterilized.

### 6.3 Instructions for Automated Cleaning with either Enzymatic Detergent or Alkaline Detergent

<b>Step</b>	<b>Procedure</b>
<b>1</b>	Ensure to remove all accessories from the rail and rail extension.
<b>2</b>	Disassemble the rail and rail extension. a) If the rail and rail extension are contained within the rail tray, remove the lid from the rail tray then remove the rail and rail extension. b) If the rail and rail extension are assembled together, pull up on the rail extension to remove from the rail (Figure 4). For the rail, ensure that the knob is completely unthreaded from the proximal clamp assembly and the threads are exposed (Figure 3).
<b>3</b>	Rinse the rail, rail extension, and rail tray under running cold utility (potable) water for 1 minute to remove gross soil. If soil is visible, use a soft-bristled brush to scrub visible soil. While rinsing the rail: a) Focus on threads, crevices, sliding surfaces, and engraved surfaces. b) Ensure that the knob remains fully unthreaded from the proximal clamp assembly threads and internal knob threads are fully exposed (Figure 3). While rinsing the rail extension: a) Focus on crevices and engraved surfaces. While rinsing the rail tray: a) Focus on crevices, latching side clips, handles, brackets, and engraved surfaces.
<b>4</b>	Place the rail and rail extension in their corresponding brackets of the rail tray (Figure 2).
<b>5</b>	Transfer the devices in the rail tray onto an appropriate rack system (e.g., 3-Level manifold rack accessory) contained inside the washer-disinfector for processing. a) Ensure the lid of the rail tray is not on the tray. Place the lid in an open spot in the washer-disinfector and at an angle to facilitate drainage.

Step	Procedure																																																								
6	<p>Select the automated wash cycle parameters as listed below for EITHER enzymatic detergent or alkaline detergent, NOT BOTH. If using enzymatic detergent, use the table of parameters for enzymatic detergent; if using alkaline detergent, use the table of parameters for alkaline detergent. Once the rail tray (containing the rail and rail extension) has been cleaned with EITHER enzymatic detergent or alkaline detergent, proceed to the next step.</p> <p>The following table is the list of parameters for automated cleaning with <b>enzymatic detergent</b>.</p> <table border="1"> <thead> <tr> <th colspan="4">AUTOMATIC CLEANING PARAMETERS FOR ENZYMATIC DETERGENT</th> </tr> <tr> <th>Phase</th> <th>Time (Minutes)</th> <th>Temperature and Water Type</th> <th>Detergent Type and Concentration</th> </tr> </thead> <tbody> <tr> <td>Pre-wash</td> <td>2:00</td> <td>Cold utility (potable) water</td> <td>N/A</td> </tr> <tr> <td>Enzyme Wash</td> <td>4:00*</td> <td>Hot utility (potable) water</td> <td>Enzymatic detergent, such as Valsure® Enzymatic Cleaner or equivalent, prepared at the detergent manufacturer's recommendation</td> </tr> <tr> <td>Rinse</td> <td>2:00</td> <td>≥ 43 °C utility (potable) water (Set Point)</td> <td>N/A</td> </tr> <tr> <td>Thermal Disinfection (optional#)</td> <td>5:00</td> <td>90 °C critical water (Set Point)</td> <td>N/A</td> </tr> <tr> <td>Drying</td> <td>7:00</td> <td>45 °C (Set Point)</td> <td>N/A</td> </tr> </tbody> </table> <p>*Or per manufacturer's recommendation, whichever is greater #Consult your facility's guidelines to determine if disinfection is required</p> <p>The following table is the list of parameters for automated cleaning with <b>alkaline detergent</b>.</p> <table border="1"> <thead> <tr> <th colspan="4">AUTOMATIC CLEANING PARAMETERS FOR ALKALINE DETERGENT</th> </tr> <tr> <th>Phase</th> <th>Time (Minutes)</th> <th>Temperature and Water Type</th> <th>Detergent Type and Concentration</th> </tr> </thead> <tbody> <tr> <td>Pre-wash</td> <td>2:00</td> <td>Cold utility (potable) water</td> <td>N/A</td> </tr> <tr> <td>Alkaline Wash</td> <td>10:00*</td> <td>Hot utility (potable) water</td> <td>Alkaline detergent, such as Neodisher® Mediclean Forte or equivalent, prepared at the detergent manufacturer's recommendation</td> </tr> <tr> <td>Rinse</td> <td>2:00</td> <td>≥ 43 °C utility (potable) water (Set Point)</td> <td>N/A</td> </tr> <tr> <td>Thermal Disinfection (optional#)</td> <td>5:00</td> <td>90 °C critical water (Set Point)</td> <td>N/A</td> </tr> <tr> <td>Drying</td> <td>7:00</td> <td>45 °C (Set Point)</td> <td>N/A</td> </tr> </tbody> </table> <p>*Or per manufacturer's recommendation, whichever is greater #Consult your facility's guidelines to determine if disinfection is required</p>	AUTOMATIC CLEANING PARAMETERS FOR ENZYMATIC DETERGENT				Phase	Time (Minutes)	Temperature and Water Type	Detergent Type and Concentration	Pre-wash	2:00	Cold utility (potable) water	N/A	Enzyme Wash	4:00*	Hot utility (potable) water	Enzymatic detergent, such as Valsure® Enzymatic Cleaner or equivalent, prepared at the detergent manufacturer's recommendation	Rinse	2:00	≥ 43 °C utility (potable) water (Set Point)	N/A	Thermal Disinfection (optional#)	5:00	90 °C critical water (Set Point)	N/A	Drying	7:00	45 °C (Set Point)	N/A	AUTOMATIC CLEANING PARAMETERS FOR ALKALINE DETERGENT				Phase	Time (Minutes)	Temperature and Water Type	Detergent Type and Concentration	Pre-wash	2:00	Cold utility (potable) water	N/A	Alkaline Wash	10:00*	Hot utility (potable) water	Alkaline detergent, such as Neodisher® Mediclean Forte or equivalent, prepared at the detergent manufacturer's recommendation	Rinse	2:00	≥ 43 °C utility (potable) water (Set Point)	N/A	Thermal Disinfection (optional#)	5:00	90 °C critical water (Set Point)	N/A	Drying	7:00	45 °C (Set Point)	N/A
AUTOMATIC CLEANING PARAMETERS FOR ENZYMATIC DETERGENT																																																									
Phase	Time (Minutes)	Temperature and Water Type	Detergent Type and Concentration																																																						
Pre-wash	2:00	Cold utility (potable) water	N/A																																																						
Enzyme Wash	4:00*	Hot utility (potable) water	Enzymatic detergent, such as Valsure® Enzymatic Cleaner or equivalent, prepared at the detergent manufacturer's recommendation																																																						
Rinse	2:00	≥ 43 °C utility (potable) water (Set Point)	N/A																																																						
Thermal Disinfection (optional#)	5:00	90 °C critical water (Set Point)	N/A																																																						
Drying	7:00	45 °C (Set Point)	N/A																																																						
AUTOMATIC CLEANING PARAMETERS FOR ALKALINE DETERGENT																																																									
Phase	Time (Minutes)	Temperature and Water Type	Detergent Type and Concentration																																																						
Pre-wash	2:00	Cold utility (potable) water	N/A																																																						
Alkaline Wash	10:00*	Hot utility (potable) water	Alkaline detergent, such as Neodisher® Mediclean Forte or equivalent, prepared at the detergent manufacturer's recommendation																																																						
Rinse	2:00	≥ 43 °C utility (potable) water (Set Point)	N/A																																																						
Thermal Disinfection (optional#)	5:00	90 °C critical water (Set Point)	N/A																																																						
Drying	7:00	45 °C (Set Point)	N/A																																																						
7	Remove the rail tray (with the rail and rail extension within) and the rail tray lid from the washer-disinfector. When unloading, check moving parts and surrounding surfaces for complete removal of visible soil. If soil is visible, repeat the automated cleaning procedure (Section 6.3) or perform a manual cleaning procedure (Section 6.1 or Section 6.2). If no soil is visible, place and latch the rail tray's lid on the rail tray and proceed to the next step.																																																								
8	After cleaning (and disinfecting), the rail tray (with the rail and rail extension inside) should be wrapped in two layers of standard sterilization wrap until sterilized.																																																								

### 7.0 Disinfection Instructions (Optional\*)

Step	Procedure								
1	While ensuring the knob threads of the rail are fully exposed (Figure 3) place the rail and rail extension in their corresponding brackets of the rail tray (Figure 2).								
2	Transfer the devices in the rail tray onto an appropriate rack system contained inside the washer-disinfector for processing. a) Ensure the lid of the rail tray is not on the tray. Place the lid in an open spot in the washer-disinfector and at an angle to facilitate drainage.								
3	Set the motor speed to HIGH, if applicable, and select the following recommended parameters:								
	<table border="1"> <thead> <tr> <th>Phase</th> <th>Time (Minutes)</th> <th>Temperature and Water Type</th> <th>Detergent Type and Concentration</th> </tr> </thead> <tbody> <tr> <td>Thermal Rinse</td> <td>5:00</td> <td>90 °C critical water (Set Point)</td> <td>N/A</td> </tr> </tbody> </table>	Phase	Time (Minutes)	Temperature and Water Type	Detergent Type and Concentration	Thermal Rinse	5:00	90 °C critical water (Set Point)	N/A
Phase	Time (Minutes)	Temperature and Water Type	Detergent Type and Concentration						
Thermal Rinse	5:00	90 °C critical water (Set Point)	N/A						

Step	Procedure
4	Remove the rail tray (with the rail and rail extension inside) and the rail tray lid from the washer-disinfector. Place the rail tray's lid on the rail tray and latch.
5	After disinfecting, the rail tray (with the rail and rail extension inside) should be wrapped in two layers of standard polypropylene sterilization wrap until sterilized.

\*Automated thermal disinfection is recommended per facility guidelines as an added safety factor for personnel

## 8.0 Sterilization Instructions

The rail set (rail, rail extension, rail tray) must be sterilized prior to each use in accordance with appropriate hospital procedures and the following recommended instructions and parameters.

**Note: When sterilizing multiple devices in one cycle, ensure devices or trays are not stacked and sterilizer max load is not exceeded.**

Step	Procedure																							
1	For the rail, ensure that the knob is completely unthreaded from the proximal clamp assembly (Figure 3).																							
2	Ensure the rail and rail extension are secured in their corresponding brackets within the rail tray and the tray lid is on and latched (Figure 2).																							
3	Sterilize using steam sterilization. Devices may be steam sterilized via pre-vacuum or gravity methods:																							
	<table border="1"> <thead> <tr> <th>Sterilizer Type</th> <th>Configuration</th> <th>Minimum Temperature</th> <th>Minimum Exposure Time</th> <th>Minimum Dry Time</th> <th>Minimum Cool-Down Time</th> </tr> </thead> <tbody> <tr> <td rowspan="2">Gravity Displacement</td> <td>Double Wrapped<sup>1</sup></td> <td>121 °C</td> <td>30 minutes</td> <td>50 minutes</td> <td>30 minutes<sup>2</sup></td> </tr> <tr> <td>Double Wrapped<sup>1</sup></td> <td>132 °C</td> <td>4 minutes</td> <td>30 minutes</td> <td>N/A</td> </tr> <tr> <td rowspan="2">Pre-vacuum</td> <td>Double Wrapped<sup>1</sup></td> <td>134 °C</td> <td>3 minutes</td> <td>30 minutes</td> <td>60 minutes<sup>2</sup></td> </tr> </tbody> </table>	Sterilizer Type	Configuration	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time	Minimum Cool-Down Time	Gravity Displacement	Double Wrapped <sup>1</sup>	121 °C	30 minutes	50 minutes	30 minutes <sup>2</sup>	Double Wrapped <sup>1</sup>	132 °C	4 minutes	30 minutes	N/A	Pre-vacuum	Double Wrapped <sup>1</sup>	134 °C	3 minutes	30 minutes	60 minutes <sup>2</sup>
	Sterilizer Type	Configuration	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time	Minimum Cool-Down Time																		
	Gravity Displacement	Double Wrapped <sup>1</sup>	121 °C	30 minutes	50 minutes	30 minutes <sup>2</sup>																		
Double Wrapped <sup>1</sup>		132 °C	4 minutes	30 minutes	N/A																			
Pre-vacuum	Double Wrapped <sup>1</sup>	134 °C	3 minutes	30 minutes	60 minutes <sup>2</sup>																			
	<sup>1</sup> Halyard Health H300 or equivalent loaded in configuration depicted in Figure 2																							
<sup>2</sup> Post steam sterilization handling may be necessary. Place outside of chamber on wire rack to cool.																								
<b>WARNING: Do not remove the rail tray from sterile packaging after steam sterilization. This may compromise sterility.</b>																								

## 9.0 Storage

Store the steam sterilized Reusable Rail Set (i.e., the rail tray with the rail and rail extension within) in double wrapped configuration in a cool, dry place.

## 10.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed.

The following table contains the Basic UDI-DIs for the Edwards Reusable Rail Set:

Product	Model	Basic UDI-DI
Edwards reusable rail	10500RL	0690103D004REU000YA
Edwards reusable rail extension	10500EX	0690103D004REU000YA
Edwards rail tray	TRAY10500RL	0690103D004REU000YA

## 11.0 Reporting

Users should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user is established.

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TMT79



Figure 2: Top-down view of rail and rail extension in rail tray with lid removed

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



**Figure 3: Knob completely unthreaded from proximal clamp assembly**

- **Figure 3 : Molette complètement déenfilée de l'ensemble de clamp proximal**
- **Abbildung 3: Knopf vollständig von der proximalen Klemmen-Einheit abgeschraubt**
- **Figura 3: Pomo completamente desenroscado del conjunto de la abrazadera proximal**
  - **Figura 3: manopola completamente libera dal gruppo clamp prossimale**
  - **Afbeelding 3: Knop volledig uit het proximale klemsamenstel gedraaid**
- **Figur 3: Knap skruet helt ud af den proksimale klemmeenhed** ■ **Figur 3: Ratt helt urskruvad från proximal klämma**
  - **Εικόνα 3: Ο περιστροφικός διακόπτης εντελώς ξεβιδωμένος από την εγγύς διάταξη σφιγκτήρα**
  - **Figura 3: Manípulo totalmente desenroscado do conjunto do grampo proximal**
- **Obrázek 3: Knoflík zcela vyšroubovaný ze sestavy proximální svorky** ■ **3. ábra: A proximális szorítóegységéről teljesen letekert gomb**
  - **Rysunek 3: Pokrętło całkowicie wykręcone z proksymalnego zespołu zacisku**
  - **Obrázok 3: Otočný ovládač je úplne odpojený od zostavy proximálnej svorky**
- **Figur 3: Knotten er helt ute av den proksimale klemmen** ■ **Kuva 3: Nuppi kokonaan kierrettynä irti proksimaalisesta puristinkokoonpanosta**
  - **Фигура 3: Копче, напълно развито от комплект на проксимална клампа**
  - **Figura 3: Buton complet desfăcut din ansamblul pensei clamp proximale**
  - **Joonis 3. Nupp on proksimaalse klambri komplektilt täielikult lahti keeratud**
- **3 pav. Rankenėlė visiškai atsukta nuo proksimalinio spaustuko mazgo** ■ **3. attēls. Poga pilnībā izskrūvēta no proksimālās spaiļes bloka**
  - **Şekil 3: Topuz yivleri proksimal klemp tertibatından tamamen çıkarılmıştır**
  - **Рис. 3. Ручка, полностью выкрученная из проксимального зажима**
- **Slika 3: Obrtno dugme potpuno odvijeno sa sklopa proksimalne kleme** ■ **Slika 3: Kotačić potpuno odvijen iz sklopa proksimalne hvataljke**



- Figure 4: Rail extension removed from the rail ■ Figure 4 : Rallonge de rail retirée du rail
- Abbildung 4: Schienenverlängerung von der Schiene entfernt ■ Figura 4: Extensi3n de riel desmontada del riel
  - Figura 4: prolunga per guida rimossa dalla guida ■ Afbeelding 4: railverlenging verwijderd van de rail
  - Figur 4: Skinneforlængerer fjernet fra skinnen ■ Figur 4: Skenf3rlængning borttagen från skenan
  - Εικόνα 4: Αφαίρεση της προέκτασης ράγας από τη ράγα ■ Figura 4: Extens3o para calha removida da calha
  - Obrázek 4: Prodloužení kolejnice odstraněné z kolejnice ■ 4. ábra: A sínr3l eltávolított sínhosszabbító
  - Rysunek 4: Przedłużenie prowadnicy wyjęte z prowadnicy ■ Obrázok 4: Nadstavec na kolajnicu odstránený z kolajnice
  - Figur 4: Skinneforlengelse fjernet fra skinnen ■ Kuva 4: Kiskon jatke irrotettuna kiskosta
  - Фигура 4: Удължител за релса, отстранен от релсата ■ Figura 4: Extensie pentru şină scoasă de pe şină
  - Joonis 4. R3öra pikendus on r3öpalt eemaldatud ■ 4 pav. Nuo b3gelio nuimtas b3gelio ilgintuvas
  - 4. att3ls. Sliedes pagarinājums noņemts no sliedes ■ Şekil 4: Ray uzatması raydan çıkarılmıştır
  - Рис. 4. Удлинитель направляющей, снятый с направляющей ■ Slika 4: Produžetak šine uklonjen sa šine
  - Slika 4: Nastavak za vodilicu uklonjen s vodilice

**Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola**

	English	Français	Deutsch	Español	Italiano
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Catalogue Number	Référence catalogue	Artikelnummer	Número de catálogo	Numero di catalogo
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Quantity	Quantité	Menge	Cantidad	Quantità
	Contents	Contenu	Inhalt	Contenido	Contenuti
	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche
	Keep away from sunlight	Tenir à l'abri de la lumière du soleil	Vor Sonneneinstrahlung schützen	Mantener alejado de la luz del sol	Tenere lontano dalla luce solare
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Consult instructions for use on the website	Consulter le mode d'emploi sur le site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile
	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union	Representante autorizado en la Comunidad Europea/Unión Europea	Rappresentante autorizzato nella Comunità Europea/Unione Europea
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione

**Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola**

	English	Français	Deutsch	Español	Italiano
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	Importer	Importateur	Importeur	Importador	Importatore

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Catalogusnummer	Katalognummer	Katalognummer	Αριθμός καταλόγου	Número de catálogo
	Lotnummer	Partinummer	Lotnummer	Αριθμός παρτίδας	Número de lote
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Inhoud	Indhold	Innehåll	Περιεχόμενο	Conteúdo
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas
	Uit het zonlicht houden	Beskyttes mod sollys	Skyddas mot solljus	Φυλάσσετε μακριά από την άμεση ηλιακή ακτινοβολία	Manter afastado da luz solar
	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό	Manter seco
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης στον ιστότοπο	Consultar as instruções de utilização no site
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Niet-steriel	Ikke-steril	Icke-sterilt	Μη αποστειρωμένο	Não esterilizado


Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/União Europeia
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Importeur	Importør	Importör	Εισαγωγέας	Importador

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
	Suomi	Български	Română	Eesti	Lietuvių
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vésioje, sausoje vietose
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба в уебсайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege kasutusjuhiseid veebisaidilt	Žr. naudojimo instrukcijas interneto svetainėje
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Epästeriili	Нестерилно	Nesteril	Mittesteriilne	Nesterilus
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Maahantuoja	Вносител	Importator	Importija	Importuotojas

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Kataloga numurs	Katalog Numarası	Номер по каталогу	Kataloški broj	Kataloški broj
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Daudzums	Miktar	Количество	Količina	Količina
	Saturs	İçindekiler	Содержимое	Sadržaj	Sadržaj
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu

**Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümboleite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola**

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Sargāt no saules gaismas	Güneş ışığından uzak tutun	Беречь от солнечных лучей	Držite podalje od sunčeve svetlosti	Držite dalje od sunčeve svjetlosti
	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Medicinska ierice	Tibbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na internet stranici	Pogledajte upute za upotrebu na internetskoj stranici
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Nesterils	Steril Değildir	Нестерильно	Nesterilno	Nesterilno
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Evropskoj zajednici / Evropskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik

**Note:** The labeling of this product may not contain every symbol depicted in this legend. ■ **Remarque :** il est possible que certains symboles de cette légende n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Die Kennzeichnung dieses Produkts enthält möglicherweise nicht jedes Symbol, das in dieser Zeichenerklärung abgebildet ist. ■ **Nota:** Es posible que no todos los símbolos mostrados en esta leyenda aparezcan en el etiquetado de este producto. ■ **Nota:** l'etichetta di questo prodotto potrebbe non contenere tutti i simboli riportati in questa legenda. ■ **Opmerking:** De etikettering van dit product bevat mogelijk niet elk symbool dat in deze legenda wordt weergegeven. ■ **Bemærk:** Produktets mærkning omfatter muligvis ikke alle symboler gengivet i denne forklaring. ■ **Obs!** Märkningen av produkten innehåller eventuellt inte alla symboler som avbildas i denna symbolförklaring. ■ **Σημείωση:** Η σήμανση του προϊόντος ενδέχεται να μην περιέχει κάθε σύμβολο που αναφέρεται στο παρόν υπόμνημα. ■ **Nota:** a rotulagem deste produto pode não conter todos os símbolos ilustrados nesta legenda. ■ **Poznámka:** Označení tohoto výrobku nemusí obsahovat všechny symboly uvedené v této legendě. ■ **Megjegyzés:** A termék címkéje nem feltétlenül tartalmazza az ezen jelmagyarázattal szereplő összes szimbólumot. ■ **Uwaga:** Etykieta tego wyrobu może nie zawierać wszystkich symboli przedstawionych w niniejszej legendzie. ■ **Poznámka:** Na štítku tohto produktu sa nemusia nachádzať všetky symboly uvedené v týchto vysvetlivkách. ■ **Merk:** Merkingen av dette produktet inneholder kanskje ikke alle symbolene som er gjengitt i denne symbolforklaringen. ■ **Huomautus:** Tämän tuotteen merkinnät eivät välttämättä sisällä kaikkia tässä selityksessä olevia symboleita. ■ **Забележка:** Етикетите по този продукт може да не съдържат всеки символ, представен в тази легенда. ■ **Notă:** este posibil ca eticheta acestui produs să nu conțină fiecare simbol descris în această legendă. ■ **Märkus.** Tootel esinev märgistus ei pruugi sisaldada kõiki tingmärkide seletuses välja toodud sümboleid. ■ **Pastaba.** Šio gaminio etiketėje gali nebūti visų šiame paaiškinime pavaizduotų simbolių. ■ **Piezīme.** Ši produkta marķējuma var nebūt visu simbolu, kuri iekļauti šajā simbolu skaidrojumā. ■ **Not:** Bu ürünün etiketi bu açıklama verilen her sembolü içermeyebilir. ■ **Примечание.** Маркировка изделия может содержать не все условные обозначения. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** oznaka ovoga proizvoda možda ne sadrži svaki simbol prikazan u ovoj legendi simbola.

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2026-03  
10060806002 A  
© Copyright 2026, Edwards Lifesciences LLC  
All rights reserved.

**Edwards Lifesciences LLC**  
One Edwards Way  
Irvine, CA 92614 USA

**Edwards Lifesciences B.V.**  
Verlengde Poolseweg 16  
4818 CL Breda, Netherlands

Telephone+1.949.250.2500  
+1.800.424.3278  
FAX +1.949.250.2525

Web IFU

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**Edwards Lifesciences GmbH**