Edwards

**Edwards EVOQUE Tricuspid Valve Replacement System**  
**Systeme de remplacement de la valve tricuspide Edwards EVOQUE**

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**Instructions for Use**

**Rx only**

**Caution: Implantation of the transcatheter tricuspid valve should only be performed by physicians who have received training for the Edwards EVOQUE tricuspid valve replacement system.**

- Do not attempt to use the Edwards EVOQUE valve (herein referred to as the EVOQUE valve), delivery system, or associated accessories before completely reading and understanding the information contained in this booklet.
- Failure to follow these instructions, warnings, and precautions may lead to device damage or patient injury. Use of the EVOQUE system should be restricted to those physicians trained to perform invasive endovascular procedures and to those physicians trained in the proper use of the system.
- Consult with authorized Edwards personnel for appropriate selection of EVOQUE valve size.
- The EVOQUE valve, delivery system, dilator kit, loading system, and stabilizer are supplied STERILE. The single-use accessories (base and plate) and reusable accessories (rail, platform, and plate) are supplied NON-STERILE.
- For the Edwards reusable platform and plate, refer to the reprocessing instructions in the Edwards Reusable Accessories IFU for wipe down steps. For the Edwards reusable rail assembly, refer to the reprocessing instructions in the Edwards Reusable Rail Set IFU. All other devices are supplied for single-use only. After use, dispose in accordance with hospital administrative and/or government policy.

**1.0 Introduction**

The Edwards EVOQUE tricuspid valve replacement system (herein referred to as the EVOQUE system) is designed to replace the native tricuspid valve in patients with tricuspid valve regurgitation without the need for conventional open-heart surgery. The EVOQUE system consists of four (4) elements which are used with either of the two (2) accessories sets, as outlined below:

Product Name	44 mm	48 mm	52 mm	56 mm
Edwards EVOQUE Valve (EVOQUE Valve)	9850EV44	9850EV48	9850EV52	9850EV56
Edwards EVOQUE Loading System	9850LS			-
	9850LSA			
Edwards EVOQUE Loading System and Trimmer	9850LSB			

Edwards, Edwards Lifesciences, the stylized E logo, Edwards EVOQUE, and EVOQUE are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Product Name	44 mm	48 mm	52 mm	56 mm
Edwards EVOQUE Tricuspid Delivery System	9850TDS			
Edwards EVOQUE Dilator Kit	9850DK			

Accessories Set 1	
Product Name	Model Number
Edwards EVOQUE Stabilizer	9850SB
Edwards EVOQUE Stabilizer Base	9850BA
Edwards EVOQUE Stabilizer Plate	9850PT

Accessories Set 2	
Product Name	Model Number
Edwards EVOQUE Stabilizer	9850SZ
Edwards Reusable Rail	10500RL
Edwards Reusable Platform	10000UP
Edwards Reusable Plate	10000PT

**Note: The EVOQUE system is designed for use with either Accessories set 1 or Accessories set 2. Accessories set 1 is considered an optional and single-use set for use with the EVOQUE system. Components from set 1 and set 2 are not interchangeable.**

**Note: The devices that comprise the EVOQUE system are intended to be used together and are not intended to be used as standalone devices. The information in this IFU is applicable to all the EVOQUE system devices.**

## 2.0 Device Description

### • Edwards EVOQUE Valve (Figure 1)

The EVOQUE valve consists of a trileaflet bovine pericardial tissue valve, nitinol frame, and fabric skirt, and is packaged and terminally sterilized in glutaraldehyde.

Valve size recommendations are based on native valve annulus size, as measured by computed tomography (CT). Patient anatomical factors and imaging modalities should be considered during valve size selection.

Device Diameter (Recommended Valve Size)	Systole		Diastole	
	Recommended Treatable Perimeter - Derived Diameter Range (mm)	Maximum Treatable Annulus Length (mm)	Recommended Treatable Perimeter - Derived Diameter Range (mm)	Maximum Treatable Annulus Length (mm)
44	36.5 - 43	45.5	39.6 - 45.5	50
48	40 - 47	49.5	43.2 - 49.5	54
52	45 - 51	53.5	46.8 - 53.5	58
56	48 - 55	57.5	50.4 - 57.5	62

### • Edwards EVOQUE Tricuspid Delivery System (Figure 2, Figure 3, and Figure 4)

The delivery system has an outer diameter of 28 F and is intended to deliver the EVOQUE valve in the crimped position via the transfemoral venous approach. The delivery system handle contains a primary flex knob, secondary flex knob, and depth knob to facilitate EVOQUE valve alignment and positioning in the native valve, and a capsule knob and release knob to control the expansion and release of the EVOQUE valve.

### • Edwards EVOQUE Dilator Kit (Figure 5)

The 24 F, 28 F, and 33 F diameter hydrophilic-coated dilators are intended to dilate the access site, facilitating delivery system insertion. All dilators accommodate a 0.035 inch (0.89 mm) guidewire and are tapered to minimize access site trauma.

### • Edwards EVOQUE Loading System (Figure 6) / Loading System and Trimmer (Figure 7)

The loading system, which consists of multiple components, is intended to facilitate loading and attachment of the EVOQUE valve onto the delivery system. The loading system assists in crimping the EVOQUE valve to the appropriate diameter, which allows the outer capsule to advance over the EVOQUE valve.

### • Edwards EVOQUE System Accessories Set 1: Stabilizer (Model 9850SB), Base, and Plate (Figure 8)

---

The single-use stabilizer with single-use base and plate are intended to secure the delivery system at an angle appropriate for the transfemoral venous approach and to enable fine adjustments of the position of the delivery system during the implantation procedure. The base is height-adjustable to accommodate patient lower extremities and is intended to provide a stable base for the stabilizer. The plate is intended to provide a stable, flat surface for the base on the operating table.

• **Edwards EVOQUE System Accessories Set 2: Stabilizer (Model 9850SZ) and Reusable Rail, Platform, and Plate (Figure 9)**

The reusable accessories compatible with the EVOQUE system include the reusable rail, platform, and plate. The reusable rail is a stable, sterilizable rail that provides a surface to secure compatible stabilizer to the reusable platform. The reusable platform allows for positioning and stabilization of Edwards transcatheter replacement delivery systems through adjustable leg height and a stable platform. The reusable plate is placed below the reusable platform, providing a flat and stable surface during use. The single-use stabilizer is intended to aid with positioning and stabilization of the EVOQUE system during implantation procedures. For the Accessories set utilizing reusable accessories (Figure 9), refer to the Edwards Reusable Accessories IFU for set up of the platform and plate. If reusable rail is non-sterile or if unpackaging for the first time, refer to the Edwards Reusable Rail Set IFU.

### **3.0 Facilities and Training Requirement**

Facilities that intend to perform an implantation procedure utilizing the EVOQUE system must have access to cine fluoroscopy and transesophageal echocardiography (TEE) throughout the procedure. In addition, the implanting physicians must have prompt access to facilities with the necessary equipment, instruments, supplies, and personnel to perform emergency tricuspid valve surgery, if required.

A comprehensive training program is provided by Edwards Lifesciences and must be completed by the implanting physicians before use of the EVOQUE system. The implanting physicians should have advanced technical knowledge and experience in related catheter-based procedures.

### **4.0 Indications for Use**

#### **4.1 Intended Use**

The EVOQUE tricuspid valve replacement system is indicated for the improvement of health status in patients with symptomatic severe tricuspid regurgitation despite being treated optimally with medical therapy for whom tricuspid valve replacement is deemed appropriate by a Heart Team.

### **5.0 Contraindications**

The EVOQUE valve is contraindicated in patients who cannot tolerate an anticoagulation/antiplatelet regimen, who have active bacterial endocarditis or other active infections, or who have untreatable hypersensitivity to nitinol alloys (nickel and titanium).

### **6.0 Warnings**

- The EVOQUE valve, delivery system, loading system, dilator kit, and stabilizer are designed, intended, and distributed for STERILE single-use only. The EVOQUE stabilizer base and EVOQUE stabilizer plate are provided nonsterile for single-use only. Do not resterilize or reuse any of the single-use devices. There are no data to support the sterility, nonpyrogenicity, or functionality of the single-use devices after reprocessing.
- Ensure proper sterile techniques are utilized during the preparation, transfer, and use of the devices.
- Do not use the valve if the tamper evident seal is broken, the storage solution does not completely cover the valve, the temperature indicator has been activated, the valve is damaged, or the expiration date has elapsed. The EVOQUE valve must remain hydrated at all times. The valve cannot be exposed to solutions, antibiotics, or chemicals other than its shipping storage solution and sterile physiologic saline solution. This will help prevent leaflet damage that may impact valve functionality. Keep the EVOQUE valve hydrated with normal saline until ready for implantation.
- Ensure the correct valve size is selected. Implantation of the improper size (i.e., undersizing or oversizing) may lead to paravalvular leak (PVL), migration, embolization, and/or annular damage.
- Patients with previously implanted devices (e.g., inferior vena cava (IVC) filter) should be carefully assessed prior to insertion of the delivery system to avoid potential damage to vasculature or a previously implanted device.
- Patients with pre-existing cardiac leads should be carefully assessed prior to implantation to avoid potential adverse interaction between devices.
- Care should be taken when implanting cardiac leads after EVOQUE valve implantation to avoid potential adverse interaction between the devices.
- Patients implanted with the EVOQUE valve should be maintained on anticoagulant/antiplatelet therapy as determined by their physicians in accordance with current guidelines, to minimize the risk of valve thrombosis or thromboembolic events.
- There are no data to support device safety or effectiveness if the patient has:
  - Echocardiographic evidence of severe right ventricular dysfunction;

- 
- Pulmonary arterial systolic pressure (PASP) > 70 mmHg by echo Doppler;
  - A trans-tricuspid pacemaker or defibrillator lead that has been implanted in the RV within the last 3 months;
  - Dependency on a trans-tricuspid pacemaker without alternative pacing options.

## 7.0 Precautions

### 7.1 Precautions Prior to Use

- The patient's eligibility depends on the anatomic conditions based on CT scan.
- It is advised that a multi-disciplinary heart team be of the opinion that EVOQUE valve implantation is preferable to alternative percutaneous device solutions, including minimally-invasive open heart surgery.
- It is advised that a multi-disciplinary heart team takes into consideration the severity of disease and the chances of reversibility of right heart failure based on a complete hemodynamic assessment.

### 7.2 Precautions

- The EVOQUE valve is to be used only with the EVOQUE delivery system and EVOQUE loading system.
- The procedure should be conducted under appropriate imaging modalities, such as transesophageal echocardiography (TEE), fluoroscopy, and/or intracardiac echocardiography (ICE).
- Glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Safety Data Sheet available from Edwards Lifesciences.
- Conduction disturbances may occur before, during, or following implantation of the EVOQUE valve, which may require continuous ECG monitoring before hospital discharge. If a patient has confirmed or suspected conduction disturbances, consider patient monitoring and/or electrophysiology evaluation. The risk of conduction disturbances may increase with the 56 mm valve size.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Long-term durability has not been established for the EVOQUE valve. Regular medical follow-up is advised to evaluate EVOQUE valve performance.
- Implantation of the EVOQUE valve should be postponed in patients with (1) a history of myocardial infarction within one month (30 days) of planned intervention, (2) pulmonary emboli within 3 months (90 days) of planned intervention, (3) cerebrovascular accident (stroke or TIA) within 3 months (90 days) of planned intervention, (4) active upper GI bleeding within 3 months (90 days) prior to procedure requiring transfusion.

## 8.0 Potential Adverse Events

Potential adverse events related to standard cardiac catheterization, use of anesthesia, the EVOQUE valve, and the implantation procedure include:

- Abnormal lab values
- Allergic reaction to anesthesia, contrast media, anti-coagulation medication, or device materials
- Anaphylactic shock
- Anemia or decreased hemoglobin (Hgb), may require transfusion
- Aneurysm or pseudoaneurysm
- Angina or chest pain
- Arrhythmia – atrial (i.e., atrial fibrillation, supraventricular tachycardia)
- Arrhythmias – ventricular (i.e., ventricular tachycardia, ventricular fibrillation)
- Arterio-venous fistula
- Bleeding
- Cardiac arrest
- Cardiac (heart) failure
- Cardiac injury, including perforation
- Cardiac tamponade / pericardial effusion
- Cardiogenic shock
- Chordal entanglement or rupture that may require intervention
- Coagulopathy, coagulation disorder, bleeding diathesis
- Conduction system injury, which may require implantation of a pacemaker (temporary or permanent)
- Conversion to open heart surgery
- Coronary artery occlusion
- Damage to or interference with function of pacemaker or implantable cardioverter defibrillator (ICD)
- Death

- 
- Edema
  - Electrolyte imbalance
  - Embolization including air, particulate, calcific material, or thrombus
  - Emergent cardiac surgery
  - Endocarditis
  - Esophageal irritation
  - Esophageal perforation or stricture
  - EVOQUE system component(s) embolization
  - Failure to retrieve any EVOQUE system components
  - Fever
  - Gastrointestinal bleeding
  - Hematoma
  - Hemodynamic compromise
  - Hemolysis / hemolytic anemia
  - Hemorrhage requiring transfusion/surgery
  - Hypertension
  - Hypotension
  - Inflammation
  - Injury to the tricuspid apparatus including chordal damage, rupture, papillary muscle damage
  - Local and systemic infection
  - Mesenteric ischemia or bowel infarction
  - Multi-system organ failure
  - Myocardial infarction
  - Nausea and/or vomiting
  - Nerve injury
  - Neurological symptoms, including dyskinesia, without diagnosis of TIA or stroke
  - Non-emergent reoperation
  - Pain
  - Pannus formation
  - Paralysis
  - Percutaneous valve intervention
  - Peripheral ischemia
  - Permanent disability
  - Pleural effusion
  - Pneumonia
  - Pulmonary edema
  - Pulmonary embolism
  - Reaction to anti-platelet or anticoagulation agents
  - Rehospitalization
  - Renal failure
  - Respiratory failure, atelectasis – may require prolonged intubation
  - Retroperitoneal bleed
  - Right ventricular outflow tract (RVOT) obstruction
  - Septicemia, sepsis
  - Skin burn, injury, or tissue changes due to exposure to ionizing radiation
  - Stroke
  - Structural deterioration (wear, fracture, calcification, leaflet tear, leaflet thickening, stenosis of implanted device, or new leaflet motion disorder)
  - Thromboembolism
  - Transient ischemic attack (TIA)
  - Valve dislodgement/embolization
  - Valve endocarditis
  - Valve explant
  - Valve leaflet entrapment
  - Valve malposition
  - Valve migration
  - Valve paravalvular leak (PVL)

- 
- Valve regurgitation (new or worsening tricuspid, aortic, mitral, pulmonary)
  - Valve thrombosis
  - Vascular injury or trauma, including dissection or occlusion
  - Vessel spasm
  - Wound dehiscence, delayed or incomplete healing

## 9.0 Additional Equipment

The implantation procedure requires additional equipment that is not supplied with the EVOQUE system. The additional equipment is provided below.

### 9.1 Equipment for EVOQUE Valve Loading

**Note: Volumes reflect adequate amount for preparing 1 implant.**

- 3500 ml (minimum) sterile physiological saline solution at ambient temperature (~ 23 °C)
- 500 ml (minimum) heparinized saline solution (2 units/ml) at ambient temperature (~ 23 °C)
- 4 sterile bowls (≥ 500 ml, ≥ 7 cm depth, plastic)
- 1 large sterile bowl (≥ 2 L, ≥ 10 cm depth, plastic)
- 1 scalpel, no. 11 scalpel blade
- 1 luer locking syringe (≥ 20 cc)
- Gauze pads
- Blunt tip forceps
- Blunt tip scissors
- Sterile towels

### 9.2 Equipment for Access, Procedure, and Monitoring

- Standard cardiac catheterization lab equipment
- Femoral vessel introducer sheath
- Fluoroscopy (fixed, mobile, or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal echocardiography capabilities
- Steerable introducer sheath
- Exchange length 0.035 inch (0.89 mm) max guidewire
- Extra small curve 0.035 inch (0.89 mm) max guidewire
- Right coronary arterial catheter and guidewire
- Sterile table for EVOQUE valve and device preparation

### 9.3 Standby Equipment

- Arterial bypass cannula (~18 F)
- Cardiopulmonary bypass machine
- Compliant balloon (> 20 mm diameter, 9 cc contrast volume)
- Diluted radiopaque contrast medium (15:85 medium to saline dilution)
- High pressure contrast injector
- Intra-aortic balloon pump and appropriately sized balloon
- Pigtail angiographic catheter
- Venous return bypass cannula (~18 F)
- Transthoracic echocardiographic (TTE) equipment
- Vascular access lubricant
- Temporary pacing equipment

## 10.0 Directions for Use

### 10.1 Prior to Use Inspection

Before using the EVOQUE system, visually examine each element and accessory for evidence of gross damage (e.g., a cracked jar or lid, leakage, broken or missing seals) that may have compromised the packaging sterility (if applicable) or functionality of the components.

**WARNING: Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging and/or sterile barriers and any components have been opened or damaged, or the expiration date has elapsed, as sterility and/or function may be compromised.**

**WARNING: Do not mishandle the EVOQUE valve or use device/container if found to be damaged, leaking, or without adequate sterilant (not fully submerged in glutaraldehyde, or missing intact seals). The EVOQUE valve must not be used for implantation, as sterility may be compromised.**

**WARNING: Do not use the EVOQUE valve if the expiration date has elapsed, as either sterility or valve function may be compromised.**

**WARNING: Do not use the EVOQUE valve if the temperature indicator has been activated, as valve function may be compromised.**

**WARNING: Do not use the EVOQUE valve if the tamper evident seal is broken, as sterility may be compromised.**

## 10.2 Accessories Setup

### 10.2.1 Accessories Set 1 Setup

Step	Procedure
1	Place the plate on the operating table beneath the patient's leg to support the base during the procedure.
2	After the patient has been positioned on the operating table, place single-use base over the patient's leg, on top of the plate, at the desired distance from the mid-sternum.
3	Establish sterile barrier.
4	After sterile draping, place single-use base adapter in-line with the access site and attached to the front of the base using a clamp (Figure 8).

**Note: In the Accessories set 1 configuration, the base adapter and clamp are both supplied with the stabilizer.**

### 10.2.2 Accessories Set 2 Setup

For the Accessories set 2 configuration utilizing reusable accessories (Figure 9), refer to the Edwards Reusable Accessories IFU for set up of the platform and plate. If reusable rail is non-sterile or if unpackaging for the first time, refer to Edwards Reusable Rail Set IFU. Once reusable rail is sterile, refer to below instructions.

Step	Procedure
1	Remove Edwards reusable rail from sterilization wrap (blue wrap) and inspect for damage.
2	Prior to sterile draping the patient, assemble and position the reusable plate and reusable platform around the legs of patient, adjusting the height and angle of the platform as needed. Refer to the Edwards Reusable Accessories IFU.
3	Establish sterile barrier.
4	After sterile draping, assemble and attach Edwards reusable rail and stabilizer to the reusable platform.
5	Following the procedure, remove the reusable rail from the reusable platform. Perform standard wipe down of reusable rail before sending for reprocessing. Refer to the Edwards Reusable Rail Set IFU.

## 10.3 Device Handling and Preparation

Follow sterile technique during device preparation and implantation. Preparation to be performed by two operators, an implant operator (operating the distal section of the delivery system) and a handle operator (operating the proximal section of the delivery system). Edwards personnel will always serve as implant operator. Handle operator can be an Edwards personnel or hospital staff. The handle operator is to follow instructions below:

**WARNING: Do not mishandle the EVOQUE valve. EVOQUE valve leaflets mishandled or damaged during any part of the loading process will require replacement of the EVOQUE valve.**

**WARNING: The EVOQUE valve should not remain fully crimped for a period of time longer than 120 minutes, as it may impair valve functionality.**

**CAUTION: Do not place the jar or the pouch of the delivery system into the sterile field. The outside of the jar and the pouch are not sterile, and the contents within the jar and pouch should be handled using standard aseptic techniques to prevent contamination.**

**CAUTION: To reduce the risk of contamination, do not open the EVOQUE valve jar until implantation is certain.**

**CAUTION: Ensure that the entire suture is removed when removing the serial number tag from the EVOQUE valve, as it may lead to emboli.**

**CAUTION: Do not allow the EVOQUE valve to come in contact with any sharp instrument, as it may impair valve function.**

**WARNING: Adequate rinsing with normal saline must be performed before implantation to reduce glutaraldehyde concentration, as it may result in glutaraldehyde toxicity.**

**CAUTION: Avoid contact of the leaflet tissue or the rinse solution with towels, linens, or other sources of lint or particulate matter that could be transferred to the leaflet tissue, as it may lead to emboli.**

### 10.3.1 Delivery System Preparation

Follow sterile technique during device preparation and implantation. Edwards personnel will direct the handle operator to do the following:

Step	Procedure
1	Wipe only the handle of the delivery system using gauze soaked in saline via syringe.
2	Remove stylets from proximal and distal end of the delivery system, then flush the guidewire lumen (Figure 4(3)) with heparinized saline as the instructed amount by Edwards personnel.
3	Flush the handle flush port (Figure 3(4)) then sheath port (Figure 3(1)) with heparinized saline as the instructed amount by Edwards personnel.
4	Retract outer capsule: Rotate the white knob (Figure 3(2)) in the direction of the arrows until it comes to a hard stop to ensure outer capsule (Figure 4(8)) is fully retracted.
5	Retract inner capsule: Rotate the blue knob (Figure 3(3)) in the direction of the arrows until it comes to a hard stop to ensure inner capsule (Figure 4(6)) is fully retracted.
6	Ensure primary and secondary knobs (Figure 3(5 and 6)) are at 0 (zero). Ensure black depth knob (Figure 3(7)) is at ½. (if not, rotate to ½ or instructed by Edwards personnel).
7	Fully advance tapered tip: Press tapered tip locking button and advance the tapered tip slider (Figure 3(8)), releasing tapered tip locking button when complete.

### 10.3.2 Implant Preparation

Authorized Edwards personnel will perform the necessary steps to ensure the implant is ready for loading.

### 10.3.3 Implant Locking and Loading

Loading system components will be attached to the implant and delivery system by authorized Edwards personnel serving as the implant operator.

Edwards personnel will direct the handle operator to:

Step	Procedure
1	Slowly turn the white knob (Figure 3(2)) counterclockwise in the opposite direction of the arrows to advance the outer capsule (Figure 4(8)), stopping as indicated by Edwards personnel.
2	Rotate white knob (Figure 3(2)) in the direction of the arrows to retract outer capsule until the knob hard-stops (Figure 4(8)).
3	Rotate the blue knob (Figure 3(3)) in the opposite direction of the arrows to advance the inner capsule (Figure 4(6)), stopping when indicated by Edwards personnel. <b>Note: Edwards personnel will remove the lock driver and prepare the implant for loading then direct the handle operator to do the following:</b>
4	Rotate the white knob (Figure 3(2)) in the opposite direction of the arrows to continuously advance the outer capsule (Figure 4(8)), stopping at several stages as indicated by Edwards personnel until the guidewire shield and the implant anchors are at the appropriate position.
5	Fully advance inner capsule (Figure 4(6)) by rotating the blue knob (Figure 3(3)) in the opposite direction of the arrows until it hard-stops.
6	Fully advance outer capsule (Figure 4(8)) by rotating the white knob (Figure 3(2)) in the opposite direction of the arrows until it hard-stops. <b>Note: Edwards personnel will remove loader funnel and direct the handle operator to do the following:</b>
7	Relieve tension in the outer and inner capsule. <b>Note: Edwards personnel will pull tension on the support tube assembly to fully seat valve, while simultaneously directing the handle operator to do the following:</b>
8	Pull tension on delivery system handle while rotating white knob (Figure 3(2)) in the opposite direction of the arrows to relieve tension. <b>Note: Edwards personnel will perform final preparations steps to trim the capsule and direct the handle operator to do the following:</b>

Step	Procedure
9	Flush the handle flush port (Figure 3(4)) with the instructed amount of saline until the capsule is filled.
10	Slowly retract the tapered tip (Figure 4(1)) by pressing tapered tip button and retracting tapered tip slider (Figure 3(8)) until it contacts capsule.
11	Rotate the black depth knob (Figure 3(7)) counterclockwise in the opposite direction of the arrows until depth knob hard-stops.
12	Flush the handle flush port (Figure 3(4)) with the instructed amount of saline.
13	Flush the sheath port (Figure 3(1)) with the instructed amount to eliminate air and advance the sheath distally until mated flush to capsule.
14	Remove stylets from ends of the delivery system.
15	Flush the handle flush port (Figure 3(4)) and guidewire lumen (Figure 4(3)) with the instructed amount of heparinized saline.
16	Wipe the coated components (sheath, capsule, and tapered tip) using heparinized saline-soaked gauze.

## 10.4 EVOQUE Valve Implant

### 10.4.1 Guidewire Placement

Prepare femoral venous access using standard interventional techniques.

**WARNING: Do not use excessive force and/or manipulation during advancement and positioning of the guidewire, as it may result in perforation/dissection of arteries, veins, and/or other cardiac structures. This may also result in cardiac arrhythmias and conduction disturbances.**

Step	Procedure
1	Advance steerable sheath into the right atrium at the exit of the IVC.
2	Insert a guidewire through the steerable sheath.
3	Advance the guidewire across the tricuspid valve. <b>Note: Other interventional devices and techniques (e.g., guide catheters) may be used to assist the guidewire in crossing the tricuspid valve.</b>
4	Establish proper guidewire path and confirm no entanglement with cardiac structures.

### 10.4.2 EVOQUE Valve Delivery

**WARNING: Avoid excessive movement of the delivery system while executing the procedure to protect vasculature or cardiac structures. Avoid excessive rotation of the delivery system to maintain functionality of the delivery system.**

**WARNING: If functionality of the delivery system is compromised (e.g. significant increase in resistance during knob rotations, lack of responsiveness of knobs), do not proceed to deploy the valve, as it may result in inability to deploy the valve or valve malposition.**

**Note: Flush the delivery system with heparinized saline as needed throughout the procedure.**

Step	Procedure
1	Ensure the hydrophilic coating on dilators and delivery system is activated prior to use. Dilate access site. An Edwards EVOQUE dilator kit should be used as needed.
2	Insert the delivery system over the guidewire.
3	Advance the delivery system until the distal end of the tapered tip is positioned at the junction between IVC and right atrium.
4	Using fluoroscopy, ensure the delivery system is oriented in the proper orientation. <b>WARNING: The primary flex on the delivery system flexes in the direction of the flush ports; care should be taken to ensure that the delivery system is orientated correctly at this point.</b>
5	Retract the sheath.
6	Flex and orient the delivery system towards the tricuspid valve.
7	Advance the delivery system to cross the tricuspid valve. <b>Note: Delivery system flex, delivery system rotation, and guidewire position may be adjusted during valve crossing to optimize crossing position.</b>

Step	Procedure
8	Using echocardiography and fluoroscopic guidance, verify that the delivery system has crossed through the tricuspid valve into the right ventricle.
9	If using the Accessories set 1 configuration, dock the stabilizer onto the base adapter and secure it to the base. If using the Accessories set 2 configuration, dock the rail onto the platform and secure it. Next, dock the stabilizer onto the rail and secure it.
10	Attach the delivery system and sheath to the stabilizer.
11	Adjust the delivery system as needed to ensure hemodynamic stability.
12	Using pre-operative CT data (if available), position the C-arm at the optimal viewing projection.
13	Position the delivery system to be coaxial to the tricuspid annulus while minimizing contact with the native anatomy.
14	Using echocardiography and fluoroscopic guidance, confirm that the EVOQUE valve is positioned at the appropriate depth and coaxiality with relation to the native valve. <b>CAUTION: Once the capsule is retracted to expose the EVOQUE valve anchors, the valve is unable to be retrieved or recaptured into the delivery system.</b> <b>WARNING: Maintain central position of the delivery system within the native valve during deployment to ensure proper positioning of the valve.</b>
15	Retract the outer capsule until the anchors are exposed.
16	Adjust the position of the EVOQUE valve so that the anchors are positioned within native leaflets as dictated by patient anatomy.
17	Retract the outer and inner capsule until the desired EVOQUE valve diameter is achieved.
18	Engage the leaflets.
19	Confirm positioning of the EVOQUE valve by using echo imaging to assess engagement of the leaflets. Adjust the EVOQUE valve positioning as needed.
20	Using echo imaging, observe the motion of the native leaflets and adjust the position of the EVOQUE valve as required to fully engage native tricuspid valve leaflets.
21	Once full engagement has been confirmed, ensure that the EVOQUE valve is perpendicular to the tricuspid annular plane.
22	Retract the tapered tip until it is positioned within the EVOQUE valve.
23	Retract the inner capsule until the EVOQUE valve is released from the delivery system. <b>CAUTION: Care should be taken during final release of the EVOQUE valve using the release knob, as premature release may impact performance of the EVOQUE valve.</b>
24	Using echo and fluoroscopic imaging, assess final position and functionality of the EVOQUE valve.

#### 10.4.3 Delivery System Removal

**WARNING: Take care to maintain position of delivery system central within the EVOQUE valve during delivery system removal, as failure to do so may impact valve functionality or lead to dislodgement of the valve.**

**Note: Delivery system can be removed from the stabilizer at any point during removal, if applicable.**

Step	Procedure
1	Fully retract the tapered tip.
2	Unflex and retract delivery system as needed until the tapered tip is above the locking tabs of the EVOQUE valve. Adjust the guidewire as needed to maintain a central position relative to the EVOQUE valve. Ensure that the locking ring is free from the EVOQUE valve.
3	Rotate the release knob so that the inner capsule is in contact with the tapered tip.
4	Unflex and retract the delivery system as needed.
5	Rotate the capsule retraction knob until the outer capsule is in contact with the inner capsule.
6	Ensure the delivery system is fully unflexed and remove the delivery system from the access site.

Step	Procedure
	<b>Note: A sheath may be used to seal the femoral vein following system removal.</b>
7	Perform femoral closure as applicable using standard interventional techniques.
8	Perform a ventriculogram if needed to assess final position of the EVOQUE valve.

## 11.0 How Supplied

### 11.1 Sterilization and Packaging

The EVOQUE valve is provided sterile by terminal liquid sterilization and is non-pyrogenic. It is packaged and sterilized in a glutaraldehyde solution inside a jar to which a seal has been applied. The outer surface of the jar is not sterile and must not be placed in the sterile field. The EVOQUE valve is supplied with a temperature indicator and should not be used if the indicator has been activated.

The delivery system, dilator kit, and loading system are supplied sterilized by ethylene oxide and provided non-pyrogenic. The components are secured on a card and packaged in a pouch and shelf box.

The stabilizer is provided ethylene oxide sterilized. The components are secured on a card and packaged in a pouch and shelf box.

The single-use base and plate or the reusable rail, platform, and plate are provided non-sterile.

The components are packaged in individual shipper boxes.

### 11.2 Storage

The EVOQUE valve should be stored between 10 °C and 25 °C (50 °F and 77 °F). Stock inspection and rotation at regular intervals are recommended so that the EVOQUE valve with an earlier expiration date is used first.

The delivery system, dilator kit, loading system, single-use stabilizer/base/plate and reusable platform/plate should be stored in a cool, dry place that is contamination free. Refer to the Edwards Reusable Rail Set IFU for storage of the reusable rail.

## 12.0 Magnetic Resonance (MR) Safety Information

Non-clinical testing has demonstrated that the Edwards EVOQUE valve is MR Conditional. A patient with the valve can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 and 3 T only
- Maximum spatial gradient magnetic field of 3000 gauss/cm (30.0 T/m) or less
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg
- Normal operating mode of the MR system for both gradients and SAR

Under the scan conditions defined above, the EVOQUE valve is expected to produce a maximum temperature rise of 4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the EVOQUE valve extends approximately 0.8 cm from the device when imaged with a gradient echo or spin echo pulse sequence and a 3 T MRI system.

## 13.0 Patient Information

Patient education brochures are provided to each site and should be given to the patient to inform them of the risks and benefits of the procedure and alternatives in adequate time before the procedure to be read and discussed with their physician. A copy of this brochure may also be obtained from Edwards Lifesciences by calling 1.888.713.1564.

A patient implant card request form is provided with each transcatheter heart valve. After implantation, all requested information should be completed on this form. The serial number may be found on the package and on the identification tag attached to the transcatheter heart valve. The original form should be returned to the Edwards Lifesciences address indicated on the form and upon receipt, Edwards Lifesciences will provide an identification card to the patient.

## 14.0 Recovered Implant and Device Disposal

Edwards Lifesciences is interested in obtaining recovered clinical specimens of the EVOQUE valve for analysis. A written report summarizing our findings will be provided upon completion of our evaluation. Please contact Edwards for return of the recovered valve.

If you do decide to return any of the devices, please follow the following instructions:

- 
- **Unopened Package with Sterile Barrier Intact:** If the pouches have not been opened, return the device in its original packaging.
  - **Package Opened but Not Implanted:** If a pouch is opened, the device is no longer sterile. Please return the device in its original packaging.
  - **Explanted Implant:** The explanted implant should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit for return to Edwards.

#### **14.1 Disposal**

Use universal precautions for biohazards and sharps to avoid user injury. Used devices (includes all devices that come in contact with patients) should be handled and disposed of in accordance with institutional guidelines on biohazardous materials and hospital waste to avoid possible cross-contamination. For disposal of the reusable platform and plate accessories, refer to the Edwards Reusable Accessories IFU.

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© Copyright 2026, Edwards Lifesciences LLC
All rights reserved.
**Edwards Lifesciences LLC**